

## Dynamic Diagnostics Ltd Strep latex test

A rapid latex test for the confirmatory identification of *Streptococci* of Lancefield Groups A, B, C, D, F and G from culture plates.

Most strains of *Streptococci* which have been isolated from human infections possess group specific antigens. Identification of the organisms includes extraction and characterisation of these antigens from organism grown in culture.

The use of specific antisera combined with an enzyme extraction method provide a speedy and reliable typing system for *Streptococci*.



### Dynamic Diagnostics Ltd Strep Kit features

- Rapid grouping of Groups A, B, C, D, F, G
- Simple enzymic extraction
- Easy-to-read agglutination
- High specificity to group D

- Easy to use
- High sensitivity and specificity
- Results in one minute after extraction



# Dynamic Diagnostics<sup>Ltd</sup>

## Strep latex test

### Ordering Information:

Code	Description	Size	Storage
DD-47	Strep Latex Kit	50 tests	2 - 8°C

### Reagents/Kit:

Reagents	Product Code	Contents	Quantity/Size
Reagent A	DD-47a	Contains rabbit <i>Strep</i> Group A antibody sensitised latex particles	2.5ml
Reagent B	DD-47b	Contains rabbit <i>Strep</i> Group B antibody sensitised latex particles	2.5ml
Reagent C	DD-47c	Contains rabbit <i>Strep</i> Group C antibody sensitised latex particles	2.5ml
Reagent D	DD-47d	Contains rabbit <i>Strep</i> Group D antibody sensitised latex particles	2.5ml
Reagent F	DD-47f	Contains rabbit <i>Strep</i> Group F antibody sensitised latex particles	2.5ml
Reagent G	DD-47g	Contains rabbit <i>Strep</i> Group G antibody sensitised latex particles	2.5ml
Positive controls	DD-47pc	Positive control, contains inactivated polyvalent antigenic extracts to groups A, B, C, D, F and G	1.0ml
Extraction enzyme	DD-47ex	Lyophilised extraction enzyme	2x10ml

---