


FRI, OCT 11, 2024



CUSTOMIZEDBONE MODIFICATIONS DURING SURGERY

As stated in the product IFU:

"Any modification to the supplied device shall be the sole and exclusive responsibility of the surgeon. If during surgery a modification of the implant is deemed necessary by the surgeon, however, it should be made with extreme caution, only by using low speed diamond drills under water irrigation."

Especially in pediatrics, surgeons are more comfortable to adjust the implant vs adjust the bone edges.

TAKE HOME MESSAGES:



Surgeon Responsibility



Extreme Caution



Diamond drill under irrigation

