

# ELISABETH

## The Mitral Valve Apparatus Visualizer

Specifically designed for mitral subvalvular apparatus exposure during valve repair. Particularly useful and helpful in minimally invasive procedures when mitral valve exposure is difficult or inadequate.

- Designed in two configurations (Left and Right) in order to have better adaptability to the anatomy of the papillary muscles. In very large anuli, they can be used together to achieve better visualization and stability.
- Stability System design assures optimal positioning of the device in the natural anulus
- Easy to use

# ELISABETH

## The Mitral Valve Apparatus Visualizer

technical  
features


### THE MAGIC TRICK





HERE YOU DON'T SEE IT...



NOW YOU SEE IT...



...AND YOUR MAGIC TRICK IS SUCCESSFUL



A = 24 mm

B = 18 mm

C = 95 mm

CODE

DESCRIPTION

UNITS PER BOX

SMA-SYST-002-M

2 Elizabeth Mitral Valve Apparatus Visualizers  
R & L in 5 sterile pouches

10 units, 5 pouches



Manufactured by:

Distributed by:



MED EUROPE s.r.l.  
Via Gramsci, 9  
40018 San Pietro in Casale BOLOGNA - ITALY  
Tel. (0039) 051 810572 - 051 817902  
051 817850  
Fax (0039) 051 817845  
E-mail: info@medeurope.eu



Medical Instruments S.p.A.  
Via Salarolo, 2/B  
40068 San Lazzaro di Savena (BO)  
Tel. +39 051 6271111  
Fax +39 051 0195378  
E-mail: info@medicalinstruments.it