# Product Manual for disposable mask





Efficient protection



Light & breathable



Dust-proof and skin-friendly



3 Layer anti crease design




Prevent pathogens

Prevent virus and bacterial

Prevent droplets and dust

#### **Features:**

- High density barrier purification air intake
- Medical material, dustproof and breathable
- Skin-friendly materials
- Three layers filter, safe to use
- 5 Efficient Protection
- 6 Solid space design , closer protection



# **Packing**

#### **Individual Package**

Product Dimension 175mm \* 95mm

Material Outer layer / New anti-sticking non-woven fabric layer

Middle layer/ Spray fleece filter material Inner layer / Superfine soft fiber material



#### **Color Box**

Four colors Printing pcs/box

190\*98\*90mm

Dimension **Master** 

Carton

Boxes/ CTN Total 2,000pcs
Dimension 505\*390\*380mm





## **Product Instruction:**

| Product Name      | Pearlteck Comfortable Mask  |
|-------------------|---|
| Product Dimension | About 175mm * 95mm  |
| Specification     | 50pcs individual package  |
| Expiration Date   | Three years   |
| Product Manual    | This product is a disposable mask, not washable;  |
|                   | suitable for adults. Wear it for daily life, work and travel.   |
| Material          | Outer layer / New anti-sticking non-woven fabric layer  |
|                   | Middle layer/ Spray fleece filter material  |
|                   | Inner layer / Superfine soft fiber material   |
|                   |   |
| Note              | Use with caution for non-manufactured cloth allergies   |
| Note              | Use with caution for non-manufactured cloth allergies and those with abnormal cardiopulmonary function.   |
| Note              | ·   |
| Note              | and those with abnormal cardiopulmonary function.   |
| Note              | and those with abnormal cardiopulmonary function.  Store in a dry, ventilated, non-corrosive environment,   |
| Note              | and those with abnormal cardiopulmonary function.  Store in a dry, ventilated, non-corrosive environment,  away from fire and flammable materials   |
|                   | and those with abnormal cardiopulmonary function.  Store in a dry, ventilated, non-corrosive environment,  away from fire and flammable materials  This product is not a substitute for anti-virus masks and medical protective masks.                              |
| Manufacturer      | and those with abnormal cardiopulmonary function.  Store in a dry, ventilated, non-corrosive environment,  away from fire and flammable materials  This product is not a substitute for anti-virus masks and medical protective masks.  SHENZHEN PEARLTECK CO.,LTD. |

### **Dimension:**

