

Modular LC Columns



Our range of un-packed HPLC and UHPLC systems are made out of the highest UHPLC grade Stainless Steel with amazingly smooth internal finishes in the 5-8 μ m range. This has proven to give better peak results, maximize your investment and be highly efficient.

- With our UHPLC hardware, you can achieve high pressures(1500+ bar) and superior peak shape.
- Every modular style column we sell has the ID size laser etched for easy identification.
- We also offer a variety of customization, ranging from end fittings for your own special look, chemically resistant hard anodized colors for part identification or both.

| 2.1mm Modular Column Hardware | |
|-------------------------------|--|
| PART # | Description |
| 21020 | 2.1mm x 20mm Tube Only |
| 21020-020 | 2.1mm x 20mm x 2.0um |
| 21020-005 | 2.1mm x 20mm x 0.5um |
| 21020-002HP | High Pressure 2.1mm x 20mm x 0.2um |
| 21020-005HP | High Pressure 2.1mm x 20mm x 0.5um |
| 21020-020HP | High Pressure 2.1mm x 20mm x 2.0um |
| 21025 | 2.1mm x 25mm Tube Only |
| 21030-020 | 2.1mm x 30mm x 2.0um |
| 21030-005 | 2.1mm x 30mm x 0.5um |
| 21030-005HP | High Pressure 2.1mm x 30mm x 0.5um |
| 21030-002HP | High Pressure 2.1mm x 30mm x 0.2um |
| 21035-005HP | High Pressure 2.1mm x 35mm x 0.5um |
| 21050-020 | 2.1mm x 50mm x 2.0um |
| 21030 | 2.1mm x 30mm Tube Only |
| 21050 | 2.1mm X 50mm Tube Only |
| 21050-005 | 2.1mm x 50mm x 0.5um |
| 21050-005HP | High Pressure 2.1mm x 50mm x 0.5um |
| 21050-005HP-W | High Pressure 2.1mm x 50mm x 0.5um - Waters Seat |
| 21050-002HP | High Pressure 2.1mm x 50mm x 0.2um |
| 21075 | 2.1mm x 75mm Tube Only |
| 21075-020 | 2.1mm x 75mm x 2.0um |
| 21075-005 | 2.1mm x 75mm x 0.5um |
| 21075-005HP | High Pressure 2.1mm x 75mm x 0.5um |
| 21075-002HP | High Pressure 2.1mm x 75mm x 0.2um |
| 21100-020 | 2.1mm x 100mm x 2.0um |
| 21100 | 2.1mm x 100mm Tube Only |
| 21100-005 | 2.1mm x 100mm x 0.5um |
| 21100-005HP | High Pressure 2.1mm x 100mm x 0.5um |
| 21100-005HPW | High Pressure 2.1mm x 100mm x 0.5um -Waters Seat |
| 21100-002HP | High Pressure 2.1mm x 100mm x 0.2um |
| 21125 | 2.1mm x 125mm Tube Only |
| 21125-020 | 2.1mm x 125mm x 2um |
| 21125-005 | 2.1mm x 125mm x 0.5um |
| 21125-005HP | High Pressure 2.1mm x 125mm x 0.5um |
| 21150 | 2.1mm x 150mm Tube Only |
| 21150-020 | 2.1mm x 150mm x 2.0um |
| 21150-005 | 2.1mm x 150mm x 0.5um |
| 21150-020HP | High Pressure 2.1mm x 150mm x 2um |
| 21150-005HP | High Pressure 2.1mm x 150mm x 0.5um |
| 21150-002HP | High Pressure 2.1mm x 150mm x 0.2um |
| 21200 | 2.1mm x 200mm Tube Only |
| 21200-020 | 2.1mm x 200mm x 2.0um |
| 21200-005 | 2.1mm x 200mm x 0.5um |
| 21250 | 2.1mm x 250mm Tube Only |
| 21250-020 | 2.1mm x 250mm x 2.0um |
| 21250-005 | 2.1mm x 250mm x 0.5um |
| 21250-005HP | High Pressure 2.1mm x 250mm x 0.5um |
| 21250-002HP | High Pressure 2.1mm x 250mm x 0.2um |
| 21300 | 2.1mm x 300mm Tube Only |
| 21300-005 | 2.1mm x 300mm x 0.5um |
| 21300-020 | 2.1mm x 300mm x 2.0um |

| 3.0mm Modular Column Hardware | |
|--------------------------------------|--|
| PART # | Description |
| 30020 | 3.0mm x 20mm Tube Only |
| 30020-020 | 3.0mm x 20mm x 2.0um |
| 30020-005 | 3.0mm x 20mm x 0.5um |
| 30020-002HP | High Pressure 3.0mm x 20mm x 0.2um |
| 30025 | 3.0mm x 25mm Tube Only |
| 30025-005HP | High Pressure 3.0mm x 25mm x 0.5um |
| 30030-020 | 3.0mm x 30mm x 2.0um |
| 30030 | 3.0mm x 30mm Tube Only |
| 30030-005 | 3.0mm x 30mm x 0.5um |
| 30030-005HP | High Pressure 3.0mm x 30mm x 0.5um |
| 30030-002HP | High Pressure 3.0mm x 30mm x 0.2um |
| 30033 | 3.0mm x 33mm Tube Only |
| 30033-020HP | High Pressure 3.0mm x 33mm x 2um |
| 30050-020 | 3.0mm x 50mm x 2.0um |
| 30050 | 3.0mm x 50mm Tube Only |
| 30050-005 | 3.0mm x 50mm x 0.5um |
| 30050-002HP | High Pressure 3.0mm x 50mm x 0.2um |
| 30050-005HP | High Pressure 3.0mm x 50mm x 0.5um |
| 30050-020-005HP | High Pressure 3.0mm x 50mm x 2um & 0.5um |
| 30060 | 3.0mm x 60mm Tube Only |
| 30060-020HP | High Pressure 3.0mm x 60mm x 2.0um |
| 30060-005HP | High Pressure 3.0mm x 60mm x 0.5um |
| 30075-020 | 3.0mm x 75mm x 2.0um |
| 30075 | 3.0mm x 75mm Tube Only |
| 30075-005 | 3.0mm x 75mm x 0.5um |
| 30075-005HP | High Pressure 3.0mm x 75mm x 0.5um |
| 30075-002HP | High Pressure 3.0mm x 75mm x 0.2um |
| 30100-020 | 3.0mm x 100mm x 2.0um |
| 30100 | 3.0mm x 100mm Tube Only |
| 30100-005 | 3.0mm x 100mm x 0.5um |
| 30100-005HP | High Pressure 3.0mm x 100mm x 0.5um |
| 30100-002HP | High Pressure 3.0mm x 100mm x 0.2um |
| 30125-020 | 3.0mm x 125mm x 2.0um |
| 30125 | 3.0mm x 125mm Tube Only |
| 30125-005 | 3.0mm x 125mm x 0.5um |
| 30125-005HP | High Pressure 3.0mm x 125mm x 0.5um |
| 30150-020 | 3.0mm x 150mm x 2.0um |
| 30150 | 3.0mm x 150mm Tube Only |
| 30150-005 | 3.0mm x 150mm x 0.5um |
| 30150-020HP | High Pressure 3.0mm x 150mm x 2um |
| 30150-005HP | High Pressure 3.0mm x 150mm x 0.5um |
| 30150-002HP | High Pressure 3.0mm x 150mm x 0.2um |
| 30200 | 3.0mm x 200mm Tube Only |
| 30250-020 | 3.0mm x 250mm x 2.0um |
| 30250 | 3.0mm x 250mm Tube Only |
| 30250-005 | 3.0mm x 250mm x 0.5um |
| 30250-005HP | High Pressure 3.0mm x 250mm x 0.5um |
| 30300 | 3.0mm x 300mm Tube Only |
| 30300-005 | 3.0mm x 300mm x 0.5um |
| 30300-020 | 3.0mm x 300mm x 2.0um |

| 4.0mm Modular Column Hardware | |
|--------------------------------------|-------------------------------------|
| PART # | Description |
| 40020 | 4.0mm x 20mm Tube Only |
| 40020-020 | 4.0mm x 20mm x 2.0um |
| 40020-005 | 4.0mm x 20mm x 0.5um |
| 40030-020 | 4.0mm x 30mm x 2.0um |
| 40030 | 4.0mm x 30mm Tube Only |
| 40030-005 | 4.0mm x 30mm x 0.5um |
| 40030-005HP | High Pressure 4.0mm x 30mm x 0.5um |
| 40050-020 | 4.0mm x 50mm x 2.0um |
| 40050 | 4.0mm x 50mm Tube Only |
| 40050-005 | 4.0mm x 50mm x 0.5um |
| 40075 | 4.0mm x 75mm Tube Only |
| 40075-020 | 4.0mm x 75mm x 2.0um |
| 40075-005 | 4.0mm x 75mm x 0.5um |
| 40080 | 4.0mm x 80mm Tube Only |
| 40100-020 | 4.0mm x 100mm x 2.0um |
| 40100 | 4.0mm x 100mm Tube Only |
| 40100-005 | 4.0mm x 100mm x 0.5um |
| 40100-005HP | High Pressure 4.0mm x 100mm x 0.5um |
| 40125-020 | 4.0mm x 125mm x 2.0um |
| 40125 | 4.0mm x 125mm Tube Only |
| 40125-005 | 4.0mm x 125mm x 0.5um |
| 40125-005HP | High Pressure 4.0mm x 125mm x 0.5um |
| 40125-020-005 | 4.0mm x 125mm x 2um,0.5um |
| 40150-020 | 4.0mm x 150mm x 2.0um |
| 40150 | 4.0mm x 150mm Tube Only |
| 40150-020-005 | 4.0mm x 150mm x 2.0um,.5um |
| 40150-005 | 4.0mm x 150mm x 0.5um |
| 40200 | 4.0mm x 200mm Tube Only |
| 40250-020-005 | 4.0mm x 250mm x 2.0um,.5um |
| 40250-020 | 4.0mm x 250mm x 2.0um |
| 40250 | 4.0mm x 250mm Tube Only |
| 40250-005 | 4.0mm x 250mm x 0.5um |
| 40250-020HP | High Pressure 4.0mm x 250mm x 2.0um |
| 40300-005 | 4.0mm x 300mm x 0.5um |
| 40300 | 4.0mm x 300mm Tube Only |
| 40300-020 | 4.0mm x 300mm x 2um |
| 40300-020HP | High Pressure 4.0mm x 300mm x 2um |

| 4.6mm Modular Column Hardware | |
|-------------------------------|-------------------------------------|
| PART # | Description |
| 46020 | 4.6mm x 20mm Tube Only |
| 46020-020 | 4.6mm x 20mm x 2.0um |
| 46020-005 | 4.6mm x 20mm x 0.5um |
| 46020-002HP | High Pressure 4.6mm x 20mm x 0.2um |
| 46025 | 4.6mm x 25mm Tube Only |
| 46025-002HP | High Pressure 4.6mm x 25mm x 0.2um |
| 46030-002HP | High Pressure 4.6mm x 30mm x 0.2um |
| 46030-020 | 4.6mm x 30mm x 2.0um |
| 46030 | 4.6mm x 30mm Tube Only |
| 46030-005 | 4.6mm x 30mm x 0.5um |
| 46033-005 | 4.6mm x 33mm x 0.5um |
| 46050-020 | 4.6mm x 50mm x 2.0um |
| 46050 | 4.6mm x 50mm Tube Only |
| 46050-020-005 | 4.6mm x 50mm x 2.0um x 0.5um |
| 46050-005 | 4.6mm x 50mm x 0.5um |
| 46050-005HP | High Pressure 4.6mm x 50mm x 0.5um |
| 46050-002HP | High Pressure 4.6mm x 50mm x 0.2um |
| 46050-020HP | High Pressure 4.6mm x 50mm x 2um |
| 46075-020 | 4.6mm x 75mm x 2.0um |
| 46075-020-005 | 4.6mm x 75mm x 2.0um & 0.5um |
| 46075 | 4.6mm x 75mm Tube Only |
| 46075-005 | 4.6mm x 75mm x 0.5um |
| 46075-005HP | High Pressure 4.6mm x 75mm x 0.5um |
| 46075-002HP | High Pressure 4.6mm x 75mm x 0.2um |
| 46075-020-005 | 4.6mm x 75mm x 2um & 0.5um |
| 46100-020 | 4.6mm x 100mm x 2.0um |
| 46100-020-005 | 4.6mm x 100mm x 2.0um & 0.5um |
| 46100 | 4.6mm x 100mm Tube Only |
| 46100-005 | 4.6mm x 100mm x 0.5um |
| 46100-020HP | High Pressure 4.6mm x 100mm x 2.0um |
| 46100-002HP | High Pressure 4.6mm x 100mm x 0.2um |
| 46100-005HP | High Pressure 4.6mm x 100mm x 0.5um |
| 46120-005 | 4.6mm x 120mm x 0.5um |
| 46125-020 | 4.6mm x 125mm x 2.0um |
| 46125 | 4.6mm x 125mm Tube Only |
| 46125-020-005 | 4.6mm x 125mm x 2.0um, 0.5um |
| 46125-005 | 4.6mm x 125mm x 0.5um |
| 46150-020 | 4.6mm x 150mm x 2.0um |
| 46150 | 4.6mm x 150mm Tube Only |
| 46150-020-005 | 4.6mm x 150mm x 2.0um x 0.5um |
| 46150-005 | 4.6mm x 150mm x 0.5um |
| 46150-020HP | High Pressure 4.6mm x 150mm x 2.0um |
| 46150-005HP | High Pressure 4.6mm x 150mm x 0.5um |
| 46150-002HP | High Pressure 4.6mm x 150mm x 0.2um |
| 46200 | 4.6mm x 200mm Tube Only |
| 46200-020 | 4.6mm x 200mm x 2.0um |
| 46220-005 | 4.6mm x 220mm x 0.5um |
| 46250-020 | 4.6mm x 250mm x 2.0um |
| 46250 | 4.6mm x 250mm Tube Only |
| 46250-020-005 | 4.6mm x 250mm x 2.0um x 0.5um |
| 46250-005 | 4.6mm x 250mm x 0.5um |
| 46250-020HP | High Pressure 4.6mm x 250mm x 2.0um |
| 46250-005HP | High Pressure 4.6mm x 250mm x 0.5um |
| 46300-020 | 4.6mm x 300mm x 2.0um |
| 46300-005 | 4.6mm x 300mm x 0.5um |
| 46300 | 4.6mm x 300mm Tube Only |