CozyZone Technical Information

USP	Description	Size (mm)
Base Receptor Sheet	 Easy to fit Easily removed for cleaning. Easy to reapply to bed after cleaning 	Cot 1200 x 600 Cot Bed 1400 x 700 Single Bed 1900 x 900 Special Sizes available on request.
Brackets	 Solid coloured brackets for easy identification. Lightweight and supplied with handy grip holes. 	SMALL - H 80 W 80 D 140 MEDIUM - H 150 W 150 D 165 LARGE - H 200 W 150 D 200
Bracket Pads	 Easy to apply to bracket and remove for cleaning. Wipe cleanable. Durable with our unique CozyTherm ACTIVE Temperature Management Technology within the PU. 	SMALL - H 80 W 80 D 140 MEDIUM - H 150 W 150 D 165 LARGE - H 200 W 150 D 200

CozyZone FLIP Overmantle CozyZone FLIP Cushions

- 'FLIP' 38mm Overmantle to offer either/or ACTIVE Temperature Management.
- White side for ACTIVE Skin Management at 33 degrees.
 CozyTherm ACTIVE Temperature management impregnated in PU surface.
- Grey side for warmth and coziness –
 ideal for winter months or for those
 who find it difficult to generate/retain
 heat.
- Additional comfort with 38mm thick and soft memory foam. Excellent body immersion qualities.
- Wipe cleanable.

Cot 1200 x 600 Cot Bed 1400 x 700 Single Bed 1900 x 900 Special Sizes available on request.







- Offering dynamic therapy through fragmentation of memory foam.
- Wipe clean PU surface.
- 'FLIP' Cushions to offer either ACTIVE Temperature Management or to support Heat generation and help maintaining heat.
- White side for ACTIVE Skin Management at 33 degrees.
 CozyTherm ACTIVE Temperature management impregnated in PU surface.
- Grey side for warmth and coziness ideal for winter months or for those

Medium Square FLIP CozyTherm Cushion PU	50 x 50
Large Square FLIP CozyTherm Cushion PU	70 x 60
Medium Rectangle FLIP CozyTherm Cushion PU	50 x 27.5
Large Rectangle FLIP CozyTherm Cushion PU	75 x 45
Small Roll CT Cushion PU	30 x Ø 8
Medium Roll CT Cushion PU	40 x Ø 13
Large Roll CT Cushion PU	60 x Ø 13.5
Extra Large Roll CT Cushion PU	100 x Ø 15

Leg Support	 who find it difficult to generate/retain heat. Fragmented 'Repurposed' memory foam from luxury furnishings that would have been sent as 'waste' material. Animal themed leg support. Encourages new Users of sleep systems with child friendly cushion. Leg support can be used in supine or side lying. Easy to roll when required. ACTIVE Temperature Management. Crib 5 Fire Retardancy Fragmented cushions for Dynamic support Prevent Internal/External rotation of legs. Symmetrical alignment. Manage Contractures/Reduce Tone 	2 Sizes - Can be used for anyone from 2 years' old to 6ft 2in tall. Small Large	
CozyZone Cushion Covers	 Super soft fabric using yarn sourced from luxury furnishings manufacturer. 	Medium Square FLIP CozyTherm Cushion PU	50 x 50
	Easy to remove and refit for	Large Square FLIP CozyTherm Cushion PU	70 x 60
	laundering.	Medium Rectangle FLIP CozyTherm Cushion PU	50 x 27.5
	4 way stretch to conform to shape of	Large Rectangle FLIP CozyTherm Cushion PU	75 x 45
	cushions.	Small Roll CT Cushion PU	30 x Ø 8
		Medium Roll CT Cushion PU	40 x Ø 13
		Large Roll CT Cushion PU	60 x Ø 13.5
		Extra Large Roll CT Cushion PU	100 x Ø 15

