




ASTM Level 2 & 3



Color	black
Filtration - BFE	>98%
Filtration - PFE	>98%
Fluid Resistance, mmHG	160
Delta P, mmH ₂ O/cm ²	<6.0
Flame Spread	Class 1
Ear loop	black elastic
Packaging	50/box - 500/cs



How to wear a mask correctly




Use elastic ear loops to place mask with the blue side as the outer layer.



Stretch mask to cover nose and mouth.



Pinch nose wire and mold to face.




Unhook elastic ear loops, discard, and wash hands.



Pieces: 50
 GW (lbs): 0.45
 Dimensions:
 SKU: PPE-M-AS-1105B



Pieces: 40
 Boxes: 8
 GW (lbs): 4.28
 Dimensions: 16x8x8
 Per Pallet: 50
 Per Layer: 10
 Layers: 5



Pieces: 2,000
 Boxes: 40
 GW (lbs): 20.23
 Dimensions: 21x16x16
 Per Pallet: 24
 Per Layer: 6
 Layers: 4

