



This 5 layer N95 respirator is produced at our Santa Ana, CA USA production facility and is being tested by NIOSH. Upon approval, these N95 respirators with yellow head straps may be used in medical environments per the U.S. Emergency Use Authorization.



US-MADE **5 LAYER N95** **RESPIRATORS**

MODEL:

Pure-5

PRODUCT NAME:

5 Layer N95 (Pure-5) Respirator

MATERIALS:

Layer 1 (Interior): Non-woven Polyester

Layer 2: Melt-blown Polypropylene

Layer 3: Melt-blown Polypropylene

Layer 4: Hot cotton fabric

Layer 5 (Exterior): Non-woven Polyester

PACKING DETAILS:

20 N95 Respirators (Individually Wrapped) Per Carton

32 Cartons Per Case (640 Respirators Per Case)

12 Cases Per Pallet, 36 Cases

PACKING SIZE SPECIFICATIONS:

Box: 5.5" x 3" x 5.5"

Carton: 11.25" x 13.4" x 22.2"

NET WEIGHT:

16.75 lbs

GROSS WEIGHT:

18.95 lbs

