



ASTM Level 3

| | |
|---------------------------------------------|------------------------------|
| Color | blue outside white inside |
| Filtration - BFE | >98% |
| Filtration - PFE | >98% |
| Fluid Resistance, mmHG | 160 |
| Delta P, mmH ₂ O/cm ² | <6.0 |
| Flame Spread | Class 1 |
| Ear loop | white elastic |
| Packaging | 50/box - 500/cs |



How to wear a mask correctly



Use elastic ear loops to place mask with the blue side as the outer layer.



Stretch mask to cover nose and mouth.



Pinch nose wire and mold to face.



Unhook elastic ear loops, discard, and wash hands.

PPE
 Get Protected With Bigfoot



Pieces: 50
 GW (lbs): 0.45
 Dimensions:
 SKU:
 PPE-M-AS-SUR-LVL3-S



Pieces: 40
 Boxes: 8
 GW (lbs): 4.28
 Dimensions: 16x8x8
 Per Pallet: 50
 Per Layer: 10
 Layers: 5



Pieces: 2,000
 Boxes: 40
 GW (lbs): 20.23
 Dimensions: 21x16x16
 Per Pallet: 24
 Per Layer: 6
 Layers: 4

