

Respokare®

Antiviral N95 NIOSH

Features

- USA Patented antiviral technology.
- Usage of Soft and comfortable material.
- Adjustable nose clip and elastic headband.
- Easy to carry around by V shape.
- Fitted to human face tightly.
- Free to breathe.
- Three sizes (L,M,S) available for option.
- NIOSH approved, at least 98% filtration efficiency.
- FDA registered for medical use.



Applications

Widely used in Hospitals, Commercial Buildings, Construction, Design & Construction, Food Processing, Food Safety, General Manufacturing, Heavy Infrastructure, Mining, Oil & Gas, Transportation.

Product Specification

Item Number	Size	Mask Size	Package	Carton Size
RK-200-3040A	S	4.0" x 5.0" x 4.5"	30 pieces/box 24 boxes/carton	585X435X445mm
RK-200-3041A	M	4.5" x 5.7" x 5.0"	30 pieces/box 24 boxes/carton	585X435X445mm
RK-200-3042A	L	5.0" x 6.5" x 5.5"	30 pieces/box 24 boxes/carton	585X435X445mm

