

Patient Benefits by product and performance level.

TruHearing®

| | Performance Level | Premium | Advanced |
|--|---|--------------------------------|--------------------------------|
|  | Hear conversations more easily with greater focus on speech and less on the challenging soundscape around you. (RIC, BTE) | ✓ | ✓ |
|  | Automatic adjustments via Motion Sensors to adapt to unique listening needs the environment changes. (RIC, BTE) | ✓ | ✓ |
|  | Natural voice and higher acceptance via Own Voice Processing to avoid having your voice sound too loud or unnatural. (RIC, BTE) | ✓ | ✓ |
|  | Hearing phone calls in both ears by using the TwinPhone feature- holding up phone to ear and having it stream automatically to the opposite ear OR take advantage of Bluetooth® (BT) streaming via a smartphone to hear the phone conversation in both ears*. (RIC, BTE, Instant-fit, Custom) | ✓ | ✓ |
|  | Enhance music with preset programs via HD Music. Included options for live music, recorded music, and personal performances. (RIC, BTE, Instant-fit, Custom) | 3 | 3 |
|  | Relief for ringing in the ears (Tinnitus) via Therapy solutions and exclusive Notch Therapy. (RIC, BTE, Custom) | ✓ | ✓ |
|  | Automatically focus on speakers in front of you via Directionality for better speech understanding in noise. (RIC, BTE, Instant-fit, Custom) | ✓ | ✓ |
|  | Customize focus in the TruHearing App via Directional Hearing to adjust as needed to better suit your personal preferences in any situation - dining in restaurants, walking outdoors, or visiting with friends. (RIC, BTE, Custom) | Front Back Right Left | Front Back Right Left |
|  | Automatic spatial awareness via Spatial SpeechFocus to hear speech and noise from all directions, like hearing those around you when in a car. (RIC, BTE, Custom) | Front Back Right Left | Front Back Surround |
|  | Reduced wind noise via eWindscreen for better clarity in speech when outdoors on a windy day or when riding in a car with the top down. (RIC, BTE, Instant-fit, Custom) | Auto | Manual |
|  | Better awareness of high-frequency sound via Extended Bandwidth to hear higher pitches with enhanced sound quality, such as a bird chirping, a small child or an alarm. (RIC, BTE, Instant-fit, Custom) | ✓ | |
|  | Increase understanding in reverberant environments via echoShield to reduce echoes in large spaces like places of worship, museums, and auditoriums. (RIC, BTE, Instant-fit, Custom) | ✓ | |

Higher performance levels include more channels for additional customization options.

Checkmarks highlight feature availability at respective levels, while numbers / text highlight added functionalities.

*Bluetooth-integrated hearing aids needed for BT streaming. The Bluetooth® word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by Sivantos is under license. TruHearing® is a registered trademark of TruHearing, Inc. Copyright © 2021 Signia, a Division of WSAUD A/S. 8/21 2.0 SIV-11315-21