

Sterile Kit Processing

Asceptic Technique for Processing PRP

prep for centrifuge

1 ☐ apply and disinfect gloves with cavicide

2 ☐ prepare tray

a ☐ double female luer lock (optional)

b ☐ PRP processing device

3 ☐ transfer blood to PRP system

process in centrifuge

4 ☐ disinfect centrifuge

5 ☐ replace and disinfect gloves

6 ☐ process blood in centrifuge

once it is completed

7 ☐ replace and disinfect gloves

8 ☐ remove PRP processing device from centrifuge

Sterile Technique for Injection

9 ☐ prepare tray

a ☐ disinfected probe

b ☐ tegaderm for probe, optional

c ☐ sterile gel

d ☐ sterile gauze soaked in stanhexidine

e ☐ treatment needle

prp treatment

10 ☐ open sterile drape

11 ☐ apply and disinfect gloves

12 ☐ position drape

13 ☐ apply sterile gloves

14 ☐ apply soaked gauze to area surrounding treatment

15 ☐ apply sterile gel to treatment area

16 ☐ complete injection

17 ☐ dry skin

18 ☐ apply bandaid or gauze