



# Revolutionizing vascular healthcare with affordable, high-quality devices

Delivering innovative vascular medical devices, with a special emphasis on Outpatient and Ambulatory Surgery Centers (OBL/ASCs).

## Falexexpand HP PTA Balloon Catheter


High-Pressure with RBP up to 24atm  
0.035" OTW, 6F, Dia:3-12mm, BL: 20-  
200mm, CL: 40cm, 75cm, and 135cm

High-pressure and semi-compliant PTA  
balloon catheters with soft and tapered  
tips for better lesion crossing with  
hydrophilic coating and small profiles for  
challenging stenosis



## Falspeed Semi- Compliance PTA Balloon Catheter

Semi-compliant with Dia: 2.0 to 7.0mm,  
BL: 20 to 300mm, CL: 150cm (OTW), Tip  
Crossing: 0.0225", 4F to 6F, NP: 6atm



## 014" Workhorse Guidewire (300cm)

Distal 40cm high-load Nitinol for fast  
recovery with Stainless Steel body for  
predictable tracking.  
Shapable tip design with 0.8g & 1.5g tip  
weight options; moderate support



## 018", 032", and 035" Workhorse Guidewire

Nitinol core wire guarantees great tip shape  
resistance, outstanding kink resistance  
and 1:1 torque control  
Straight, Angled, and Shapeable tips,  
Standard and Stiff shaft



## Balancing Quality & Affordability

Contact Sepehr Fariabi (+1-949-726-2846) or visit [www.worldwidevascular.com](http://www.worldwidevascular.com) to learn more