

Mobile LTE

ON-THE-GO SOLUTION

1-888-311-3120

<https://advancedmedicalalertllc.com/>



Equipment included:
All-in-one neck button, charging cradle, and interchangeable neck or wrist button

What is it?

Designed to work on the AT&T® cellular network, this system integrates a built-in two-way speaker and an interchangeable, lightweight neck or wrist button. This simple mobile solution includes an advanced GPS locator* and clear voice communication, so that users have immediate access to help anytime, anywhere. It is compatible with our optional fall detection pendant which can call for help if a fall is detected, even if a user cannot push the help button.**

How does it work?



When you need help, just press your button.



Your information is sent to the emergency response center. An emergency response operator connects you through the speaker on your device and decides on the proper action.



The operator sends the appropriate responder based on your needs - emergency services, caregiver, family member or neighbor.

Who is it for?

This simple safety solution is ideal for anyone looking for an easy way to maintain their independence wherever they go.

It works for users who do not have a cellular phone or landline, so anyone can use it. It can be worn in the shower or bath where accidents often occur. Plus, it has a battery life of up to five days, depending on your activity, which means it's ready to go when you are.

This system is also particularly beneficial for anyone who lives or travels alone, and those with medical conditions.

Mobile LTE

ON-THE-GO SOLUTION

Key Features

- No landline or cellular phone required
- Works on the AT&T® cellular network
- No cellular contract required
- Advanced GPS technology*
- One button access to emergency help through built-in two-way speaker
- Optional fall detection pendant can call for help if a fall is detected **
- Interchangeable, lightweight neck or wrist button can be worn in the shower or bath where accidents are common
- Up to 5 days battery depending on activity and usage

Specifications

- | | |
|-----------------------------|--|
| • Dimensions | 1.75 " x 2.25 " (W x L) |
| • Weight | 1.9 ounces |
| • Installation Requirements | Power outlet, AT&T® cellular coverage† |
| • Maintenance Requirements | Test monthly, charge daily |
| • Battery | Up to 5 days |



* Location accuracy may vary depending on environmental conditions.

** Fall detection feature does not detect 100% of falls. If able, users should always press their button when they need help.

† Cellular strength may vary depending on environmental conditions.