

# Mobile

ON-THE-GO SOLUTION

1-888-311-3120

<https://advancedmedicalalertllc.com/>



SOLUTION SHEET

## What is it?

Designed to work on the AT&T® cellular network, this system integrates a built-in two-way speaker with a lightweight, all-in-one help button that can be worn around the neck. This simple mobile solution includes an advanced GPS\* locator with clear voice communication, so that users have immediate access to help anytime, anywhere. Optional fall detection is available, which can call for help if a fall is detected and a user cannot push the help button.\*\*



Equipment included:  
All-in-one neck button with two-way speaker and charging cradle

## How does it work?



When you need help, just press your button.



Your information is sent to the emergency response center. An emergency response operator connects you through the speaker on your device and decides on the proper action.



The operator sends the appropriate responder based on your needs - emergency services, caregiver, family member or neighbor.

## Who is it for?

This simple safety solution is ideal for anyone looking for an easy way to maintain their independence wherever they go.

It works for users who do not have a cellular phone or landline, so anyone can use it. It can be worn in the shower or bath where accidents often occur. Plus, it has a battery life of up to five days, depending on your activity, which means it's ready to go when you are.

This system is also particularly beneficial for anyone who lives or travels alone, and those with medical conditions.

# Mobile

## ON-THE-GO SOLUTION

### Key Features

- No landline or cellular phone required
- Works on the AT&T® cellular network
- No cellular contract required
- Advanced GPS technology\*
- One-button access to emergency help through built-in two-way speaker
- Optional built-in fall detection can call for help if a fall is detected\*\*
- Can be worn in the shower or bath where accidents are common
- Up to 5 days battery depending on activity and usage

### Specifications

- |                             |  |
|-----------------------------|--|
| • Dimensions                | 1.75 " x 2.25 " (W x L)                |
| • Weight                    | 1.9 ounces                             |
| • Installation Requirements | Power outlet, AT&T® cellular coverage† |
| • Maintenance Requirements  | Test monthly, charge daily             |
| • Battery                   | Up to 5 days                           |



\* Location accuracy may vary depending on environmental conditions.

\*\* Fall detection feature does not detect 100% of falls. If able, users should always press their button when they need help

† Cellular strength may vary depending on environmental conditions.