



PROUDLY MADE

BY AMERICANS
with Domestic & Imported Materials

Golden Lift Chairs

AMERICA'S #1 SELLING BRAND



Our power lift recliners come in a variety of shapes and sizes for the perfect fit. Our premium fabrics match nearly every décor and our patented technologies give you the most comfortable seating experience.



 GOLDEN

MaxiComfort SERIES



TWILIGHT



Discover Weightlessness in the MaxiComfort Cloud with Twilight positioning. Our patented, three-motor technology offers an array of rejuvenating options to cradle the body and ease the mind. Achieve extreme Zero-gravity and TV positions while enjoying the world's finest rejuvenating seating experience!

MaxiComfort Cloud with Twilight Positioning

PR514 shown in Brisa® CoffeeBean*




Quick Ship Fabrics for PR514

Hazelnut, Sterling, Copper, Calypso, Shiraz, Brisa® Coffee Bean*, Brisa® Buckskin* & Brisa® Bridle*

MaxiComfort Twilight PR514 Fabrics




Hazelnut



Sterling



Calypso



Shiraz




Copper



Brisa® Buckskin*



Brisa® Bridle*



Brisa® Coffee Bean*


Please view fabric colors in person for accuracy.

*Upcharge applies



TWILIGHT

Patented 3-Motor Seating Technology



Up


TV



TV in Twilight



Lounge in Twilight



Lounge Reclined
in Twilight



Zero Gravity



Zero Gravity
in Twilight



Trendelenburg



Motor #1
Adjusts back
positioning

Motor #2
Controls the footrest
movement and the
patented movement of the
rising seat box

Motor #3
Activates Twilight technology, a
tilting movement that provides
a weightless experience



PR514

SCAN WITH SMART
PHONE CAMERA FOR
"HOW-TO" VIDEO



MaxiComfort SERIES



The Cloud™ power lift recliners offer our unique, ergonomic seating system for long-term comfort and support. They feature our custom bucket seat and plush backrest. Available in two styles!



MaxiComfort Cloud
PR510-SME shown in Calypso



Quick Ship Fabrics for PR510-SME Hazelnut, Sterling, Copper, Shiraz, Calypso & Brisa® Coffee Bean*



AutoFlex optional.

MaxiComfort Cloud
PR510-MLA shown in Hazelnut



Quick Ship Fabrics for PR510-MLA Hazelnut, Sterling, Calypso, Shiraz, Copper, Brisa® Buckskin*, Brisa® Coffee Bean* & Brisa® Bridle*



The Cloud™ PR512 features our exclusive PowerPillow™ articulating headrest and a modern back design.

MaxiComfort Power Cloud
PR512 shown in Brisa® Coffee Bean*

PowerPillow™ STANDARD



Quick Ship Fabrics for PR512 Hazelnut, Sterling, Copper, Brisa® Buckskin*, Brisa® Coffee Bean* & Brisa® Bridle*

Shown with PowerPillow™ Option



SCAN FOR
MAXICOMFORT
PATENTED
POSITIONING
VIDEO

Cloud PR510, PR512 and Relaxer PR756 Fabrics



DayDreamer™ PR632 Fabrics



The DayDreamer™ is a medium sized power lift recliner offering our exclusive PowerPillow™ articulating headrest and The Cloud™ seating system, featuring our custom bucket seat design. Recline, relax, and adjust your head position with the touch of a button.



DayDreamer™

PR632 shown in Brisa® Coffee Bean*

PowerPillow™ STANDARD



SCAN FOR PR632 VIDEO

The Pub chair is a medium sized power lift recliner with a classic look, featuring a two pillow waterfall back, nail head trim on the arms and wings, and cut back arms with a full chaise pad and pocketed coil seat.



MaxiComfort Pub Chair

PR713 shown in Brisa® Bridle*



Quick Ship for PR632M Hazelnut, Calypso, Brisa® Coffee Bean* & Brisa® Bridle*



Quick Ship Fabrics for PR713 Mahogany, Geranium & Brisa® Bridle*



Relaxer

PR756MC shown in Sterling



Relaxer

PR756L shown in Copper

The Relaxer, the original MaxiComfort power lift recliner is offered in two sizes, featuring four individual zipper backed pillows for truly customizable comfort.



Quick Ship Fabrics for PR756MC & PR756L: Hazelnut, Sterling, Copper & Brisa® Coffee Bean*. Additional Quick Ship Fabric for PR756L: Brisa® Bridle*

Pub Chair PR713 Fabrics



Mahogany



Geranium



Brisa®



Bridle* = PowerPillow™ option

= pocketed coil spring seat std.

Please view fabric colors in person for accuracy.

*Upcharge applies

MaxiComforters

The MaxiComforters power lift recliners are now offered in SIX SIZES, ensuring a perfect fit for just about anyone!



MaxiComforter Jr. Petite
PR505JP shown in Evergreen



Quick Ship for PR505JPT: Pearl



MaxiComforter Small
PR505S shown in Palomino

PowerPillow™ Option Available



Quick Ship for PR505S: Palomino, Admiral, Evergreen, Cabernet & Pearl



MaxiComforter Medium Extra Wide

PR505M26 shown in Admiral

(AutoDrive not available on this model.)



Quick Ship for PR505M26: Admiral & Palomino

MaxiComforter Fabrics



Palomino



Admiral




Evergreen



Cabernet



Pearl



Coffee Bean*

Please view fabric colors in person for accuracy.

*Upcharge applies



MaxiComforter Medium

PR505M shown in Palomino

Shown with PowerPillow™



MaxiComforter Large

PR505L shown in Cabernet

PowerPillow™ Option Available



Quick Ship Fabrics for PR505M: Palomino, Admiral, Evergreen, Cabernet, Pearl & Brisa® Coffee Bean*



Quick Ship for PR505L: Palomino, Admiral, Cabernet & Pearl



MaxiComforter Tall PR505T

shown in Palomino

PowerPillow™ Option Available



Quick Ship for PR505T: Admiral, Palomino & Cabernet



= PowerPillow™ option



= pocketed coil spring seat std.

Please view fabric colors in person for accuracy.

*Upcharge applies



Comforter SERIES (a perfect fit)

The Comforter series was designed to fit individuals small to tall with weight capacities varying from 300 to 375 lbs. These chairs give the perfect sizes for the perfect fit.



Jr. Petite
PR501JP shown in Admiral



Small
PR501S shown in Pearl



Medium
PR501M shown in Palomino



Quick Ship for PR501JP & PR501S:
Palomino, Admiral, Cabernet & Pearl



Quick Ship Fabrics for PR501M: Palomino, Admiral, Evergreen, Cabernet & Pearl

Comforter Wides

The Comforter Wides are perfect for those who need more space between the armrests. Constructed using the highest quality components for reliable, steady lifting and reclining, these chairs offer weight capacities varying from 375 to 700 lbs.



Small 23
PR501S-23 shown in Cabernet



Medium 26
PR501M-26 shown in Pearl



Large 26
PR501L-26 shown in Evergreen
(Dual Motor)



Quick Ship for PR501S-23:
Palomino & Admiral



Quick Ship for PR501M-26:
Palomino, Cabernet & Admiral



Quick Ship for PR501L-26
Palomino

Comforter & Comforter Wide Fabrics



Palomino Admiral Evergreen Cabernet Pearl

Please view fabric colors in person for accuracy.



Did you know that Golden power lift recliners can be custom tailored to fit you perfectly? Golden Flagship Dealers are specially trained to customize a Comforter to fit you perfectly. So whether you are small, tall, or up to 700 lbs., the perfect fit is only a phone call away.



Large
PR501L shown in Evergreen



Tall
PR501T shown in Cabernet



Quick Ship Fabrics for PR501L:
Palomino, Admiral, Evergreen,
Cabernet & Pearl



Quick Ship for PR501T: Palomino



Tall 28
PR501T-28 shown in Admiral
(Dual Motor)



Super 33
PR502 shown in Palomino
(Triple Motor)
(AutoDrive not available on this model.)
(Freight Supplement applies.)



Quick Ship Fabrics for PR501T-28 & PR502:
Palomino



 = pocketed
coil spring
seat std.

SIGNATURE SERIES
(multi-functional)



**USB Port
Standard!**



Introducing the world's first three-position recliner with Twilight Technology! Our dual-motor system comes standard and allows for a deeper recline to achieve rejuvenating positions no other three-position chair on the market offers.




Orion with Twilight Positioning

PR405-MLA *shown in Sandstorm*



Quick Ship for PRPR405-MLA: Brisa® Coffee Bean*, Brisa® Buckskin*, Anchor, Sandstorm & Bittersweet

Orion Twilight PR405 Fabrics



Brisa®
Coffee Bean*



Brisa®
Buckskin*



Anchor



Sandstorm



Bittersweet

Please view fabric colors in person for accuracy.



= pocketed coil spring seat std.



*Upcharge applies

INTRODUCING!



TWILIGHT

Patented 3-POSITION Seating Technology



Motor #1
Controls the movement of both the backrest and footrest.

Motor #2
Activates Twilight technology, a tilting movement that provides a weightless experience.

SCAN FOR
PR405 VIDEO



SIGNATURE SERIES

(multi-functional)



New compact 6-button remote



America's favorite lift and recline chair has a powered articulating headrest and power lumbar system to truly customize to your comfort from top to bottom!

Imperial PR404

shown in Brisa® Coffee Bean*

PowerPillow™ & Power Lumbar STANDARD



Quick Ship Fabrics for PR404:
Hazelnut, Sterling, Brisa® Buckskin* & Brisa® Coffee Bean*



Designed for the Spine

Articulating headrest adjusts for personalized comfort



Power Lumbar



Imperial PR404 Fabrics



Please view fabric colors in person for accuracy.



= PowerPillow™ option



= pocketed coil spring seat std.

*Upcharge applies



Cambridge
PR401SME
shown in Sterling



Cambridge
PR401MLA
shown in Hazelnut



Quick Ship Fabrics for PR401-SME & PR401-MLA: Hazelnut, Sterling, Copper, Calypso, Shiraz & Brisa® Coffee Bean*. **Additional Quick Ship Fabric for PR401-MLA:** Bamboo*



Available in three sizes, this chair glides forward and fits perfectly in almost any room.



Space Saver
PR931S shown in Sandstorm
Needs only 4" from wall to recline.



Space Saver
PR931M shown in Anchor
Needs only 5" from wall to recline.



Space Saver
PR931L shown in Bittersweet
Needs only 6" from wall to recline.



Quick Ship Fabrics for PR931 S, M, L: Anchor, Bittersweet & Sandstorm

SCAN
FOR PR931
VIDEO



Cambridge PR401 Fabrics



Space Saver PR931 Fabrics



Please view fabric colors in person for accuracy.

*Upcharge applies

SIESTA SERIES

Siesta

Enjoy infinite positioning with our shear reducing motion seat! Power lumbar and a power headrest come standard with a USB port on the hand control.

Luxurious fabrics with designer contrast stitching will stand out in any décor!



Siesta
PR445* shown in Hickory



Power Headrest Standard



Power Lumbar Standard



Luxurious Fabric with Contemporary Contrast Stitching



Unique 2-way Motion Seat Reduces Shear



Soft-feel Hand Control with USB Charging Port

Siesta PR445 Fabrics



Ash



Timber



Hickory

Please view fabric colors in person for accuracy.

*This chair is designed by Golden and manufactured by our offshore partner.

IMPORTED PRODUCT

VALUE SERIES

(cost-effective)

This series combines ultra-soft lumbar pillow support with our unique "grip-rite" arm design for easy transfer.



Capri
PR200 shown in Cornflower



Monarch Medium
PR355M shown in Autumn



Monarch Large
PR355L shown in Spring



Quick Ship Fabrics for PR200, PR355M & PR355L:
Autumn, Cornflower, Spring, Rosewood & Chestnut

Value Fabrics



Autumn Cornflower Spring Rosewood Chestnut

Please view fabric colors in person for accuracy.

**Upcharge applies.*



Features that make Golden Power Lift Recliners America's #1 brand!

Golden Technologies is the world's leading manufacturer of power lift recliners. Founded in 1985, Golden has become America's most trusted and reliable brand. All chairs are handcrafted by American workers with domestic and imported materials (with the exception of the PR445 Siesta Series chair, which is designed by Golden and manufactured by our offshore partner). Each chair undergoes over 42 quality control inspections before the product is approved for shipping. Golden is the pioneer of innovative products and features.



(Excluding PR445 Siesta Series. IMPORTED PRODUCT)



Self Diagnostic System



Intuitive Hand Control



Lumbar Relief Technology



Relaxing Heat and Massage



Brisa



49 Pocketed Coil Seat Springs



GOLDEN
AutoDrive™ 3.0
CUSTOMIZED COMFORT

TOTAL CONTROL

Operates the basic movement of the chair.

POSITION CONTROL

Programmable positioning buttons for individualized comfort.

CUSTOM CONTROL

Fine tunes any position by moving the backrest and footrest independently.



AutoDrive™ is the exclusive hand control system for MaxiComfort power lift recliners. This smart system has taken the guesswork out of getting into or out of your favorite positions. The AutoDrive control assists in finding favorite comfort positions and it can memorize up to four positions.



AUTODRIVE
SCAN FOR
VIDEO

* Not available on models PR505M26 & PR502.

THE BRISA® COLLECTION

GOLDEN OFFERS THE ONLY LIFETIME* WARRANTY for all Brisa® breathable fabrics on lift chairs in the industry!



BRISA®
GOLDEN
SCAN FOR
VIDEO



** Limited Lifetime Warranty, as long as normal use is maintained. Restrictions apply. See Owner's Manual for details.*

NOTE: Crossing fabrics & series is available. Charges will apply and lead times will change.


NEW! HIGH PERFORMANCE FABRIC



SEE IT IN
ACTION!



Bamboo
(DBA)



Morning Glory
(DMG)



Eucalyptus
(DEU)



www.goldentech.com



Complies with
UL 73 and
CSA C22.2 No. 68

YouTube

[www.YouTube.com/
GoldenTechChannel](http://www.YouTube.com/GoldenTechChannel)

f Find us on
Facebook

www.facebook.com/goldentechnologies

PRODUCT SPECIFICATIONS

Name	Model	Weight Capacity	Overall Width	Overall Height	Seat to Top of Back	Width Between Arms	Seat Depth	Floor to Top of Seat	Extension - While Reclined from Back of Seat to End of Footrest	Floor to Top of Footrest - Reclined	Heat & Massage Option * Available	Footrest Extension Option Available	Straight Lift Option Available*		
MaxiComfort®	Cloud w/ Twilight	PR514	375 lbs.	39.5"	46.5"	32.5"	20"	20"	20.5"	43"	25"	No	No	No	
	Pub Chair	PR713	375 lbs.	35"	45"	26.5"	21"	18"	20"	43.5"	23"	Yes	No	No	
	DayDreamer	PR632	375 lbs.	38.5"	45"	31"	22.5"	20.5"	20"	42.5"	24"	Yes	No	No	
	Cloud SME	PR510-SME	375 lbs.	33.5"	41"	27"	18.5"	19.5"	19"	43"	23"	Yes	No	No	
	Cloud MLA	PR510-MLA	375 lbs.	38"	45"	29"	20"	21.5"	20"	44"	24.5"	Yes	No	No	
	Power Cloud	PR512-MLA	375 lbs.	38"	45"	30"	20"	21"	20"	44"	24.5"	Yes	No	No	
	Relaxer	PR756MC	375 lbs.	36.5"	44"	29"	20"	21"	19"	43"	23"	Yes	Yes	No	
	Relaxer	PR756L	375 lbs.	39"	46"	30"	22"	22"	20.5"	44"	24"	Yes	Yes	No	
	Comforter JP	PR505JP	300 lbs.	30.5"	39"	25"	19"	18.5"	18"	38"	22"	Yes	No	No	
	Comforter Sm.	PR505S	300 lbs.	31.5"	40"	25"	20"	19.5"	19"	39"	23"	Yes	Yes	No	
	Comforter Md.	PR505M	375 lbs.	32.5"	42"	27"	21"	21"	19.5"	42"	23"	Yes	Yes	No	
	Comforter Lg.	PR505L	375 lbs.	34.5"	44"	28.5"	23"	22"	20.5"	43"	23"	Yes	Yes	No	
	Comforter Tall	PR505T	375 lbs.	34.5"	47"	30"	23"	23"	22.5"	45"	23.5"	Yes	Yes	No	
Comforter M-26	PR505M-26	500 lbs.	38"	42"	28"	26"	20.5"	20"	42"	23"	Yes	Yes	No		
Comforter	Junior Petite	PR501JP	300 lbs.	30.5"	39"	25"	19"	18.5"	18"	38"	22"	Yes	No	Yes	
	Small	PR501S	300 lbs.	31.5"	40"	25"	20"	19.5"	19"	39"	23"	Yes	Yes	Yes	
	Medium	PR501M	375 lbs.	32.5"	42"	27"	21"	21"	19.5"	42"	23"	Yes	Yes	Yes	
	Large	PR501L	375 lbs.	34.5"	44"	28.5"	23"	22"	20.5"	43"	23"	Yes	Yes	Yes	
	Tall	PR501T	375 lbs.	34.5"	47"	30"	23"	24"	22.5"	45"	23.5"	Yes	Yes	Yes	
	Comforter Wides	Small X-wide	PR501S-23	375 lbs.	34.5"	40"	25"	23"	19.5"	18"	39"	23"	Yes	Yes	Yes*
		Medium-26 Dbl.	PR501M-26D	500 lbs.	38.5"	42"	27"	27.5"	21.5"	19.5"	42"	23"	Yes	Yes	Yes*
		Large-26 Dbl.	PR501L-26D	500 lbs.	38.5"	44"	28.5"	27.5"	22.5"	20.5"	43"	23"	Yes	Yes	Yes*
		Tall-28 Dbl.	PR501T-28D	500 lbs.	41"	48"	30"	29.5"	24"	22.5"	45"	23"	Yes	Yes	Yes*
		Super 33	PR502	700 lbs.	45"	47"	29"	33"	22"	23"	45"	21"	Yes	No	No
Signature	Orion w/ Twilight	PR405	375 lbs.	40.5"	46.5"	30.5"	23"	21"	21"	43"	25"	No	No	No	
	Imperial	PR404	375 lbs.	32.5"	43"	28.5"	19"	22"	19"	44"	23"	No	No	Yes	
	Cambridge	PR401-SME	375 lbs.	33"	41.5"	27"	20"	19"	19"	41.5"	22"	Yes	Yes	Yes	
	Cambridge	PR401-MLA	375 lbs.	35"	44.5"	31"	22'	22"	22"	44"	23"	Yes	Yes	Yes	
	SpaceSaver Sm.	PR931S	350 lbs.	30.5"	39.5"	26.5"	17.5"	20"	18"	41"	22"	Yes	No	No	
	Space Saver Md.	PR931M	350 lbs.	32.5"	41.5"	28.5"	20"	21"	19"	42.5"	22"	Yes	No	No	
	Space Saver Lg.	PR931L	350 lbs.	35"	43.5"	30.5"	22.5"	22.5"	20"	44"	22"	Yes	No	No	
Siesta	PR445	375 lbs.	34"	45"	30.5"	21.5"	19.5"	20.5"	42"	18.5"	No	No	No		
Value	Capri	PR200	325 lbs.	31"	41"	27"	20"	19"	18"	40"	23"	No	No	Yes	
	Monarch Md.	PR355M	375 lbs.	32"	42"	28"	21"	20.5"	18.5"	43"	23"	Yes	No	Yes	
	Monarch Lg.	PR355L	375 lbs.	34"	44"	30"	23"	21.5"	19.5"	43"	22"	Yes	No	Yes	

*Must be set at factory

**All chairs feature a battery backup system for power outages (9 volt batteries not included).



PROUDLY MADE



BY AMERICANS

with Domestic & Imported Materials



**Proudly Made by American
Workers in Old Forge, Pa.**



 **GOLDEN**