



# Golden Scooters & Power Wheelchairs

SUPERIOR DESIGNS | FIRST CLASS LUXURY





## Relax & Ride in Style!

**Our scooters and power wheelchairs offer excellent maneuverability and performance to get you where you want to go in comfort and style. Explore our line of mobility products and choose the one that's perfect for you! Plus, all models are backed by the best warranties in the industry!**

Highly Dependable • Maximum Maneuverability  
Unparalleled Stability • Greatest Functionality  
Superior Designs • First Class Luxury



## Luxury Meets Rugged in TWO NEW Powerful Scooters!

**INTRODUCING America's ONLY Full-time Luxury Travel Scooter! The GB119 3-wheel and GB149 4-wheel Buzzaround LX models offer comfort and convenience in a stylish new design.**

Enjoy upgraded features in a whole new class of scooters  
in our Buzzaround Series!

# BUZZAROUND LX

LUXURY

Full-size comfort comes with the convenience of true portability on the Buzzaround LX, along with standard features that include a stylish captain's seat, an infinitely adjustable tiller at your fingertips, and a USB charging port right under the dashboard. With a 375 lb. weight capacity, comfort suspension springs, and up to 18 mile range on two U1 batteries, the Buzzaround LX offers ample foot room to take you where you want to go in comfort and style!

**BEST IN CLASS  
RANGE!**

**BEST IN CLASS  
WEIGHT CAPACITY!**



Buzzaround LX  
GB119

**WARRANTIES:**

**FRAME:** Lifetime Warranty

**DRIVE TRAIN:** 2 Year Warranty

**ELECTRONICS:** 2 Year Warranty

**BATTERIES:** 1 Year Warranty



**BUZZAROUND LX  
FEATURES:**

- Stylish Captain's Seat
- 375 lb. Weight Capacity
- Infinitely Adjustable Tiller with EZ Reach Adjustment Lever
- Up to 18 Mile Range with two U1 Batteries
- Ample Foot Room
- USB Charging Port
- Comfort Spring Suspension for a Smoother Ride
- Red & Blue Color Shrouds Standard



Buzzaround LX  
GB149

# BUZZ<sup>TM</sup> AROUND EX EXTREME

## FULL-SIZE PORTABLE

The all new Golden BuzzAround Extreme Scooter features the portability and easy disassembly of a lightweight scooter with the amazing benefits of a full size scooter. Both the three and four-wheel BuzzAround Extreme Scooters come standard with 2-U1 batteries – industry leading and the largest battery in this category – for maximum operating range up to 18 miles. The Extreme also features a brand new ultrabright LED headlight and taillight, and comfort suspension springs for the smoothest ride.

Go the extra mile on a BuzzAround Extreme Scooter!

### BUZZAROUND EX FEATURES:

- 350 lb. weight capacity
- 4" ground clearance
- Full size & portable, yet easy to disassemble
- Larger 90 amp controller
- Larger 420 watt motor
- Larger 9" tires
- 2-U1 batteries for maximum range
- Ultrabright LED headlight and taillight
- Comfort suspension springs
- Operating range up to 18 miles per charge

Optional color  
shrouds available.  
\*Upcharge applies.



**BEST IN CLASS  
RANGE!**





GB118

**WARRANTIES:**  
**FRAME:** Lifetime Warranty  
**DRIVE TRAIN:** 2 Year Warranty  
**ELECTRONICS:** 2 Year Warranty  
**BATTERIES:** 1 Year Warranty



### OPTIONAL BATTERY DOCKING STATION

The Battery Docking Station for models GB118, GB148, GB119 and GB149 is portable. It connects and charges the U1 batteries. Perfect for charging batteries off the scooter, especially when visiting family or friends, or at a vacation destination. Enjoy being on the go!



GB148



### OPTIONAL BUMPER CONVERSION KIT

Can be ordered as a KIT or pre-installed on models GB118, GB148, GB119 & GB149.





Optional color shrouds available.  
\*Upcharge applies.



**BUZZYL  
AROUND XL™**

GB117D

Optional color shrouds available.  
\*Upcharge applies.



**BUZZYL  
AROUND XL™**



GB147D

## BUZZAROUND XL FEATURES:

- 300 lb. weight capacity
- Ultra Bright LED headlight standard
- Charging port below the control panel
- 3 inch ground clearance
- Red and Blue Shroud panels included free
- Disassembles quickly and easily for transport

## WARRANTIES:

**FRAME:** Lifetime Warranty  
**DRIVE TRAIN:** 2 Year Warranty  
**ELECTRONICS:** 2 Year Warranty  
**BATTERIES:** 1 Year Warranty



Optional color shrouds available.  
\*Upcharge applies.



**BUZZVI  
AROUND  
AL  
HD**



GB117H

Optional color shrouds available.  
\*Upcharge applies.



**BUZZVI  
AROUND  
AL  
HD**



GB147H

## **BUZZAROUND XLHD FEATURES:**

- Increased weight capacities up to 325 lbs.
- Ultra Bright LED headlight standard
- Charging port below the control panel
- 3 inch ground clearance
- Red & Blue Shroud panels included free
- Disassembles quickly and easily for transport
- 22Amp batteries standard
- Operating range up to 14 miles per charge

## **WARRANTIES:**

**FRAME:** Lifetime Warranty  
**DRIVE TRAIN:** 2 Year Warranty  
**ELECTRONICS:** 2 Year Warranty  
**BATTERIES:** 1 Year Warranty



**BUZZYLS  
AROUNDXLS™**  
Comfort Spring Suspension

Optional color  
shrouds available.  
\*Upcharge applies.



GB117S

**BUZZYLS  
AROUNDXLS™**  
Comfort Spring Suspension

Optional color  
shrouds available.  
\*Upcharge applies.



GB147S

## BUZZAROUND XLS FEATURES:

- 300 lb. weight capacity
- Ultra Bright LED headlight standard
- Charging port below the control panel
- 3 inch ground clearance
- Red and Blue Shroud panels included free
- Disassembles quickly and easily for transport
- NEW Comfort Spring Suspension for a smoother ride

## WARRANTIES:

**FRAME:** Lifetime Warranty  
**DRIVE TRAIN:** 2 Year Warranty  
**ELECTRONICS:** 2 Year Warranty  
**BATTERIES:** 1 Year Warranty





# BUZZAROUND XL<sup>SM</sup> HD

Optional color shrouds available.  
\*Upcharge applies.



GB117Z



# BUZZAROUND XL<sup>SM</sup> HD

Optional color shrouds available.  
\*Upcharge applies.



GB147Z

## BUZZAROUND XL<sup>SM</sup> HD FEATURES:

- Increased weight capacity up to 325 lbs.
- Ultra Bright LED headlight standard
- Charging port below the control panel
- 3 inch ground clearance
- Red and Blue Shroud panels included free
- Disassembles quickly and easily for transport
- 22 Amp batteries standard
- Operating range up to 14 miles per charge
- NEW Comfort Spring Suspension for a smoother ride

## WARRANTIES:

- FRAME:** Lifetime Warranty
- DRIVE TRAIN:** 2 Year Warranty
- ELECTRONICS:** 2 Year Warranty
- BATTERIES:** 1 Year Warranty



Optional color shrouds available.  
\*Upcharge applies.



GB107



**BUZZAROUND LT FEATURES:**

- 300 lb. weight capacity
- Increased foot room
- Color through chip resistant plastic (red only)
- Charger port in the tiller
- Quickly disassembles into five, lightweight pieces for easy transport
- Comes standard in matte red
- Pink, yellow, silver and camouflage color shrouds available
- Accepts new interchangeable color shrouds

**WARRANTIES:**

**FRAME:** Lifetime Warranty

**DRIVE TRAIN:** 2 Year Warranty

**ELECTRONICS:** 2 Year Warranty

**BATTERIES:** 1 Year Warranty



(Buzz XL model shown)

# LiteRider<sup>GOLDEN</sup>

## LARGE PORTABLE

The unique tubular design of the LiteRider scooters gets you wherever you need to go in a style all their own. The only scooters to come with two cargo baskets standard, the LiteRiders are easy to disassemble, and maneuver smoothly indoors and out.



LiteRider<sup>GOLDEN</sup>  
LiteRider GL111

### **WARRANTIES:**

**FRAME:** Lifetime Warranty

**DRIVE TRAIN:** 2 Year Warranty

**ELECTRONICS:** 2 Year Warranty

**BATTERIES:** 1 Year Warranty



LiteRider<sup>GOLDEN</sup>  
LiteRider GL141

### **GL111 & GL141 FEATURES:**

- 300 lb. weight capacity
- Wraparound Delta Tiller standard
- Easy to read control panel with backlit battery gauge
- Two storage baskets standard
- Comfortable 17"x16" stadium style seat standard
- Adjustable armrests
- Optional 20" x 18" stadium style seat
- 2-22 AH Batteries; increased operating range up to 12 miles
- Completely wireless disassembly
- 9" front and rear tires
- New Ultrabright LED adjustable headlight and rugged plastic basket standard!

# GOLDEN COMPANION

## FULL-SIZE LUXURY

**ALL NEW Rugged Steel Frame  
with Weight Capacity up to 400 lbs.!**

The Golden Companion scooters are your ticket to ride first class every day! Allow the luxurious, full-size Companions to enhance your life and take you to the places you've always had trouble getting to. Shop in the mall, attend sporting events and enjoy the great outdoors with the unmatched comfort of your favorite new Companion! All Companion models feature:



Companion GC340



*Companion  
GC340 shown  
with optional  
power lift seat.*

### **WARRANTIES:**

**FRAME:** Lifetime Warranty  
**DRIVE TRAIN:** 2 Year Warranty  
**ELECTRONICS:** 2 Year Warranty  
**BATTERIES:** 1 Year Warranty

### **COMPANIONS FEATURE:**

- All new design for easy disassembly
- Wraparound Delta Tiller standard
- Adjustable LED headlight
- EZ Reach tiller adjustment lever
- Auto lock tiller
- Easy to read control panel
- Easy access charging port
- Black, non-marking, low profile tires on sturdy steel rims
- LED battery gauge
- Comfortable, fully adjustable padded armrests
- Under floor mat troubleshooting guide
- Removable shopping basket with handle
- Four mobility transportation tie downs



Companion GC240



Companion GC440

# golden **Avenger™**



Avenger GA541  
(Freight Supplement applies.)

When you want rugged, heavy duty scooters to tackle the challenging terrain of the great outdoors, look no further than the Golden Avenger™ and Patriot™. Both feature oversized tires with sporty mag wheels, a full-size captain's seat for all day riding comfort, and an impressive operating range of up to 18 miles.

## **AVENGER FEATURES:**

- 500 lb. weight capacity
- Adjustable tiller for riding comfort
- Deluxe, reclining high back captain's seat
- Large, 13" front and rear tires
- Adjustable front headlight and four-way flashers
- Large front cargo basket standard
- Operating range up to 18 miles per charge

# Golden **Patriot™**

## **PATRIOT FEATURES:**

- Full front and rear lighting package with tail lights, headlights and flashers
- Four-wheel independent suspension for smoothly traversing outdoor terrain
- Rugged rubber mat and anti-skid material on the floor deck
- Convenient water bottle holders for two cold beverage bottles
- Full control panel featuring safety horn, flashers, lights, turn signals, speed knob and backlit battery gauge
- Black, pneumatic tires with brushed alloy wheels
- 13-Inch Pneumatic Tire and Rims

## **WARRANTIES:**

**FRAME:** 5 Year Warranty  
**DRIVE TRAIN:** 2 Year Warranty  
**ELECTRONICS:** 2 Year Warranty  
**BATTERIES:** 1 Year Warranty

## **WARRANTIES:**

**FRAME:** Lifetime Warranty  
**DRIVE TRAIN:** 2 Year Warranty  
**ELECTRONICS:** 2 Year Warranty  
**BATTERIES:** 1 Year Warranty



Patriot GR575  
(Freight Supplement applies.)

# Live the Golden Life

Test drive the industry's most comprehensive line of consumer power wheelchairs today!

## THE NEW PORTABLE LITERIDER ENVY'S DESIGN MAKES OUR POWER WHEELCHAIRS THE EASIEST TO TRANSPORT!

- Lightweight & Convenient
- Disassembles into 4 pieces
- 22 Amp Hour battery for all-day range
- Off board battery charging for on the go

### **WARRANTIES:**

**FRAME:** 5 Year Warranty  
**DRIVE TRAIN:** 13 Month Warranty  
**ELECTRONICS:** 13 Month Warranty  
**BATTERIES:** 12 Month Warranty

**ONLY  
36 lbs.!**





LiteRider GP162

# LiteRider<sup>®</sup> GOLDEN ENVY

Available in 6 Vibrant Colors with a Larger Foot Plate and New Shroud Design!

## Golden's RETAIL CASH POWER WHEELCHAIR features:

- 300 lb. weight capacity
- Operating range up to 15.5 miles
- Turning radius: 28.5"
- NEW premium 17" Seat Option

### *SOLAR FLARE ORANGE*

Shown with NEW high back premium seat option

### *ENVY GREEN*



### *CHERRY RED*



### *PARADISE BLUE*



### *SUNBURST YELLOW*



### *SATIN SILVER*



# Experience Constant Contact



Inclines and thresholds are no match for the Golden Compass Sport's true center wheel drive design with our patented springless articulating chassis!

This unique suspension system keeps all six wheels on the ground simultaneously, eliminating the rocking sensation when traversing inclines, giving both Compass™ models unparalleled curb and threshold climbing ability. No other power base can match the constant wheel to ground contact of a Compass!



golden  
**COMPASS™  
SPORT**

Model GP605

## **COMPASS SPORT FEATURES:**

- NEW All Black Seat
- NEW aluminum wheels provide increased durability
- Medicare Code: K0823
- 300 lb. Weight Capacity
- Turning Radius: 19.5"
- Batteries 2-U1
- Single Post Seat Mount
- 40 Amp Controller Standard

## **WARRANTIES:**

**FRAME:** Lifetime Warranty  
**DRIVE TRAIN:** 13 Month Warranty  
**ELECTRONICS:** 13 Month Warranty  
**BATTERIES:** 12 Month Warranty

## **THE COMPASS HD HAS THE MOST STABLE BASE IN ITS CLASS!**

Compass™ models provide superior forward stability making them the safest power wheelchairs to mount or dismount. The Compass line has fewer moving parts than a power base with springs for the best forward stability and no tipping!



golden  
**COMPASS™  
HD**

Model GP620

## **COMPASS HD FEATURES:**

- NEW All Black Seat
- Medicare Code: K0825
- 450 lb. Weight Capacity
- Turning Radius: 20"
- Batteries 2-NF22
- Four Post Seat Mount with Slide Seat Bracket
- 75 Amp Controller Standard



*You can't stand like this on any other power wheelchair base.*

## **WARRANTIES:**

**FRAME:** 5 Year Warranty  
**DRIVE TRAIN:** 13 Month Warranty  
**ELECTRONICS:** 13 Month Warranty  
**BATTERIES:** 12 Month Warranty



# Alanté SPORT

- Low-maintenance Color Through Shroud Design!
- Larger, Lightweight Composite Footplate!
- Generous High Back Captain's Seat



## Alanté SPORT

### ALANTÉ SPORT FEATURES:

- Low center of gravity
- Excellent maneuverability
- 28" turning radius
- 300 lb weight capacity
- Upgraded user friendly motor connections
- Auto shut-off charger system
- Electromagnetic braking system
- Padded flip up arms standard
- 9" driving wheels
- Larger 6" rear caster wheels

**ALL  
NEW!**

Model GP208



### WARRANTIES:

**FRAME:** 5 Year Warranty  
**DRIVE TRAIN:** 13 Month Warranty  
**ELECTRONICS:** 13 Month Warranty  
**BATTERIES:** 12 Month Warranty

# AVAILABLE ACCESSORIES

\*Only available on Companion models.



Cane Holder



Scooter Cover



Quad Cane Holder



Rear View Mirror



Crutch Holder



Rear Basket



Oxygen Tank Holder



Power Elevating Seat\*



Walker Holder



Safety Flag

*(Included with most models. Not available on GB119 & GB149.)*



Cup Holder



Pack'n'Go Armbag



Cell Phone Holder



Deluxe Pack'n'Go

*(Fits seats 20" wide or larger and high back Captain's Seat.)*



Optional BuzzAround Shroud Cover Colors:  
Pink, Silver, Yellow & Camouflage

*(Red & Blue come standard. GB107 standard in red only.)*



Back Pack

# SPECIFICATIONS\*

Item	Description	Length			Width			Front Tire	Rear Tire	Center Tire	Casters	Weight Capacity (lbs.)	Maximum Speed (mph)	Ground Clearance	Operating Range (miles)	Turning Radius	Batteries (2)	Weight w/ Batteries	Front Half Weight	Rear Half Weight	Seat Weight w/ Arms	Heaviest Piece	Base Weight
GB107	Buzz LT (3-wheel)	40.5"	21.5"	8"	8"	N/A	N/A	300	4	3"	8	37"	12AH	97	33	25	20	33	N/A				
GB117D	Buzz XL (3-wheel)	42.5"	21.5"	8"	8"	N/A	N/A	300	4	3"	8	37"	12AH	102	35	25	23	35	N/A				
GB117S	Buzz XLS (3-wheel)	42.5"	21.5"	8"	8"	N/A	N/A	300	4	3"	8	37"	12AH	105	35	28	23	35	N/A				
GB117H	Buzz XLHD (3-wheel)	42.5"	21.5"	8"	8"	N/A	N/A	325	4	3"	14	37"	22AH	112	35	25	23	35	N/A				
GB117Z	Buzz XLSHD (3-wheel)	42.5"	21.5"	8"	8"	N/A	N/A	325	4	3"	14	37"	22AH	115	35	28	23	35	N/A				
GB147D	Buzz XL (4-wheel)	41.5"	21.5"	8"	8"	N/A	N/A	300	4	3"	8	47"	12AH	110	43	25	23	43	N/A				
GB147S	Buzz XLS (4-wheel)	41.5"	21.5"	8"	8"	N/A	N/A	300	4	3"	8	47"	12AH	113	43	28	23	43	N/A				
GB147H	Buzz XLHD (4-wheel)	41.5"	21.5"	8"	8"	N/A	N/A	325	4	3"	14	47"	22AH	120	43	25	23	43	N/A				
GB147Z	Buzz XLSHD (4-wheel)	41.5"	21.5"	8"	8"	N/A	N/A	325	4	3"	14	47"	22AH	123	43	28	23	43	N/A				
GB118	Buzz EX (3-wheel)	45"	22"	9"	9"	N/A	N/A	350	5	4"	18	40.5"	U1	154	46	36	23	46	N/A				
GB148	Buzz EX (4-wheel)	44"	22"	9"	9"	N/A	N/A	350	5	4"	18	56"	U1	161	53	36	23	53	N/A				
GB119	Buzz LX (3-wheel)	45"	22"	9"	9"	N/A	N/A	375	5	4"	18	40.5"	U1	168	50	36	32	50	N/A				
GB149	Buzz LX (4-wheel)	44"	22"	9"	9"	N/A	N/A	375	5	4"	18	54"	U1	174	57	36	32	57	N/A				
GL111	LiteRider (3-wheel)	41.5"	22"	9"	9"	N/A	N/A	300	5	3.75"	12	39"	22AH	125	44	28	21	44	N/A				
GL141	LiteRider (4-wheel)	43.25"	22"	9"	9"	N/A	N/A	300	5	3.75"	12.5	48"	22AH	133	52	28	21	52	N/A				
GC240	Companion (3-wheel)	42.25"	23.5"	9"	10.75"	N/A	N/A	350	5.3	3.75"	19.9	37"	U1	163	45	48	26	48	N/A				
GC340	Companion (3-wheel)	45.5"	23.5"	9"	10.75"	N/A	N/A	400	5.3	4.5"	20.9	41"	U1	176	56	48	26	56	N/A				
GC440	Companion (4-wheel)	47.5"	23.5"	10.75"	10.75"	N/A	N/A	400	5.3	4.5"	18.9	56.5"	U1	194	74	48	26	74	N/A				
GA541	Avenger	57.5"	24.5"	13"	13"	N/A	N/A	500	7	4"	18 <sup>2</sup>	67"	NF22	265	68	68	55	68	N/A				
GR575	Patriot	57"	28"	13"	13"	N/A	N/A	400	8	6.75"	24	76.5"	75AH	302	N/A	N/A	44	152	152				
GP162	LiteRider Envy	36"	22.5"	N/A	9"	N/A	N/A	300	3.5	2.5"	15.5	28.5"	22AH	123	37	36	23	37	N/A				
GP208	Alanté Sport	41.5"	25"	9"	N/A	N/A	N/A	300	4.3	2.25"	15 <sup>2</sup>	28"	U1	146.5	N/A	N/A	37	66	66				
GP605	Compass Sport	39"	24"	N/A	N/A	10"	N/A	300	4	3.25"	24	19.5"	U1	185	N/A	N/A	46	92	92				
GP620	Compass HD	40"	25"	N/A	N/A	10"	N/A	450	4	3.5"	12 <sup>2</sup>	20"	NF22	226	N/A	N/A	51	101	101				

\* All weight shown in lbs. All products feature Freewheel Mode, Electronic Speed Control and Off-Board Charger.

NOTE: Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

<sup>1</sup>. Battery range at 200 lbs., and will vary due to rider weight, drive surface, terrain and battery type.  
<sup>2</sup>. Battery range at maximum weight capacity, and will vary due to rider weight, drive surface, terrain and battery type.

**Brakes:** All Golden Technologies scooters come standard with Dynamic/Electromechanical braking systems.  
**Rear wheel drive:** All Golden Technologies scooters come standard with a sealed transaxle, 24 volt motor.

# Choose from a Complete Line of Lifestyle Solution Products by Golden



## SCOOTERS

- BuzzAround Series
- Lite Rider Series
- Companion Series
- Heavy Duty Series



## POWER WHEELCHAIRS

- Compass Sport
- Compass HD
- LiteRider Envy
- Alanté Sport



## POWER LIFT RECLINERS

- Value Series
- Signature Series
- Comforter Series
- Comforter Wide Series
- MaxiComforter Series
- **PATENTED** Twilight Technology Positioning Features!

Visit Your Local Authorized Golden Dealer.



[www.goldentech.com](http://www.goldentech.com)