

PRIDE®

MOBILITY
SCOOTERS

Get Out There!



Don't let a change in mobility affect your lifestyle. Live your best with Pride® Mobility Scooters and enjoy sleek, sporty style in a high-performance package.

With an array of models and colors to choose from, Pride Mobility Scooters make it easier for you to get out there and enjoy life.



What's Inside

Travel Mobility	3-7	Victory® 9.....	9
Jazzy® Zero Turn	3	Victory® 10.....	9
Jazzy® Passport	3	Victory® 10 LX w/CTS	10
Go Chair®	4	Victory® LX Sport.....	10
Revo™ 2.0	4	Full-Size Luxury	11-12
Go-Go® ES 2.....	5	Pursuit®	11
Go-Go® Ultra X.....	5	Pursuit® XL.....	11
Go-Go Elite Traveller®.....	5	Raptor	12
Go-Go Elite Traveller® Plus.....	6	Heavy Duty	12
Go-Go® LX with CTS.....	6	Maxima	12
Go-Go® Sport	6	Accessories	13
Go-Go® Folding Scooter.....	7	Innovative Technology	14
Mid-Size Luxury	8-10	Specifications	15-18
Victory® 9.2.....	8	Travel Mobility	15-16
Victory® 10.2.....	8	Full/Mid-Size and Heavy Duty	17-18



Choosing the right Pride® Mobility Scooter

When choosing a Pride® Mobility Scooter, specific requirements and conditions can affect your process.

Indoor maneuverability

- Travel mobility to mid-size luxury scooters offer enhanced maneuverability

Stability

- 4-wheel scooters provide increased stability

Indoor and outdoor versatility

- Mid-size and full-size luxury scooters are great for both indoor and outdoor use

Outdoor capabilities

- Mid-size and full-size luxury scooters offer excellent outdoor capabilities

Increased weight capacity

- Mid-size and full-size luxury scooters provide increased weight capacity

Easy portability

- Travel mobility scooters are easy to disassemble, transport and store

Speed and performance

- Mid-size and full-size luxury scooters offer increased speed and performance

Performance and Convenience on the Go



Jazzy® ZERO TURN

4-WHEEL SCOOTER

Features

- 325 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into five Patent Pending, super lightweight pieces for convenient transport
- iTurn Technology™ provides 4-wheel stability with 3-wheel maneuverability
- 38.25" turning radius navigates tight corners and small spaces effortlessly
- Dual motors for added traction in rough terrain
- Standard under seat storage
- Delta tiller with wraparound handles
- Charger port in tiller
- LED lighting package
- USB charger built into the tiller for convenient charging of smart phone and portable devices

Featuring intelligent-turning



Color Choices

INTERCHANGEABLE SHROUDS⁷



Fire Opal Green Moonstone Onyx Black Pewter Sapphire Blue Pearl White



Jazzy® PASSPORT

Features

- 250 lb. weight capacity
- Maximum speed up to 3.6 mph
- Up to 9.6 miles per charge range at 250 lbs., 11 miles per charge range at 200 lbs.
- A compact and lightweight design makes traveling a breeze
- Back seat pocket for extra storage while on-the-go
- Folds in just a few simple steps for effortless transport
- Includes convenient, under seat mesh storage bag, mesh cup holder and 40" lap belt
- Large foot platform to accommodate the user's comfort
- Specialty foam seat design for extreme comfort

Now available with optional airline-approved travel battery

Note: Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Jazzy® Passport and its lithium-ion battery are permitted.

Color Choice



Pewter



Feather-touch disassembly

GO CHAIR®

Features

- #1 for portability
- 300 lb. weight capacity
- One hand feather-touch disassembly
- Up to 13.2 miles per charge range at 200 lbs., up to 8.7 miles per charge range at 300 lbs.
- 18" x 17", 5-point swivel seat with latch release
- Faster speed up to 3.7 mph
- Compact turning radius 25.5"
- Dual swing-away, under seat storage bins
- Height and width adjustable armrests
- Heaviest piece is only 36 lbs.

Color Choices

INTERCHANGEABLE SHROUDS⁷



Ruby Red Sapphire Blue Rose Quartz Citrine Yellow Amber Orange Pearl White



3-wheel



4-wheel

REVO™ 2.0

3 AND 4-WHEEL SCOOTER

Features

- 375 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into five super lightweight pieces for convenient transport
- CTS Suspension (Comfort-Trac front and rear suspension)
- Up to 12 miles per charge range at 375 lbs., 17.4 miles per charge range at 200 lbs. (3W); Up to 13 miles per charge range at 375 lbs., 19.8 miles per charge range at 200 lbs. (4W)
- Standard under seat storage
- Delta tiller with wraparound handles
- Charger port in tiller
- LED headlight
- Easily charge the batteries on or off the unit
- USB mobile device charging port

Color Choices



Grey Street True Blue

GO-GO[®] es2

3-WHEEL SCOOTER



Features

- 250 lb. weight capacity
- Two speed settings high and low (indoor/outdoor), maximum speed up to 4 mph
- Up to 6.8 miles per charge range
- Lightweight, one-piece frame for added stability and easy portability

Color Choices



Black

GO-GO[®] ULTRAX

3 AND 4-WHEEL SCOOTER



3-wheel



4-wheel

Features

- 260 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into five super lightweight pieces for convenient transport
- Maximum speed up to 4 mph
- Up to 6.9 miles per charge range (3W), 7.2 miles per charge range (4W)

Color Choice



Red

GO-GO[®] ELITE TRAVELLER[®]

3 & 4-WHEEL SCOOTER



3-wheel



4-wheel

Features

- 300 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into five super lightweight pieces for convenient transport
- Maximum speed up to 4 mph
- Up to 12.4 miles per charge range (3W), 12 miles per charge range (4W)
- Charger port in tiller

Color Choices



Red



Blue

INTERCHANGEABLE SHROUDS⁷



3-wheel



4-wheel

GOOGO ELITE TRAVELLER[®] PLUS

3 AND 4-WHEEL SCOOTER

Features

- 300 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into five super lightweight pieces for convenient transport
- Maximum speed up to 4.5 mph
- Up to 14.5 miles per charge range
- Delta tiller with wraparound handles
- Charger port in tiller
- LED headlight

Color Choices

INTERCHANGEABLE SHROUDS⁷



Red Blue



3-wheel



4-wheel

GOOGO[®] LX WITH CTS SUSPENSION

3 AND 4-WHEEL SCOOTER

Features

- 300 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into five super lightweight pieces for convenient transport
- CTS Suspension (front and rear suspension)
- Maximum speed up to 4.4 mph
- Up to 11.7 miles per charge range (3W), 12.6 miles per charge range (4W)
- Delta tiller with wraparound handles
- Charger port in tiller
- LED headlight

Color Choices

INTERCHANGEABLE SHROUDS⁷



Red Blue **Optional** Orion Yellow Onyx Black Pearl Pink Sunburst Orange



3-wheel



4-wheel

GOOGO[®] SPORT

3 AND 4-WHEEL SCOOTER

Features

- 325 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into five super lightweight pieces for convenient transport
- Maximum speed up to 4.7 mph
- Up to 15.9 miles per charge range (3W), 14.5 miles per charge range (4W)
- Delta tiller with wraparound handles
- Charger port in tiller
- LED headlight

Color Choices

INTERCHANGEABLE SHROUDS⁷



Red Blue

GO-GO[®] FOLDING SCOOTER

4-WHEEL SCOOTER



Features

- 250 lb. weight capacity
- Maximum speed up to 4 mph
- Up to 13 miles per charge range (lithium-ion), 9.3 miles per charge range (sealed lead acid)
- Height-adjustable tiller
- Lighting package
- Optional lithium-ion battery
- Optional armrests

Color Choice



White



Note: Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Go-Go[®] Folding Scooter and its lithium-ion battery is permitted.



The Ultimate in Style and Performance®

VICTORY® SERIES

VICTORY® 9.2

3 AND 4-WHEEL SCOOTER



3-wheel



4-wheel

Features

- 300 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand
- Maximum speeds up to 5.6 mph
- Up to 17.5 miles per charge
- Delta tiller with wraparound handles
- Frosted LED headlights and low-mounted LED pathway light brightens your route
- Under-tiller puddle light improves visibility to USB charging port and tiller-mounted storage
- Integrated USB mobile device charging port in tiller
- Integrated storage/cup holder in tiller

Color Choices



Candy Apple Red Ocean Blue

VICTORY® 10.2

3 AND 4-WHEEL SCOOTER



3-wheel



4-wheel

Features

- 400 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand
- Maximum speeds up to 5.6 mph
- Up to 13.5 miles per charge
- Delta tiller with wraparound handles
- Frosted LED headlights and low-mounted LED pathway light brightens your route
- Under-tiller puddle light improves visibility to USB charging port and tiller-mounted storage
- Integrated USB mobile device charging port in tiller
- Integrated storage/cup holder in tiller

Color Choices



Candy Apple Red Ocean Blue



3-wheel



4-wheel

VICTORY[®] 9

3 AND 4-WHEEL SCOOTER

Features

- 300 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand
- Maximum speed up to 5.3 mph
- Up to 13 miles per charge range
- Delta tiller with wraparound handles
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Easy access tiller-mounted charger port

Color Choices



Candy Apple Red Viper Blue



3-wheel



4-wheel

VICTORY[®] 10

3 AND 4-WHEEL SCOOTER

Features

- 400 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand
- Maximum speed up to 5.3 mph
- Up to 16 miles per charge range
- Delta tiller with wraparound handles
- Long lasting, high intensity LED headlight offers optimal pathway illumination
- Easy access tiller-mounted charger port

Color Choices



Candy Apple Red Viper Blue

VICTORY[®] 10 LX WITH CTS SUSPENSION

4-WHEEL SCOOTER



Features

- 400 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand
- CTS Suspension (Comfort-Trac front and rear suspension)
- Maximum speed up to 5.6 mph
- Up to 16 miles per charge range
- Delta tiller with wraparound handles
- Black vinyl reclining high-back seat with sliders, height adjustable headrest, flip-up, height and width adjustable armrests
- Full LED lighting package that includes curb and rear lights and directional signals
- Easy access tiller-mounted charger port

Color Choices



Candy Apple Red True Blue

VICTORY[®] LX^{SPORT}

4-WHEEL SCOOTER



Features

- 400 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand
- CTS Suspension (Comfort-Trac front and rear suspension)
- Maximum speed up to 8 mph
- Up to 19.7 miles per charge range (200 lbs.),
Up to 13 miles per charge range (400 lbs.)
- Delta tiller with wraparound handles
- Black vinyl reclining high-back seat with sliders, height adjustable headrest, flip-up, height and width adjustable armrests
- Full LED lighting package includes curb and rear lights, directional signals and courtesy light under tiller
- Easy access tiller-mounted charger port

Color Choices



Candy Apple Red True Blue

Outstanding Outdoor Performance with Luxurious Features



4-WHEEL SCOOTER



Features

- 400 lb. weight capacity
- Infinitely-adjustable tiller angle with conveniently located handle
- Maximum speed up to 8 mph
- Front and rear suspension
- Up to 14 miles per charge range (NF-22), 19 miles per charge range (G24)
- Feather-touch disassembly
- Delta tiller with wraparound handles
- Full lighting package that includes directional signals and lower headlight
- Reclining high-back seat with sliders and headrest

Color Choices



Candy Apple Red Silver



4-WHEEL SCOOTER



Features

- 400 lb. weight capacity
- Infinitely-adjustable tiller angle with conveniently located handle
- Maximum speed up to 7.5 mph
- Up to 30 miles per charge range (G27), 22 miles per charge range (G24)
- Front and rear suspension
- Delta tiller with wraparound handles
- Deluxe reclining high-back seat with sliders and headrest
- 13" low-profile, pneumatic tires for excellent outdoor performance

Color Choices



Candy Apple Red Silver

FULL-SIZE LUXURY

RAPTOR

3-WHEEL SCOOTER



**This product is not a medical device and is not intended to assist, treat, diagnose, or alleviate any medical condition or disability.*

Features

- 400 lb. weight capacity
- Front and rear suspension
- Maximum speed up to 14 mph
- Full LED lighting package with directional signals and brake lights that offer full visibility for increased safety
- Hand brake with regenerative braking system
- Secure, lockable under seat storage compartment with seat support rod
- Easy-to-read LCD console which shows estimated speed, mileage and battery display
- Limited recline fold-flat seat with manual adjustment and angle/width adjustable armrests
- Sleek and solid design for enhanced stability
- Front storage compartment located on the tiller
- Tubeless pneumatic tires for a smooth ride

Color Choice



Silver

HEAVY-DUTY

Durable Performance at an Exceptional Value

MAXIMA

3 AND 4-WHEEL SCOOTER



3-wheel

Features

- 500 lb. weight capacity
- Ultra heavy-duty drivetrain with 110 amp controller
- Maximum speed up to 5.3 mph
- Easily adjustable tiller
- Up to 18 miles per charge range (3W, NF-22), 9.6 miles per charge range (power elevating seat, U-1), 16 miles per charge range (4W)
- Durable and stylish rear polyurethane bumper
- User-friendly console

Color Choices



Candy Apple Red



Viper Blue



4-wheel

Accessories

A wide selection of Pride Mobility Scooter accessories are available to personalize your Pride Mobility Scooter to better meet your needs and make it uniquely yours.



Cane clip attachment



Cup holder



Single cane holder



Safety flag



Double cane holder



Under seat storage bins



Rear basket



XLR USB charger



Oxygen holder



RAM® X-Grip®
cell phone holder



Wishbone crutch holder



Additional 12-amp
battery box



Lapbelt



Saddlebag



Walker holder



Additional 18-amp
battery box



Rearview mirror



Weather breaker



Weather cover

Military Patches Available



DISCLAIMER: Not all accessories are available for all products. Please speak with your Provider for details and availability.

Innovative Technology



Our patent-pending, intelligent-turning **iTurn Technology™** ensures effortless navigation in small spaces and around corners. Experience the stability of a 4-wheel scooter with the tight turning radius of a 3-wheel scooter. iTurn Technology allows you to make a U-turn in a small space without ever having to put the scooter in reverse! From elevators to narrow grocery store aisles, iTurn Technology can handle the tightest spaces!



Comfort-Trac Suspension (CTS)



Front Suspension



Rear Suspension

Comfort-Trac Suspension (CTS) is the next generation in mobility scooter ride comfort. Pride's **industry-first** technology incorporates advanced suspension components for improved absorption over rough surfaces, and a smoother ride experience over various terrains.

Feather-Touch Disassembly



Scooter/travel mobility disassembly



Go Chair® disassembly

Most Go-Go® Travel Mobility and Pride® mid-size scooters feature **Feather-Touch**—one-hand disassembly for hassle-free transport and storage.

Travel Mobility Specifications

MODEL NAME	Jazzy® Zero Turn	Revo™ 2.0	Jazzy® Passport	Go Chair®	Go-Go® ES 2
MODEL NUMBER	JAZZYZTS	S66 (3W) S67 (4W)	FJPT1700	Go Chair	S81
FRONT TIRES	2" x 7" solid	3" x 8" solid (3W) 3" x 9" solid (4W)	8" solid (caster)	6" solid (caster)	2" x 7.5" solid
REAR TIRES	2.5" x 9" solid	3" x 9" solid	12" (Pneumatic drive)	9" solid (drive)	3" x 7.5" solid
MAXIMUM SPEED ²	Up to 4 MPH	Up to 5 MPH	Up to 3.6 MPH	Up to 3.7 MPH	Up to 4 MPH
GROUND CLEARANCE ^{4,11,14}	1.4" at motor 1.5" at anti-tip 3.5" at side deck	1.25" at transaxle 3.25" at mid-frame (3W) 3.5" at mid-frame (4W)	3" at motor	1.8" at motor 1.7" at anti-tip 2.7" at foot plate bracket	1.5" at motor 0.75" at anti-tip
3-W TURNING RADIUS ⁴	N/A	40"	N/A	N/A	31"
4-W TURNING RADIUS ⁴	38.25"	52"	24.25"	25.5" with foot platform	N/A
3-W LENGTH/WIDTH ^{4,10}	N/A	46"/20.75"	N/A	N/A	36.75"/19.5"
4-W LENGTH/WIDTH ^{4,10}	45.75"/20"	47"/20.75"	40"/23.5"	33"/25.5"	N/A
SEAT-TO-GROUND HEIGHT RANGE	21.5"-23.5"	21"-23.5"	19"	22"-23"	19"-20"
SEAT TYPE (STANDARD)	Compact, foldable (17"x17")	Compact, foldable (20"x18")	(17.5"x16.5")	Compact, foldable (18"x17")	Compact, foldable (17"x15")
(3/4-WHEEL) WEIGHT ^{3,4,12}	120 lbs.	120 lbs./132 lbs.	54 lbs.	67 lbs.	58 lbs.
SEAT-TO-DECK HEIGHT RANGE	16"-18"	15"-18"	14.5"	16.25"-17.25"	16.25"-17.25"
BATTERY WEIGHT ⁹	13.5 lbs. each	25 lbs. each	7 lbs. each	14 lbs. each	9 lbs. each
BATTERY REQUIREMENTS ⁶	(2) 22AH	(2) U1	(1) 24V, 18AH (Lithium-ion)	(2) 18AH	(2) 12AH
PER CHARGE RANGE (MAX.) ^{2,8}	Up to 13.4 miles at 200 lbs., Up to 11.6 miles at 325 lbs.	3W: Up to 12 miles at 375 lbs., Up to 13.2 miles at 200 lbs., 4W: Up to 12.4 miles at 375 lbs., Up to 17.8 miles at 200 lbs.	Up to 9.6 miles at 250 lbs., Up to 11 miles at 200 lbs. Up to 5.7 miles (Optional travel battery)	Up to 8.7 miles at 375 lbs., Up to 11 miles at 200 lbs.	Up to 6.8 miles
BATTERY CHARGER	Off-board, 3.5A	Off-board, 3.5A	Off-board, 2A	Off-board, 2A	Off-board, 2A
SUSPENSION	Front & rear	Front & rear	Limited	Limited	N/A
MAXIMUM WEIGHT CAPACITY	325 lbs.	375 lbs.	250 lbs.	300 lbs.	250 lbs.
STANDARD BODY COLORS	Fire Opal, Green Moonstone, Onyx Black, Pewter, Pearl White, Sapphire Blue	Grey Street, True Blue	Pewter	Ruby Red, Sapphire Blue, Citrine Yellow, Amber Orange, Pearl White, Rose Quartz	Black
WEIGHT OF HEAVIEST PIECE ⁴	48 lbs. (rear section)	48 lbs. (3-wheel) (rear section) 57 lbs. (4-wheel) (front section)	N/A	36 lbs. (rear section)	41 lbs. (3W) (scooter base)
HCPCS	K0899*	K0899*	K0899*	K0899*	K0800 — Group 1 standard (3W)

*K0899 — power mobility device, not coded by dme pdac or does not meet criteria.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement.

DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.

Go-Go® Ultra X	Go-Go Elite Traveller®	Go-Go Elite Traveller® Plus	Go-Go® Sport	Go-Go® LX w/ CTS Suspension	Go-Go® Folding Scooter
SC40X (3W) SC44X (4W)	SC40E (3W) SC44E (4W)	SC53 (3W) SC54 (4W)	S73 (3W) S74 (4W)	SC50LX (3W) SC54LX (4W)	S19WH1001 (SLA) S19WH1002 (Lithium-ion)
2" x 8" solid (3W) 2" x 7" solid (4W)	2" x 8" solid (3W) 2" x 7" solid (4W)	2" x 8" solid	3" x 8" solid (3W) 3" x 9" solid (4W)	2" x 8" solid (3W) 2" x 7" solid (4W)	1.5" x 7" solid
2.5" x 8" solid	2.5" x 8" solid	3" x 9" solid	3" x 9" solid	2.5" x 8" solid	2" x 8" solid
Up to 4 MPH	Up to 4 MPH	Up to 4.5 MPH	Up to 4.7 MPH	Up to 4.4 MPH	Up to 4 MPH
1.3" at transaxle 0.67" at anti-tip 3.1" at deck (3W) 3" at deck (4W)	1.38" 0.75" at anti-tip 2.25" at mid-frame	1.7" at transaxle (3W) 0.83" at anti-tip (3W) 2.6" at side deck (3W)	1.5" at motor 0.75" at anti-tip 1.25" at rear frame lock (3W) 3" at controller (4W)	1.5" at motor 0.75" at anti-tip 2.5" at frame	2.3"
33"	33"	34.5"	36"	33"	32"
44"	44"	45.75"	50.5"	44.75"	N/A
37"/19.5"	37"/19.5"	40"/21.25"	42"/21.25"	39.75"/20.5"	37"/19"
39.5"/19.5"	39.5"/19.5"	40.25"/21.25"	43"/21.25"	41"/20.5"	N/A
18.25"-20.25"	18.25"-20.25"	20.5"-22"	20.5"-22"	20.75"-22.25"	20.5"
Compact, foldable (17"x17")	Compact, foldable (17"x17")	Compact, foldable (17"x17")	Compact, foldable (17"x17")	Compact, foldable (17"x17")	Compact, foldable (16"x13")
76 lbs./80 lbs.	76 lbs./80 lbs.	86 lbs./87 lbs.	93 lbs./99 lbs.	85 lbs./87 lbs.	51 lbs.
15.75"-17.25"	15.75"-17.25"	16.25"-17.75"	17"-18"	16.25"-17.25"	13"
9 lbs. each ⁵	9 lbs. each ⁵ / 14 lbs. each ⁷	14 lbs. each ⁷	14 lbs. each ⁷	9 lbs. each ⁵ / 14 lbs. each ⁷	9 lbs. each (Sealed Lead Acid / 5 lbs. each (Lithium-ion))
(2) 12AH	(2) 12AH or (2) 18AH	(2) 18AH	(2) 18AH	(2) 12AH or (2) 18AH	(2) 12AH or (1) 12AH (Lithium-ion)
3W: Up to 6.9 miles 4W: Up to 7.2 miles	3W: Up to 8.5 miles (12AH) Up to 12.4 miles (18AH) 4W: Up to 8.2 miles (12AH) Up to 12 miles (18AH)	3W: Up to 14.5 miles 4W: Up to 14.5 miles	3W: Up to 15.9 miles 4W: Up to 14.5 miles	3W: Up to 8 miles (12AH) Up to 11.7 miles (18AH) 4W: Up to 8.6 miles (12AH) Up to 12.6 miles (18AH)	Up to 13 miles (Lithium-Ion) Up to 9.3 miles (Sealed Lead Acid)
Off-board, 2A	Off-board, 2A (12AH) Off-board, 3.5A (18AH)	Off-board, 3.5A	Off-board, 3.5A	Off-board, 2A (12AH) Off-board, 3.5 A (18AH)	Off-board, 2A (sealed lead acid) Off-board, 3.5A (Lithium-ion)
N/A	N	N	N	Front & rear	N
260 lbs.	300 lbs.	300 lbs.	325 lbs.	300 lbs.	250 lbs.
Red	Red, Blue ¹	Red, Blue ¹	Red, Blue ¹	Red, Blue ¹	White
29 lbs. (3-wheel) (front section) 34 lbs. (4-wheel) (front section)	29 lbs. (3-wheel) (front section) 34 lbs. (4-wheel) (front section)	32 lbs. (3-wheel) (front section) 34 lbs. (4-wheel) (front section)	35 lbs. (3-wheel) (front section) 44 lbs. (4-wheel) (front section)	34 lbs. (3-wheel) (front section) 37 lbs. (4-wheel) (front section)	51 lbs. (unit does not disassemble — unit weight without battery pack)
K0899*	K0800 — Group 1 standard	K0806 — Group 2 standard ¹³	K0801 — Group 1 heavy duty	K0899*	K0899*

1. All shroud panels included.

2. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.

3. Battery weight not included.

4. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.

5. Battery pack is 18 lbs. with 12AH batteries.

6. Deep cycle (AGM or gel cell type recommended).

7. 18 AH battery pack weighs 29.5 lbs.

8. Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards.

Results derived from theoretical calculations used on battery specifications and drive system performance. Test conducted at maximum weight capacity.

9. Battery weight may vary based on manufacturer.

10. Width of base (not including seat).

11. Does not account for anti-tip clearance.

12. Includes standard seat weight.

13. Must be billed to Medicare as either a beneficiary requested upgrade with a properly executed ABN, or as an upgrade at no additional charge to the beneficiary.

14. At midframe.

15. Light package.

16. This product is not a medical device and is not intended to assist, treat, diagnose, or alleviate any medical condition or disability.

Mid/Full-Size, Heavy-Duty

MODEL NAME	Victory® 9	Victory® 10	Victory® 9.2	Victory® 10.2	Victory 10® LX w/ CTS Suspension ¹⁵
MODEL NUMBER	SC609 (3W) SC709 (4W)	SC610 (3W) SC710 (4W)	S6092 (3W) S7092 (4W)	S6102 (3W) S7102 (4W)	S710LX
FRONT TIRES	3W: 3" x 8" solid 4W: 3" x 9" solid	3.5" x 10" solid	3W: 3" x 8" solid 4W: 3" x 9" solid	3.5" x 10" solid	3.5" x 10" solid
REAR TIRES	3" x 9" solid	3.5" x 10" solid	3" x 9" solid	3.5" x 10" solid	3.5" x 10" solid
MAXIMUM SPEED ²	Up to 5.3 MPH	Up to 5.3 MPH	Up to 5.6 MPH	Up to 5.2 MPH	Up to 5.6 MPH
GROUND CLEARANCE ^{4,11, 14}	1.8" at transaxle 1.4" at anti-tip 2.7" at side deck (3W) 2.1" at motor bracket bolt (4W)	2.2" at motor (3W) 1.26" at anti-tip (3W) 3.4" at side deck (3W) 2" at motor (4W) 1.5" at anti-tip (4W) 3.1 at side deck (4W)	1.8" at transaxle 1.4" at anti-tip 2.7" at side deck	2.4" at transaxle 2.1" at anti-tip 3.0" at side deck	2.4" at transaxle 0.9" at anti-tip 2.9" at front frame
3-W TURNING RADIUS ⁴	40.25"	42"	40.25"	42"	N/A
4-W TURNING RADIUS ⁴	51.5"	54"	52"	52.25"	60"
3-W LENGTH/WIDTH ^{4,10}	43.5"/22.25"	45.25"/22.25"	43.5"/22.25"	45.25"/22.25"	N/A
4-W LENGTH/WIDTH ^{4,10}	45.5"/22.5"	47"/22.25"	46"/22.25"	47"/22.25"	47"/22.25"
SEAT-TO-GROUND HEIGHT RANGE	19.25"-20.75"	20.5"-22"	21.5"-24.5"	22"-25"	23.5"- 25.5"
SEAT TYPE (STANDARD)	Compact, foldable (18"x17")	Compact, foldable (18"x17")	Compact, foldable (18"x17")	Compact, foldable (18"x17")	Limited reclining H/B w/sliders & headrest (18"x18")
(3/4-WHEEL) WEIGHT ^{3,4,12}	117 lbs./131 lbs.	133 lbs./137 lbs.	125.75 lbs./136 lbs.	130.25 lbs./139 lbs.	179 lbs.
SEAT-TO-DECK HEIGHT RANGE	17.25"-18.75"	17.5"-19"	15.25"-18.25"	15.25" – 18.25"	16.5"-18.5"
BATTERY WEIGHT ⁹	24 lbs. each	24 lbs. each/ 32 lbs. each	24.5 lbs. each	24 lbs. each 32 lbs. each	32 lbs. each
BATTERY REQUIREMENTS ⁶	(2) U1	(2) U1 or (2) 40AH	(2) U1	(2) U1 or (2) 40AH	(2) 40AH
PER CHARGE RANGE (MAX.) ^{2,8}	Up to 13 miles	Up to 16 miles	3W: Up to 17.5 miles 4W: Up to 16.5 miles	3W/4W: Up to 11 miles (U-1) 3W: Up to 13 miles (40AH) 4W: Up to 13.5 miles (40AH)	Up to 16 miles
BATTERY CHARGER	Off-board, 3.5A (U-1)	Off-board, 3.5A (U-1)/ Off-board, 5A (40AH)	Off-board, 3.5A (U-1)/ Off-board, 5A (40AH)	Off-board, 3.5A (U-1)/ Off-board, 5A (40AH)	Off-board, 5A (40AH)
SUSPENSION	N/A	N/A	N/A	N/A	Front & rear
MAXIMUM WEIGHT CAPACITY	300 lbs.	400 lbs.	300 lbs.	400 lbs.	400 lbs.
STANDARD BODY COLORS	Candy Apple Red, Viper blue	Candy Apple Red, Viper blue	Candy Apple Red, Ocean Blue	Candy Apple Red, Ocean Blue	Candy Apple Red, True Blue
WEIGHT OF HEAVIEST PIECE ⁴	48 lbs. (3-wheel) (front section) 59 lbs. (4-wheel) (front section)	54 lbs. (3-wheel) (front section) 61 lbs. (4-wheel) (front section)	51.75 lbs. (3-wheel) (front section) 61.5 lbs. (4-wheel) (front section)	56.5 lbs. (3-wheel) (front section) 64.5 lbs. (4-wheel) (front section)	78 lbs. (4-wheel) (front section)
HCPCS	K0800 — Group 1 standard	K0801 — Group 1 heavy duty	K0800 — Group 1 standard	K0801 — Group 1 heavy duty	K0899*

*K0899 — Power mobility device, not coded by dme pdac or does not meet criteria.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement.

DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %.

Specifications

Victory® LX Sport ¹⁵	Maxima	Pursuit®15	Pursuit® XL15	Raptor ¹⁵
S710LXW	SC900 (3W) SC940 (4W)	SC713	SC714	R3-1700
3.5" x 10" Solid	3" x 10" Solid	3" x 13" Solid	13" Pneumatic	3" x 16" Pneumatic
3.5" x 10" Solid	4" x 10" Solid	13" Solid	13" Pneumatic	3" x 16" Pneumatic
Up to 8 MPH	Up to 5.3 MPH	Up to 8 MPH	Up to 7.5 MPH	Up to 14 MPH
2.2" at motor 3.15" at rear mid-frame	2.5" at motor 2" at anti-tip (3w) 4.5" at side deck (3w) 1.5" at anti-tip (4w) 3.25" (4w)	3.7" at motor 1.8" at anti-tip 4.5" at side deck	2.75" at motor 5.25" at mid-frame	3.5" at center of frame 1.37" at anti-tip 6.18"
N/A	50.5"	N/A	N/A	60"
61.25"	54"	67.75" Left 66.75" Right	82.5"	N/A
N/A	46"/25.5"	N/A	N/A	59"/24.75"
47"/25.62"	47"/25.5"	53.25"/25.25"	56"/27.25"	N/A
23"-25"	22.5"-24.5"	24.75"-26.25"	26.25"-27.75"	28"
Limited reclining H/B w/sliders & headrest (18"x18")	Deluxe contoured Medium-Back (22"x18")	Limited reclining H/B w/sliders & headrest (18"x18")	Limited reclining H/B pillow top w/sliders & headrest (19.5"x16")	High Back (20" x 18") (100° limited recline)
160 lbs.	171 lbs./176 lbs.	225 lbs.	254 lbs.	208 lbs.
16.5"-18.5"	3W: 16"-18" 4W: 16"-18"	17"-18.75"	21"-22.5"	19"
30.2 lbs. each	38 lbs. each (NF-22 or 55AH) 24 lbs. each (U-1)	38 lbs. each (NF-22 or 55AH) 53 lbs. each (75AH)	53 lbs. each (75AH) 60 lbs. each (100AH)	53 lbs. each
(2) 50AH	(2) NF-22 (2) U-1 (power elevating seat)	(2) NF-22 or (2) Group 24	(2) Group 27 or (2) Group 24	(2) 75AH
Up to 19.7 miles (200 lbs.) Up to 13 miles (400 lbs.)	3W: Up to 18 miles, U-1 (power elevating seat) Up to 9.6 miles 4W: Up to 16 miles	G24: Up to 19 miles NF-22: Up to 14 miles	G27: Up to 30 miles G24: Up to 22 miles	19 miles (up to 400 lbs. and full speed) 31 miles (up to 220 lbs. and at low speed)
Off-board, 5A	Off-board, 5A	Off-board, 5A (NF-22) Off-board, 8A (Group 24)	Off-board, 8A	Off-board, 8A
Front & rear	N/A	Front & rear	Front & rear	Front & rear
400 lbs.	500 lbs.	400 lbs.	400 lbs.	400 lbs.
Candy Apple Red, True Blue	Candy Apple Red, Viper Blue	Candy Apple Red, Silver	Candy Apple Red, Silver	Silver
77 lbs. (4W) (front section)	76 lbs. (3W) (rear section) 76 lbs. (4W) (rear section)	95 lbs. (front section)	205 lbs. (main frame)	N/A (unit does not disassemble)
K0899*	K0802 — Group 1 very heavy duty (3W) K0808 — Group 2 ¹⁴ very heavy duty (4W)	K0899*	K0899*	NM ¹⁶

1. All shroud panels included.

2. Varies with user weight, terrain type, battery amp- hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.

3. Battery weight not included.

4. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.

5. Battery pack is 18 lbs. with 12AH batteries.

6. Deep cycle (AGM or gel cell type recommended).

7. 18 AH battery pack weighs 29.5 lbs.

8. Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations used on battery specifications and drive system performance. Test conducted at maximum weight capacity.

9. Battery weight may vary based on manufacturer.

10. Width of base (not including seat).

11. Does not account for anti-tip clearance.

12. Includes standard seat weight.

13. Must be billed to Medicare as either a beneficiary requested upgrade with a properly executed ABN, or as an upgrade at no additional charge to the beneficiary.

14. At midframe.

15. Light package.

16. This product is not a medical device and is not intended to assist, treat, diagnose, or alleviate any medical condition or disability.

PRIDE[®]

MOBILITY
SCOOTERS

pridemobility.com |    

