

SIGVARIS

The Worldwide Leader
in Compression Therapy



COMPRESSION SOLUTIONS 2018

GRADUATED COMPRESSION SOCKS & HOSIERY
FOR CHRONIC VENOUS DISEASE

HOW TO ORDER

HOW TO ORDER MEDICAL

Series Number 23 = COTTON <i>(See individual pages for series number)</i>	Compression Level 1 = 15–20mmHg 2 = 20–30mmHg* 3 = 30–40mmHg 4 = 40–50mmHg** 5 = 50–60mmHg**	Product Style C = Knee-high or Calf <i>(See individual pages for available styles)</i>	Size ML = Medium Long <i>(See individual pages for sizing)</i>	Gender M = Men W = Women O = Open toe	Color 99 = Black <i>(See individual pages for available colors)</i>
23 2 C ML M 99					
* 18–25mmHg available in DIABETIC COMPRESSION SOCKS only					
* 40–50mmHg and 50–60mmHg available in NATURAL RUBBER only					
To order Plus Calf with Grip-top or Pantyhose: add “/P” after the color in the ordering code (<i>860P for women only</i>)					
To order Calf Style with Grip-top: add “/S” after the color in the ordering code (<i>see product pages for availability</i>)					
To order Thigh with Waist Attachment: add “/L” to indicate the Left leg or “/R” to indicate the Right leg after the color in the ordering code (<i>500W and 860W only</i>)					
To order TRAVERSE SOCKS and MERINO WOOL: The digit for Gender is not used in the ordering code. Example: 412CSM59 = TRAVERSE SOCKS, 20–30mmHg, Closed toe Calf, Size Small Medium, Color Red					

HOW TO ORDER WELL BEING

Series Number 120 = SHEER FASHION <i>(See individual pages for series number)</i>	Product Style N = Thigh-high <i>(See individual pages for available styles)</i>	Size A = Small <i>(See individual pages for sizing)</i>	Color 33 = Natural <i>(See individual pages for available colors)</i>
120 N A 33			
To order ATHLETIC RECOVERY SOCKS, PERFORMANCE SLEEVES, and MERINO OUTDOOR SOCKS: Size is only one letter and the digit for Gender is not used in the ordering code. Example: 421CL31 = MERINO OUTDOOR SOCK, 15–20mmHg, Closed toe Calf, Size Large, Color Olive			

TO ORDER:
 1.800.322.7744, ext. 1
Online Dealer Portal:
sigvaris-online.com

SIGVARIS DEALER PORTAL

The SIGVARIS Dealer Portal is your online location for securely ordering SIGVARIS products and managing your account. Place and track your orders, review order history, contact Customer Care, and much more at sigvaris-online.com.

DEALER PORTAL BENEFITS

24/7 access-order any time!

Three ways to securely place an order online:
Catalog, Quick Entry or Order Upload

Access to order history and tracking

Shipping options including drop ship and easy access to an address table for quick shipment of reorders

Check your account balance and print invoices



Email sigvariscares@sigvaris.com and you will receive your active account within one business day!

SIGVARIS EDUCATION

SIGVARIS has several new and exciting educational opportunities for 2018 to help you get a leg up on leg health knowledge. These opportunities include online learning, webinars, and certified fitter in-services, as well as a basic or advanced course, in a city near you. Learn more at sigvariseducation.com.

BECOME A SIGVARIS CERTIFIED FITTER!

Contact your Territory Manager

or call Customer Care at 1-800-322-7744, ext. 1

or email sigvariscares@sigvaris.com



HOW TO MEASURE

Please follow these instructions carefully to determine the correct size for different product styles.

Measurements should always be taken in the following order, and *as soon after waking as possible when swelling is minimal.*

CALF STYLE


Open/Closed Toe

1	Measure the ankle circumference at the narrowest part of the ankle, above the ankle bone. Use the appropriate series sizing chart to determine if ankle is SMALL, MEDIUM, LARGE or X-LARGE. Note: This should be your first measurement.
2	Measure the calf circumference at the fullest part of the calf.
3	Determine the calf length by measuring from the floor up to the bend of the knee. This measurement is taken on the outside of the leg.

PANTYHOSE AND THIGH-HIGH STYLE

Open/Closed Toe

1	Measure the ankle circumference at the narrowest part of the ankle, above the ankle bone. Use the appropriate series sizing chart to determine if ankle is SMALL, MEDIUM or LARGE. Note: This should be your first measurement.
2	Measure the calf circumference at the fullest part of the calf.
3	Skip Step 3 (but refer to General Guidelines #6 below).
4	Measure the thigh at its widest circumference.
5	Determine the leg length by measuring from the floor up to the gluteal fold.
6	Determine the hip circumference by measuring the widest part of the hips.



GENERAL GUIDELINES

1	Measurements should always be taken as soon after waking as possible, when swelling is minimal.
2	Measure against bare skin, not over clothing.
3	The measuring system is based on having the foot flat against the floor and the calf or leg at a right angle to the foot (see above images for examples).
4	Measurements should be taken by pulling the tape around the leg so that it is snug but not constricting.
5	Length measurements should be taken from the ground up.
6	Take all measurements, even if the patient only requires calf-length products, as the patient might decide at a later date to purchase additional styles.

WHY SIGVARIS?

Our objective is to produce the best compression garments to enhance our customers' health and quality of life.

We are dedicated to helping legs look and feel their best. Using the world's premium yarns, we have spun ourselves into the worldwide leader in graduated compression socks and hosiery by embracing our core beliefs of **QUALITY**, **COMFORT**, **FIT** and **INNOVATION**. Staying true to these values, we are continuing our efforts by expanding our business to help more patients with various edemas.

QUALITY

More than 20 people handle our products from start to finish, ensuring every garment is of the highest quality. We want to make sure that when patients choose SIGVARIS, they can count on our quality.



COMFORT

We offer the most comfortable fabric options available, including the finest cottons, wools and microfibers. We also double-cover the yarns on our premium compression stockings for wearing ease.



FIT

We offer premium products designed for all shapes and sizes, from petite to plus, including compression for a broad spectrum of edemas. All our products are designed using the most recent anthropometric studies and thorough wear tests to ensure a proper fit for exceptional comfort and efficacy.



INNOVATION

With more than 60 years in the compression market, we continue to invest in medical research and product development to improve the quality of life for patients suffering from venous and lymphatic disorders. To do this, we develop medically efficacious circular-knit products using the finest yarns and the most sophisticated knitting machines.



SIGVARIS, PFS², and Sensinnov are registered trademarks of SIGVARIS AG, CH-9014 St.Gallen/Switzerland, in many countries worldwide. © Copyright 2018 SIGVARIS, Inc.

It Stays is a trademark of KIOUS, O.H. Lycra is a registered trademark of INVISTA. Supima is a registered trademark of the Supima Association of America Corporation. NILIT is a registered trademark of Nilit. All rights reserved. *drirelease* and FreshGuard are copyright ©2018 *drirelease*. *drirelease* is an Optimizer Brands company. All Rights Reserved. Sim-Slide and Magride are registered trademarks of Arion Holding B.V.

SIGVARIS MEDICAL

	SERIES	PAGE	COMPRESSION LEVELS	STYLES
NATURAL FIBERS	COTTON (230)	6	20–30mmHg 30–40mmHg	closed toe: calf, calf w/grip-top, thigh-high w/grip-top open toe: calf, thigh-high w/grip-top
	CUSHIONED COTTON (360)	8	20–30mmHg	closed toe: calf
	SEA ISLAND COTTON (220) <small>ZÜRICH COLLECTION</small>	10	20–30mmHg	closed toe: calf
	ALL-SEASON MERINO WOOL (240) <small>ZÜRICH COLLECTION</small>	12	20–30mmHg	closed toe: calf
SHEER	ALLURE (710) <i>(for women only)</i>	14	15–20mmHg 20–30mmHg	closed toe: thigh-high w/grip-top, pantyhose
	EVERSHEER (780) <i>(for women only)</i>	16	15–20mmHg 20–30mmHg 30–40mmHg	closed & open toe: calf, thigh-high w/grip-top, pantyhose, maternity pantyhose
	MIDSHEER (750) <i>(for women only)</i>	18	20–30mmHg	closed toe: calf, thigh-high w/grip-top, pantyhose
OPAQUE	SOFT OPAQUE (840) <i>(for women only)</i>	20	15–20mmHg 20–30mmHg 30–40mmHg	closed toe: calf, calf w/grip-top, thigh-high w/grip-top, pantyhose, maternity pantyhose open toe: calf, thigh-high w/grip-top, pantyhose
	SELECT COMFORT (860)	22	20–30mmHg 30–40mmHg	closed toe: calf, petite calf, calf plus, calf w/grip-top, thigh w/grip-top, petite thigh-high w/grip-top, extra-large thigh w/grip-top, pantyhose, petite pantyhose, plus-sized pantyhose open toe: calf, calf plus, calf w/grip-top, thigh w/grip-top, extra-large thigh w/grip-top, thigh-high w/waist attachment
	MIDTOWN MICROFIBER (820) <i>(for men only)</i>	24	15–20mmHg 20–30mmHg 30–40mmHg	closed toe: calf, calf w/grip-top, thigh w/grip-top, extended foot sizes
	MICROFIBER SHADES (830)	26	20–30mmHg	closed toe: calf
ACTIVE	DAILY COMFORT (850)	28	20–30mmHg	closed toe: calf
	TRAVERSE SOCKS (412)	30	20–30mmHg	closed toe: calf
	MERINO WOOL (422)	32	20–30mmHg	closed toe: calf
SPECIALTY	DIABETIC COMPRESSION SOCKS (602)	34	18–25mmHg	closed toe: calf
	ANTI-EMBOLISM STOCKINGS (930)	36	18–23mmHg	open toe: calf, thigh-high
	NATURAL RUBBER (500)	38	30–40mmHg 40–50mmHg 50–60mmHg	open toe: calf, thigh-high, thigh w/grip-top, thigh-high w/waist attachment, pantyhose
	ADVANCE ARMSLEEVE (910)	40	15–20mmHg 20–30mmHg 30–40mmHg	grip-top with gauntlet grip-top without gauntlet

SIGVARIS WELL BEING & CONSUMER HEALTH

SERIES	PAGE	COMPRESSION LEVEL	STYLES
CASUAL COTTON (146/186)	44	15–20mmHg	closed toe: calf
CUSHIONED COTTON (142/182)	46	15–20mmHg	closed toe: calf
SEA ISLAND COTTON (151/191) <small>ZURICH COLLECTION</small>	48	15–20mmHg	closed toe: calf
ALL-SEASON MERINO WOOL (152/192) <small>ZURICH COLLECTION</small>	50	15–20mmHg	closed toe: calf
MICROFIBER SHADES (143/183)	52	15–20mmHg	closed toe: calf
SHEER FASHION (120) <i>(for women only)</i>	54	15–20mmHg	closed toe: calf, thigh-high w/grip-top, pantyhose, maternity pantyhose open toe: calf
SOFT SILHOUETTE LEGGINGS (170) <i>(for women only)</i>	56	15–20mmHg	open foot: leggings
BUSINESS CASUAL (189) <i>(for men only)</i>	58	15–20mmHg	closed toe: calf
ATHLETIC RECOVERY SOCKS (401)	60	15–20mmHg	closed toe: calf
PERFORMANCE SLEEVES (412)	62	20–30mmHg	open foot: calf sleeve
MERINO OUTDOOR SOCKS (421)	64	15–20mmHg	closed toe: calf
EVERSOFT DIABETIC SOCKS (160)	66	8–15mmHg	closed toe: calf

Custom sizes are available for a cotton product, a microfiber product, and our NATURAL RUBBER series. For more information, please speak with your SIGVARIS Sales Representative or Customer Care: 1.800.322.7744, ext. 1.

SIGVARIS IS THE GLOBAL LEADER IN COMPRESSION THERAPY

From our comfortable, patented top bands to our double-covered yarns that allow our products to glide on and off with ease, SIGVARIS designs garments to meet patients' needs. We use the finest textiles that provide optimal comfort and improved patient compliance.



As you look through the pages of our catalog, you'll discover products that not only deliver medically sound therapeutic results, but also provide the ultimate product assortment designed to fit any lifestyle. From petite to plus, and everything in between, legs love living in SIGVARIS!



COTTON

for women & men



FEATURES

Only Supima Cotton touches the skin to ensure superior comfort

Breathable and durable

Classic ribbed look

Ideal for everyday wear and for sensitive skin

BEST SUITED FOR

Containment of various types of edema, including venous edema, phlebolympheidema, and early stages of lymphedema

All-day cotton comfort in all types of climates

Ideal for dress or casual shoes

Ideal for long-distance travel

FIBER CONTENT

50% Nylon


25% Spandex

25% Supima® Cotton

(content can vary slightly by style)

Latex-Free

SUPIMA®
WORLD'S FINEST COTTONS





20–30mmHg Closed Toe Calf
Medium Short in Black Mist
232CMSW14



Silicone band
grip-top

20–30mmHg Closed Toe Thigh-high
Medium Long in Black
232NMLM99

AVAILABLE STYLES AND COLORS

Closed Toe = ●	White (00)	Navy (10)	Black Mist (14)	Crispa (66)	Chocolate (88)*	Black (99)
Open Toe = ○						
CALF (C)	●	●	●	●●	●	●
CALF WITH GRIP-TOP (C/S)				●*		●**
THIGH WITH GRIP-TOP (N)				●●		●

* 20–30mmHg ONLY

** 30–40mmHg CALF WITH GRIP-TOP ONLY AVAILABLE IN BLACK

30–40mmHg EXTRA LARGE CALF ONLY AVAILABLE IN BLACK WITH GRIP-TOP

COMPRESSION LEVELS

20–30mmHg	30–40mmHg
-----------	-----------

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"–8.5" (18–21.5cm)		MEDIUM 8.5"–10" (21.5–25.5cm)		LARGE 10"–11.5" (25.5–29cm)		EXTRA LARGE 11.5"–14" (29–36cm)	
CALF CIRCUMFERENCE	SMALL 11"–15.5" (28–39cm)		MEDIUM 12.5"–17.5" (31.5–45cm)		LARGE 14"–20" (35.5–51cm)		EXTRA LARGE 16"–24" (41–61cm)	
CALF LENGTH	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)
THIGH CIRCUMFERENCE	SMALL 18"–26" (45.5–66cm)		MEDIUM 19"–28.5" (48–73cm)		LARGE 20.5"–33" (52–84cm)			
LEG LENGTH: WOMEN & OPEN TOE	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)	XS & XL not available in Thigh style	
LEG LENGTH: MEN	SHORT ≤30" (76cm)	LONG >30" (76cm)	SHORT ≤30" (76cm)	LONG >30" (76cm)	SHORT ≤30" (76cm)	LONG >30" (76cm)		
SIZE	SS	SL	MS	ML	LS	LL	XS	XL

MEN'S SHOE SIZES: SHORT LEG 7–11, LONG LEG 9.5–14 **WOMEN'S SHOE SIZES:** SHORT LEG 5.5–9.5, LONG LEG 8–12

Product Code Example: For product at the top left of page: 232CMSW14	Series	mmHg Level	Product Style	Size	Gender/ Open Toe	Color
	23	2	C	MS	W	14

Product Code Example: For 20–30mmHg Closed Toe Calf w/ Grip-top Large Long for Men in Crispa: 232CLLM66/S	Series	mmHg Level	Product Style	Size	Gender/ Open Toe	Color	Grip-top
	23	2	C	LL	M	66	/S

CUSHIONED COTTON

for women & men



FEATURES

Thick, padded sole, heel and toe to protect feet

Flat, low-friction toe seam

Soft, stretchable knit-in band prevents pinching

No compression in heel and instep for ease of donning

Ideal for everyday wear and for sensitive skin

Recommended for patients with diabetes

BEST SUITED FOR

Management of CVD symptoms

People who stand for prolonged periods of time

Diabetic friendly

Patients with sensitive toes

Work shoes and boots

Casual shoes and sports shoes

FIBER CONTENT

64% Supima Cotton

28% Nylon

8% Spandex

Latex-Free

SUPIMA®
WORLD'S FINEST COTTONS

WHITE
(00)

BLACK
(99)



20–30mmHg Closed Toe Calf
Small Long in White
362CSLW00



20–30mmHg Closed Toe Calf
Medium Long in Black
362CMLM99

AVAILABLE STYLE AND COLORS

Closed Toe = ●	White (00)	Black (99)
CALF (C)	●	●

COMPRESSION LEVEL

20–30mmHg
ALSO AVAILABLE IN 15–20mmHg, SEE PAGE 46

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"–8.5" (18–21.5cm)		MEDIUM 8.5"–10" (21.5–25.5cm)		LARGE 10"–11.5" (25.5–29cm)		EXTRA LARGE 11.5"–14" (29–36cm)	
CALF CIRCUMFERENCE	SMALL 11"–15.5" (28–39cm)		MEDIUM 12.5"–17.5" (31.5–45cm)		LARGE 14"–20" (35.5–51cm)		EXTRA LARGE 16"–24" (41–61cm)	
CALF LENGTH	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)
SIZE	SS	SL	MS	ML	LS	LL	XS	XL

MEN'S SHOE SIZES: SHORT LEG 7–11, LONG LEG 9.5–14 **WOMEN'S SHOE SIZES:** SHORT LEG 5.5–9.5, LONG LEG 8–12

Product Code Example: For product at the top left of page: 362CSLW00	Series	mmHg Level	Product Style	Size	Gender/ Open Toe	Color
	36	2	C	SL	W	00

SEA ISLAND COTTON

for women & men

FEATURES

Knitted from Sea Island Cotton, the world's rarest cotton, for a fashionable look and a soft, luxurious feel

Only cotton touches the skin to ensure comfort

Breathable and comfortable for sensitive skin

BEST SUITED FOR

Management of symptoms related to chronic venous disease

Patients who prefer natural fibers for everyday wear

Long-distance travel

FIBER CONTENT

43% Sea Island Cotton


39% Nylon

11% Spandex

7% Supima Cotton

Latex-Free

ZURICH
COLLECTION





20–30mmHg Closed Toe Calf
Small Long in Brown
222CSLW11



20–30mmHg Closed Toe Calf
Medium Long in Black
222CMLM99

AVAILABLE STYLE AND COLORS

Closed Toe = ●	Navy (10)	Brown (11)	Black (99)
CALF (C)	●	●	●

COMPRESSION LEVEL

20–30mmHg
ALSO AVAILABLE IN 15–20mmHg, SEE PAGE 48

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"–8.5" (18–21.5cm)		MEDIUM 8.5"–10" (21.5–25.5cm)		LARGE 10"–11.5" (25.5–29cm)		EXTRA LARGE 11.5"–14" (29–36cm)	
CALF CIRCUMFERENCE	SMALL 11"–15.5" (28–39cm)		MEDIUM 12.5"–17.5" (31.5–45cm)		LARGE 14"–20" (35.5–51cm)		EXTRA LARGE 16"–24" (41–61cm)	
CALF LENGTH	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)
SIZE	SS	SL	MS	ML	LS	LL	XS	XL

MEN'S SHOE SIZES: SHORT LEG 7–11, LONG LEG 9.5–14 **WOMEN'S SHOE SIZES:** SHORT LEG 5.5–9.5, LONG LEG 8–12

Product Code Example: For product at the top left of page: 222CSLW11	Series	mmHg Level	Product Style	Size	Gender	Color
	22	2	C	SL	W	11

ALL-SEASON MERINO WOOL

for women & men



FEATURES

Soft, comfortable, and itch-free

Made from the highest quality extra-fine Australian Merino Wool

Natural properties of wool make these socks odor-resistant

Thermo-regulating for year-round wear and keeps feet dry

BEST SUITED FOR

Management of symptoms related to chronic venous disease

Patients who prefer fashionable and premium fabric options

Long-distance travel

FIBER CONTENT

53% Merino Wool

38% Nylon

9% Spandex

Latex-Free

NAVY
(10)

BROWN
(11)

BLACK
(99)

ZURICH
COLLECTION



20–30mmHg Closed Toe Calf
Medium Long in Brown
242CMLM10



20–30mmHg Closed Toe Calf
Small Long in Navy
242CSLW10

AVAILABLE STYLE AND COLORS

Closed Toe = ●	Navy (10)	Brown (11)	Black (99)
CALF (C)	●	●	●

COMPRESSION LEVEL

20–30mmHg

ALSO AVAILABLE IN 15–20mmHg,
SEE PAGE 50

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"–8.5" (18–21.5cm)		MEDIUM 8.5"–10" (21.5–25.5cm)		LARGE 10"–11.5" (25.5–29cm)		EXTRA LARGE 11.5"–14" (29–36cm)	
CALF CIRCUMFERENCE	SMALL 11"–15.5" (28–39cm)		MEDIUM 12.5"–17.5" (31.5–45cm)		LARGE 14"–20" (35.5–51cm)		EXTRA LARGE 16"–24" (41–61cm)	
CALF LENGTH	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)
SIZE	SS	SL	MS	ML	LS	LL	XS	XL

MEN'S SHOE SIZES: SHORT LEG 7–11, LONG LEG 9.5–14 **WOMEN'S SHOE SIZES:** SHORT LEG 5.5–9.5, LONG LEG 8–12

Product Code Example: For product at the top left of page: 242CMLM10	Series	mmHg Level	Product Style	Size	Gender	Color
	24	2	C	ML	M	10

ALLURE

for women



FEATURES

Stylish diamond-dot pattern for fashionable daily wear

Durable construction with reinforced heel

Excellent breathability and moisture management for added comfort

Cotton-lined gusset on pantyhose style

Sensinnov® band grip-top on thigh style

BEST SUITED FOR

Management of early stages of chronic venous disease

Pre- and post-vein treatment with minimal edema

Patients who prefer fashionable options

Ideal for business or dress

FIBER CONTENT

64% Nylon

36% Spandex

(can vary slightly by style)

Latex-Free





20–30mmHg Closed Toe Pantyhose
Small Long in Graphite
712PSLW91



15–20mmHg Closed Toe Thigh-high
Medium Long in Black
711NMLW99



Sensinnov
band grip-top

AVAILABLE STYLES AND COLORS

Closed Toe = ●	Graphite (91)	Black (99)
THIGH W/GRIP-TOP (N)	●	●
PANTYHOSE (P)	●	●

COMPRESSION LEVELS

15–20mmHg	20–30mmHg
-----------	-----------

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"–8.5" (18–21.5cm)		MEDIUM 8.5"–10" (21.5–25.5cm)		LARGE 10"–11.5" (25.5–29cm)	
CALF CIRCUMFERENCE	SMALL 11"–15.5" (28–39cm)		MEDIUM 12.5"–17.5" (31.5–45cm)		LARGE 14"–20" (35.5–51cm)	
THIGH CIRCUMFERENCE	SMALL 18"–26" (45.5–66cm)		MEDIUM 19"–28.5" (48–73cm)		LARGE 20.5"–33" (52–84cm)	
LEG LENGTH: WOMEN	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)
GUIDELINES FOR HIP CIRCUMFERENCE	SMALL 32.5"–44.5" (83–112.5cm)		MEDIUM 34.5"–50" (87–127cm)		LARGE 38.5"–59.5" (98–150.5cm)	
SIZE	SS	SL	MS	ML	LS	LL

WOMEN'S SHOE SIZES: SHORT LEG 5.5–9.5, LONG LEG 8–12

Product Code Example: For product at the top left of page: 712PSLW91	Series	mmHg Level	Product Style	Size	Gender	Color
	71	2	P	SL	W	91

EVERSHEER

for women



FEATURES

- Fashionable, sheer medical graduated compression hosiery
- Constructed from double-covered inlay yarns for softness against the skin
- Excellent breathability and moisture management for added comfort
- Cotton-lined gusset on pantyhose styles
- Sensinno band grip-top on thigh style

BEST SUITED FOR

- Management of early stages of chronic venous disease
- Pre- and post-vein treatment with minimal edema
- Patients who prefer fashionable options
- Ideal for business or dress

FIBER CONTENT

- 64% Nylon
- 36% Spandex
- (content can vary slightly by style)
- Latex-Free*

DARK NAVY (08)	NATURAL (33)	
SUNTAN (36)	CAFÉ (73)	MOCHA (85)
	NIGHTSHADE (94)	BLACK (99)



Sensinno band grip-top



AVAILABLE STYLES AND COLORS

Closed Toe = ●	Dark Navy (08)	Natural (33)	Suntan (36)	Café (73)	Mocha (85)	Nightshade (94)	Black (99)
Open Toe = ○							
CALF (C)	●	●○	●○	●○	●○	●○	●○
THIGH W/ GRIP-TOP (N)	●	●○	●○	●○	●○	●○	●○
PANTYHOSE (P)	●	●○	●○	●○	●○	●○	●○
MATERNITY PANTYHOSE (M)*		●					●

* 20–30mmHg ONLY

COMPRESSION LEVELS

15–20mmHg	20–30mmHg	30–40mmHg
-----------	-----------	-----------

15–20mmHg Closed Toe Thigh-high Medium Long in Café
781NMLW73

20–30mmHg Closed Toe Maternity Pantyhose Small Long in Natural
782MSLW33

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"–8.5" (18–21.5cm)		MEDIUM 8.5"–10" (21.5–25.5cm)		LARGE 10"–11.5" (25.5–29cm)	
CALF CIRCUMFERENCE	SMALL 11"–15.5" (28–39cm)		MEDIUM 12.5"–17.5" (31.5–45cm)		LARGE 14"–20" (35.5–51cm)	
CALF LENGTH	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)
THIGH CIRCUMFERENCE	SMALL 18"–26" (45.5–66cm)		MEDIUM 19"–28.5" (48–73cm)		LARGE 20.5"–33" (52–84cm)	
LEG LENGTH: WOMEN	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)
GUIDELINES FOR HIP CIRCUMFERENCE	SMALL 32.5"–44.5" (83–112.5cm)		MEDIUM 34.5"–50" (87–127cm)		LARGE 38.5"–59.5" (98–150.5cm)	
SIZE	SS	SL	MS	ML	LS	LL

WOMEN'S SHOE SIZES: SHORT LEG 5.5–9.5, LONG LEG 8–12

MATERNITY HIP CIRCUMFERENCE GUIDELINES: SMALL: <56" (142cm), MEDIUM: <62" (157.5cm), LARGE: <70" (178cm)

Product Code Example: For product at the top left of page: 781NMLW73	Series	mmHg Level	Product Style	Size	Gender/ Open Toe	Color
	78	1	N	ML	W	73

MIDSHEER

for women



FEATURES

Elegant and durable fabric for every woman, every day

Consistent sheerness from ankle to thigh

Excellent breathability and moisture management

Comfortable silicone dot band on thigh-high style to prevent rolling

Extra-wide calf band for comfort

Cotton-lined gusset on pantyhose style

BEST SUITED FOR

Management of early stages of chronic venous disease

Pre- and post-vein treatment with minimal edema

Daily wear for business or casual

Active professionals needing a durable sheer




FIBER CONTENT

68% Nylon

32% Spandex


(content can vary slightly by style)

Latex-Free





20–30mmHg Closed Toe Pantyhose
Small Long in Dark Navy
752PSLW08



20–30mmHg Closed Toe Thigh-high
Medium Long in Black
752NMLW99

AVAILABLE STYLES AND COLORS

	Dark Navy (08)	Natural (33)	Nightshade (94)	Black (99)
Closed Toe = ●	●	●	●	●
CALF (C)	●	●	●	●
THIGH W/ GRIP-TOP (N)	●	●	●	●
PANTYHOSE (P)	●	●	●	●

COMPRESSION LEVEL

20–30mmHg

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"–8.5" (18–21.5cm)		MEDIUM 8.5"–10" (21.5–25.5cm)		LARGE 10"–11.5" (25.5–29cm)	
CALF CIRCUMFERENCE	SMALL 11"–15.5" (28–39cm)		MEDIUM 12.5"–17.5" (31.5–45cm)		LARGE 14"–20" (35.5–51cm)	
CALF LENGTH	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)
THIGH CIRCUMFERENCE	SMALL 18"–26" (45.5–66cm)		MEDIUM 19"–28.5" (48–73cm)		LARGE 20.5"–33" (52–84cm)	
LEG LENGTH: WOMEN	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)
GUIDELINES FOR HIP CIRCUMFERENCE	SMALL 32.5"–44.5" (83–112.5cm)		MEDIUM 34.5"–50" (87–127cm)		LARGE 38.5"–59.5" (98–150.5cm)	
SIZE	SS	SL	MS	ML	LS	LL

WOMEN'S SHOE SIZES: SHORT LEG 5.5–9.5, LONG LEG 8–12

Product Code Example: For product at the top left of page: 752PSLW08	Series	mmHg Level	Product Style	Size	Gender	Color
	75	2	P	SL	W	08

SOFT OPAQUE

for women



FEATURES

- Comfortable tights with lasting microfiber softness
- Advanced moisture management that keeps legs cool and dry
- Comfortable wide top band prevents knee-highs from rolling down
- Cotton-lined gusset on pantyhose styles
- Sensinno band grip-top on thigh style

BEST SUITED FOR

- Management of early stages of chronic venous disease
- Pre- and post-vein treatment with minimal edema
- Patients who prefer fashionable options
- Pregnancy
- Long-distance travel


FIBER CONTENT

- 79% Nylon
- 21% Spandex
- (content can vary slightly by style)
- Latex-Free*

MIDNIGHT BLUE (09)	NUDE (35)	
PECAN (40)	MULBERRY (60)	ESPRESSO (89)
GRAPHITE (91)	BLACK (99)	



20–30mmHg Closed Toe Calf
Medium Long in Mulberry
842CMLW60



15–20mmHg Closed Toe Thigh-high
Medium Long in Black
841NMLW99



Sensinno
band grip-top

AVAILABLE STYLES AND COLORS

Closed Toe = ●	Midnight Blue (09)	Nude (35)	Pecan (40)	Mulberry (60)*	Espresso (89)**	Graphite (91)	Black (99)
Open Toe = ○							
CALF (C)	●	○	●	●	●	●	●○
CALF W/ BEADED GRIP-TOP (C/S)**							●
THIGH W/ GRIP-TOP (N)	●	○	●	●	●	●	●○
PANTYHOSE (P)	●	○	●	●	●	●	●○
MATERNITY PANTYHOSE (M)*							●

* 15–20 & 20–30mmHg ONLY ** 20–30mmHg ONLY

COMPRESSION LEVELS

15–20mmHg	20–30mmHg	30–40mmHg
-----------	-----------	-----------

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"–8.5" (18–21.5cm)		MEDIUM 8.5"–10" (21.5–25.5cm)		LARGE 10"–11.5" (25.5–29cm)	
CALF CIRCUMFERENCE	SMALL 11"–15.5" (28–39cm)		MEDIUM 12.5"–17.5" (31.5–45cm)		LARGE 14"–20" (35.5–51cm)	
CALF LENGTH	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)
THIGH CIRCUMFERENCE	SMALL 18"–26" (45.5–66cm)		MEDIUM 19"–28.5" (48–73cm)		LARGE 20.5"–33" (52–84cm)	
LEG LENGTH: WOMEN	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)
GUIDELINES FOR HIP CIRCUMFERENCE	SMALL 32.5"–44.5" (83–112.5cm)		MEDIUM 34.5"–50" (87–127cm)		LARGE 38.5"–59.5" (98–150.5cm)	
SIZE	SS	SL	MS	ML	LS	LL

WOMEN'S SHOE SIZES: SHORT LEG 5.5–9.5, LONG LEG 8–12

MATERNITY HIP CIRCUMFERENCE GUIDELINES: SMALL: <56" (142cm), MEDIUM: <62" (157.5cm), LARGE: <70" (178cm)

Product Code Example: For product at the top left of page: 842CMLW60	Series	mmHg Level	Product Style	Size	Gender/ Open Toe	Color
	84	2	C	ML	W	60

Product Code Example: For 20–30mmHg Closed Toe Calf w/ Grip-top Small Short in Black: 842CSSW99/S	Series	mmHg Level	Product Style	Size	Gender/ Open Toe	Color	Grip-top
	84	2	C	SS	W	99	/S

SELECT COMFORT

*closed toe for women and
open toe for women & men*



FEATURES

From petite to plus, SELECT COMFORT offers a complete line of sizes designed to fit more than 99% of the population

Durable construction combined with a smooth, opaque look

Comfortable, stretchable material constructed from double-covered yarn easily glides on and off

Plus-size pantyhose is a great alternative for maternity

Cotton-lined gusset on pantyhose styles

BEST SUITED FOR

Containment of various types of edema including venous edema, phlebolympheidema, and early stages of lymphedema

Post-vein treatment including thermal and non-thermal techniques

Patients looking for petite and plus-sized options


FIBER CONTENT

64% Nylon

36% Spandex

(content can vary slightly by style)

Latex-Free



Product Code Example: For product at the top left of facing page: 863CSLW99	Series	mmHg Level	Product Style	Size	Gender/ Open Toe	Color
	86	3	C	SL	W	99

Product Code Example: For 20–30mmHg Plus Size Closed Toe Calf Large Long for Women in Crispa: 862CLLW66/P	Series	mmHg Level	Product Style	Size	Gender/ Open Toe	Color	Plus Size
	86	2	C	LL	W	66	/P



30–40mmHg Closed Toe Calf
Small Long in Black
863CSLW99



20–30mmHg Open Toe Thigh with
Grip-top Medium Long in Crispa
862NMLO66

AVAILABLE STYLES AND COLORS

	Natural (33)	Suntan (36)	Crispa (66)	Black (99)
Closed Toe = ● Open Toe = ○				
CALF (C)	●	●	○	●○
CALF PLUS W/GRIP-TOP (C/P)			●○	
CALF W/GRIP-TOP (C/S)*			●○	●
THIGH W/GRIP-TOP (N)	●	●	○	●
XS & XL THIGH W/GRIP-TOP (N)			●○	
PANTYHOSE (P)	●	●		●
PANTYHOSE PLUS (P/P)**			●	●
THIGH WITH WAIST ATTACHMENT (W) (INDICATE LEG WITH /L OR /R)			○	

NOTE: 860 CLOSED TOE PRODUCTS ARE SIZED FOR WOMEN, OPEN TOE PRODUCTS ARE UNISEX
* CALF PLUS STYLE AVAILABLE ONLY WITH GRIP-TOP
** PANTYHOSE PLUS IS AN ALTERNATIVE FOR MATERNITY

COMPRESSION LEVELS

20–30mmHg	30–40mmHg
-----------	-----------

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"–8.5" (18–21.5cm)			MEDIUM 8.5"–10" (21.5–25.5cm)			LARGE 10"–11.5" (25.5–29cm)		EXTRA LARGE 11.5"–14" (29–36cm)	
CALF CIRCUMFERENCE	SMALL 11"–15.5" (28–39cm)			MEDIUM 12.5"–17.5" (31.5–45cm)			LARGE 14"–20" (35.5–51cm)		EXTRA LARGE 16"–24" (41–61cm)	
PLUS SIZE CALF CIRCUMFERENCE	Not available			MEDIUM 15.5"–21" (39.5–53cm)			LARGE 18"–24" (45.5–61cm)		EXTRA LARGE 20"–28" (51–71cm)	
CALF LENGTH	PETITE** ≤14" (35.5cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	PETITE** ≤14" (35.5cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)
THIGH CIRCUMFERENCE	SMALL 18"–26" (45.5–66cm)			MEDIUM 19"–28.5" (48–73cm)			LARGE 20.5"–33" (52–84cm)		EXTRA LARGE 25"–37.5" (64–95cm)	
LEG LENGTH: WOMEN & OPEN TOE	PETITE** ≤25.5" (65cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)	PETITE** ≤25.5" (65cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)
GUIDELINES FOR HIP CIRCUMFERENCE	SMALL 32.5"–44.5" (83–112.5cm)			MEDIUM 34.5"–50" (87–127cm)			LARGE 38.5"–59.5" (98–150.5cm)		Not available in Full Pantyhose or Thigh w/ Waist Attachment	
SIZE	SP**	SS	SL	MP**	MS	ML	LS	LL	XS	XL
PLUS SIZE CALF	Not available				MS/P	ML/P	LS/P	LL/P	XS/P	XL/P

WOMEN'S SHOE SIZES: PETITE LEG 4–8, SHORT LEG 5.5–9.5, LONG LEG 8–12

PANTYHOSE PLUS SIZES HIP CIRCUMFERENCE GUIDELINES: SMALL: <56" (142cm), MEDIUM: <62" (157.5cm), LARGE: <70" (178cm)

** PETITE SIZING (SP & MP): AVAILABLE IN CLOSED TOE ONLY, WITHOUT GRIP-TOP, 20–30mmHg ONLY, UNAVAILABLE IN PLUS SIZES

MIDTOWN MICROFIBER

for men



FEATURES

100% of microfiber rests against the skin to ensure comfort

Excellent breathability and moisture management for added comfort

Easy to put on and take off with a non-binding top band

Soft and durable sock for everyday wear

Sensinnov band grip-top on thigh style

BEST SUITED FOR

Management of symptoms related to chronic venous disease

Patients who prefer fashionable options

Men with larger feet

Long-distance travel

FIBER CONTENT

87% Nylon

13% Spandex

(content can vary slightly by style)

Latex-Free




NAVY
(10)




TAN KHAKI
(32)



STEEL GREY
(93)



BLACK
(99)



20–30mmHg Closed Toe Calf with Grip-top Medium Long in Black
822CMLM99/S



20–30mmHg Closed Toe Thigh-high with Grip-top Large Long in Tan Khaki
822NLLM32

AVAILABLE STYLES AND COLORS

Closed Toe = ●	Navy (10)**	Tan Khaki (32)	Steel Grey (93)**	Black (99)
CALF (C)	●	●	●	●
CALF W/GRIP-TOP (C/S)*		●		●
THIGH W/GRIP-TOP (N)		●		●

* CALF W/GRIP-TOP NOT AVAILABLE IN 15–20mmHg

** NAVY AND STEEL GREY ONLY AVAILABLE IN 15–20 & 20–30mmHg

*Looking for a men's open toe option?
See SELECT COMFORT, page 23!*

COMPRESSION LEVELS

15–20mmHg	20–30mmHg	30–40mmHg
-----------	-----------	-----------

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"–8.5" (18–21.5cm)		MEDIUM 8.5"–10" (21.5–25.5cm)		LARGE 10"–11.5" (25.5–29cm)		EXTRA LARGE 11.5"–14" (29–36cm)	
CALF CIRCUMFERENCE	SMALL 11"–15.5" (28–39cm)		MEDIUM 12.5"–17.5" (31.5–45cm)		LARGE 14"–20" (35.5–51cm)		EXTRA LARGE 16"–24" (41–61cm)	
CALF LENGTH	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)
THIGH CIRCUMFERENCE	SMALL 18"–26" (45.5–66cm)		MEDIUM 19"–28.5" (48–73cm)		LARGE 20.5"–33" (52–84cm)		XS & XL only available in Calf styles	
LEG LENGTH: MEN	SHORT ≤30" (76cm)	LONG >30" (76cm)	SHORT ≤30" (76cm)	LONG >30" (76cm)	SHORT ≤30" (76cm)	LONG >30" (76cm)		
SOCK SIZE / SHOE SIZE	SS 7–11	SL 9.5–14	MS 7–11	ML 9.5–14	LS 7–11	LL 9.5–14	XS 7–11	XL 9.5–14
				MX*** 12.5–17.5		LX*** 12.5–17.5		XX*** 12.5–17.5

MEN'S SHOE SIZES: SHORT LEG 7–11, LONG LEG 9.5–14 (FOR EXTENDED FOOT SIZES: LONG LEG 12.5–17.5)

*** MX, LX, AND XX EXTENDED FOOT SIZES ONLY AVAILABLE IN CALF STYLE IN BLACK & TAN KHAKI IN 15–20 & 20–30mmHg ONLY.

EXTENDED FOOT SIZES NOT AVAILABLE WITH GRIP TOP

Product Code Example: For product at the top left of page: 822CMLM99/S	Series	mmHg Level	Product Style	Size	Gender/ Open Toe	Color	Grip-top
	82	2	C	ML	M	99	/S

MICROFIBER SHADES

for women & men

FEATURES

- Advanced microfiber technology for ultimate comfort
- Glides on and off with ease
- Excellent breathability and moisture management
- Soft and durable sock for everyday wear

BEST SUITED FOR

- Management of symptoms related to chronic venous disease
- Patients who prefer fashionable options
- Long-distance travel

FIBER CONTENT

- 87% Nylon
- 13% Spandex
- Latex-Free*



Actual color may vary



20–30mmHg Closed Toe Calf
Medium Short in Pink Stripe
832CMSW15



20–30mmHg Closed Toe Calf
Large Long in Graphite Heather
832CLLM27

AVAILABLE STYLE AND COLORS

Closed Toe = ●	Dark Navy Mini-Stripe (03)	Graphite Mini-Stripe (07)	Dark Navy Stripe (13)	Pink Stripe (15)	Graphite Stripe (17)	Onyx Stripe (19)
CALF (C)	●	●	●	●	●	●
Closed Toe = ●	Navy Heather (23)	Graphite Heather (27)	Pink Argyle (45)	Purple Argyle (46)	Graphite Argyle (47)	Royal Blue Argyle (48)
CALF (C)	●	●	●	●	●	●

COMPRESSION LEVEL

20–30mmHg

ALSO AVAILABLE IN 15–20mmHg,
SEE PAGE 52

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"–8.5" (18–21.5cm)		MEDIUM 8.5"–10" (21.5–25.5cm)		LARGE 10"–11.5" (25.5–29cm)	
CALF CIRCUMFERENCE	SMALL 11"–15.5" (28–39cm)		MEDIUM 12.5"–17.5" (31.5–45cm)		LARGE 14"–20" (35.5–51cm)	
CALF LENGTH	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)
SIZE	SS	SL	MS	ML	LS	LL

MEN'S SHOE SIZES: SHORT LEG 7–11, LONG LEG 9.5–14 **WOMEN'S SHOE SIZES:** SHORT LEG 5.5–9.5, LONG LEG 8–12

Product Code Example: For product at the top left of page: 832CMSW15	Series	mmHg Level	Product Style	Size	Gender/ Open Toe	Color
	83	2	C	MS	W	15

DAILY COMFORT

for women & men

FEATURES

NILIT® Aquarius fibers provide thermal control to keep feet cool, dry and comfortable

Cushioned soles to absorb impact and reduce blisters

Anatomically-contoured Achilles tendon padding

Smooth, flat toe seam

Additional compression in calf improves circulation

BEST SUITED FOR

Management of symptoms related to chronic venous disease

Healthcare professionals

Business travelers

Individuals working in physically-demanding occupations

Active individuals

FIBER CONTENT

86% Nylon

14% Spandex

Latex-Free

BLACK
(99)



20–30mmHg Closed Toe Calf
Large Long in Black
852CLLM99



20–30mmHg Closed Toe Calf
Medium Long in Black
852CMLW99

AVAILABLE STYLE AND COLOR

Closed Toe = ●	Black (99)
CALF (C)	●

COMPRESSION LEVEL

20–30mmHg

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"–8.5" (18–21.5cm)		MEDIUM 8.5"–10" (21.5–25.5cm)		LARGE 10"–11.5" (25.5–29cm)	
CALF CIRCUMFERENCE	SMALL 11"–15.5" (28–39cm)		MEDIUM 12.5"–17.5" (31.5–45cm)		LARGE 14"–20" (35.5–51cm)	
CALF LENGTH	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)
SIZE	SS	SL	MS	ML	LS	LL

MEN'S SHOE SIZES: SHORT LEG 7–11, LONG LEG 9.5–14 **WOMEN'S SHOE SIZES:** SHORT LEG 5.5–9.5, LONG LEG 8–12

Product Code Example: For product at the top left of page: 852CLLM99	Series	mmHg Level	Product Style	Size	Gender	Color
	85	2	C	LL	M	99

TRAVERSE SOCKS

for women & men

FEATURES

Provides thermal control to keep feet cool, dry and comfortable

Anatomical shape of the foot bed and toe cushion zone help prevent blisters

Additional compression in calf improves circulation

BEST SUITED FOR

Management of symptoms related to chronic venous disease

People who stand on their feet all day

Walking and everyday activities

Casual and occupational footwear

Active individuals


FIBRE CONTENT

64% Nylon

17% Polyester

19% Spandex

Latex-Free





20–30mmHg Closed Toe Calf
Small Medium in Pink
412CSM56



20–30mmHg Closed Toe Calf
Medium Large in Steel Blue
412CML62

AVAILABLE STYLE AND COLORS

Closed Toe = ●	White (00)	Blue (50)	Lime (54)	Pink (56)	Red (59)
CALF (C)	●	●	●	●	●
Closed Toe = ●	Steel Blue (62)	Blood Orange (63)	Limeade (64)	Black (99)	
CALF (C)	●	●	●	●	

COMPRESSION LEVEL

20–30mmHg

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"–8.5" (18–21.5cm)			MEDIUM 8.5"–10" (21.5–25.5cm)			LARGE 10"–11.5" (25.5–29cm)			EXTRA LARGE 11.5"–13" (29–33cm)		
CALF CIRCUMFERENCE	SMALL 11"–15.5" (28–39cm)			MEDIUM 12.5"–17.5" (31.5–45cm)			LARGE 14"–20" (35.5–51cm)			EXTRA LARGE 16"–23" (40.5–58cm)		
MEN'S SHOE SIZE	≤6.5	7–9.5	10–12.5	≤6.5	7–9.5	10–12.5	7–9.5	10–12.5	13–15.5	7–9.5	10–12.5	13–15.5
WOMEN'S SHOE SIZE	4.5–7.5	8–10.5	11–13.5	4.5–7.5	8–10.5	11–13.5	8–10.5	11–13.5	>14	8–10.5	11–13.5	>14
SIZE	SS	SM	SL	MS	MM	ML	LS	LM	LL	XS	XM	XL

Product Code Example: For product at the top left of page: 412CSM56	Series	mmHg Level	Product Style	Size	Color
	41	2	C	SM	56

MERINO WOOL

for women & men

FEATURES

- Soft, comfortable and itch-free
- Natural properties of wool make these socks odor-resistant
- Mesh flex zones for comfort and ventilation
- Cushioned, padded soles protect against blisters
- No compression in heel and instep for ease of donning

BEST SUITED FOR

- Management of symptoms related to chronic venous disease
- People who stand on their feet all day
- Work boots, steel-toed shoes, and most occupational footwear
- Outdoor work environments

FIBER CONTENT

67% Fine Merino Wool

26% Nylon

7% Spandex

Latex-Free

CHARCOAL
(12)

OLIVE
(31)



20–30mmHg Calf
Large Large in Charcoal
422CLL12



20–30mmHg Calf
Small Large in Olive
422CSL31

AVAILABLE STYLE AND COLORS

Closed Toe = ●	Charcoal (12)	Olive (31)
CALF (C)	●	●

COMPRESSION LEVEL

20–30mmHg

ALSO AVAILABLE IN 15–20mmHg,
SEE PAGE 64

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"–8.5" (18–21.5cm)			MEDIUM 8.5"–10" (21.5–25.5cm)			LARGE 10"–11.5" (25.5–29cm)		
CALF CIRCUMFERENCE	SMALL 11"–15.5" (28–39cm)			MEDIUM 12.5"–17.5" (31.5–45cm)			LARGE 14"–20" (35.5–51cm)		
MEN'S SHOE SIZE	≤6.5	7–9.5	10–12.5	≤6.5	7–9.5	10–12.5	7–9.5	10–12.5	13–15.5
WOMEN'S SHOE SIZE	4.5–7.5	8–10.5	11–13.5	4.5–7.5	8–10.5	11–13.5	8–10.5	11–13.5	>14
SIZE	SS	SM	SL	MS	MM	ML	LS	LM	LL

Product Code Example: For product at the top left of page: 422CLL12	Series	mmHg Level	Product Style	Size	Color
	42	2	C	LL	12

DIABETIC COMPRESSION SOCKS

for women & men



FEATURES

The only clinically-proven and patented graduated compression sock for individuals with diabetes

Compression from ankle to knee, compression-free foot

Non-irritating toe seam and non-binding top band

drirelease provides moisture wicking properties and odor control

Thick, padded soles protect feet

BEST SUITED FOR

Individuals with diabetes and symptoms related to chronic venous disease and leg edema

Individuals with peripheral neuropathy and small vessel disease

FIBER CONTENT

76% *drirelease*® (65% Polyester, 11% Cotton)

19% Nylon

5% Spandex

Latex-Free

WHITE
(00)

ACCEPTED





18–25mmHg Closed Toe Calf
Medium Long in White
602CMLM00



18–25mmHg Closed Toe Calf
Small Short in White
602CSSW00

AVAILABLE STYLE AND COLOR

Closed Toe = ●	White (00)
CALF (C)	

COMPRESSION LEVEL

18-25mmHg
ALSO AVAILABLE: EVERSOFT DIABETIC SOCKS IN 8–15mmHg, SEE PAGE 66

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"–8.5" (18–21.5cm)		MEDIUM 8.5"–10" (21.5–25.5cm)		LARGE 10"–11.5" (25.5–29cm)		EXTRA LARGE 11.5"–14" (29–36cm)	
CALF CIRCUMFERENCE	SMALL 11"–15.5" (28–39cm)		MEDIUM 12.5"–17.5" (31.5–45cm)		LARGE 14"–20" (35.5–51cm)		EXTRA LARGE 16"–24" (41–61cm)	
CALF LENGTH	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)
SIZE	SS	SL	MS	ML	LS	LL	XS*	XL*

MEN'S SHOE SIZES: SHORT LEG 7–11, LONG LEG 9.5–14 **WOMEN'S SHOE SIZES:** SHORT LEG 5.5–9.5, LONG LEG 8–12

*XS & XL ONLY AVAILABLE IN MEN'S SIZING

Product Code Example: For product at the top left of page: 602CMLM00	Series	mmHg Level	Product Style	Size	Gender/ Open Toe	Color
	60	2	C	ML	M	00

ANTI-EMBOLISM STOCKINGS

for women & men



FEATURES

Unique open-toe design allows easy inspection of feet

Anatomical foot design makes it easy to wear from bed to physical therapy

Wide top-band designed to stay comfortably in place

BEST SUITED FOR

Individuals needing

Post-operative compression

Bed rest

Temporary use when transitioning from recovery to early rehabilitation

FIBER CONTENT

71% Nylon

29% Spandex

(content can vary slightly by style)

Latex-Free

WHITE
(00)



18–23mmHg Open Toe Calf
Medium in White
932CM00



18–23mmHg Open Toe Thigh
Small in White
932NS00

AVAILABLE STYLE AND COLOR

Open Toe = <input checked="" type="radio"/>	White (00)
CALF (C) <input checked="" type="radio"/>	
THIGH (N) <input checked="" type="radio"/>	

COMPRESSION LEVEL

18-23mmHg

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"–8.5" <i>(18–22cm)</i>	MEDIUM 8.5"–10" <i>(22–26cm)</i>	LARGE 10"–11.5" <i>(26–30cm)</i>	X LARGE 11.5"–13" <i>(30–33cm)</i>	XX LARGE 11.5"–13" <i>(30–33cm)</i>
CALF CIRCUMFERENCE	SMALL 11.5"–15" <i>(30–38cm)</i>	MEDIUM 13.5"–16.5" <i>(34–42cm)</i>	LARGE 15.5"–18.5" <i>(39–47cm)</i>	X LARGE 17.5"–21" <i>(45–53cm)</i>	XX LARGE 19.5"–23" <i>(50–58cm)</i>
THIGH CIRCUMFERENCE	SMALL 18"–23.5" <i>(46–60cm)</i>	MEDIUM 21"–26.5" <i>(53–68cm)</i>	LARGE 23.5"–30" <i>(60–76cm)</i>	X LARGE 26.5"–33" <i>(68–84cm)</i>	XX LARGE 30"–36" <i>(76–92cm)</i>
SIZE	S	M	L	X	2

CALF LENGTH: <18" (46cm), **THIGH LENGTH:** <31.5" (80cm)

Product Code Example: For product at the top left of page: 932CM00	Series	mmHg Level	Product Style	Size	Open Toe	Color
	93	2	C	M	O	00

NATURAL RUBBER

for women & men



FEATURES

- Fits a broad range of medical patients, including those with extra-large ankle and calf circumferences
- Durable, therapeutic solution for those with long-term compression needs
- Minimal fiber fatigue throughout the day
- This product can be layered with our Cotton Underliner (see Accessories on page 69) to increase compression yet glide on and off easily

BEST SUITED FOR

- Severe varicosities
- Chronic edema post DVT
- Severe venous stasis
- Prevention of ulcer recurrence
- Early stages of lymphedema
- Active venous ulcers
- C3–C6 (CEAP)

FIBER CONTENT

- 70% Nylon
- 30% Natural Latex Rubber
- (content can vary slightly by style)
- (Caution: this product contains Natural Rubber Latex which may cause allergic reactions)





Silicone band grip-top

30–40mmHg Open Toe Thigh-high w/Grip-top Small Average Long in Beige
503NS2077



40–50mmHg Open Toe Calf Medium Full Long in Beige
504CM4077

AVAILABLE STYLES AND COLOR

Open Toe = ○	BEIGE (77)	
CALF (C)	○	CALF ONLY AVAILABLE IN S, M, AND L FULL CIRCUMFERENCE FOR 50–60mmHg
THIGH (T) WITHOUT GRIP-TOP	○	THIGH ALSO AVAILABLE IN X-SMALL (XS)
THIGH (N) WITH GRIP-TOP	○	NOT AVAILABLE IN 50–60mmHg
PANTYHOSE (P)	○	PANTYHOSE ONLY AVAILABLE IN 30–40mmHg AND ONLY IN FULL THIGH CIRCUMFERENCE
THIGH (W) WITH WAIST ATTACHMENT (INDICATE LEG WITH /L OR /R)	○	THIGH WITH WAIST ATTACHMENT ONLY AVAILABLE IN 30–40mmHg AND 40–50mmHg AND ONLY AVAILABLE IN FULL THIGH CIRCUMFERENCE

*CAUTION: THIS PRODUCT CONTAINS NATURAL RUBBER LATEX WHICH MAY CAUSE ALLERGIC REACTIONS. FOR ALL INDICATIONS & CONTRAINDICATIONS, PLEASE SEE SIGVARIS PRESCRIPTION PAD. SIGVARIS ACCEPTS NO LIABILITY FOR NON-OBSERVANCE OF CONTRAINDICATIONS OR CAUTIONS.

COMPRESSION LEVELS

30–40mmHg	40–50mmHg	50–60mmHg
-----------	-----------	-----------

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL				MEDIUM				LARGE				X-LARGE (CALF ONLY)			
	7.5"–9" (19–23cm)				9"–12" (23–30cm)				12"–15.5" (30–40cm)				15.5"–18" (40–46cm)			
CALF CIRCUMFERENCE	AVERAGE		FULL		AVERAGE		FULL		AVERAGE		FULL		AVERAGE		FULL	
	11"–15" (28–38cm)		≤17" (≤43cm)		12.5"–17" (32–43cm)		≤19" (≤48cm)		14.5"–19" (36–48cm)		≤22.5" (≤57cm)		18"–22.5" (45–57cm)		≤27" (≤68cm)	
CALF CIRCUMFERENCE FOR 50S—FULL ONLY	505C: FULL CALF ONLY		12.5"–17" (32–43cm)		505C: FULL CALF ONLY		14.5"–19" (36–48cm)		505C: FULL CALF ONLY		18"–22.5" (45–57cm)		505C NOT AVAILABLE IN XL		505C NOT AVAILABLE IN XL	
CALF LENGTH	SHORT	LONG	SHORT	LONG	SHORT	LONG	SHORT	LONG	SHORT	LONG	SHORT	LONG	SHORT	LONG	SHORT	LONG
	≤15.5" (40cm)	>15.5" (40cm)	≤15.5" (40cm)	>15.5" (40cm)	≤15.5" (40cm)	>15.5" (40cm)	≤15.5" (40cm)	>15.5" (40cm)	≤15.5" (40cm)	>15.5" (40cm)	≤15.5" (40cm)	>15.5" (40cm)	≤15.5" (40cm)	>15.5" (40cm)	≤15.5" (40cm)	>15.5" (40cm)
THIGH CIRCUMFERENCE	AVERAGE		FULL		AVERAGE		FULL		AVERAGE		FULL		CUSTOM SIZES AVAILABLE. PLEASE CONTACT CUSTOMER CARE FOR DETAILS			
	16"–21" (41–54cm)		≤24" (≤61cm)		19"–24" (48–61cm)		≤28" (≤71cm)		23"–28" (58–71cm)		≤32" (≤81cm)					
THIGH CIRCUMFERENCE FOR 50S—FULL ONLY	505T: FULL THIGH ONLY		≤24" (≤61cm)		505T: FULL THIGH ONLY		≤28" (≤71cm)		505T: FULL THIGH ONLY		≤32" (≤81cm)					
LEG LENGTH	SHORT	LONG	SHORT	LONG	SHORT	LONG	SHORT	LONG	SHORT	LONG	SHORT	LONG				
	≤28" (71cm)	>28" (71cm)	≤28" (71cm)	>28" (71cm)	≤28" (71cm)	>28" (71cm)	≤28" (71cm)	>28" (71cm)	≤28" (71cm)	>28" (71cm)	≤28" (71cm)	>28" (71cm)				

SIZE S1 S2 S3 S4 M1 M2 M3 M4 L1 L2 L3 L4 X1 X2 X3 X4

GUIDELINES FOR HIP CIRCUMFERENCES: SMALL FULL <48" (<122cm), MEDIUM FULL <54" (<137cm), LARGE FULL <61" (<155cm)

Product Code Example: For product at the top left of page: 503NS2077	Series	mmHg Level	Product Style	Size	Gender/ Open Toe	Color
	50	3	N	S2	O	77

ADVANCE ARMSLEEVE



FEATURES

Developed in collaboration with patients and clinicians to create a comfortable and effective armsleeve

Odor control and silver ion properties maintain skin hygiene

Moisture-wicking properties to keep skin dry

Sensinno silicone top-band keeps armband in place all day

BEST SUITED FOR

Patients with arm lymphedema

Post-mastectomy care

Patients with previous armsleeve compliance issues due to lack of comfort

FIBER CONTENT

50% Nylon

28% Rayon

22% Spandex

(content can vary slightly by style)

Latex-Free

BEIGE
(77)



Silicone band grip-top

15–20mmHg Armsleeve Without Gauntlet
Small Long in Beige
911ASLO77/S



Silicone band grip-top

20–30mmHg Armsleeve With Gauntlet
Small Long in Beige
912ASLG77/S

AVAILABLE STYLES AND COLOR


	Beige (77)
WITH GAUNTLET (G)	●
WITHOUT GAUNTLET (O)	●

COMPRESSION LEVELS

15–20mmHg*	20–30mmHg	30–40mmHg
------------	-----------	-----------

*15–20mmHg NOT AVAILABLE WITH GAUNTLET AND NOT AVAILABLE IN PLUS SIZING

SIZING MEASUREMENTS



SIZING CHART

PALM (a)	SMALL 6.5"–8" (17–21cm)		MEDIUM 7"–8.5" (18–22cm)				LARGE 7.5"–9" (19–23cm)				EXTRA-LARGE 8"–9.5" (20–24cm)			
WRIST (c)	SMALL 5.5"–6.5" (14–16cm)		MEDIUM 6.5"–7" (16–18cm)				LARGE 7"–8" (18–20cm)				EXTRA-LARGE 8"–9" (20–23cm)			
UPPER MID-ARM (f)	REGULAR 10"–13.5" (26–34cm)		REGULAR 11.5"–14.5" (29–37cm)		PLUS 14.5"–17" (37–43cm)		REGULAR 12.5"–15.5" (32–40cm)		PLUS 15.5"–18" (40–46cm)		REGULAR 13.5"–16.5" (34–42cm)		PLUS 16.5"–19" (42–48cm)	
LENGTH (C–G)	REGULAR ≤17" (≤43cm)	LONG >17" (>43cm)	REG. ≤17" (≤43cm)	LONG >17" (>43cm)	REG. ≤17" (≤43cm)	LONG >17" (>43cm)	REG. ≤17" (≤43cm)	LONG >17" (>43cm)	REG. ≤17" (≤43cm)	LONG >17" (>43cm)	REG. ≤17" (≤43cm)	LONG >17" (>43cm)	REG. ≤17" (≤43cm)	LONG >17" (>43cm)
SIZE	SR	SL	MR	ML	MR/P	ML/P	LR	LL	LR/P	LL/P	XR	XL	XR/P	XL/P

Product Code Example: For product at the top left of page: 911ASLO77/S	Series	mmHg Level	Product Style	Size	Gauntlet	Color	Grip-top
	91	1	A	SL	O	77	/S

Product Code Example: For 20–30mmHg Armsleeve with Gauntlet Large Regular Plus Size in Beige: 912ALRG77/PS	Series	mmHg Level	Product Style	Size	Gauntlet	Color	Plus Size	Grip-top
	91	2	A	LR	G	77	/P	S

Note: all armsleeves come with Grip-top (/S)



IMPROVE LEG HEALTH WITH A SIMPLE CHANGE OF SOCKS


For many people, experiencing tired or achy legs is the first step in learning more about the health benefits of graduated compression.

The SIGVARIS WELL BEING
& CONSUMER HEALTH
line introduces many of
the therapeutic benefits of
graduated compression with
stylish options in a simplified
sizing system.



CASUAL COTTON

for women & men



FEATURES

Made with premium Supima Cotton for optimum softness and durability

Flat, low-friction toe seam to prevent irritation

Soft, stretchable knit-in band prevents pinching

BEST SUITED FOR

Tired, heavy, restless, aching legs

Travelers, occupational leg health risks, expecting mothers, family history of chronic vein disease

FIBER CONTENT


66% Supima Cotton

26% Nylon

8% Spandex

Latex-Free

SUPIMA®
WORLD'S FINEST COTTONS



146: WOMEN'S CALF STYLE



15–20mmHg Closed Toe Calf
Size B in Khaki
146CB30

186: MEN'S CALF STYLE



15–20mmHg Closed Toe Calf
Size C in Black
186CC99

AVAILABLE STYLE AND COLORS

Closed Toe = ●	White (00)	Navy (10)	Brown (11)	Khaki (30)	Black (99)
CALF (C)	●	●	●	●	●

COMPRESSION LEVEL

15–20mmHg

SIZING CHART

WOMEN'S SHOE	5–7	7.5–9.5	10–12
ANKLE CIRCUMFERENCE	7"–9" <i>(18–23cm)</i>	7.5"–10" <i>(19–25cm)</i>	8"–10.5" <i>(20–27cm)</i>
MEN'S SHOE	6.5–8.5	9–11	11.5–14
ANKLE CIRCUMFERENCE	7.5"–10" <i>(19–25cm)</i>	8"–10.5" <i>(20–27cm)</i>	8.5"–11.5" <i>(21–29cm)</i>
SIZE	A	B	C

Product Code Example: For product at the top left of page: 146CB30	Series	Product Style	Size	Color
	146	C	B	30

CUSHIONED COTTON

for women & men

FEATURES

- Cushioned sole, heel and toe area to protect feet
- Flat, low-friction toe seam
- Soft, stretchable knit-in band prevents pinching
- Ideal for everyday wear and for sensitive skin
- Suitable for patients with diabetes
- No compression in heel and instep for ease of donning

BEST SUITED FOR

- Tired, heavy, restless, aching legs
- Travelers, occupational leg health risks, expecting mothers, family history of chronic vein disease

FIBER CONTENT

- 68% Supima Cotton
- 25% Nylon
- 7% Spandex
- Latex-Free*

SUPIMA®
WORLD'S FINEST COTTONS



WHITE
(00)

BLACK
(99)

SEA ISLAND COTTON

for women & men

FEATURES

Knitted from the world's rarest cotton

Sea Island Cotton yarns are lightweight and have great strength, a silky luster, and an incredibly soft feel

Breathable and comfortable for sensitive skin

Ideal for travel and everyday wear

BEST SUITED FOR

Tired, heavy, restless, aching legs

Travelers, occupational leg health risks, expecting mothers, family history of chronic vein disease

FIBER CONTENT

43% Sea Island Cotton

39% Nylon

11% Spandex

7% Supima Cotton

Latex-Free

NAVY
(10)

BROWN
(11)

BLACK
(99)

ZURICH
COLLECTION

151: WOMEN'S CALF STYLE



15–20mmHg Closed Toe Calf
Size A in Brown
151CA11

191: MEN'S CALF STYLE



15–20mmHg Closed Toe Calf
Size B in Navy
191CB10

AVAILABLE STYLE AND COLORS

Closed Toe = ●	Navy (10)	Brown (11)	Black (99)
CALF (C)	●	●	●

COMPRESSION LEVEL

15–20mmHg
ALSO AVAILABLE IN 20–30mmHg, SEE PAGE 11

SIZING CHART

WOMEN'S SHOE	5–7	7.5–9.5	10–12
ANKLE CIRCUMFERENCE	7"–9" <i>(18–23cm)</i>	7.5"–10" <i>(19–25cm)</i>	8"–10.5" <i>(20–27cm)</i>
MEN'S SHOE	6.5–8.5	9–11	11.5–14
ANKLE CIRCUMFERENCE	7.5"–10" <i>(19–25cm)</i>	8"–10.5" <i>(20–27cm)</i>	8.5"–11.5" <i>(21–29cm)</i>
SIZE	A	B	C

Product Code Example: For product at the top left of page: 151CA11	Series	Product Style	Size	Color
	151	C	A	11

ALL-SEASON MERINO WOOL

for women & men



FEATURES

Soft, comfortable and itch-free

Natural properties of wool make these socks odor resistant

Thermoregulating for year-round wear and keeps feet dry

BEST SUITED FOR

Tired, heavy, restless, aching legs

Travelers, occupational leg health risks, expecting mothers, family history of chronic vein disease

FIBER CONTENT

53% Extra-fine Merino Wool

38% Nylon

9% Spandex

Latex-Free

NAVY
(10)

BROWN
(11)

BLACK
(99)

ZURICH
COLLECTION

152: WOMEN'S CALF STYLE



15–20mmHg Closed Toe Calf
Size A in Brown
152CA11

192: MEN'S CALF STYLE



15–20mmHg Closed Toe Calf
Size C in Navy
192CC10

AVAILABLE STYLE AND COLORS

Closed Toe = ●	Navy (10)	Brown (11)	Black (99)
CALF (C)	●	●	●

COMPRESSION LEVEL

15–20mmHg
ALSO AVAILABLE IN 20–30mmHg, SEE PAGE 13

SIZING CHART

WOMEN'S SHOE	5–7	7.5–9.5	10–12
ANKLE CIRCUMFERENCE	7"–9" <i>(18–23cm)</i>	7.5"–10" <i>(19–25cm)</i>	8"–10.5" <i>(20–27cm)</i>
MEN'S SHOE	6.5–8.5	9–11	11.5–14
ANKLE CIRCUMFERENCE	7.5"–10" <i>(19–25cm)</i>	8"–10.5" <i>(20–27cm)</i>	8.5"–11.5" <i>(21–29cm)</i>
SIZE	A	B	C

Product Code Example: For product at the top left of page: 152CA11	Series	Product Style	Size	Color
	152	C	A	11

MICROFIBER SHADES

for women & men

FEATURES

- Advanced microfiber technology for ultimate comfort
- Glides on and off with ease
- Excellent breathability and moisture management
- Soft and durable sock for everyday wear

BEST SUITED FOR

- Tired, heavy, restless, aching legs
- Travelers, occupational leg health risks, expecting mothers, family history of chronic vein disease

FIBER CONTENT

- 87% Nylon
- 13% Spandex
- Latex-Free*



Actual color may vary

143: WOMEN'S CALF STYLE

183: MEN'S CALF STYLE



15–20mmHg Closed Toe Calf
Size B in Purple Argyle
143CB46



15–20mmHg Closed Toe Calf
Size C in Graphite Stripe
183CC17

AVAILABLE STYLE AND COLORS

Closed Toe = ●	Dark Navy Mini-Stripe (03)	Graphite Mini-Stripe (07)	Dark Navy Stripe (13)	Pink Stripe (15)	Graphite Stripe (17)	Onyx Stripe (19)
CALF (C)	●	●	●	●	●	●
Closed Toe = ●	Navy Heather (23)	Graphite Heather (27)	Pink Argyle (45)	Purple Argyle (46)	Graphite Argyle (47)	Royal Blue Argyle (48)
CALF (C)	●	●	●	●	●	●

COMPRESSION LEVEL

15–20mmHg

ALSO AVAILABLE IN 20–30mmHg,
SEE PAGE 26

SIZING CHART

WOMEN'S SHOE	5–7	7.5–9.5	10–12
ANKLE CIRCUMFERENCE	7"–9" <i>(18–23cm)</i>	7.5"–10" <i>(19–25cm)</i>	8"–10.5" <i>(20–27cm)</i>
MEN'S SHOE	6.5–8.5	9–11	11.5–14
ANKLE CIRCUMFERENCE	7.5"–10" <i>(19–25cm)</i>	8"–10.5" <i>(20–27cm)</i>	8.5"–11.5" <i>(21–29cm)</i>
SIZE	A	B	C

<p>Product Code Example: For product at the top left of page: 143CB46</p>	Series	Product Style	Size	Color
	143	C	B	46

SHEER FASHION

for women



FEATURES

- Fabric feels soft against the skin
- Incredible sheer design in a variety of fashionable colors
- Sheer toe and heel area
- Comfortable stretch is excellent for everyday wear
- Cotton-lined gusset on pantyhose styles

BEST SUITED FOR

- Tired, heavy, restless, aching legs
- Travelers, occupational leg health risks, expecting mothers, family history of chronic vein disease

FIBER CONTENT

- 79% Nylon
- 21% Spandex
- (content can vary slightly by style)
- Latex-Free*



CALF AND THIGH SIZING CHART

WOMEN'S SHOE	5-7	7.5-9.5	10-12
ANKLE CIRCUMFERENCE	7.5"-9" (19-23cm)	8"-10" (20.5-25cm)	8.5"-10.5" (22-27cm)
SIZE	A	B	C

WEIGHT GUIDELINES FOR THIGH STYLE:
 SIZE A: <135LBS., SIZE B: <180LBS., SIZE C: <230LBS.



15–20mmHg Closed Toe Calf
Size B in Taupe
120CB29



15–20mmHg Closed Toe Pantyhose
Size C in Black
120PC99

AVAILABLE STYLES AND COLORS

	Charcoal (12)	Taupe (29)	Natural (33)	Suntan (36)	Café (73)	Black (99)
Closed Toe = ●						
Open Toe = ○						
CALF (C)	●	○	○	○	●	●
THIGH (N)	●	○	○	○	●	●
PANTYHOSE (P)	●	○	○	○	●	●
MATERNITY PANTYHOSE (M)			○			●

COMPRESSION LEVEL

15–20mmHg

PANTYHOSE AND MATERNITY SIZING CHART

Height	Lbs/Kg	100 / 45	105 / 48	110 / 50	115 / 52	120 / 54	125 / 57	130 / 59	135 / 61	140 / 64	145 / 66	150 / 68	155 / 70	160 / 73	165 / 75	170 / 77	175 / 80	180 / 82	185 / 84	190 / 86	195 / 89	200 / 91	205 / 93	210 / 95	215 / 98	220 / 100	225 / 102	230 / 105	
4'9" / 145cm																													
4'10" / 147cm																													
4'11" / 150cm																													
5'0" / 152cm																													
5'1" / 155cm																													
5'2" / 157cm																													
5'3" / 160cm																													
5'4" / 162cm																													
5'5" / 165cm																													
5'6" / 167cm																													
5'7" / 170cm																													
5'8" / 173cm																													
5'9" / 175cm																													
5'10" / 178cm																													
5'11" / 180cm																													
6' / 183cm																													

NOTE: PLEASE USE PRE-PREGNANCY WEIGHT FOR MATERNITY PANTYHOSE.

* MATERNITY PANTYHOSE NOT AVAILABLE IN SIZE F.

Product Code Example: For product at the top left of page: 120CB29	Series	Product Style	Size	Color
	120	C	B	29

Product Code Example: For 15–20mmHg Open Toe Calf Size C in Suntan: 120CC036	Series	Product Style	Size	Open Toe	Color
	120	C	C	O	36

SOFT SILHOUETTE LEGGINGS

for women



FEATURES

Stylish, subtle ribbed design looks chic and pairs perfectly with dresses, blouses, tunics and shirts

Foldable, comfort-stretch waist band provides a seamless look for all-day wearing comfort

Advanced soft stretch microfiber fabric shapes legs and smooths skin

BEST SUITED FOR

Tired, heavy, restless, aching legs

Occupational leg health risks

Family history of chronic vein disease

FIBER CONTENT

83% Nylon

17% Spandex

Latex-Free

BLACK
(99)



15–20mmHg Legging
Size B in Black
170LB99



Non-binding ankle cuff stays
comfortably in place for all day wear


**AVAILABLE STYLE
AND COLOR**

Open Foot = ●	Black (99)
LEGGING (L)	●

COMPRESSION LEVEL

15–20mmHg

SIZING CHART



GUIDELINES FOR ANKLE CIRCUMFERENCE
SIZE A: 7"–9", SIZE B & C: 8"–10", SIZE D: 8.5"–11.5"

Product Code Example: For product at the top left of page: 170LB99	Series	Product Style	Size	Color
	170	L	B	99

BUSINESS CASUAL

for men



FEATURES

Stretchable toe fabric and heel design increase comfort and durability

Classic ribbed look for everyday wear

Wide, non-binding top band for comfortable fit

Soft, microfiber inlay yarns

BEST SUITED FOR

Tired, heavy, restless, aching legs

Travelers, occupational leg health risks, family history of chronic vein disease

FIBER CONTENT

88% Nylon

12% Spandex

Latex-Free

NAVY
(10)

BROWN
(11)

CHARCOAL
(12)

BLACK
(99)



15–20mmHg Closed Toe Calf
Size B in Charcoal
189CB12



15–20mmHg Closed Toe Calf
Size C in Brown
189CC11

AVAILABLE STYLE AND COLORS

Closed Toe = ●	Navy (10)	Brown (11)	Charcoal (12)	Black (99)
CALF (C)	●	●	●	●

COMPRESSION LEVEL

15–20mmHg

SIZING CHART

MEN'S SHOE	6.5–8.5	9–11	11.5–14
ANKLE CIRCUMFERENCE	7.5"–10" <i>(19–25cm)</i>	8"–10.5" <i>(20–27cm)</i>	8.5"–11.5" <i>(21–29cm)</i>
SIZE	A	B	C

Product Code Example: For product at the top left of page: 189CB12	Series	Product Style	Size	Color
	189	C	B	12

ATHLETIC RECOVERY SOCKS

for women & men

FEATURES

- Helps reduce lactic acid buildup
- Comfortable and easy to put on and take off
- drirelease* yarns provides moisture wicking and odor control
- No compression in heel and instep for ease of donning

BEST SUITED FOR

- Tired, heavy, restless, aching legs
- Post-workout recovery
- Mild athletic activity
- Walking

FIBER CONTENT

- 67% *drirelease* Polyester
(57% Polyester, 10% Cotton)
- 26% Nylon
- 7% Spandex
- Latex-Free*





15–20mmHg Calf
Size XX-Large in Black
401C299



15–20mmHg Calf
Size Medium in White
401CM00

AVAILABLE STYLE AND COLORS

Closed Toe = ●	White (00)	Black (99)
CALF (C)	●	●

COMPRESSION LEVEL

15–20mmHg

SIZING CHART

WOMEN'S SHOE SIZE	<7.5	7.5–9.5	10–12	12.5–15	≥15.5
MEN'S SHOE SIZE	<6.5	6.5–8.5	9–11	11.5–14	14.5–17
ANKLE CIRCUMFERENCE	7"–9.5" (18–24cm)	7.5"–10" (19–26cm)	8"–10.5" (20–27cm)	8.5"–11.5" (21–29cm)	8.5"–11.5" (21–29cm)
SIZE	S	M	L	X	2

Product Code Example: For product at the top left of page: 401C299	Series	mmHg Level	Product Style	Size	Color
	40	1	C	2	99

PERFORMANCE SLEEVES

for women & men



FEATURES

Fabric resists deterioration in salt, chlorine or fresh water

Stabilizes and diminishes vibrations of active muscles

Quick drying

High-visibility reflective transfer

BEST SUITED FOR

Running, triathlons, endurance sports and gym workouts

Team sports

FIBER CONTENT

75% Nylon

25% Spandex

Latex-Free





20–30mmHg Calf Sleeve
Size Large in Blood Orange
412VL63



20–30mmHg Calf Sleeve
Size Medium in Steel Blue
412VM62

AVAILABLE STYLE AND COLORS

	White (00)	Blue (50)	Lime (54)	Pink (56)*	Purple (58)*
CALF SLEEVE (V)	●	●	●	●	●
	Red (59)	Steel Blue (62)	Blood Orange (63)	Limeade (64)	Black (99)
CALF SLEEVE (V)	●	●	●	●	●

* PINK AND PURPLE NOT AVAILABLE IN X-LARGE

COMPRESSION LEVEL

20–30mmHg

SIZING CHART

ANKLE CIRCUMFERENCE	7"–8.5" (18–21.5cm)	8.5"–10" (21.5–25.5cm)	10"–11.5" (25.5–29cm)	11.5"–13" (29–33cm)
CALF CIRCUMFERENCE	11"–15.5" (28–39cm)	12.5"–17.5" (31.5–45cm)	14"–20" (35.5–51cm)	16"–23" (40.5–58cm)
SIZE	S	M	L	X

Product Code Example: For product at the top left of page: 412VL63	Series	mmHg Level	Product Style	Size	Color
	41	2	V	L	63

MERINO OUTDOOR SOCKS

for women & men



FEATURES

Soft, comfortable and itch-free

Natural properties of wool make these socks odor-resistant

Mesh flex zones for comfort and ventilation

Cushioned, padded soles protect against blisters

No compression in heel and instep for ease of donning

BEST SUITED FOR

Tired, heavy, restless, aching legs

Hiking

Fishing & hunting

Post-workout recovery

Outdoor sports and walking

FIBER CONTENT

67% Fine Merino Wool

26% Nylon

7% Spandex

Latex-Free

CHARCOAL
(12)

OLIVE
(31)



15–20mmHg Calf
Size Large in Olive
421CL31



15–20mmHg Calf
Size Medium in Charcoal
421CM12

AVAILABLE STYLE AND COLORS

Closed Toe = ●	Charcoal (12)	Olive (31)
CALF (C)	●	●

COMPRESSION LEVEL

15–20mmHg
ALSO AVAILABLE IN 20–30mmHg, SEE PAGE 32

SIZING CHART

MEN'S SHOE SIZE	<6.5	6.5–8.5	9–11	11.5–14
WOMEN'S SHOE SIZE	<7.5	7.5–9.5	10–12	12.5–15
ANKLE CIRCUMFERENCE	7"–9.5" (18–24cm)	7.5"–10" (19–26cm)	8"–10.5" (20–27cm)	8.5"–11.5" (21–29cm)
SIZE	S	M	L	X

Product Code Example: For product at the top left of page: 421CL31	Series	mmHg Level	Product Style	Size	Color
	42	1	C	L	31

EVERSOFT DIABETIC SOCKS

for women & men



FEATURES

Thick, padded soles, heels, and toes for comfort and protection

Flat, low friction toe seam

Soft, stretchable knit-in band prevents pinching

Compression from ankle to knee, compression-free foot

BEST SUITED FOR

Individuals with diabetes

Individuals with limited dexterity or lower tolerance for compression

Individuals with sensitive skin

Individuals wanting less compression in the foot

FIBER CONTENT

67% drirelease (57% Polyester, 10% Cotton)

26% Nylon

7% Spandex

Latex-Free

WHITE
(00)



8–15mmHg Closed Toe Calf
Size Small in White
160CS00



8–15mmHg Closed Toe Calf
Size Large in White
160CL00

AVAILABLE STYLE AND COLOR

Closed Toe = ●	White (00)
	CALF (C) ●

COMPRESSION LEVEL

8–15mmHg
ALSO AVAILABLE: DIABETIC COMPRESSION SOCKS IN 18–25mmHg, SEE PAGE 34

SIZING CHART

WOMEN'S SHOE	<7½	7.5–9.5	10–12	12.5–15
MEN'S SHOE	<6½	6.5–8.5	9–11	11.5–14
ANKLE CIRCUMFERENCE	7"–9½" (18–24cm)	7.5"–10" (19–26cm)	8"–10.5" (20–27cm)	8.5"–11.5" (21–29cm)
SIZE	S	M	L	X

Product Code Example: For product at the top left of page: 160CS00	Series	Product Style	Size	Color
	160	C	S	00

ACCESSORIES

for women & men



Rubber Gloves – Ridged

The practical and protective way to apply and remove your compression stockings.

Package: 12 per pack

SIZE	PRODUCT CODE
SMALL	591R400S
MEDIUM	591R400M
LARGE	591R400L
X-LARGE	591R400X

Package: Single pair

SIZE	PRODUCT CODE
SMALL	591RPRS
MEDIUM	591RPRM
LARGE	591RPRL
X-LARGE	591RPRX



Latex-Free Gloves

Gripping donning gloves recommended for use in a latex-free environment.

Package: 12 per pack

SIZE	PRODUCT CODE
SMALL	592R400S
MEDIUM	592R400M
LARGE	592R400L
X-LARGE	592R400X

Package: Single pair

SIZE	PRODUCT CODE
SMALL	592RPRS
MEDIUM	592RPRM
LARGE	592RPRL
X-LARGE	592RPRX

GLOVE SIZING CHART

MEASURE CIRCUMFERENCE OF THE HAND AT THE PALM

SMALL	MEDIUM	LARGE	X-LARGE
6.5"–7" (16.5–18cm)	7.5"–8" (19–20.5cm)	8.5"–9" (21.5–23cm)	9.5"–10" (24–25.5cm)



ARION Sim-Slide®

SIZE	PRODUCT CODE
SMALL	587SLIDESM
MEDIUM	587SLIDEMD
LARGE	587SLIDELG
X-LARGE	587SLIDEXL

SIZING CHART

FOR USE WITH OPEN-TOE STYLES

WOMEN'S SHOE SIZE	≤7	7½–9½	10–13	≥13½
MEN'S SHOE SIZE	≤5½	6–8	8½–11½	≥12
ARION SIM-SLIDE SIZE	S	M	L	XL



Video: How to Use the ARION Sim-Slide

To watch the video go to: bit.ly/SIGssyt



ARION Magnide On/Off®

SIZE	PRODUCT CODE
MEDIUM	587MAGMD
LARGE	587MAGLG
X-LARGE	587MAGXL

SIZING CHART

FOR USE WITH CLOSED-TOE STYLES

WOMEN'S SHOE SIZE	≤8½	9–12	≥12½
MEN'S SHOE SIZE	≤7.5	8–11	≥11½
ARION MAGNIDE SIZE	M	L	XL*

* XL ALSO SUITABLE FOR EXTRA WIDE FEET

Video: How to Use the ARION Magnide On/Off
To watch the video go to: bit.ly/SIGmofyt



For more SIGVARIS videos, including donning tips and product demonstrations, visit our YouTube channel: bit.ly/SIGytch



Washing Solution

Gentle washing solution specifically formulated to help maintain the quality of SIGVARIS compression socks and stockings. Also great for fine washables!

Package Options

CONTAINS	PRODUCT CODE
SIX 16 OZ. BOTTLES (473mL)	586W900
ONE 16 OZ. BOTTLE (473mL)	586W950
TWELVE 2 OZ. BOTTLES (60mL)	586W975
ONE 2 OZ. BOTTLE (60mL)	586W985



It Stays®

Roll-on water-soluble liquid adhesive to keep stockings in place.

Package Options

CONTAINS	PRODUCT CODE
TWELVE 2 OZ. BOTTLES (60mL)	581I300
ONE 2 OZ. BOTTLE (60mL)	581I350



SIMON

A donning and doffing aid so versatile that it can easily be stowed to accompany you wherever you go.

SIZE	CALF CIRCUMFERENCE	PRODUCT CODE
X-SMALL	25–29cm	589SXS
SMALL	30–38cm	589SSM
MEDIUM	38–46cm	589SMD
LARGE	46–54cm	589SLG



Video: *How to Use the SIMON*

To watch the video go to: bit.ly/SIGSIMONC



Cotton Underliner

10–15mmHg stocking. Can be used as an under layer to protect fragile skin, secure wound dressings, add 10–15mmHg of compression, or create an easy glide surface to aide in putting on and taking off compression stockings.

For sizing, please refer to the NATURAL RUBBER Calf sizing on page 39.

Package: 2 per pack

SIZE	PRODUCT CODE
SMALL/SHORT	588USS
SMALL/LONG	588USL
MEDIUM / SHORT	588UMS
MEDIUM / LONG	588UML
LARGE-XL / SHORT	588ULS
LARGE-XL / LONG	588ULL

ACCESSORIES

for women & men



SOS: Slip-On-SIGVARIS

The easy way to put on compression stockings.

SIZE	PRODUCT CODE
ONE SIZE	585S500



SIGVARIS Doff N' Donner

A soft device that allows you to quickly guide a compression sock or stocking up and down your leg or arm, making it simple to put on and remove all styles of graduated compression. It is perfect for both self-donning and caregiver use.

SIZE	PRODUCT CODE
ONE SIZE	583D100



SIGVARIS Doff N' Donner Cone

This loading device makes it easy to roll stockings onto the SIGVARIS Doff N' Donner. The bottom of the base firmly attaches to most clean, flat surfaces via suction.

SIZE	PRODUCT CODE
ONE SIZE	583D200




Foot Slip

Comfortable applicator for use with open toe stockings.

Package: 6 per pack

SIZE	PRODUCT CODE
ONE SIZE	580F200



Video: How to Use the SIGVARIS Doff N' Donner
To watch the video go to:
bit.ly/SIGdndyt

SYMPTOMS & INDICATIONS

SYMPTOMS / INDICATIONS	COMPRESSION LEVEL AT THE ANKLE
<ul style="list-style-type: none"> ▪ Tired, achy legs ▪ Occupational or evening edema ▪ Leg discomfort from long hours of standing or sitting ▪ Prevent varicosities & swelling during pregnancy ▪ Preventative for legs predisposed to risk ▪ Reduce swelling during travel ▪ Prevent DVT during travel 	15–20mmHg
<ul style="list-style-type: none"> ▪ Diabetics who experience edema (see Cautions below) 	18–25mmHg
<ul style="list-style-type: none"> ▪ Heavy, fatigued, aching legs ▪ Mild edema in lower extremities ▪ Mild varicosities with minimal edema ▪ Mild varicosities during pregnancy ▪ Lymphedema upper extremities ▪ Post-procedure of small veins ▪ Prevention or management of DVT 	20–30mmHg
<ul style="list-style-type: none"> ▪ Moderate varicosities with mild to moderate edema (also during pregnancy) ▪ Post-procedure of larger veins to maintain treatment success ▪ CEAP C3 (edema), C4 (skin changes without ulceration), C5 (skin changes with healed ulceration), C6 (skin changes with active ulcerations) ▪ Management and treatment of DVT or Post-Thrombotic Syndrome ▪ Superficial phlebitis ▪ Lymphedema after decongestant therapy to maintain reduction ▪ Orthostatic / postural hypotension ▪ Post-fracture, post-traumatic edema <p data-bbox="165 1172 677 1225"><i>Note: A prescription is recommended for stockings with the indications of 30mmHg or higher</i></p>	30–40mmHg
<ul style="list-style-type: none"> ▪ Lymphedema after decongestant therapy to maintain reduction ▪ Severe Post-Thrombotic Syndrome ▪ Severe varicosities and/or edema 	40mmHg OR HIGHER
<ul style="list-style-type: none"> ▪ Severe arterial insufficiency, intermittent claudication, ischemia ▪ Uncontrolled congestive heart failure 	CONTRAINDICATIONS
<ul style="list-style-type: none"> ▪ Skin sensitivities or allergies ▪ Advanced neuropathy with sensory loss in leg or foot ▪ Diabetes with arterial compromise ▪ Confinement to bed or non-ambulatory use unless otherwise prescribed by the physician ▪ Eczema, hypodermatitis, dermatosclerosis, dermatitis 	CAUTIONS

The therapeutic effect of graduated compression stockings is based on the compression level. With the proper CEAP classification, the correct compression level at the ankle can be chosen. For more severe symptoms, consult with a vascular specialist.

No liability accepted for non-observance of contraindications and cautions.

Fronek HS, Bergan JJ, et al. The Fundamentals of Phlebology: Venous Disease for Clinicians. 2004. pg 151. Blaettner, et al. Leg Symptoms of Healthy People and Their Treatment with Compression Hosiery. Phlebology 2008. Compression therapy of the extremities. Partsch, Rabe, Stemmer. ISBN: 2.85480.770.7.


GLOSSARY OF TERMS

TERM	DEFINITION	SIGVARIS USAGE
CHRONIC VENOUS DISEASE (CVD)	A disease state comprised of disorders ranging from tired, achy legs, varicose veins, skin issues to venous leg ulcers (open and healed).	SIGVARIS manufactures MEDICAL and WELL BEING products to help reduce the risk of, improve symptoms of, and provide treatment for CVD.
DOUBLE-COVERED YARN	Spandex fibers are covered twice with yarn to prevent bare spandex from touching the skin for improved wearing comfort.	Used in nearly all SIGVARIS graduated compression socks and hosiery to create products that are easy to put on and take off, while also providing greater comfort.
DEEP VEIN THROMBOSIS (DVT)	A blood clot in a deep vein, most commonly in the lower leg or thigh.	SIGVARIS cares about DVT prevention. Clinical evidence supports the use of graduated compression to help prevent DVTs and the more serious complication of PE (Pulmonary Embolism), which can be life-threatening.
GRADUATED COMPRESSION VS. ANTI-EMBOLISM	Graduated compression is for ambulatory patients, whereas anti-embolism stockings are for non-mobile patients.	SIGVARIS offers true graduated compression, meaning products are tightest at the ankle and decrease in pressure going up the leg to help improve circulation.
HYPOALLERGENIC	Designed to minimize the chance of an allergic reaction.	SIGVARIS offers latex-free products and cotton fabrics for those with sensitive skin and allergies. Also, many SIGVARIS accessories are designed for consumers with sensitive skin.
ISO CERTIFIED	SIGVARIS is ISO (International Organization for Standardization) certified, indicating we have a quality system in place that meets the requirements to manufacture Class II medical devices.	SIGVARIS follows ISO Standard 13485:2003 for our products manufactured in North and South America.
MICROFIBER	An extremely fine synthetic fiber knitted into textiles with the texture and drape of natural-fiber fabric. This fabric has excellent breathability.	Used to create products with a soft and comfortable feel. Find this fiber in SOFT OPAQUE, MICROFIBER SHADES, MIDTOWN MICROFIBER and BUSINESS CASUAL.
mmHg (MILLIMETERS OF MERCURY)	Defined as the extra pressure generated by a column of mercury one millimeter high; the same measurement used for blood pressure.	SIGVARIS products are offered in several ranges of compression levels based on this measurement.
MERINO WOOL	Merino is the softest and finest wool available in the world and is inherently itch-free.	SIGVARIS uses this premium wool yarn to create soft products that are thermoregulating for use in all seasons. Used in MERINO OUTDOOR SOCKS, MERINO OUTDOOR PERFORMANCE SOCKS and ALL-SEASON MERINO WOOL.
SEA ISLAND COTTON	Sea Island Cotton is the rarest cotton in the world and represents the finest in cotton softness.	SIGVARIS uses the highest-quality yarns in the industry to create products that are best in class. Look for this fine quality cotton in our SEA ISLAND COTTON SOCKS.
SENSINNOV®	A patented SIGVARIS silicone top band that gently grips the skin and keeps the band in place all day.	This top band is featured in EVERSHEER, SOFT OPAQUE, ALLURE and MIDTOWN MICROFIBER to help keep thigh-high stockings in place all day long.
Supima COTTON SUPIMA WORLD'S FINEST COTTONS	Supima accounts for only about three percent of annual cotton production in the United States. Its fineness and longer staple length make it a premium cotton fiber.	SIGVARIS uses this premium cotton yarn in the COTTON, SEA ISLAND COTTON, CASUAL COTTON and CUSHIONED COTTON products, making them comfortable and ideal for those with sensitive skin.



The Worldwide Leader
in Compression Therapy

Please turn this Catalog over to see our
Inelastic Wraps and Bandaging Systems
for Lymphedema, Wound Care and
Complex Edemas



Please turn this Catalog over to see our
Graduated Compression Socks & Hosiery
for Chronic Venous Disease



The Worldwide Leader
in Compression Therapy

SIGVARIS COMPRESSION THERAPY FOR ANY CONDITION



Lymphedema, venous disease, and phlebolymphe­dema are just a few of the conditions that SIGVARIS compression garments can help manage.

Clinicians, dealers and patients have all contributed to the design and features of every SIGVARIS product.

The result is a range of products to address many types of edema, but with a single purpose: to achieve an improved quality of life.



SIGVARIS

SIGVARIS, INC. 103055 Riley Street Suite 30 Holland, MI 49424

SIGVARIS, COMPREFLEX, and COMPREFIT are registered trademarks of SIGVARIS AG, CH-9014 St.Gallen/Switzerland, in many countries worldwide. © Copyright 2018 SIGVARIS, Inc. Breathe-O-Prene is a registered trademark of Accumed Innovative Technologies, LLC. DURAPORE is a registered trademark of 3M. All rights reserved. drirelease is copyright ©2018 drirelease. drirelease is an Optimizer Brands company. All Rights Reserved.

UPPER EXTREMITY

	SERIES	PAGE	STAGE 0-I*	STAGE II*	STAGE III*
ARM	CHIPSLEEVE	6	●	●	
	COMPRESLEEVE	8	●	●	
	COMPREFLEX	10	●	●	
	MEDAFIT	12		●	●

	SERIES	PAGE	STAGE 0-I*	STAGE II*
FACE	NECK/MANDIBLE PAD	18	●	●
	COMPVEST	20	●	●
TORSO	CHIPVEST CUSTOM	22	●	●
	COMPRESHORTS	24	●	●

HAND	SERIES	PAGE
	MEDAHAND	14
	MEDAGLOVE	14
	DORSAL POCKET GLOVE	14

LOWER EXTREMITY

	SERIES	PAGE	STAGE 0-I*	STAGE II*	STAGE III*	WOUND CARE	CVI
LEG	CHIPSLEEVE BK/TH	28/30	●	●			
	COMPRESLEEVE BK	32	●	●			
	COMPREFIT BK	34	●	●		●	●
	COMPREFIT NF	36	●	●		●	●
	COMPREFIT PLUS	38	●	●		●	●
	COMPREFIT TC	40	●	●			●
	COMPREFLEX BK	42	●	●		●	●
	COMPREFLEX NF	44	●	●		●	●
	COMPREFLEX LITE TRANSITION	46	●	●			●
	COMPREFLEX LITE COMPLETE	48	●	●			●
	COMPREFLEX TC	50	●	●			●
	COOLFLEX NF	52	●	●			●
	MEDAFIT BK/TH	54/56		●	●		

	SERIES	PAGE
FOOT	COMPBOOT LITE	58
	COMPBOOT	58
	COMPBOOT PLUS	58

BK = BELOW KNEE
 TC = THIGH COMPONENT
 TH = THIGH HIGH
 *LYMPHEDEMA

SIGVARIS ACCESSORIES & THERAPY SUPPLIES

ACCESSORY	PAGE
SOCKS & LINERS	62
STRAP EXTENDERS	63
COMPRESKNEE	63
CHIP PADS	64
GENIFIT CHIP PAD	65

THERAPY SUPPLIES	PAGE
COMPREFLEX REDUCE BK	66
COMPREFLEX REDUCE BOOT	66
COMPREFLEX REDUCE TC	68
COMPREFLEX REDUCE KNEE	68
BANDAGE KITS	70
BANDAGES	72
GAUZE & TAPE	72
STOCKINETTE	72
PADDING & FOAM	73

SIGVARIS CUSTOM

SERIES	PAGE
CUSTOM CHIPSLEEVE ARM	76
ARMASSIST	78
MEDABOOT	81
LEGASSIST BELOW KNEE	82
LEGASSIST BELOW KNEE WITH SHELF	84
LEGASSIST LOBULE COMPRESSION SYSTEM	87
LEGASSIST THIGH HIGH	90



COMFORTABLE SOLUTIONS FOR BOTH DAY & NIGHT

SIGVARIS arm compression garments are ideal for treating all stages of lymphedema. The garments are comfortable and easy to use, promoting patient compliance and improved quality of life.



SIGVARIS hand garments effectively treat stage I and II lymphedema, post-operative edema and other swelling issues.

The DORSAL POCKET GLOVE treats areas that are difficult to compress, including the full length of individual fingers, and features finger loops for doffing assistance.

CHIPSLEEVE ARM



Shown with PINK OVERSLEEVE

FEATURES

15–25mmHg of mild, gradient compression from hand to axilla when used with the OVERSLEEVE (included)

Open-cell foam chips massage fibrotic tissue, stimulating fluid movement and helping to soften the superficial tissue

Use as a bandage liner or independent compression garment

Four-way stretch fabric contains moisture-wicking, odor control properties

Safe for day or nighttime use

Custom size available (see page 76)

Easy care - machine wash and dry

Latex-free

INDICATIONS

Stage I or II Lymphedema

Mild fibrosis

Post-surgical edema

INCLUDES

CHIPSLEEVE Arm

OVERSLEEVE (5–10mmHg)

Pair cotton liners




FOAM CHIPS

Open-cell foam chips massage superficial lymphatics

DIRECTIONAL CHANNELS


Channel formations facilitate efficient lymph drainage



LEOPARD



BLACK



TIE DYE



PINK



BLUE HORIZON



NAVY



INCLUDES OVERSLEEVE

Six different OVERSLEEVE patterns available! The CHIPSLEEVE provides 15–25mmHg of mild, gradient compression when used with the OVERSLEEVE

P.O. #	Account #	Contact	Date
--------	-----------	---------	------

Phone	Patient Name	Age	Height	Weight
-------	--------------	-----	--------	--------

Company Name

Dropship Name & Address

PRODUCT INFORMATION


Arm	Size	Length	Item #
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT			

OVERSLEEVE Color

BLACK NAVY PINK TIE DYE BLUE HORIZON LEOPARD **NOTE:** If no color is selected, BLACK will be included.

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

Measuring instructions on page 93



CHIPSLEEVE ARM SIZING CHART (A–C)

	SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
C	23–32cm	28–37cm	33–43cm	39–49cm	44–54cm
B1	22–30cm	26–34cm	30–39cm	35–44cm	39–49cm
B	20–27cm	24–31cm	28–35cm	32–39cm	36–43cm
A	14–17cm	15–18cm	16–19cm	18–20cm	19–22cm
LEFT					
SHORT	2631-ARS-L	2632-ARS-L	2633-ARS-L	2634-ARS-L	2635-ARS-L
REGULAR	2631-ARR-L	2632-ARR-L	2633-ARR-L	2634-ARR-L	2635-ARR-L
LONG	2631-ARL-L	2632-ARL-L	2633-ARL-L	2634-ARL-L	2635-ARL-L
RIGHT					
SHORT	2631-ARS-R	2632-ARS-R	2633-ARS-R	2634-ARS-R	2635-ARS-R
REGULAR	2631-ARR-R	2632-ARR-R	2633-ARR-R	2634-ARR-R	2635-ARR-R
LONG	2631-ARL-R	2632-ARL-R	2633-ARL-R	2634-ARL-R	2635-ARL-R

Custom sizes and configurations are available for CHIPSLEEVE. For more information, please see page 76 or speak with your SIGVARIS Sales Representative or Customer Care: 1-800-322-7744.

COMPRESLEEVE ARM



FEATURES

- 20–30mmHg of moderate, gradient compression from hand to axilla
- Designed to replace day or nighttime bandaging
- Unique lycra / cotton blend creates comfortable, short-stretch compression
- WAVEFOAM creates pressure ridges to facilitate fluid drainage
- Attached hand piece provides additional palm compression
- Easy care - machine wash and dry
- Latex-free

INDICATIONS

- Stage I or II Lymphedema
- Mild fibrosis
- Post-surgical edema

INCLUDES

- COMPRESLEEVE Arm
- Donning Loop
- Pair cotton liners



HIGHLY ADJUSTABLE
Garment secures with adjustable tabs that create a snug, custom-like fit



WAVEFOAM
The inner WAVEFOAM creates high / low pressure channels to promote fluid movement



ATTACHED HAND PIECE
Features an attached, adjustable hand piece for added palm compression


P.O. #	Account #	Contact	Date		
Phone	Patient Name		Age	Height	Weight
Company Name					
Dropship Name & Address					

PRODUCT INFORMATION

Arm	Size	Length	Item #
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT			

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

Measuring instructions on page 93



COMPRESLEEVE ARM SIZING CHART (A-C)

	SMALL	MEDIUM	LARGE	XLARGE
C	23–32cm	28–37cm	33–43cm	39–49cm
B1	22–30cm	26–34cm	30–39cm	35–44cm
B	20–27cm	24–31cm	28–35cm	32–39cm
A	14–17cm	15–18cm	16–19cm	18–20cm
LEFT				
SHORT	1501-ARS-L	1502-ARS-L	1503-ARS-L	1504-ARS-L
REGULAR	1501-ARR-L	1502-ARR-L	1503-ARR-L	1504-ARR-L
LONG	1501-ARL-L	1502-ARL-L	1503-ARL-L	1504-ARL-L
RIGHT				
SHORT	1501-ARS-R	1502-ARS-R	1503-ARS-R	1504-ARS-R
REGULAR	1501-ARR-R	1502-ARR-R	1503-ARR-R	1504-ARR-R
LONG	1501-ARL-R	1502-ARL-R	1503-ARL-R	1504-ARL-R

COMPREFLEX ARM



FEATURES

- Provides mild/moderate compression
- Unique elbow design for 90 degree range of motion
- Soft, conforming Breathe-O-Prene® provides a comfortable fit and wicks away moisture
- Easy to adjust straps
- Available in black or beige
- Easy care - machine wash and dry
- Latex-free

INDICATIONS


- Stage I or II Lymphedema
- Post-surgical edema

INCLUDES

- COMPREFLEX® Arm
- Pair cotton liners

OPTIONAL HAND PIECE

DORSAL POCKET GLOVE may be paired with COMPREFLEX Arm. See page 14 for additional information and order form



COMPREFLEX Arm features an inner sleeve for easy donning and doffing


P.O. #	Account #	Contact	Date		
Phone	Patient Name		Age	Height	Weight
Company Name					
Dropship Name & Address					

PRODUCT INFORMATION

Arm	Size	Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

Measuring instructions on page 93



COMPREFLEX ARM SIZING CHART (A–C)

	SMALL	MEDIUM	LARGE	XLARGE
C	23–32cm	28–37cm	33–43cm	39–49cm
B1	22–30cm	26–34cm	30–39cm	35–44cm
B	20–27cm	24–31cm	28–35cm	32–39cm
A	14–19cm	15–20cm	16–21cm	18–23cm
BLACK				
SHORT	1401-ARS	1402-ARS	1403-ARS	1404-ARS
REGULAR	1401-ARR	1402-ARR	1403-ARR	1404-ARR
LONG	1401-ARL	1402-ARL	1403-ARL	1404-ARL
BEIGE				
SHORT	1411-ARS	1412-ARS	1413-ARS	1414-ARS
REGULAR	1411-ARR	1412-ARR	1413-ARR	1414-ARR
LONG	1411-ARL	1412-ARL	1413-ARL	1414-ARL

MEDAFIT ARM



FEATURES

30–40mmHg of firm, gradient compression from hand to axilla

The inner WAVEFOAM liner creates high/low pressure ridges to promote lymph drainage

Outer shell fabric is very durable, and provides firm compression

Replaceable straps and foam liner extend the garment life

Separate MEDAHAND garment constructed with chipped foam

UNISTRAP can be wrapped around the hand piece for additional compression

Safe for day or night use

Easy care - machine wash and dry

Latex-free

INDICATIONS

Stage II or III lymphedema

Limb fills quickly when compression is removed

Fibrotic tissue is present

Higher compression levels needed

INCLUDES

MEDAFIT Arm UniStrap

MEDAHAND Pair cotton liners

Donning Loop

WAVEFOAM LINER

The inner foam liner has high / low pressure ridges help facilitate lymph movement



INCLUDES MEDAHAND

The MEDAHAND provides comfortable, moderate compression

HAND GARMENT UPGRADES

Upgraded hand compression options are available at an additional cost (see page 14 for additional information)



MEDAGLOVE



DORSAL POCKET GLOVE


P.O. #	Account #	Contact	Date		
Phone	Patient Name		Age	Height	Weight
Company Name					
Dropship Name & Address					

PRODUCT INFORMATION

Arm <input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT	Hand Options
Size	<input type="checkbox"/> MEDAHAND (Included) Size
Length	<input type="checkbox"/> MEDAGLOVE (Upcharge) Size
Item #	<input type="checkbox"/> DORSAL POCKET GLOVE (Upcharge) Size

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

Measuring instructions on page 93



MEDAFIT ARM SIZING CHART (A-C)

	SMALL	MEDIUM	LARGE	XLARGE
C	23-32cm	28-37cm	33-43cm	39-49cm
B1	22-30cm	26-34cm	30-39cm	35-44cm
B	20-27cm	24-31cm	28-35cm	32-39cm
A	14-17cm	15-18cm	16-19cm	18-20cm

LEFT

SHORT	1201-ARS-L	1202-ARS-L	1203-ARS-L	1204-ARS-L
REGULAR	1201-ARR-L	1202-ARR-L	1203-ARR-L	1204-ARR-L
LONG	1201-ARL-L	1202-ARL-L	1203-ARL-L	1204-ARL-L

RIGHT

SHORT	1201-ARS-R	1202-ARS-R	1203-ARS-R	1204-ARS-R
REGULAR	1201-ARR-R	1202-ARR-R	1203-ARR-R	1204-ARR-R
LONG	1201-ARL-R	1202-ARL-R	1203-ARL-R	1204-ARL-R

HAND SIZING CHART (A-E)

	SMALL	MEDIUM	LARGE	X LARGE
E	15-19cm	19-22cm	22-26cm	26-29cm
A	14-20cm	16-22cm	18-24cm	20-26cm

MEDAHAND

RIGHT	1201-HDR-R	1202-HDR-R	1203-HDR-R	1204-HDR-R
LEFT	1201-HDR-L	1202-HDR-L	1203-HDR-L	1204-HDR-L

HAND COMPRESSION

MEDAHAND



Dorsal view

FEATURES

- Breathable / moisture wicking / anti-odor and soil-release fabric
- Open-cell foam chips gently massage tissue
- Conforming design ensures proper fit and gradient compression

MEDAGLOVE



Anterior view

FEATURES

- Breathable / moisture wicking / anti-odor and soil-release fabric
- Open-cell foam chips gently massage tissue
- Conforming design with finger separation ensures proper fit
- Effective as a compression garment or bandage liner
- Gradient compression from fingers to wrist

DORSAL POCKET GLOVE



Dorsal view

FEATURES

- Dorsal pocket holds CHIP PAD Dorsum or other therapeutic insert
- Full finger compression with doffing loops for easy use
- Palm pad provides soft, mild compression
- Ultra-thin profile for full hand function during use



P.O. #	Account #	Contact	Date		
Phone	Patient Name		Age	Height	Weight

Company Name

Dropship Name & Address

PRODUCT INFORMATION

Hand Style


<input type="checkbox"/> MEDAHAND	Size	Item #
<input type="checkbox"/> MEDAGLOVE	Size	Item #
<input type="checkbox"/> DORSAL POCKET GLOVE	Size	Item #

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

CIRCUMFERENCES

Measure circumference at Palm and record on line E.

Measure at Wrist, and record on line A.



HAND SIZING CHART (A-E)

	SMALL	MEDIUM	LARGE	X LARGE
E	15-19cm	19-22cm	22-26cm	26-29cm
A	14-20cm	16-22cm	18-24cm	20-26cm

MEDAHAND

RIGHT	1201-HDR-R	1202-HDR-R	1203-HDR-R	1204-HDR-R
LEFT	1201-HDR-L	1202-HDR-L	1203-HDR-L	1204-HDR-L

MEDAGLOVE

RIGHT	1201-GLR-R	1202-GLR-R	1203-GLR-R	1204-GLR-R
LEFT	1201-GLR-L	1202-GLR-L	1203-GLR-L	1204-GLR-L

DORSAL POCKET GLOVE

RIGHT	1201-DGR-R	1202-DGR-R	1203-DGR-R	1204-DGR-R
LEFT	1201-DGR-L	1202-DGR-L	1203-DGR-L	1204-DGR-L

SIGVARIS TORSO PRODUCTS



THERAPEUTIC RESULTS PATIENTS CAN SEE & FEEL

SIGVARIS torso compression garments are designed to treat mild edema and stage I or stage II lymphedema. Swelling in the lower torso, chest and neck is often difficult to safely and effectively address. SIGVARIS offers customizable options, such as the addition of drain tube pockets and holes for a trach tube or a catheter to address a wide variety of patient needs.

NECK/MANDIBLE PAD



FEATURES

- Provides 15–25mmHg of mild compression for the neck and lower face
- Conforming design and materials ensure a comfortable fit
- Straps are easily adjusted and feature durable, medical grade material
- Easy care - machine wash and dry
- Four-way stretch fabric contains moisture-wicking, odor control properties
- Latex-free

INDICATIONS

- Stage I or II Lymphedema
- Post-surgical edema

INCLUDES

- 1 NECK/MANDIBLE PAD

ONE PIECE GARMENT

The NECK/MANDIBLE PAD is a one piece garment that provides soft, mild compression to the lower face and neck.

Straps secure at the top of the head and behind the neck to ensure a comfortable fit. The strap fabric is high quality, medical grade, and allows for some adjustability.



DIRECTIONAL CHANNELS

Channel formations facilitate efficient lymph drainage

FOAM CHIPS

Open-cell foam chips massage superficial lymphatics

P.O. #	Account #	Contact	Date		
Phone	Patient Name		Age	Height	Weight


Company Name

Dropship Name & Address


PRODUCT INFORMATION

Size Item #
 SMALL MEDIUM LARGE X-LARGE

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms



HEAD - measure circumference starting at the top of the head down to under the chin, back to top of head.



NECK - measure circumference at largest part of neck.

NECK/MANDIBLE PAD SIZING CHART

	SMALL	MEDIUM	LARGE	XLARGE
HEAD	46-56cm	52-62cm	56-66cm	58-68cm
NECK	28-33cm	33-38cm	38-43cm	43-58cm
	2231-NM	2232-NM	2233-NM	2234-NM



FEATURES

- 15–25mmHg of mild, supportive compression
- Features anterior hook adjustment for optimal fit
- Straps are easily adjusted and feature durable, medical grade material
- Interior pockets hold CHIP PAD Bra or prosthetic insert
- Easy care - machine wash and dry

Caution: elastic band contains Natural Rubber Latex which may cause allergic reactions

INDICATIONS

- Stage I or II Lymphedema
- Pitting edema, non-fibrotic
- Post-surgical edema from mastectomy, lumpectomy or reconstructive procedures

INCLUDES

- 1 COMPREVEST

ACCESSORIES (see page 64)



CHIP PADS


CHIP PAD Bra is easily inserted into the interior side pocket of the COMPREVEST

DRAIN TUBE POCKET

This accessory attaches to the inside of the waist band and holds a drain tube and pouch



Multiple size options for hook enclosures



Adjustable shoulder straps



P.O. #	Account #	Contact	Date
--------	-----------	---------	------

Phone	Patient Name	Age	Height	Weight
-------	--------------	-----	--------	--------

Company Name

Dropship Name & Address

PRODUCT INFORMATION


Color <input type="checkbox"/> BLACK <input type="checkbox"/> WHITE	Accessories (Additional cost)	
Size	<input type="checkbox"/> CHIP PAD HALF (breast to under arm, see page 52)	Color <input type="checkbox"/> BLACK <input type="checkbox"/> WHITE
Item #	<input type="checkbox"/> CHIP PAD FULL (breast to spine, see page 52)	Color <input type="checkbox"/> BLACK <input type="checkbox"/> WHITE
	<input type="checkbox"/> DRAIN TUBE POCKET (see page 52)	Color <input type="checkbox"/> BLACK <input type="checkbox"/> WHITE

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

STEP ONE

Measure the circumference at the fullest part of bust.


BUST SIZE



STEP TWO

Measure the circumference in inches directly under the bust (where a regular bra sits).

BAND SIZE



STEP THREE

Determine the **Cup Size**: subtract the **Band** from the **Bust**.

CUP SIZE

2" or less = **A/B**

2" - 4" = **C/D**

4" - 6" = **DD/DDD**

COMPREVEST SIZING CHART

	SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
BAND SIZE	28-30in	32-34in	36-38in	40-42in	44-46in

BLACK

A/B	2201-BRR-A	2202-BRR-A	2203-BRR-A	2204-BRR-A	2205-BRR-A
C/D	2201-BRR-C	2202-BRR-C	2203-BRR-C	2204-BRR-C	2205-BRR-C
DD/DDD	2201-BRR-D	2202-BRR-D	2203-BRR-D	2204-BRR-D	2205-BRR-D

WHITE

A/B	2281-BRR-A	2282-BRR-A	2283-BRR-A	2284-BRR-A	2285-BRR-A
C/D	2281-BRR-C	2282-BRR-C	2283-BRR-C	2284-BRR-C	2285-BRR-C
DD/DDD	2281-BRR-D	2282-BRR-D	2283-BRR-D	2284-BRR-D	2285-BRR-D

CHIPVEST CUSTOM



Bilateral channeling shown



Shoulder straps are adjustable

FEATURES

- Provides mild/moderate compression
- Open-cell foam chips help soften fibrotic areas
- Straps are easily adjusted and feature durable, medical grade material
- Four-way stretch fabric contains moisture-wicking, odor control properties
- Back drawstring and corset design allows patient to adjust garment easily
- Easy care - machine wash and dry

Caution: back draw string contains Natural Rubber Latex which may cause allergic reactions

INDICATIONS

- Stage I or II Lymphedema
- Mild fibrosis
- Post-surgical edema

INCLUDES

- 1 CHIPVEST


DIRECTIONAL CHANNELS

Channel formations facilitate efficient lymph drainage

FOAM CHIPS

Open-cell foam chips massage superficial lymphatics

The drawstring may be adjusted to ensure a comfortable fit



P.O. #	Account #	Contact	Date
--------	-----------	---------	------

Phone	Patient Name	Age	Height	Weight
-------	--------------	-----	--------	--------

Company Name

Dropship Name & Address

PRODUCT INFORMATION

Channeling
 FULL (Bilateral) LEFT SIDE (Unilateral) RIGHT SIDE (Unilateral)

Item # **2239-VS**


A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

STEP ONE

Measure the **Length** in centimeters from the Suprasternal Notch to the waist.

STEP TWO

Measure **Circumferences** in centimeters.



COMPRESHORT



FEATURES

10–15mmHg of supportive, mild compression for truncal / abdominal lymphedema

Soft, breathable fabric reduces moisture containment

Reinforced seams for durability

Machine wash and dry

May be paired with the GENIFIT CHIP PAD for additional genital support (see page 65)

Caution: elastic band on sizes 2XLRG and above contain Natural Rubber Latex which may cause allergic reactions

INDICATIONS

Stage I or II Lymphedema


Mild fibrosis

Post-surgical edema

INCLUDES

1 COMPRESHORTS

STANDARD SIZE



COMPRESHORTS SIZING CHART (WAIST, HIP, MIDTHIGH)

	WAIST	HIP	MIDTHIGH	ITEM NUMBER
SMALL	69–76cm	84–91cm	41–46cm	1301-CS
SMALL MAX	76–86cm	91–102cm	41–46cm	1301-CM
MEDIUM	76–86cm	91–102cm	46–53cm	1302-CS
MED. MAX	86–97cm	102–112cm	46–53cm	1302-CM
LARGE	86–97cm	102–112cm	53–61cm	1303-CS
LRG MAX	97–112cm	112–127cm	53–61cm	1303-CM
XLARGE	97–112cm	112–127cm	61–71cm	1304-CS
XLRG MAX	112–132cm	127–147cm	61–71cm	1304-CM
2XLARGE	112–132cm	127–147cm	71–81cm	1305-CS
2XLRG MAX	132–152cm	147–168cm	71–81cm	1305-CM
3XLARGE	132–152cm	147–168cm	81–91cm	1306-CS
3XLRG MAX	152–172cm	168–188cm	81–91cm	1306-CM
4XLARGE	152–172cm	168–188cm	91–101cm	1307-CS
4XLRG MAX	172–193cm	188–208cm	91–101cm	1307-CM
5XLARGE	172–193cm	188–208cm	101–112cm	1308-CS

P.O. #	Account #	Contact	Date		
Phone	Patient Name		Age	Height	Weight

Company Name

Dropship Name & Address

PRODUCT INFORMATION		
Size <input type="checkbox"/> STANDARD	Size <input type="checkbox"/> CUSTOM	
Standard Size	Custom Waist Measurement	Custom Hip Measurement
Item #	Custom Midthigh Measurement	

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

CUSTOM SIZE

Excellent choice for patients with lipedema, or who are outside of standard size ranges.

STEP ONE

Record the **WAIST**, **HIP**, and **MIDTHIGH** measurements


STEP TWO

Measure from the groin down and indicate the garment inseam length, record in the corresponding box

STEP THREE

Mark the end of the inseam and take a circumference measurement at that point, record in the corresponding box

After recording these measurements, please call Customer Care: 1-800-322-7744.





SIGVARIS below knee garments feature high quality, bacteriostatic liners. The moisture-wicking, odor control properties allow the patient to feel comfortable while wearing their compression garment.

SOLUTIONS FOR PATIENTS WITH CVI, ULCERS AND OTHER COMPLEX EDEMAS

SIGVARIS leg garments are comfortable and easy to don and doff, promoting patient compliance and improved quality of life.

The variety of the SIGVARIS compression line ranges from mild compression with chipped foam, to firm, non-elastic compression.



CHIPSLEEVE BELOW KNEE



FEATURES

15–25mmHg of mild, gradient compression from foot to below knee when used with the OVERSLEEVE (included)

Open-cell foam chips massage fibrotic tissue, stimulating fluid movement and helping to soften the superficial tissue

Use as a bandage liner or independent compression garment

Four-way stretch fabric contains moisture-wicking, odor control properties

Safe for day or nighttime use

Easy care - machine wash and dry

Latex-free

INDICATIONS

Stage I or II Lymphedema

Mild fibrosis

Post-surgical edema

INCLUDES

1 CHIPSLEEVE BK

1 OVERSLEEVE (5–10mmHg)

Pair of non-compressive BASIC liners



FOAM CHIPS

Open-cell foam chips massage superficial lymphatics

DIRECTIONAL CHANNELS

Channel formations facilitate efficient lymph drainage



BASIC LINER


Non-compressive

- Easy to apply knee-high liners
- *direrelease*® moisture wicking and odor control
- Flat seams for added comfort
- Silver ionized treatment for bacteriostatic properties



LEOPARD

BLACK



TIE DYE

PINK



BLUE HORIZON

NAVY



Shown with BLACK OVERSLEEVE

INCLUDES OVERSLEEVE

Six different OVERSLEEVE patterns available! The CHIPSLEEVE provides 15–25mmHg of mild, gradient compression when used with the OVERSLEEVE

P.O. #	Account #	Contact	Date		
Phone	Patient Name		Age	Height	Weight

Company Name

Dropship Name & Address

PRODUCT INFORMATION

Leg LEFT RIGHT Size Length Item #

Leg LEFT RIGHT Size Length Item #

OVERSLEEVE Color

BLACK NAVY PINK TIE DYE BLUE HORIZON LEOPARD **NOTE: If no color is selected, BLACK will be included.**


A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

Measuring instructions on page 93

CHIPSLEEVE - BK SIZING CHART (A-C1, NO FOOT)

	SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
C1	32-42cm	37-47cm	42-52cm	51-61cm	58-68cm
C	29-39cm	34-44cm	39-49cm	48-58cm	55-65cm
B	24-34cm	29-39cm	33-43cm	41-51cm	44-55cm
A	20-29cm	21-30cm	25-36cm	32-42cm	33-43cm
REGULAR	2631-BKR	2632-BKR	2633-BKR	2634-BKR	2635-BKR
TALL	2631-BKT	2632-BKT	2633-BKT	2634-BKT	2635-BKT

CIRCUMFERENCE



CHIPSLEEVE THIGH HIGH



FEATURES

15–25mmHg of mild, gradient compression from foot to gluteal fold when used with the OVERSLEEVE (included)

Open-cell foam chips massage fibrotic tissue, stimulating fluid movement and helping to soften the superficial tissue

Use as a bandage liner or independent compression garment

Four-way stretch fabric contains moisture-wicking, odor control properties

Safe for day or nighttime use

Easy care - machine wash and dry

Latex-free

INDICATIONS

Stage I or II Lymphedema

Mild fibrosis

Post-surgical edema

INCLUDES

1 CHIPSLEEVE TH

1 OVERSLEEVE (5–10mmHg)

Pair of non-compressive cotton liners



FOAM CHIPS

Open-cell foam chips massage superficial lymphatics

DIRECTIONAL CHANNELS

Channel formations facilitate efficient lymph drainage



NAVY



PINK



BLACK



BLUE HORIZON



TIE DYE



LEOPARD

INCLUDES OVERSLEEVE

Six different OVERSLEEVE patterns available! The CHIPSLEEVE provides 15–25mmHg of mild, gradient compression when used with the OVERSLEEVE

P.O. #	Account #	Contact	Date		
Phone	Patient Name		Age	Height	Weight

Company Name

Dropship Name & Address

PRODUCT INFORMATION

Leg	Size	Length	Item #
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT			

Leg	Size	Length	Item #
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT			

OVERSLEEVE Color

BLACK NAVY PINK TIE DYE BLUE HORIZON LEOPARD **NOTE: If no color is selected, BLACK will be included.**

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

Measuring instructions on page 93

CHIPSLEEVE - TH SIZING CHART (A-F, NO FOOT)


	SMALL	MEDIUM	LARGE	XLARGE
F	48-58cm	56-66cm	64-74cm	74-84cm
E	43-53cm	51-61cm	58-68cm	68-78cm
D	38-48cm	46-56cm	53-63cm	63-73cm
C1	32-42cm	37-47cm	42-52cm	51-61cm
C	29-39cm	34-44cm	39-49cm	48-58cm
B	24-34cm	29-39cm	33-43cm	41-51cm
A	20-29cm	21-30cm	25-36cm	32-42cm

LEFT

SHORT	2631-THS-L	2632-THS-L	2633-THS-L	2634-THS-L
REGULAR	2631-THR-L	2632-THR-L	2633-THR-L	2634-THR-L
TALL	2631-THT-L	2632-THT-L	2633-THT-L	2634-THT-L

RIGHT

SHORT	2631-THS-R	2632-THS-R	2633-THS-R	2634-THS-R
REGULAR	2631-THR-R	2632-THR-R	2633-THR-R	2634-THR-R
TALL	2631-THT-R	2632-THT-R	2633-THT-R	2634-THT-R



COMPRESLEEVE BELOW KNEE



FEATURES

- 20–30mmHg of moderate, gradient compression from foot to below knee
- Designed to replace day or nighttime bandaging
- Unique lycra / cotton blend creates comfortable, short-stretch compression
- WAVEFOAM creates pressure ridges to facilitate fluid drainage
- Easy care - machine wash and dry
- Latex-free

INDICATIONS

- Stage I or II Lymphedema
- Mild fibrosis
- Post-surgical edema

INCLUDES

- 1 COMPRESLEEVE BK
- 1 COMPREBOOT (30–40mmHg)
- Pair of non-compressive BASIC liners



WAVEFOAM

The inner WAVEFOAM creates high / low pressure channels to promote fluid movement



COMPREBOOT

Includes a separate 30–40mmHg foot piece



HIGHLY ADJUSTABLE

Garment features a zipper panel and adjustable tabs that make donning and doffing easy

P.O. # Account # Contact Date

Phone Patient Name Age Height Weight

Company Name

Dropship Name & Address

PRODUCT INFORMATION

Leg Size Length Item # CompreBoot Size

LEFT RIGHT

Leg Size Length Item # CompreBoot Size

LEFT RIGHT

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms


Measuring instructions on page 93

COMPRESLEEVE - BK SIZING CHART (A-C1)

	SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
C1	32-42cm	38-48cm	42-52cm	51-61cm	58-68cm
C	29-39cm	34-44cm	39-49cm	48-58cm	55-65cm
B	24-34cm	29-39cm	33-43cm	41-51cm	44-55cm
A	20-29cm	21-30cm	25-36cm	32-42cm	33-43cm
SHORT	1501-BKS	1502-BKS	1503-BKS	1504-BKS	1505-BKS
REGULAR	1501-BKR	1502-BKR	1503-BKR	1504-BKR	1505-BKR
TALL	1501-BKT	1502-BKT	1503-BKT	1504-BKT	1505-BKT

COMPREBOOT SIZING CHART (I-K)

	SMALL		MED / LARGE		XL / XXL	
	REGULAR	LONG	REGULAR	LONG	REGULAR	LONG
I	28-36cm	28-36cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22-26cm	22-26cm	≤30cm	≤30cm	≤33cm	≤33cm
K	14-18cm	19-23cm	18-20cm	23-25cm	20-22cm	25-27cm



COMPREFIT BELOW KNEE



FEATURES

30–40mmHg of active, inelastic, gradient compression for reduction of edema

Medicare approved for ulcers under Code A6545-AW

Ideal for treating chronic venous insufficiency or mild/moderate lymphedema

Strap extenders may be applied if the patient exceeds standard size range (see page 63)

Combine with COMPREFIT® TC for full leg compression

Easy care - machine wash and dry

Latex-free

INDICATIONS

Stage I or II Lymphedema

Chronic Venous Insufficiency

Post-surgical Edema

INCLUDES

1 COMPREFIT BK garment

1 COMPREBOOT

Pair of BASIC liners



STRAP EXTENDERS

For patients exceeding standard size ranges, use a Strap Extender to add 10cm of circumferential range to a strap

BASIC LINER

Non-compressive



- Easy to apply knee-high liners
- *drirelease* moisture wicking and odor control
- Flat seams for added comfort
- Silver ionized treatment for bacteriostatic properties

COMFORTABLE FIT

Select the open slot that provides the most conforming fit

EASY TO USE

Strap ends fold back on themselves to prevent tangling



COMPREBOOT

Provides 30–40mmHg in the foot and ankle



P.O. #	Account #	Contact	Date
Phone	Patient Name	Age	Height Weight
Company Name			
Dropship Name & Address			

PRODUCT INFORMATION

Leg	Size	Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE
Leg	Size	Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

Measuring instructions on page 93

COMPREFIT - BK REGULAR SIZING CHART (A-C1)

	SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
C1	29-39cm	34-44cm	39-49cm	48-58cm	55-65cm
C	29-39cm	34-44cm	39-49cm	48-58cm	55-65cm
B	24-34cm	29-39cm	33-43cm	41-51cm	44-55cm
A	20-29cm	21-30cm	25-36cm	32-42cm	33-43cm

BLACK

REGULAR	1101-BKR	1102-BKR	1103-BKR	1104-BKR	1105-BKR
---------	----------	----------	----------	----------	----------

BEIGE

REGULAR	1111-BKR	1112-BKR	1113-BKR	1114-BKR	1115-BKR
---------	----------	----------	----------	----------	----------

COMPREFIT - BK TALL SIZING CHART (A-C1)


	SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
C1	29-39cm	34-44cm	39-49cm	48-58cm	55-65cm
C	26-36cm	31-41cm	35-45cm	44-54cm	50-60cm
B	21-31cm	25-35cm	30-40cm	36-46cm	40-50cm
A	20-29cm	21-30cm	25-36cm	32-42cm	33-43cm

BLACK

TALL	1101-BKT	1102-BKT	1103-BKT	1104-BKT	1105-BKT
------	----------	----------	----------	----------	----------

BEIGE

TALL	1111-BKT	1112-BKT	1113-BKT	1114-BKT	1115-BKT
------	----------	----------	----------	----------	----------



COMPREBOOT SIZING CHART (I-K)

	SMALL		MED / LARGE		XL / XXL	
	REGULAR	LONG	REGULAR	LONG	REGULAR	LONG
I	28-36cm	28-36cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22-26cm	22-26cm	≤30cm	≤30cm	≤33cm	≤33cm
K	14-18cm	19-23cm	18-20cm	23-25cm	20-22cm	25-27cm

NOTE: If foot size doesn't match COMPREFIT size, please order the COMPREFIT NF (pg 37) with separate foot piece (pg 58).

COMPREFIT NO FOOT



FEATURES

30–40mmHg of active, inelastic, gradient compression for reduction of edema

Ideal for treating chronic venous insufficiency or mild/moderate lymphedema

Good option for patients in need of a different foot garment (size or style)

Strap extenders may be applied if the patient exceeds standard size range (see page 63)

May be combined with COMPREFIT TC

Easy care - machine wash and dry

Latex-free

INDICATIONS

Stage I or II Lymphedema

Chronic Venous Insufficiency

Post-surgical Edema

INCLUDES

1 COMPREFIT BK garment

Pair of BASIC liners

BASIC LINER Non-compressive



- Easy to apply knee-high liners
- *drirelease* moisture wicking and odor control
- Flat seams for added comfort
- Silver ionized treatment for bacteriostatic properties



STRAP EXTENDERS

For patients exceeding standard size ranges, use a Strap Extender to add 10cm of circumferential range to a strap



COMFORTABLE FIT
Select the open slot that provides the most conforming fit

EASY TO USE
Strap ends fold back on themselves to prevent tangling

P.O. #	Account #	Contact	Date		
Phone	Patient Name		Age	Height	Weight
Company Name					
Dropship Name & Address					

PRODUCT INFORMATION

Leg	Size	Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

Leg	Size	Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

NOTE: If foot garment is desired, please see page 58 for options and order form.

Measuring instructions on page 93

COMPREFIT - NF REGULAR SIZING CHART (A-C1)

	SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
C1	29-39cm	34-44cm	39-49cm	48-58cm	55-65cm
C	29-39cm	34-44cm	39-49cm	48-58cm	55-65cm
B	24-34cm	29-39cm	33-43cm	41-51cm	44-55cm
A	20-29cm	21-30cm	25-36cm	32-42cm	33-43cm

BLACK

REGULAR	1101-NF-BKR	1102-NF-BKR	1103-NF-BKR	1104-NF-BKR	1105-NF-BKR
---------	-------------	-------------	-------------	-------------	-------------

BEIGE

REGULAR	1111-NF-BKR	1112-NF-BKR	1113-NF-BKR	1114-NF-BKR	1115-NF-BKR
---------	-------------	-------------	-------------	-------------	-------------

COMPREFIT - NF TALL SIZING CHART (A-C1)


	SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
C1	29-39cm	34-44cm	39-49cm	48-58cm	55-65cm
C	26-36cm	31-41cm	35-45cm	44-54cm	50-60cm
B	21-31cm	25-35cm	30-40cm	36-46cm	40-50cm
A	20-29cm	21-30cm	25-36cm	32-42cm	33-43cm

BLACK

TALL	1101-NF-BKT	1102-NF-BKT	1103-NF-BKT	1104-NF-BKT	1105-NF-BKT
------	-------------	-------------	-------------	-------------	-------------

BEIGE

TALL	1111-NF-BKT	1112-NF-BKT	1113-NF-BKT	1114-NF-BKT	1115-NF-BKT
------	-------------	-------------	-------------	-------------	-------------



COMPREFIT PLUS



FEATURES

30–40mmHg of active, inelastic, gradient compression for reduction of edema

Used to protect healing wounds, delicate granulation tissue, and treating chronic venous insufficiency or mild/moderate lymphedema

Easy care - machine wash and dry

Latex-free

INDICATIONS

Stage I or II Lymphedema

Chronic Venous Insufficiency

Post-surgical edema

INCLUDES

1 COMPREFIT PLUS

1 Flat foam liner

1 COMPREBOOT

Pair of BASIC liners

BASIC LINER

Non-compressive



- Easy to apply knee-high liners
- *drirelease* moisture wicking and odor control
- Flat seams for added comfort
- Silver ionized treatment for bacteriostatic properties

FLAT FOAM LINER

The 1/4" dense, Open-cell flat foam liner provides a soft interface under the COMPREFIT to protect fragile skin or recently healed wounds



COMPREBOOT

Provides 30–40mmHg in the foot and ankle



P.O. # Account # Contact Date

Phone Patient Name Age Height Weight

Company Name

Dropship Name & Address

PRODUCT INFORMATION

Leg Size Length Item # CompreBoot Size
 LEFT RIGHT

Leg Size Length Item # CompreBoot Size
 LEFT RIGHT

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

Measuring instructions on page 93

COMPREFIT PLUS REGULAR SIZING CHART (A-C1)

	SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
C1	24-34cm	29-39cm	34-44cm	43-53cm	50-60cm
C	24-34cm	29-39cm	34-44cm	43-53cm	50-60cm
B	19-29cm	24-34cm	28-38cm	36-46cm	39-50cm
A	15-24cm	16-25cm	20-31cm	27-37cm	28-38cm

BLACK


REGULAR	1101-BPR	1102-BPR	1103-BPR	1104-BPR	1105-BPR
---------	----------	----------	----------	----------	----------

COMPREFIT PLUS TALL SIZING CHART (A-C1)

	SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
C1	24-34cm	29-39cm	34-44cm	43-53cm	50-60cm
C	21-31cm	26-36cm	30-40cm	39-49cm	45-55cm
B	16-26cm	20-30cm	25-35cm	31-41cm	35-45cm
A	15-24cm	16-25cm	20-31cm	27-37cm	28-38cm

BLACK

TALL	1101-BPT	1102-BPT	1103-BPT	1104-BPT	1105-BPT
------	----------	----------	----------	----------	----------



COMPREBOOT SIZING CHART (I-K)

	SMALL		MED / LARGE		XL / XXL	
	REGULAR	LONG	REGULAR	LONG	REGULAR	LONG
I	28-36cm	28-36cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22-26cm	22-26cm	≤30cm	≤30cm	≤33cm	≤33cm
K	14-18cm	19-23cm	18-20cm	23-25cm	20-22cm	25-27cm

COMPREFIT THIGH COMPONENT



FEATURES


- 30–40mmHg of active, inelastic, gradient compression for reduction of edema
- Ideal for treating chronic venous insufficiency or mild/moderate lymphedema
- 3 KNEE options available
- Strap extenders may be applied if the patient exceeds standard size range (see page 63)
- Combine with COMPREFIT BK for full leg compression
- Easy care - machine wash and dry
- Latex-free

INDICATIONS

- Stage I or II Lymphedema
- Chronic Venous Insufficiency
- Post-surgical Edema

INCLUDES

- 1 COMPREFIT TC
- 1 BELT & HIP ATTACHMENT
- Pair of Cotton thigh liners



BELT & HIP ATTACHMENT

These items are included to ensure a comfortable fit, and help secure garment in place

HIGHLY ADJUSTABLE

Straps may be trimmed as needed

KNEE OPTIONS *(additional charge)*



COMPREFIT KNEE

Closed Patella

Soft, conforming compression with straps that wrap around the back of the knee



COMPREFIT KNEE

Open Patella

Soft, conforming compression with straps that wrap around the back of the knee



COMPREFLEX KNEE

Easy, comfortable compression with straps that wrap around the front of the knee. Features inner stretch panel for simple donning and doffing

P.O. #	Account #	Contact	Date		
Phone	Patient Name		Age	Height	Weight

Company Name

Dropship Name & Address

PRODUCT INFORMATION

Leg	Size/Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT			<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE
Leg	Size/Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT			<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE
COMPREFLEX STYLE	Item #	Size	Color
<input type="checkbox"/> OPEN PATELLA <input type="checkbox"/> CLOSED PATELLA <input type="checkbox"/> COMPREFLEX KNEE			<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

Measuring instructions on page 93

COMPREFIT - TC SIZING CHART (D-F)

	SMALL	MEDIUM	LARGE	XLARGE
F	48-58cm	56-66cm	64-74cm	74-84cm
E	43-53cm	51-61cm	58-68cm	68-78cm
D	38-48cm	46-56cm	53-63cm	63-73cm
BLACK	1101-TC	1102-TC	1103-TC	1104-TC
BEIGE	1111-TC	1112-TC	1113-TC	1114-TC


COMPREFLEX SIZING CHART

	SMALL	MEDIUM	LRG / X-LRG
D	38-48	46-56	53-73
C1	29-39	34-44	39-65
CLOSED PATELLA			
BLACK	1101-KP	1102-KP	1103-KP
BEIGE	1111-KP	1112-KP	1113-KP
OPEN PATELLA			
BLACK	1101-KP-OP	1102-KP-OP	1103-KP-OP
BEIGE	1111-KP-OP	1112-KP-OP	1113-KP-OP

COMPREFLEX KNEE SIZING CHART

	SMALL	MEDIUM	LARGE	X LRG
D	38-48	46-56	53-63	58-73
C1	29-39	34-44	39-49	44-65
BLACK	1401-KP	1401-KP	1403-KP	1404-KP
BEIGE	1411-KP	1412-KP	1413-KP	1414-KP

CIRCUMFERENCE



COMPREFLEX BELOW KNEE



FEATURES

20–50mmHg of gradient, short-stretch compression

Soft, conforming Breathe-O-Prene provides a comfortable fit and wicks away moisture

Strap extenders may be applied if the patient exceeds standard size ranges (see page 63)

Combine with COMREFLEX THIGH COMPONENT for full leg compression

Available in black or beige

Easy care - machine wash and dry

Latex-free

INDICATIONS

Stage I or II Lymphedema

Chronic Venous Insufficiency

Post-surgical Edema

Lipodermatosclerosis

INCLUDES

1 COMPREFLEX BK

1 COMPREBOOT

Set of ACCUTAB devices

Pair of BASIC liners

BASIC LINER Non-compressive




- Easy to apply knee-high liners
- *drirelease* moisture wicking and odor control
- Flat seams for added comfort
- Silver ionized treatment for bacteriostatic properties



INCLUDES COMPREBOOT

- Ideal for mild-moderate edema/lymphedema
- Features neoprene material with hook and loop tabs
- 4-way stretch for firm compression around ankle
- Provides moderate, therapeutic compression



ACCUTAB

Each COMPREFLEX lower leg wrap includes a set of ACCUTAB compression measuring devices

P.O. #	Account #	Contact	Date
--------	-----------	---------	------

Phone	Patient Name	Age	Height	Weight
-------	--------------	-----	--------	--------

Company Name

Dropship Name & Address

PRODUCT INFORMATION

Leg	Size	Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

Leg	Size	Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

Measuring instructions on page 93

COMPREFLEX - BK SIZING CHART (A-C1)

	SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
C1	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm
C	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm
B	24-34cm	29-39cm	34-44cm	39-49cm	44-55cm
A	19-26cm	21-30cm	26-36cm	31-41cm	36-46cm

BLACK

REGULAR	1401-BKR	1402-BKR	1403-BKR	1404-BKR	1405-BKR
TALL	1401-BKT	1402-BKT	1403-BKT	1404-BKT	1405-BKT


BEIGE

REGULAR	1411-BKR	1412-BKR	1413-BKR	1414-BKR	1415-BKR
TALL	1411-BKT	1412-BKT	1413-BKT	1414-BKT	1415-BKT

COMPREFLEX SIZING CHART (I-K)

	SMALL		MED / LARGE		XL / XXL	
	REGULAR	LONG	REGULAR	LONG	REGULAR	LONG
I	28-36cm	28-36cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22-26cm	22-26cm	≤30cm	≤30cm	≤33cm	≤33cm
K	14-18cm	19-23cm	18-20cm	23-25cm	20-22cm	25-27cm

NOTE: If foot size doesn't match COMPREFLEX size, please order the COMPREFLEX NF (pg 45) with separate foot piece (pg 58).



COMPREFLEX NO FOOT



FEATURES

- 20–50mmHg of gradient, short-stretch compression
- Soft, conforming Breathe-O-Prene provides a comfortable fit and wicks away moisture
- Strap extenders may be applied if the patient exceeds standard size ranges (see page 63)
- Good option for patients in need of a different foot option (size or boot)
- Easily worn under most standard footwear
- Available in black or beige
- Easy care - machine wash and dry
- Latex-free

INDICATIONS

- Stage I or II Lymphedema
- Chronic Venous Insufficiency
- Post-surgical Edema

INCLUDES

- 1 COMPREFLEX BK
- Set of ACCUTAB devices
- Pair of BASIC liners

BASIC LINER Non-compressive




- Easy to apply knee-high liners
- *drirelease* moisture wicking and odor control
- Flat seams for added comfort
- Silver ionized treatment for bacteriostatic properties

CONTOURED DESIGN

Garment features a curved design to provide a conforming fit

STRETCH PANEL

The front stretch panel secures garment in place for easy application and high compliance



ACCUTAB

Each COMPREFLEX lower leg wrap includes a set of ACCUTAB compression measuring devices



EASY TO USE

Straps rolls back and attach to garment to prevent tangling

P.O. # Account # Contact Date

Phone Patient Name Age Height Weight

Company Name

Dropship Name & Address

PRODUCT INFORMATION

Leg Size Length Item # Color
 LEFT RIGHT BLACK BEIGE

Leg Size Length Item # Color
 LEFT RIGHT BLACK BEIGE


A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

NOTE: If foot garment is desired, please see page 58 for options and order form.

Measuring instructions on page 93

COMPREFLEX - NF SIZING CHART (A-C1)

	SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
C1	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm
C	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm
B	24-34cm	29-39cm	34-44cm	39-49cm	44-55cm
A	19-26cm	21-30cm	26-36cm	31-41cm	36-46cm
BLACK					
REGULAR	1401-NF-BKR	1402-NF-BKR	1403-NF-BKR	1404-NF-BKR	1405-NF-BKR
TALL	1401-NF-BKT	1402-NF-BKT	1403-NF-BKT	1404-NF-BKT	1405-NF-BKT
BEIGE					
REGULAR	1411-NF-BKR	1412-NF-BKR	1413-NF-BKR	1414-NF-BKR	1415-NF-BKR
TALL	1411-NF-BKT	1412-NF-BKT	1413-NF-BKT	1414-NF-BKT	1415-NF-BKT



COMPREFLEX LITE TRANSITION



FEATURES

20–50mmHg of gradient, short-stretch compression

Medicare approved for ulcers under Code A6545-AW

Soft, conforming Breathe-O-Prene provides a comfortable fit and wicks away moisture

Strap extenders may be applied if the patient exceeds standard size ranges (see page 63)

Easily worn under most standard footwear

Available in black or beige

Easy care - machine wash and dry

Latex-free

INDICATIONS

Stage I or II Lymphedema

Chronic Venous Insufficiency

Post-surgical Edema

INCLUDES

1 COMPREFLEX BK

Set of ACCUTAB devices

Pair of TRANSITION liners

TRANSITION LINER

10–15mmHg in Foot and Ankle




- Soft, flat seam knee-high liners
- *drirelease* moisture wicking and odor control
- Extra thick sole from toe to heel for added comfort
- Silver ionized treatment for bacteriostatic properties

CONTOURED DESIGN

Garment features a curved design to provide a conforming fit

STRETCH PANEL

The front stretch panel secures garment in place for easy application and high compliance



ACCUTAB

Each COMPREFLEX lower leg wrap includes a set of ACCUTAB compression measuring devices



EASY TO USE

Straps roll back and attach to garment to prevent tangling

P.O. #	Account #	Contact	Date
--------	-----------	---------	------

Phone	Patient Name	Age	Height	Weight
-------	--------------	-----	--------	--------

Company Name

Dropship Name & Address

PRODUCT INFORMATION

Leg	Size	Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

Leg	Size	Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

Measuring instructions on page 93

COMPREFLEX LITE TRANSITION SIZING CHART (A-C1)


	SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
C1	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm
C	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm
B	24-34cm	29-39cm	34-44cm	39-49cm	44-55cm
A	19-26cm	21-30cm	26-36cm	31-41cm	36-46cm

BLACK

REGULAR	1401-UC-BKR	1402-UC-BKR	1403-UC-BKR	1404-UC-BKR	1405-UC-BKR
TALL	1401-UC-BKT	1402-UC-BKT	1403-UC-BKT	1404-UC-BKT	1405-UC-BKT

BEIGE

REGULAR	1411-UC-BKR	1412-UC-BKR	1413-UC-BKR	1414-UC-BKR	1415-UC-BKR
TALL	1411-UC-BKT	1412-UC-BKT	1413-UC-BKT	1414-UC-BKT	1415-UC-BKT



TRANSITION LINER SIZING CHART

SIZE	ANKLE CIRC. (A)	CALF CIRC. (C)	WOMEN'S SHOE	MEN'S SHOE	PAIR (BLACK)	PAIR (BEIGE)
SMALL	19-26cm	29-39cm	<7½	<6½	1601-CK	1611-CK
MEDIUM	21-30cm	34-44cm	7.5-9.5	6.5-8.5	1602-CK	1612-CK
LARGE	26-36cm	40-50cm	10-12	9-11	1603-CK	1613-CK
X-LARGE	31-41cm	46-56cm	12.5-15	11.5-14	1604-CK	1614-CK
XX-LARGE	36-46cm	54-64cm	12.5-15	11.5-14	1605-CK	1615-CK

COMPREFLEX LITE COMPLETE



FEATURES

20–50mmHg of gradient, short-stretch compression

Medicare approved for ulcers under Code A6545-AW

Soft, conforming Breathe-O-Prene provides a comfortable fit and wicks away moisture

Strap extenders may be applied if the patient exceeds standard size ranges (see page 63)

Easily worn under most standard footwear

Available in black or beige

Easy care - machine wash and dry

Latex-free

INDICATIONS

Stage I or II Lymphedema

Chronic Venous Insufficiency

Post-surgical Edema

INCLUDES

1 COMPREFLEX BK

Set of ACCUTAB devices

Pair of COMPLETE liners

COMPLETE LINER

10–15mmHg in Foot, Ankle and Calf



- Comfortable knee-high liners
- *drirelease* moisture wicking and odor control
- Extra thick sole from toe to heel for added comfort
- Silver ionized treatment for bacteriostatic properties

CONTOURED DESIGN


Garment features a curved design to provide a conforming fit

EASY TO USE

Straps rolls back and attach to garment to prevent tangling

STRETCH PANEL

The front stretch panel secures garment in place for easy application and high compliance



ACCUTAB

Each COMPREFLEX lower leg wrap includes a set of ACCUTAB compression measuring devices



P.O. # Account # Contact Date

Phone Patient Name Age Height Weight

Company Name

Dropship Name & Address

PRODUCT INFORMATION

Leg Size Length Item # Color
 LEFT RIGHT BLACK BEIGE

Leg Size Length Item # Color
 LEFT RIGHT BLACK BEIGE

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

Measuring instructions on page 93

COMPREFLEX LITE COMPLETE SIZING CHART (A-C1)

	SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
C1	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm
C	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm
B	24-34cm	29-39cm	34-44cm	39-49cm	44-55cm
A	19-26cm	21-30cm	26-36cm	31-41cm	36-46cm

BLACK


REGULAR	1401-FC-BKR	1402-FC-BKR	1403-FC-BKR	1404-FC-BKR	1405-FC-BKR
TALL	1401-FC-BKT	1402-FC-BKT	1403-FC-BKT	1404-FC-BKT	1405-FC-BKT

BEIGE

REGULAR	1411-FC-BKR	1412-FC-BKR	1413-FC-BKR	1414-FC-BKR	1415-FC-BKR
TALL	1411-FC-BKT	1412-FC-BKT	1413-FC-BKT	1414-FC-BKT	1415-FC-BKT

CIRCUMFERENCE

LEFT RIGHT



COMPLETE LINER SIZING CHART

SIZE	ANKLE CIRC. (A)	CALF CIRC. (C)	WOMEN'S SHOE	MEN'S SHOE	PAIR (BLACK)	PAIR (BEIGE)
SMALL	19-26cm	29-39cm	<7½	<6½	1601-CK-F1	1611-CK-F1
MEDIUM	21-30cm	34-44cm	7.5-9.5	6.5-8.5	1602-CK-F1	1612-CK-F1
LARGE	26-36cm	40-50cm	10-12	9-11	1603-CK-F1	1613-CK-F1
X-LARGE	31-41cm	46-56cm	12.5-15	11.5-14	1604-CK-F1	1614-CK-F1
XX-LARGE	36-46cm	54-64cm	12.5-15	11.5-14	1605-CK-F1	1615-CK-F1

COMPREFLEX THIGH COMPONENT



FEATURES

- Includes ACCUTAB for easy, accurate compression
- Soft, conforming Breathe-O-Prene provides a comfortable fit and wicks away moisture
- Strap extenders may be applied if the patient exceeds standard size ranges (see page 63)
- Combine with COMPREFLEX BK for full leg compression
- 3 KNEE options available
- Available in black or beige
- Easy care - machine wash and dry
- Latex-free

INDICATIONS

- Stage I or II Lymphedema
- Chronic Venous Insufficiency
- Post-surgical Edema

INCLUDES

- 1 COMPREFLEX TC
- Set of ACCUTAB devices
- Belt & Hip Attachment
- Pair of non-compressive cotton liners



BELT & HIP ATTACHMENT

These items are included to ensure a comfortable fit, and help secure garment in place

EASY TO USE

Straps rolls back and attach to garment to prevent tangling

STRETCH PANEL

The front stretch panel secures garment in place for easy application and high compliance

KNEE OPTIONS *(additional charge)*



COMPREFLEX KNEE

Closed Patella

Soft, conforming compression with straps that wrap around the back of the knee



COMPREFLEX KNEE

Open Patella

Soft, conforming compression with straps that wrap around the back of the knee



COMPREFLEX KNEE

Easy, comfortable compression with straps that wrap around the front of the knee. Features inner stretch panel for simple donning and doffing

P.O. #	Account #	Contact	Date
Phone	Patient Name	Age	Height Weight

Company Name

Dropship Name & Address

PRODUCT INFORMATION

Leg	Size/Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT			<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE
Leg	Size/Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT			<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE
KNEE STYLE	Item #	Size	Color
<input type="checkbox"/> OPEN PATELLA <input type="checkbox"/> CLOSED PATELLA <input type="checkbox"/> COMPREFLEX KNEE			<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

Measuring instructions on page 93

COMPREFLEX - TC SIZING CHART (D-F)

	SMALL	MEDIUM	LARGE	XLARGE
F	48-58cm	56-66cm	64-74cm	74-84cm
E	43-53cm	51-61cm	58-68cm	68-78cm
D	38-48cm	46-56cm	53-63cm	63-73cm

LEFT

BLACK	REGULAR	1401-TCR-L	1402-TCR-L	1403-TCR-L	1404-TCR-L
	TALL	1401-TCT-L	1402-TCT-L	1403-TCT-L	1404-TCT-L

RIGHT

BEIGE	REGULAR	1411-TCR-R	1412-TCR-R	1413-TCR-R	1414-TCR-R
	TALL	1411-TCT-R	1412-TCT-R	1413-TCT-R	1414-TCT-R

COMPREFLEX KNEE SIZING CHART


	SMALL	MEDIUM	LRG / X-LRG
D	38-48	46-56	53-73
C1	29-39	34-44	39-65

CLOSED PATELLA

BLACK	1101-KP	1102-KP	1103-KP
BEIGE	1111-KP	1112-KP	1113-KP

OPEN PATELLA

BLACK	1101-KP-OP	1102-KP-OP	1103-KP-OP
BEIGE	1111-KP-OP	1112-KP-OP	1113-KP-OP



COMPREFLEX KNEE SIZING CHART

	SMALL	MEDIUM	LARGE	X LRG
D	38-48	46-56	53-63	58-73
C1	29-39	34-44	39-49	44-65
BLACK	1401-KP	1402-KP	1403-KP	1404-KP
BEIGE	1411-KP	1412-KP	1413-KP	1414-KP

COOLFLEX NO FOOT



FEATURES

- 20–40mmHg of gradient compression
- Mesh fabric for heat and moisture transfer
- Pressure measurements on straps
- Hook closures for maximum strength and durability
- One-handed side-closure system
- Latex-free

INDICATIONS

- Chronic Venous Insufficiency
- Venous Leg Ulcers
- Post-thrombotic Syndrome
- Phlebolymphe^dema
- Lymphedema
- Dependent Edema

INCLUDES

- COOLFLEX NF
- Pair of BASIC LINERS



BASIC LINER Non-compressive

- Easy to apply knee-high liners
- *drirelease* moisture wicking and odor control
- Flat seams for added comfort
- Silver ionized treatment for bacteriostatic properties

FOOT COMPRESSION OPTIONS



Foot garment options on page 60

Additional liner options are on page 64



EASY DONNING

Front stretch panel secures garment in place for easy one-handed application

BREATHABLE MESH

Mesh fabric has superior breathability for a cool and comfortable fit

ACCURATE PRESSURE

Straps feature easy-to-use pressure indicators for 20–30mmHg or 30–40mmHg of compression

P.O. #	Account #	Contact	Date		
Phone	Patient Name		Age	Height	Weight
Company Name					
Dropship Name & Address					

PRODUCT INFORMATION

Leg	Size	Length	Item #
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT			


Leg	Size	Length	Item #
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT			

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

Measuring instructions on page 93

COOLFLEX - NF REGULAR SIZING CHART (A-C1)

	SMALL	SMALL MAX	MEDIUM / LARGE	MED / LRG MAX	XLARGE / XXLARGE
C1	28-39cm	35-50cm	34-50cm	45-64cm	45-64cm
C	28-39cm	35-50cm	34-50cm	45-64cm	45-64cm
B	23-33cm	27-38cm	29-43cm	35-49cm	37-54cm
A	16-26cm	16-26cm	23-33cm	23-33cm	30-42cm
LEFT					
REGULAR	2101-NF-BKR-L	2101-NF-BMR-L	2103-NF-BKR-L	2103-NF-BMR-L	2105-NF-BKR-L
TALL	2101-NF-BKT-L	2101-NF-BMT-L	2103-NF-BKT-L	2103-NF-BMT-L	2105-NF-BKT-L
RIGHT					
REGULAR	2101-NF-BKR-R	2101-NF-BMR-R	2103-NF-BKR-R	2103-NF-BMR-R	2105-NF-BKR-R
TALL	2101-NF-BKT-R	2101-NF-BMT-R	2103-NF-BKT-R	2103-NF-BMT-R	2105-NF-BKT-R



MEDAFIT BELOW KNEE



FEATURES

- 30–40mmHg of gradient, inelastic compression
- Highly adjustable compression system
- Outer shell fabric is very durable, and provides inelastic compression
- Inner WAVEFOAM liner creates high-low pressure ridges to promote lymph drainage
- COMPREBOOT PLUS features chipped foam liner
- Replaceable straps (call Customer Care for details)
- Easy care - machine wash and dry
- 1 year manufacturers warranty
- Latex-free

INDICATIONS

- Stage II or III chronic lymphedema
- Chronic Venous Insufficiency
- Post-surgical edema

INCLUDES

- 1 MEDAFIT BK
- 1 COMPREBOOT PLUS
- Pair of BASIC liners



WAVEFOAM LINER

The inner foam liner is constructed with WAVEFOAM, the high / low pressure ridges help facilitate lymph movement

BASIC LINER Non-compressive



- Easy to apply knee-high liners
- *drirelease* moisture wicking and odor control
- Flat seams for added comfort
- Silver ionized treatment for bacteriostatic properties



INCLUDES COMPREBOOT PLUS

Firm compression with inner chipped foam liner, adjustable straps and gripper patches for added traction

P.O. #	Account #	Contact	Date		
--------	-----------	---------	------	--	--

Phone	Patient Name	Age	Height	Weight
-------	--------------	-----	--------	--------

Company Name

Dropship Name & Address

PRODUCT INFORMATION

Leg	Size	Length	Item #	COMPREBOOT PLUS Size	Item #
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT					

Leg	Size	Length	Item #	COMPREBOOT PLUS Size	Item #
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT					

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

Measuring instructions on page 93

MEDAFIT - BK SIZING CHART (A-C1)

	SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
C1	32-42cm	38-48cm	42-52cm	51-61cm	58-68cm
C	29-39cm	34-44cm	39-49cm	48-58cm	55-65cm
B	24-34cm	29-39cm	33-43cm	41-51cm	44-55cm
A	20-29cm	21-30cm	25-36cm	32-42cm	33-43cm

LEFT

SHORT	1201-BKS-L	1202-BKS-L	1203-BKS-L	1204-BKS-L	1205-BKS-L
REGULAR	1201-BKR-L	1202-BKR-L	1203-BKR-L	1204-BKR-L	1205-BKR-L
TALL	1201-BKT-L	1202-BKT-L	1203-BKT-L	1204-BKT-L	1205-BKT-L

RIGHT

SHORT	1201-BKS-R	1202-BKS-R	1203-BKS-R	1204-BKS-R	1205-BKS-R
REGULAR	1201-BKR-R	1202-BKR-R	1203-BKR-R	1204-BKR-R	1205-BKR-R
TALL	1201-BKT-R	1202-BKT-R	1203-BKT-R	1204-BKT-R	1205-BKT-R

COMPREBOOT PLUS SIZING CHART (I-K)

	SMALL		MED / LARGE		XL / XXL	
	REGULAR	LONG	REGULAR	LONG	REGULAR	LONG
I	28-36cm	28-36cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22-26cm	22-26cm	≤30cm	≤30cm	≤33cm	≤33cm
K	14-18cm	19-23cm	18-20cm	23-25cm	20-22cm	25-27cm
BLACK	1801-BPR	1801-BPL	1803-BPR	1803-BPL	1805-BPR	1805-BPL

CIRCUMFERENCE

LEFT RIGHT

_____ **C1**

_____ **C**

_____ **B**


_____ **A**

_____ **Ø POINT**

_____ **I**

_____ **J**

_____ **K:** Measure **length** on medial side of foot from heel to 1st metatarsal head



MEDAFIT THIGH HIGH



FEATURES

- 30–40mmHg of gradient, inelastic compression
- Highly adjustable compression system
- Outer shell fabric is very durable, and provides inelastic compression
- Inner WAVEFOAM liner creates high-low pressure ridges to promote lymph drainage
- COMPREBOOT PLUS features chipped foam liner
- Replaceable straps (call Customer Care for details)
- Easy care - machine wash and dry
- 1 year manufacturers warranty
- Latex-free

INDICATIONS

- Stage II or III chronic lymphedema
- Chronic Venous Insufficiency
- Post-surgical edema

INCLUDES

- 1 MEDAFIT TH
- 1 COMPREBOOT PLUS
- Pair of cotton liners



INCLUDES COMPREBOOT PLUS

Firm compression with inner chipped foam liner, adjustable straps and gripper patches for added traction



WAVEFOAM LINER

The inner foam liner is constructed with WaveFoam. High / low pressure ridges help facilitate lymph movement

P.O. #	Account #	Contact	Date		
--------	-----------	---------	------	--	--

Phone	Patient Name	Age	Height	Weight
-------	--------------	-----	--------	--------

Company Name

Dropship Name & Address

PRODUCT INFORMATION

Leg	Size	Length	Item #	COMPREBOOT PLUS Size	Item #
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT					

Leg	Size	Length	Item #	COMPREBOOT PLUS Size	Item #
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT					

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

Measuring instructions on page 93

MEDAFIT - TH SIZING CHART (A-F)

	SMALL	SM. MAX	MEDIUM	MED. MAX	LARGE	LRG. MAX	XLRG
F	48-58cm	56-66cm	56-66cm	64-74cm	64-74cm	74-84cm	74-84cm
E	43-53cm	51-61cm	51-61cm	58-68cm	58-68cm	68-78cm	68-78cm
D	38-48cm	46-56cm	46-56cm	53-63cm	53-63cm	63-73cm	63-73cm
C1	32-42cm		37-47cm		42-52cm		51-61cm
C	29-39cm		34-44cm		39-49cm		48-58cm
B	24-34cm		29-39cm		33-43cm		41-51cm
A	20-29cm		21-30cm		25-36cm		32-42cm

LEFT


SHORT	1201-THS-L	1201-TMS-L	1202-THS-L	1202-TMS-L	1203-THS-L	1203-TMS-L	1204-THS-L
REGULAR	1201-THR-L	1201-TMR-L	1202-THR-L	1202-TMR-L	1203-THR-L	1203-TMR-L	1204-THR-L
TALL	1201-THT-L	1201-TMT-L	1202-THT-L	1202-TMT-L	1203-THT-L	1203-TMT-L	1204-THT-L

RIGHT

SHORT	1201-THS-R	1201-TMS-R	1202-THS-R	1202-TMS-R	1203-THS-R	1203-TMS-R	1204-THS-R
REGULAR	1201-THR-R	1201-TMR-R	1202-THR-R	1202-TMR-R	1203-THR-R	1203-TMR-R	1204-THR-R
TALL	1201-THT-R	1201-TMT-R	1202-THT-R	1202-TMT-R	1203-THT-R	1203-TMT-R	1204-THT-R

COMPREBOOT PLUS SIZING CHART (I-K)

	SMALL		MED / LARGE		XL / XXL	
	REGULAR	LONG	REGULAR	LONG	REGULAR	LONG
I	28-36cm	28-36cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22-26cm	22-26cm	≤30cm	≤30cm	≤33cm	≤33cm
K	14-18cm	19-23cm	18-20cm	23-25cm	20-22cm	25-27cm
BLACK	1801-BPR	1801-BPL	1803-BPR	1803-BPL	1805-BPR	1805-BPL



FOOT COMPRESSION

COMPRESSBOOT LITE



10–15mmHg

- Ideal for mild edema/lymphedema
- Ultra-thin profile that easily fits into most footwear
- Highly adjustable design (may be trimmed to fit)
- Breathable stretch fabric conforms to ankle
- Provides mild, therapeutic compression

COMPRESSBOOT LITE PRODUCT CODES

	SMALL/ MEDIUM / LARGE	XLRG / XX LRG
REGULAR	1803-BLR	1805-BLR
LONG	1803-BLL	1805-BLL

COMPRESSBOOT



30–40mmHg

- Ideal for mild-moderate edema/lymphedema
- Features neoprene material with hook and loop tabs
- 4-way stretch for firm compression around ankle
- Provides moderate, therapeutic compression

COMPRESSBOOT PRODUCT CODES

	SMALL	MEDIUM / LARGE	XLRG / XX LRG
BLACK			
REGULAR	1801-BTR	1803-BTR	1805-BTR
LONG	1801-BTL	1803-BTL	1805-BTL
BEIGE			
REGULAR	1811-BTR	1813-BTR	1815-BTR
LONG	1811-BTL	1813-BTL	1815-BTL

COMPRESSBOOT PLUS



30–40mmHg

- Ideal for chronic/acute edema/lymphedema
- Chipped foam lining
- Removable gripper patches provide traction
- Straps are removable for custom positioning or replacement
- Provides moderate, therapeutic compression

COMPRESSBOOT PLUS PRODUCT CODES

	SMALL	MEDIUM / LARGE	XLRG / XX LRG
REGULAR	1801-BPR	1803-BPR	1805-BPR
LONG	1801-BPL	1803-BPL	1805-BPL

P.O. #	Account #	Contact	Date
--------	-----------	---------	------

Phone	Patient Name	Age	Height	Weight
-------	--------------	-----	--------	--------

Company Name

Dropship Name & Address

PRODUCT INFORMATION

Foot <input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT	
Size	<input type="checkbox"/> COMPREBOOT LITE
Length	<input type="checkbox"/> COMPREBOOT
Item #	<input type="checkbox"/> COMPREBOOT PLUS


A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

MEASURING INSTRUCTIONS

- I:** Measure circumference of ankle bend and heel.
- J:** Measure circumference across metatarsal heads.
- K:** Measure length on the medial side of the foot from the heel to the 1st metatarsal head.

CIRCUMFERENCE

LEFT	RIGHT	
_____	_____	I
_____	_____	J
_____	_____	K



FOOT SIZING CHART (I-K)

	SMALL		MED / LARGE		XL / XXL	
	REGULAR	LONG	REGULAR	LONG	REGULAR	LONG
I	28-36cm	28-36cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22-26cm	22-26cm	≤30cm	≤30cm	≤33cm	≤33cm
K	14-18cm	19-23cm	18-20cm	23-25cm	20-22cm	25-27cm

PATIENT OPTIONS FOR EASE OF USE



SIGVARIS Accessories such as strap extenders, gripper patches, knee pieces and other accessories allow patients to maximize the value of their garments.

The CHIP PAD line features open-cell chipped foam sewn in channels, making the pads ideal for massaging fibrotic tissue and breaking up pockets of lymph fluid.



SIGVARIS therapy supplies include BIAFORM, a premium, 100% cotton short-stretch bandage with superior conforming characteristics, ring-spun cotton tubular stockinette for exceptional softness and stretch, and a range of padding options including BIASOFT™ fleece padding for comfort and durability.



ACCESSORIES & LINERS



COTTON SOCKS & LINERS

Soft, seamless cotton liners are non-compressive and can be worn under all compression products.

STYLE	SIZE	1 PACK	3 PACK	10 PACK
BK	S-XXL	1795-BK1 (1 Pair)	1795-BK3 (3 Pairs)	1795-BKB (10 Pairs)

STYLE	SIZE	
TH	S / M / L	1793-TH (One liner per pack)
TH	XL	1794-TH (One liner per pack)
ARM	S-XL	1794-AR (One liner per pack)



BASIC LINER

Non-compressive

- Easy to apply knee-high liners
- *drirelease* moisture wicking and odor control
- Flat seams for added comfort
- Silver ionized treatment for bacteriostatic properties

Size: One size fits all

PAIR	10 PAIR PACK
1700-NK1	1700-NKC



TRANSITION LINER

10–15mmHg in Foot & Ankle

- Soft, flat seam knee-high liners
- *drirelease* moisture wicking and odor control
- Extra thick sole from toe to heel for added comfort
- Silver ionized treatment for bacteriostatic properties

TRANSITION LINER SIZING CHART

SIZE	ANKLE CIRC. (A)	CALF CIRC. (C)	WOMEN'S SHOE	MEN'S SHOE	PAIR (BLACK)	PAIR (BEIGE)
SMALL	19–26cm	29–39cm	<7½	<6½	1601-CK	1611-CK
MEDIUM	21–30cm	34–44cm	7.5–9.5	6.5–8.5	1602-CK	1612-CK
LARGE	26–36cm	40–50cm	10–12	9–11	1603-CK	1613-CK
X-LARGE	31–41cm	46–56cm	12.5–15	11.5–14	1604-CK	1614-CK
XX-LARGE	36–46cm	54–64cm	12.5–15	11.5–14	1605-CK	1615-CK



COMPLETE LINER

10–15mmHg in Foot, Ankle & Calf

- Comfortable knee-high liners
- *drirelease* moisture wicking and odor control
- Extra thick sole from toe to heel for added comfort
- Silver ionized treatment for bacteriostatic properties



ANKLE PAD

Moisture wicking mesh lining. Open-cell foam chips.

SIZE	EACH
S / M / L	2233-AP
XL / XXL	2235-AP

COMPLETE LINER SIZING CHART

SIZE	ANKLE CIRC. (A)	CALF CIRC. (C)	WOMEN'S SHOE	MEN'S SHOE	PAIR (BLACK)	PAIR (BEIGE)
SMALL	19–26cm	29–39cm	<7½	<6½	1601-CK-F1	1611-CK-F1
MEDIUM	21–30cm	34–44cm	7.5–9.5	6.5–8.5	1602-CK-F1	1612-CK-F1
LARGE	26–36cm	40–50cm	10–12	9–11	1603-CK-F1	1613-CK-F1
X-LARGE	31–41cm	46–56cm	12.5–15	11.5–14	1604-CK-F1	1614-CK-F1
XX-LARGE	36–46cm	54–64cm	12.5–15	11.5–14	1605-CK-F1	1615-CK-F1



UNISTRAP

Neoprene wrap for added compression around hand.

EACH
1902



COMPREKNEE

Compression garment with either open or closed patella.

	SMALL	MEDIUM	LRG / X-LRG
D	38-48	46-56	53-73
C1	29-39	34-44	39-65
CLOSED PATELLA			
BLACK	1101-KP	1102-KP	1103-KP
BEIGE	1111-KP	1112-KP	1113-KP
OPEN PATELLA			
BLACK	1101-KP	1102-KP	1103-KP
BEIGE	1111-KP	1112-KP	1113-KP



COMPREFLEX KNEE

Garment with straps that wrap in front of the patella. Stretch panel secures garment in place.

	SMALL	MEDIUM	LARGE	X LRG
D	38-48	46-56	53-63	58-73
C1	29-39	34-44	39-49	44-65
BLACK	1401-KP	1402-KP	1403-KP	1404-KP
BEIGE	1411-KP	1412-KP	1413-KP	1414-KP



COMPREFIT STRAP EXTENDERS

Adds 10 cm of adjustable range to a single strap.
1 set extends 1 strap.

	1 SET	10 SET PACK
BLACK	1100-SE1	1100-SEB
BEIGE	1110-SE1	1110-SEB



GRIPPER PATCHES

Rubber material attaches to foot piece and provides added traction.

1 SET

1800-GP



COMPREFLEX STRAP EXTENDER

Adds 10 cm of adjustable range to a single strap.

	1 SET	10 SET PACK
BLACK	1401-SE1	1401-SEB
BEIGE	1411-SE1	1411-SEB



ACCUTAB

Compression measuring device design to use with COMPREFLEX leg wraps.

SIZE 1 SET

SMALL	1491-TBI
MEDIUM	1492-TBI
LARGE	1493-TBI
X-LARGE	1494-TBI
XX-LARGE	1495-TBI



FOAM LINER

High density, open-cell foam with breathable, moisture-wicking mesh lining. *Liner not included.*

	SIZE	REGULAR	TALL
SMALL		1101-FLR	1101-FLT
MEDIUM		1102-FLR	1102-FLT
LARGE		1103-FLR	1103-FLT
X-LARGE		1104-FLR	1104-FLT
XX-LARGE		1105-FLR	1105-FLT



COMPREVEST DRAIN TUBE POCKET

This accessory attaches to the inside of the waist band and holds a drain tube and pouch.

COLOR 2 PER PK.

BLACK	2200-DP
WHITE	2280-DP

ACCESSORIES & CHIP PADS



CHIP PAD DORSUM

Use on dorsum of hand or foot. Helps with softening scars from carpal tunnel surgery, or underneath compression bandages over pockets of lymphedema. Fits into DORSAL POCKET GLOVE.

SIZE	SINGLE
SMALL	2231-DO (8cm x 9cm)
LARGE	2233-DO (9cm x 14cm)



CHIP PAD AXILLA

Use under arm for treatment of lymphedema, especially over areas of post-operative scarring.

SIZE	SINGLE
SMALL	2231-AX (10cm x 19cm)
LARGE	2233-AX (18cm x 27cm)



CHIP PAD RADIAL / HALF

Radial pathways creates low pressure channels to promote fluid drainage.

SIZE	SINGLE
ONE SIZE FITS ALL	2230-RH (15cm x 25cm)



CHIP PAD RADIAL / FULL

Versatile, effective design is ideal for post-surgical scarring and/or pockets of lymphedema.

SIZE	SINGLE
SMALL	2231-RF (18cm diameter)
LARGE	2233-RF (23cm diameter)



CHIP PAD BRA

Ideal for post mastectomy/lumpectomy edema. Pair with COMPREVEST.

SIZE	SINGLE
BLACK	
HALF	2230-BRH
FULL	2230-BRF
WHITE	
HALF	2280-BRH
FULL	2280-BRF

CHIP PAD GENIFIT

Comfortable and mild standard compression wear for genital edema


The GENIFIT may be paired with the COMPRESHORTS, or other compression wear. Open-cell foam chips in channel formation assist in massaging and facilitating drainage of interstitial fluid. Machine wash / dry.



MALE TALL
For genital and abdominal swelling



MALE SHORT
For genital swelling



MALE GENIFIT


	SMALL	MEDIUM	LARGE	X-LARGE	XX-LARGE
WIDTH	28cm	32cm	36cm	40cm	44cm
SCROTAL CIRCUM.	26cm	32cm	38cm	44cm	51cm
RISE - SHORT	14cm	16cm	18cm	20cm	22cm
RISE - TALL	32cm	35cm	38cm	41cm	44cm
SHORT	2231-GMS	2232-GMS	2233-GMS	2234-GMS	2235-GMS
TALL	2231-GMT	2232-GMT	2233-GMT	2234-GMT	2235-GMT



FEMALE TALL
For genital and abdominal swelling



FEMALE SHORT
For genital swelling



FEMALE GENIFIT

	SMALL	MEDIUM	LARGE	X-LARGE
WIDTH	28cm	32cm	36cm	40cm
RISE - SHORT	14cm	16cm	18cm	20cm
RISE - TALL	32cm	35cm	38cm	41cm
SHORT	2231-GFS	2232-GFS	2233-GFS	2234-GFS
TALL	2231-GFT	2232-GFT	2233-GFT	2234-GFT

COMPREFLEX REDUCE BELOW KNEE & BOOT



FEATURES

Lower leg options are the BELOW KNEE and BOOT (each sold separately)

Short-stretch compression mimics bandage wrapping

Soft, conforming Breathe-O-Prene provides a comfortable fit and wicks away moisture

Replaces bandage wrapping during decongestive therapy

May be combined with other COMPREFLEX REDUCE systems

Available in black or beige

Easy care - machine wash and dry

Latex-free

INDICATIONS

Decongestive phase of lymphedema treatment

INCLUDES (each product sold separately)

1 COMPREFLEX REDUCE system

1 pair of cotton liners

The new SIGVARIS COMPREFLEX REDUCE is a simple and effective system used during early decongestive phase of lymphedema, as it replaces the need for traditional bandage wrapping.

Each COMPREFLEX REDUCE garment is easily cut to fit the patient. The zipper creates simple garment assembly.



P.O. #	Account #	Contact	Date
--------	-----------	---------	------

Phone	Patient Name	Age	Height	Weight
-------	--------------	-----	--------	--------

Company Name

Dropship Name & Address

PRODUCT INFORMATION

COMPREFLEX REDUCE BK	Size	Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

COMPREFLEX REDUCE BOOT	Size	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT			<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

MEASURING INSTRUCTIONS

COMPREFLEX REDUCE BELOW KNEE

Measure circumferences where indicated by the black dots, beginning with the dot 5cm above the **Lateral Malleolus** (outer ankle bone), continue to the 25cm dot. Record in the corresponding lines A and C (Left/Right). Measure length from the **Lateral Malleolus** (outer ankle bone) to lateral aspect of the **Fibular Head**, and record in box G.

COMPREFLEX REDUCE BOOT


Measure circumference of ankle bend and heel and record on line I. Measure circumference across metatarsal heads and record on line J. Measure length on medial side of foot from heel to 1st metatarsal head and record on line K.

COMPREFLEX REDUCE BK SIZING CHART

	RANGE 1	RANGE 3
C	<65cm	<90cm
A	<55cm	<75cm
BLACK		
REGULAR	2401-BKR	2403-BKR
TALL	2401-BKT	2403-BKT
BEIGE		
REGULAR	2411-BKR	2413-BKR
TALL	2411-BKT	2413-BKT

COMPREFLEX REDUCE BOOT SIZING CHART

	RANGE 1	RANGE 3
I	30-44cm	39-54cm
J	24-33cm	30-43cm
K	16-25cm	20-27cm
BLACK	2801-BTL	2803-BTL
BEIGE	2811-BTL	2813-BTL



COMPREFLEX REDUCE THIGH COMPONENT & KNEE



FEATURES

- Upper leg options are the THIGH COMPONENT and KNEE PIECE (each sold separately)
- Short-stretch compression mimics bandage wrapping
- Soft, conforming Breathe-O-Prene provides a comfortable fit and wicks away moisture
- Replaces bandage wrapping during decongestive therapy
- May be combined with other COMPREFLEX REDUCE systems
- Available in black or beige
- Easy care - machine wash and dry
- Latex-free

INDICATIONS

- Decongestive phase of lymphedema treatment

INCLUDES (each product sold separately)

- 1 COMPREFLEX REDUCE system
- 1 pair of cotton liners

The new SIGVARIS COMPREFLEX REDUCE is a simple and effective system used during early decongestive phase of lymphedema, as it replaces the need for traditional bandage wrapping.

Each COMPREFLEX REDUCE garment is easily cut to fit the patient. The zipper creates simple garment assembly.



P.O. #	Account #	Contact	Date
--------	-----------	---------	------

Phone	Patient Name	Age	Height	Weight
-------	--------------	-----	--------	--------

Company Name

Dropship Name & Address

PRODUCT INFORMATION

COMPREFLEX REDUCE TC	Size	Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

COMPREFLEX REDUCE KNEE	Size	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT			<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

MEASURING INSTRUCTIONS

COMPREFLEX REDUCE THIGH COMPONENT

Measure circumferences where indicated by the black dots, beginning with the dot 5cm above the lateral aspect **Top of Patella** and continue to the 23cm dot. Record in the corresponding lines D and F (Left/Right). Measure length from the lateral aspect **Top of Patella** to lateral aspect **Gluteal Fold**, and record in box G1.

COMPREFLEX REDUCE KNEE


Measure circumferences where indicated by the black dots, beginning with the dot 5cm below the lateral aspect **Fibular Head**. Continue with the dot 5cm above the lateral aspect **Top of Patella**. Record in the corresponding lines C1 and D (Left/Right).

COMPREFLEX REDUCE TC SIZING CHART

	RANGE 1	RANGE 3
F	<80cm	<125cm
D	<65cm	<110cm

BLACK	LEFT	
	REGULAR	2401-TCR-L 2403-TCR-L
	TALL	2401-TCT-L 2403-TCT-L
	RIGHT	
BEIGE	LEFT	
	REGULAR	2411-TCR-L 2413-TCR-L
	TALL	2411-TCT-L 2413-TCT-L
	RIGHT	
BLACK	LEFT	
	REGULAR	2411-TCR-R 2413-TCR-R
	TALL	2411-TCT-R 2413-TCT-R
	RIGHT	
BEIGE	LEFT	
	REGULAR	2411-TCR-L 2413-TCR-L
	TALL	2411-TCT-L 2413-TCT-L
	RIGHT	
BLACK	LEFT	
	REGULAR	2411-TCR-R 2413-TCR-R
	TALL	2411-TCT-R 2413-TCT-R
	RIGHT	
BEIGE	LEFT	
	REGULAR	2411-TCR-L 2413-TCR-L
	TALL	2411-TCT-L 2413-TCT-L
	RIGHT	

CIRCUMFERENCE



COMPREFLEX REDUCE KNEE SIZING CHART

	RANGE 1	RANGE 3
D	<75cm	<115cm
C1	<65cm	<90cm
BLACK	2401-KP	2403-KP
BEIGE	2411-KP	2413-KP

BANDAGE KITS

STANDARD BANDAGE KITS

Each kit contains supplies for ONE limb

STANDARD BK	
PRODUCT CODE	5901-BK
CONTENTS	
1	BIAFORM Short Stretch Bandage 8cm x 5m
2	BIAFORM Short Stretch Bandage 10cm x 5m
1	BIAFORM Short Stretch Bandage 12cm x 5m
2	Non-Woven Padding 4" x 3.8yds / 10cm x 3.5M
6	ELASTOSOFT elastic gauze 4" X 4.1yds (6 rolls)
6	ELASTOSOFT elastic gauze 6" X 4.1yds
1	Cotton Tubular Stockinette 4" x 10 yds
1	DURAPORE Surgical Tape - Hypoallergenic 1" x 10 yds

STANDARD TH	
PRODUCT CODE	5901-TH
CONTENTS	
1	BIAFORM Short Stretch Bandage 8cm x 5m
3	BIAFORM Short Stretch Bandage 10cm x 5m
3	BIAFORM Short Stretch Bandage 12cm x 5m
2	Non-Woven Padding 4" x 3.8yds / 10cm x 3.5M
2	Non-Woven Padding 6" x 3.8yds / 15cm x 3.5M
8	ELASTOSOFT elastic gauze 4" X 4.1yds (6 rolls)
8	ELASTOSOFT elastic gauze 6" X 4.1yds
1	Cotton Tubular Stockinette 4" x 10 yds
1	DURAPORE Surgical Tape - Hypoallergenic 1" x 10yd

STANDARD ARM	
PRODUCT CODE	5901-ARM
CONTENTS	
1	BIAFORM Short Stretch Bandage 6cm x 5m
1	BIAFORM Short Stretch Bandage 8cm x 5m
1	BIAFORM Short Stretch Bandage 10cm x 5m
2	Non-Woven Padding 4" x 3.8yds / 10cm x 3.5M
6	ELASTOSOFT elastic gauze 4" X 4.1yds (6 rolls)
6	ELASTOSOFT elastic gauze 6" X 4.1yds
1	Cotton Tubular Stockinette 3" x 10 yds
1	DURAPORE Surgical Tape - Hypoallergenic 1" x 10yd



PREMIUM BANDAGE KITS

Each kit contains supplies for TWO limbs

PREMIUM BK	
PRODUCT CODE	5902-BK
CONTENTS	
2	BIAFORM Short Stretch Bandage 8cm x 5m
4	BIAFORM Short Stretch Bandage 10cm x 5m
2	BIAFORM Short Stretch Bandage 12cm x 5m
4	Non-Woven Padding 4" x 3.8yds / 10cm x 3.5M
12	ELASTOSOFT elastic gauze 4" X 4.1yds (6 rolls)
12	ELASTOSOFT elastic gauze 6" X 4.1yds
1	Cotton Tubular Stockinette 4" x 10 yds
1	DURAPORE Surgical Tape - Hypoallergenic 1" x 10yd

PREMIUM TH	
PRODUCT CODE	5902-TH
CONTENTS	
2	BIAFORM Short Stretch Bandage 8cm x 5m
6	BIAFORM Short Stretch Bandage 10cm x 5m
6	BIAFORM Short Stretch Bandage 12cm x 5m
4	Non-Woven Padding 4" x 3.8yds / 10cm x 3.5M
4	Non-Woven Padding 6" x 3.8yds / 15cm x 3.5M
16	ELASTOSOFT elastic gauze 4" X 4.1yds (6 rolls)
16	ELASTOSOFT elastic gauze 6" X 4.1yds
1	Cotton Tubular Stockinette 4" x 10 yds
1	DURAPORE Surgical Tape - Hypoallergenic 1" x 10yd

PREMIUM ARM	
PRODUCT CODE	5902-ARM
CONTENTS	
2	BIAFORM Short Stretch Bandage 6cm x 5m
2	BIAFORM Short Stretch Bandage 8cm x 5m
2	BIAFORM Short Stretch Bandage 10cm x 5m
4	Non-Woven Padding 4" x 3.8yds / 10cm x 3.5M
12	ELASTOSOFT elastic gauze 4" X 4.1yds (6 rolls)
12	ELASTOSOFT elastic gauze 6" X 4.1yds
1	Cotton Tubular Stockinette 3" x 10 yds
1	DURAPORE Surgical Tape - Hypoallergenic 1" x 10yd



SHORT STRETCH BANDAGES



BIAFORM

Premium short stretch bandage. Superior stretch and conforming abilities. Exceptional durability.

Case holds 20 rolls

SIZE	QUANTITY	PRODUCT CODE
6cm x 5m	1 CASE	5052-C
8cm x 5m	1 CASE	5053-C
10cm x 5m	1 CASE	5054-C
12cm x 5m	1 CASE	5055-C
10cm x 10m	1 CASE	5056-C
12cm x 10m	1 CASE	5057-C



ISOFORM

Short-stretch warp, perfect for lymphedema in the thigh or trunk areas.

Case holds 10 rolls

SIZE	QUANTITY	PRODUCT CODE
15cm x 5m	1 CASE	5061-C
20cm x 5m	1 CASE	5062-C

STOCKINETTE & DRESSING RETAINERS



Cotton Stockinette

Non-compressive, 100% cotton. Ideal for use with nighttime bandaging or compression garments. Available in a range of sizes to fit almost any application.

LATEX FREE

SIZE	QUANTITY	PRODUCT CODE
2" x 25yds	1 ROLL	5402
3" x 25yds	1 ROLL	5403
4" x 25yds	1 ROLL	5404
5" x 25yds	1 ROLL	5405
6" x 25yds	1 ROLL	5406
8" x 25yds	1 ROLL	5408



BIAGRIP

Compressive stockinette. Can be washed and reused to augment compression garments. Cut to size.

LATEX FREE

SIZE	QUANTITY	PRODUCT CODE
A (1.5" x 11yds)	1 ROLL	5430
B (2.5" x 11yds)	1 ROLL	5431
C (2.75" x 11yds)	1 ROLL	5432
D (3" x 11yds)	1 ROLL	5433
E (3.5" x 11yds)	1 ROLL	5434
F (4" x 11yds)	1 ROLL	5435
G (5" x 11yds)	1 ROLL	5436
J (7" x 11yds)	1 ROLL	5437
K (8" x 11yds)	1 ROLL	5438
L (12" x 11yds)	1 ROLL	5439

GAUZE & TAPE



ELASTOSOFT

Soft cotton elastic gauze. Ideal for wrapping fingers & toes.

SIZE	QUANTITY	PRODUCT CODE
1" x 4.1yds	BAG OF 24 ROLLS	5601
2" x 4.1yds	BAG OF 12 ROLLS	5602
3" x 4.1yds	BAG OF 12 ROLLS	5603
4" x 4.1yds	BAG OF 12 ROLLS	5604
6" x 4.1yds	BAG OF 6 ROLLS	5606



3M DURAPORE™ SILK SURGICAL TAPE

Exceptionally durable, easy to cut to size. Optimal for securing bandages.

SIZE	QUANTITY	PRODUCT CODE
1" x 10yds	1 BOX (12 Rolls)	5691-B
2" x 10yds	1 BOX (6 Rolls)	5692-B

FOAM



LYMPHESOFT

Soft, open-cell rolled foam provides even compression distribution when bandaging. Use with BIAFORM short stretch bandages. Air permeable. Tear resistant.


SIZE	QUANTITY	PRODUCT CODE
10cm x 2.5m	1 ROLL	5101
10cm x 2.5m	1 CASE	5101-C
12cm x 2.5m	1 ROLL	5102
12cm x 2.5m	1 CASE	5102-C
15cm x 2.5m	1 ROLL	5103
15cm x 2.5m	1 CASE	5103-C



WAVEFOAM (1 sided)

High-pressure ridges gently massage thickened fibrotic areas. Low-pressure channels allow pathways for fluid evacuation.

SIZE	QUANTITY	PRODUCT CODE
1.2cm x 69cm x 90cm	1 SHEET	5202



WAVEFOAM II (2 sided)

High-pressure ridges gently massage thickened fibrotic areas. Low-pressure channels allow pathways for fluid evacuation.

SIZE	QUANTITY	PRODUCT CODE
1.2cm x 69cm x 90cm	1 SHEET	5212



GRAYFOAM

Moderately rigid open-cell foam helps soften fibrotic areas. Provides firm, stable support under BIAFORM short stretch bandages.

SIZE	QUANTITY	PRODUCT CODE
1/4" x 18" x 72"	1 SHEET	5301
1/2" x 18" x 72"	1 SHEET	5302

PADDING



BIASOFT

Soft, durable, and seamless fleece padding. Use under short stretch bandages. Can be machine washed and dried for repeated use. 10 rolls per case.

SIZE	QUANTITY	PRODUCT CODE
6cm x 3m	1 ROLL	5251
6cm x 3m	1 CASE	5251-C
10cm x 3m	1 ROLL	5252
10cm x 3m	1 CASE	5252-C
15cm x 3m	1 ROLL	5253
15cm x 3m	1 CASE	5253-C



Non-Woven Soft, Synthetic Padding

Optimal for use under short stretch bandages.

SIZE	QUANTITY	PRODUCT CODE
10cm x 3.5m (4" x 3.8 yds)	1 CASE (30 ROLLS)	5661-C

SOLUTIONS

FOR PATIENTS WITH CHRONIC, ACUTE AND SEVERE LYMPHEDEMA



SIGVARIS custom compression garments are ideal for treating stage III / chronic lymphedema.

For 20 years, the MEDASSIST line has been the industry standard for custom lymphedema compression wear.

Our custom CHIPSLEEVE ARM presents an affordable option for chipped foam nighttime garments.

The medical grade textiles used in SIGVARIS MEDASSIST products is high quality and extremely durable.

HOW TO ORDER CUSTOM COMPRESSION WEAR

Measuring instructions and order forms for all custom compression wear are provided in the following pages, and may also be found online at **bit.ly/SIGmccCustom**. Our sales support team is available to assist you through the process of ordering a custom garment. From determining if your patient is a good candidate for a custom garment to providing measuring guidance, we are here to help.

We require pictures of the patient after measurements have been taken, with the measurement marks still present*. These photos may be emailed to the address below.

Contact us:
866.931.0876
us_orders@sigvaris.com

All custom products in this catalog are handmade and patterned according to the measurements provided at our manufacturing facility in Holland, Michigan.

* SIGVARIS, Inc abides by HIPPA guidelines and protects all patient-sensitive information. SIGVARIS does not disclose patient sensitive information to third parties.

CUSTOM CHIPSLEEVE ARM

INDICATIONS

- Stage I or II lymphedema
- Mild fibrosis
- Post-surgical edema
- Patient measurements are outside of standard size ranges

FEATURES

- 15–25mmHg of mild, gradient compression from hand to axilla when used with the OVERSLEEVE
- Only 7 easy measurements required
- Open-cell foam chips massage fibrotic tissue, stimulating fluid movement
- Use as a bandage liner or independent compression garment
- 1 year manufacturer warranty (visit sigvarisusa.com for details)
- Four-way stretch fabric contains moisture-wicking, odor control properties
- Safe for day or nighttime use
- Easy care - machine wash and dry
- Latex-free




ARMASSIST

INDICATIONS

- Stage II or III lymphedema
- Patient measurements are outside of standard size ranges

FEATURES

- 30–40mmHg firm, gradient compression
- Constructed with Flat Foam or WaveFoam
- Exceptionally thin and lightweight design
- The medical grade fabric features non-static, non-absorbent, and non-stretch qualities
- Palm pad for added dorsum compression





CUSTOM CHIPSLEEVE ARM

P.O. # Account # Contact Date

Phone Patient Name Age Height Weight

Company Name

Dropship Name & Address

PRODUCT INFORMATION

Arm

LEFT RIGHT

OVERSLEEVE Color

BLACK NAVY PINK TIE DYE BLUE HORIZON LEOPARD

NOTE: If no color is selected, BLACK will be included.

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

MEASURING INSTRUCTIONS

With patient seated, place the arm extended and elbow bent slightly, with the palm down, on a flat surface.

CIRCUMFERENCES

Measure circumference at Palm and record on line E.

Measure circumference at Wrist, and record on line A.
Mark dorsal aspect at distal edge of tape.

Measure circumference at Forearm, and record on line A1.

Measure circumference at Elbow, and record on line B.

Measure circumference at mid Bicep, and record on line B1.


Measure circumference at Axilla, and record on line C.

Mark dorsal aspect at proximal edge of tape.

LENGTH

Measure length at dorsal aspect from mark at Point A to mark at Point C and record in box D.

NOTES:





P.O. #	Account #	Contact	Date
--------	-----------	---------	------

Phone	Patient Name	Age	Height	Weight
-------	--------------	-----	--------	--------


Company Name

Dropship Name & Address

PRODUCT INFORMATION

Arm LEFT RIGHT Foam FLAT WAVEFOAM

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms



ARMASSIST MEASURING INSTRUCTIONS

IMPORTANT

Measuring for this garment requires special measurements critical to proper fit. Before you begin, attend live training at a SIGVARIS Certified Compression Specialist course near you (info at sigvariseducation.com). Alternatively, call Customer Care at 800-322-7744, or e-mail us_orders@sigvaris.com, to receive a remote consultation/training.

SUPPLIES NEEDED

1. Cell phone with camera. Photos of the arm(s) with measurement markings MUST be e-mailed to: us_orders@sigvaris.com.
2. Measuring Instructions and Forms.
3. SIGVARIS Measuring tape and Body Pen (or eyeliner pencil).
4. Signed Custom Order Terms & Conditions Form (include with order).

INSTRUCTIONS

All measurements should be recorded in centimeters. Apply slight tension to hold the tape measure in place.

STEP 1: LOCATE LANDMARKS

1. Encircle arm with tape measure located as high as possible at Axilla. Mark on the proximal edge of tape at the Anterior, Dorsal, and Posterior aspects of arm, and record the circumference on the line labeled Posterior Axilla.
2. Encircle elbow with the distal edge of a tape measure aligned at the elbow crease and Olecranon Process (elbow). Mark the dorsal aspect of the elbow along this edge, and record the circumference on the line labeled Elbow (Ø Point).
3. Encircle wrist with the distal edge of a tape measure aligned with the ulnar styloid process. Mark the dorsal aspect of the wrist along this edge, and record the circumference on the line labeled Ulnar Styloid.

STEP 2: LOCATE CIRCUMFERENCE MEASUREMENT POSITIONS

1. Starting with "0" of the tape measure at your mark on the dorsal aspect of the elbow, mark the dorsal aspect of the upper arm in 5cm increments up to the mark made from the Axilla on the dorsal aspect of the arm.
2. Starting with the "0" of the tape measure at your mark on the dorsal aspect of the elbow, mark the dorsal aspect of the forearm in 5cm increments down to the mark made for the Ulnar Styloid on the dorsal aspect of the arm.

STEP 3: MEASURE LENGTHS

1. Measure upper arm from anterior axilla (inside armpit) to the elbow crease and record length on Box A.
2. Measure upper arm from posterior axilla (outside armpit) to the Olecranon Process (elbow) and record length on Box B.
3. Measure the forearm along dorsal aspect from elbow (as marked at the mid-point between Olecranon Process and elbow crease) to the Ulnar Styloid as indicated on diagram and record length on Box C.
4. Measure hand from your mark on dorsal aspect of wrist to the third metacarpal head and record length on Box D.
5. Measure width of hand across dorsal metacarpal heads and record width on Box E.

STEP 4: MEASURE CIRCUMFERENCES

1. Encircle arm with a tape measure at each mark located between the elbow and axilla, and record each circumference on the line corresponding to each position.
2. Encircle arm with a tape measure at each mark located between the elbow and Ulnar Styloid, and record each circumference on the line corresponding to each position.

TAKE DIGITAL PHOTOS (REQUIRED) AFTER MEASUREMENT MARKINGS ARE MADE. EMAIL TO: us_orders@sigvaris.com

NOTE: ARMASSIST is available in: 7099-AR FlatFoam, 7199-AR WAVEFOAM, 7299-AR Super FlatFoam, 7399-AR Super WAVEFOAM.

LEGASSIST

INDICATIONS

Stage II or III lymphedema

Patient measurements are outside of standard size ranges

Lobule is present (use LCS or BK with SHELF)

FEATURES

30–40mmHg firm, gradient compression

Constructed with Flat Foam or WAVEFOAM

Exceptionally thin and lightweight design

The medical grade fabric features non-static, non-absorbent and non-stretch qualities

Zipper panel for easy donning / doffing

COMPREBOOT PLUS is included (not applicable with LCS)

OPTIONS

MEDABOOT (page 81)

LEGASSIST BELOW KNEE (page 82)

LEGASSIST BELOW KNEE WITH SHELF (page 84)

LEGASSIST LOBULE COMPRESSION SYSTEM (page 87)

LEGASSIST THIGH HIGH (page 90)



LEGASSIST
BELOW KNEE (BK)



LEGASSIST
THIGH HIGH (TH)



LEGASSIST
LOBULE COMPRESSION SYSTEM (LCS)

P.O. #	Account #	Contact	Date		
--------	-----------	---------	------	--	--

Phone	Patient Name	Age	Height	Weight
-------	--------------	-----	--------	--------

Company Name

Dropship Name & Address

PRODUCT INFORMATION

Foot

LEFT RIGHT

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

IMPORTANT

Measuring for this garment requires special measurements critical to proper fit. Before you begin, attend live training at a SIGVARIS Certified Compression Specialist course near you (info at sigvariseducation.com). Alternatively, call Customer Care at 800-322-7744, or e-mail us_orders@sigvaris.com, to receive a remote consultation/training.

SUPPLIES NEEDED

1. Cell phone with camera. Photos of the foot with measurement markings MUST be e-mailed to: us_orders@sigvaris.com.
2. Measuring Instructions and Forms.
3. SIGVARIS Measuring tape and Body Pen (or eyeliner pencil).
4. Signed Custom Order Terms & Conditions Form (include with order).

INSTRUCTIONS

All measurements should be recorded in centimeters. Apply slight tension to hold the tape measure in place.

MEASURE FOOT

1. Measure foot from heel to 1st metatarsal head (or desired boot length) and record length on Line K.
2. Measure Top of Foot from the 3rd Metatarsal Head to Ankle Bend and record length on Line M.
3. Encircle the Ankle Bend and Heel with a tape measure and record the circumference on Line I.
4. Encircle the foot across the metatarsal heads and record the circumference on Line J.
5. Measure the ankle circumference and record on Line L.

CIRCUMFERENCE

LEFT RIGHT

_____ I

_____ J


_____ L

LENGTH

LEFT RIGHT

_____ K

_____ M





P.O. #	Account #	Contact	Date
--------	-----------	---------	------

Phone	Patient Name	Age	Height	Weight
-------	--------------	-----	--------	--------

Company Name


Dropship Name & Address

PRODUCT INFORMATION

Leg LEFT RIGHT Foam FLAT WAVEFOAM Foot Options COMPREBOOT PLUS CUSTOM MEDABOOT (additional charge)

Size* REGULAR SUPER ***NOTE:** If the greatest circumference measurement is >60cm, order a Super.

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms



LEGASSIST BK MEASURING INSTRUCTIONS

IMPORTANT

Measuring for this garment requires special measurements critical to proper fit. Before you begin, attend live training at a SIGVARIS Certified Compression Specialist course near you (info at sigvariseducation.com). Alternatively, call Customer Care at 800-322-7744, or e-mail us_orders@sigvaris.com, to receive a remote consultation/training.

SUPPLIES NEEDED

1. Cell phone with camera. Photos of the leg(s) with measurement markings MUST be e-mailed to: us_orders@sigvaris.com.
2. Measuring Instructions and Forms.
3. SIGVARIS Measuring tape and Body Pen (or eyeliner pencil).
4. Signed Custom Order Terms & Conditions Form (include with order).

INSTRUCTIONS

All measurements should be recorded in centimeters. Apply slight tension to hold the tape measure in place.

STEP 1: LOCATE LANDMARKS

1. The ideal placement for the top edge of the garment is two finger widths below the popliteal crease, and across the fibular head and tibial tuberosity. However, the position of the Top of Garment can be customized, if needed, to overlap better with a LEGASSIST LCS thigh piece. Encircle leg with top edge of tape measure at the desired location. Using the SIGVARIS body pen, trace a line around the leg along the top edge of the tape and measure and record the circumference on Top of Garment line.
2. Encircle leg with bottom edge of tape measure aligned with the ankle bend. Mark leg on this edge of tape at the anterior, medial, posterior and lateral aspects, and record the circumference on the line labeled Ankle Bend (Ø Point).

STEP 2: LOCATE CIRCUMFERENCE MEASUREMENT POSITIONS

1. Starting with "0" of the tape measure at your mark on the lateral aspect of the ankle bend, mark the leg in 5cm increments up to your mark on the lateral aspect of the Top of Garment.

STEP 3: MEASURE LENGTHS

1. Measure leg from your mark on the anterior aspect of the Ankle Bend to the anterior aspect of Top of Garment and record length on line marked A.
2. Measure leg from your mark on the medial aspect of the Ankle Bend to the medial aspect of Top of Garment and record length on line marked B.
3. Measure leg from your mark on the posterior aspect of the Ankle Bend to the posterior aspect of Top of Garment and record length on line marked C.
4. Measure leg from your mark on the lateral aspect of the Ankle Bend to the lateral aspect of Top of Garment and record length on line marked D.

STEP 4: MEASURE CIRCUMFERENCES

1. Encircle leg with a tape measure at each mark located between the ankle bend and Top of Garment, and record each circumference on the line corresponding to each position.

STEP 5: MEASURE FOOT

1. Measure foot from heel to 1st metatarsal head (or desired boot length) and record length on Line K.
2. Measure Top of Foot from the 3rd Metatarsal Head to Ankle Bend and record length on Line M.
3. Encircle the instep of foot and apex of Heel with a tape measure and record the circumference on Line I.
4. Encircle the foot across the metatarsal heads and record the circumference on Line J.

TAKE DIGITAL PHOTOS (REQUIRED) AFTER MEASUREMENT MARKINGS ARE MADE. EMAIL TO: us_orders@sigvaris.com

NOTE: LEGASSIST BELOW KEE is available in:
7099-BK Flat Foam, 7199-BK WAVEFOAM, 7299-BK Super, Flat Foam, 7399-BK Super, WAVEFOAM.



P.O. # _____ Account # _____ Contact _____ Date _____

Phone _____ Patient Name _____ Age _____ Height _____ Weight _____

Company Name _____


Dropship Name & Address _____

PRODUCT INFORMATION

Leg LEFT RIGHT Foam FLAT WAVEFOAM Foot Options COMPREBOOT PLUS CUSTOM MEDABOOT (additional charge)

Size* REGULAR SUPER ***NOTE:** If the greatest circumference measurement is >60cm, order a Super.

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms



IMPORTANT

Measuring for this garment requires special measurements critical to proper fit. Before you begin, attend live training at a SIGVARIS Certified Compression Specialist course near you (info at sigvariseducation.com). Alternatively, call Customer Care at 800-322-7744, or e-mail us_orders@sigvaris.com, to receive a remote consultation/training.

SUPPLIES NEEDED

1. Cell phone with camera. Photos of the leg(s) with measurement markings MUST be e-mailed to: us_orders@sigvaris.com.
2. Measuring Instructions and Forms.
3. SIGVARIS Measuring tape and Body Pen (or eyeliner pencil).
4. Signed Custom Order Terms & Conditions Form (include with order).

INSTRUCTIONS

All measurements should be recorded in centimeters. Apply slight tension to hold the tape measure in place.

STEP 1: CHOOSE THE ZERO POINT (Ø POINT)

The Zero Point is unique to the LEGASSIST BK WITH SHELF garment. It is the reference point from which all other measurements are taken, in order to determine the lengths needed for the custom garment to fit the patient's individual ankle shape. Encircle leg with bottom edge of tape measure immediately above the shelf, at the lowest aspect of the curve of the ankle lobule where a tape measure can be applied without sliding off. Trace a line around the leg along the bottom edge of the tape using the SIGVARIS body pen, and record the circumference on Ø Point line. Document your zero point location by measuring the distance from the floor.

STEP 2: CHOOSE WHERE GARMENT WILL END DISTALLY - "BOTTOM OF GARMENT"

Encircle leg with bottom edge of tape measure where the ankle lobule touches the ankle and/or foot. Do not pull tape under shelf. Mark leg on this edge of tape at the anterior, medial, posterior and lateral aspects, and record the circumference on the line labeled Bottom of Garment. The smallest ankle measurement is not relevant for this product. The product is designed to hug and lift the ankle lobule.

STEP 3: CHOOSE WHERE THE GARMENT WILL END PROXIMALLY - "TOP OF GARMENT"

The ideal placement for the top edge of the garment is two finger widths below the popliteal crease, and across the fibular head and tibial tuberosity. However, the position of the Top of Garment can be customized, if needed, to overlap better with a LEGASSIST LCS thigh piece. Encircle leg with top edge of tape measure at the desired location. Using the SIGVARIS body pen, trace a line around the leg along the top edge of the tape and measure and record the circumference on Top of Garment line.

STEP 4: MEASURE AND RECORD LENGTHS

Measure each of the lengths below from the Zero Point, following the contour of the leg and ankle, but not extending the tape measure into folds. Lengths measured distal to the Zero Point should end where the ankle lobule touches the ankle or foot at the Bottom of Garment versus extending the tape measure up under the lobule. Record lengths as follows:

LATERAL (A) LENGTHS:

- A1: Measure & record lateral length from Ø Point to Top of Garment
- A2: Measure & record lateral length from Ø Point to Bottom of Garment

MEDIAL (B) LENGTHS:

- B1: Measure & record medial length from Ø Point to Top of Garment
- B2: Measure & record medial length from Ø Point to Bottom of Garment

POSTERIOR (C) LENGTHS:

- C1: Measure & record posterior length from Ø Point to Top of Garment
- C2: Measure & record posterior length from Ø Point to Bottom of Garment

ANTERIOR (D) LENGTHS:

- D1: Measure & record anterior length from Ø Point to Top of Garment
- D2: Measure & record anterior length from Ø Point to Bottom of Garment

LEGASSIST BK W/ SHELF MEASURING INSTRUCTIONS

STEP 5: MEASURE AND RECORD CIRCUMFERENCES

On the LATERAL aspect of the leg ONLY, mark every 5cm moving proximally from the Zero Point. Next, record the circumferences at each point, making sure that the measuring tape is parallel to the floor.

STEP 6: MEASURE FOOT

1. Measure the circumference from where the shelf touches instep of the foot to the apex of the heel and record on line I.
2. Measure the foot circumference at the metatarsal heads and record on line J.
3. Measure the length from the 1st metatarsal head to the posterior aspect of the heel, and record on line K.
4. Measure the length from the dorsal aspect of the third metatarsal head to where the shelf touches the instep of the foot and record on line M.

STEP 7: PICTURES

Now that you have completed your measurements, take pictures of the leg (front, back, side) that show all of the measurement markings, as well as the overall shape of the leg. This will assist our expert design team to further customize the garment according to your patient's individual needs. E-mail the pictures with your completed forms to: us_orders@sigvaris.com.

Don't forget to record the patient's age, height, and weight, as this also helps with customizing the garment.

STEP 8: SELECT ADDITIONAL OPTIONS

You are ready to complete your forms and place the order.

Choose *Right* or *Left* leg

Choose *Size*: If the greatest leg circumference is >60cm, order a Super.

Choose *Flat* or *WAVEFOAM* (recommended)

Choose COMPREBOOT Plus (included) or Custom MEDABOOT (additional charge)

TAKE DIGITAL PHOTOS (REQUIRED) AFTER MEASUREMENT MARKINGS ARE MADE. EMAIL TO: us_orders@sigvaris.com

NOTE: LEGASSIST BELOW KNEE WITH SHELF is available in:

7099-BK Flat Foam, 7199-BK WAVEFOAM, 7299-BK Super, Flat Foam, 7399-BK Super, WAVEFOAM.

LEGASSIST LOBULE COMPRESSION SYSTEM

P.O. #	Account #	Contact	Date
--------	-----------	---------	------

Phone	Patient Name	Age	Height	Weight
-------	--------------	-----	--------	--------

Company Name


Dropship Name & Address

PRODUCT INFORMATION

Leg LEFT RIGHT Foam FLAT WAVEFOAM Optional HIP ATTACHMENT

Size* REGULAR SUPER ***NOTE:** If the greatest circumference measurement is >90cm, order a Super.

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms



LEGASSIST LCS MEASURING INSTRUCTIONS

IMPORTANT

Measuring for this garment requires special measurements critical to proper fit. Before you begin, attend live training at a SIGVARIS Certified Compression Specialist course near you (info at sigvariseducation.com). Alternatively, call Customer Care at 800-322-7744, or e-mail us_orders@sigvaris.com, to receive a remote consultation/training.

SUPPLIES NEEDED

1. Cell phone with camera. Photos of the leg(s) with measurement markings MUST be e-mailed to: us_orders@sigvaris.com.
2. Measuring Instructions and Forms.
3. SIGVARIS Measuring tape and Body Pen (or eyeliner pencil).
4. Signed Custom Order Terms & Conditions Form (include with order).

INSTRUCTIONS

All measurements should be recorded in centimeters. Apply slight tension to hold the tape measure in place.

STEP 1: CHOOSE THE ZERO POINT (Ø POINT)

The Zero Point is unique to the LEGASSIST LCS garment. It is the reference point from which all other measurements are taken, in order to determine the lengths needed for the custom garment to fit the patient's individual leg shape. Choose a location proximal to the knee joint where the leg shape is the most cylindrical and the tape measure can be positioned circumferentially around the leg and exactly parallel to the floor. Trace a line around the leg at this Zero Point using the SIGVARIS body pen and record the circumference on Ø Point line. If possible, document the location of the Zero Point by measuring the distance from the head of the fibula to the Zero Point.

STEP 2: CHOOSE WHERE THE GARMENT WILL END DISTALLY - "BOTTOM OF GARMENT"

The length of the LEGASSIST LCS will depend on how low the medial knee lobule hangs down. Start by marking the lowest aspect of the curve of the lobule where a tape measure could be applied without sliding off (lobule apex). Next, choose and mark a point on the lateral aspect of the leg where the garment will overlap the lower leg piece by at least 3 centimeters and create no more than a 30degree downward angle or no more than 10degree upward angle to the lobule apex. Encircle the leg from this point to the lobule apex, with the tape measure positioned just above your marking points. Trace this circumference with your body pen along the bottom edge of the tape and record the circumference measurement on Bottom of Garment line.

PROBLEM SOLVING:

If the Bottom of Garment is proximal to the knee, due to a lobule positioned very high on the thigh, you will need to measure for a customized LEGASSIST TH system that overlaps with the LEGASSIST LCS.

If the Bottom of Garment extends down to the mid-lower leg or ankle, due to a lobule positioned very low on the leg, you will need to measure for a customized LEGASSIST BK that is shortened to best overlap with the LEGASSIST LCS.

STEP 3: CHOOSE WHERE THE GARMENT WILL END PROXIMALLY - "TOP OF GARMENT"

You may choose to end the garment circumferentially at the gluteal crease, or extend the garment proximally on the lateral aspect of the thigh to the greater trochanter. Mark your preferred position for the top of the garment along the medial thigh, at least 4cm distal to the groin. Next, mark your preferred position for the top of the garment along the lateral thigh, anywhere from the level of the gluteal fold to the lateral greater trochanter. Encircle the leg with the proximal edge of the tape measure positioned just below your markings, trace a line around the leg along the top edge of the tape with your body pen, and record the circumference on Top of Garment line.

STEP 4: MARK THE MIDLINE OF GARMENT AND RECORD LENGTHS

Step back to visualize the leg and your marking points to determine where the midline of the garment would be, medial to lateral, encompassing the lobules. Mark the midline of the garment anteriorly (D) and posteriorly (C) at the Top of Garment, Zero Point, and Bottom of Garment lines. Measure each of the lengths from the Zero Point, following the contour of the leg but not extending the tape measure into folds. Record lengths as follows:

LATERAL (A) LENGTHS:

A1: Measure & record lateral length from Zero Point to Top of Garment

A2: Measure & record lateral length from Zero Point to Bottom of Garment

MEDIAL (B) LENGTHS:

B1: Measure & record medial length from Zero Point to Top of Garment

B2: Measure & record medial length from Zero Point to Bottom of Garment

LEGASSIST LCS MEASURING INSTRUCTIONS

POSTERIOR (C) LENGTHS:

- C1: Measure & record posterior length from Zero Point to Top of Garment
- C2: Measure & record posterior length from Zero Point to Bottom of Garment

ANTERIOR (D) LENGTHS:

- D1: Measure & record anterior length from Zero Point to Top of Garment
- D2: Measure & record anterior length from Zero Point to Bottom of Garment

STEP 5: MEASURE AND RECORD CIRCUMFERENCES

Starting from the Zero Point, on the LATERAL aspect of the leg ONLY, mark every 5cm proximally from the zero point and mark every 5 cm distally from the zero point. Record the circumferences at each point, making sure that the measuring tape is parallel to the floor (Note that some of the spaces to record circumferences may not be used, or you may need to add additional spaces, depending on the length of the limb).

STEP 6: MEASURE LOBULE

Measure the size of the lobule for the part of the LCS system that will lift and compress the lobule, called the "Lobule Sling". Starting at the medial aspect of the lobule, measure the distance from the Bottom of Garment line to the medial aspect of the leg. Follow the curve of the lobule only to where it touches the leg. Do not extend the tape measure up into the crevices between the lobule and the leg. Record this length as E. At the midline of the lobule, measure from the anterior line of the Bottom of Garment, underneath the lobule and up to the posterior line of the Bottom of Garment. Record this length as F.

HIP ATTACHMENT (OPTIONAL)

1. Measure from lateral aspect of Gluteal Fold to lateral aspect of hip at the desired location for the bottom of the belt, and record length on Line H.
2. Measure circumference of torso at the desired location for the bottom of the belt, and record circumference on Line G.

STEP 8: PICTURES

Now that you have completed your measurements, take pictures of the leg (front, back, side) that show all of the measurement markings, as well as the overall shape of the leg. This will assist our expert design team to further customize the garment according to your patient's individual needs. E-mail the pictures with your completed forms to: us_orders@sigvaris.com. Don't forget to record the patient's age, height, and weight, as this also helps with customizing the garment.

STEP 9: SELECT ADDITIONAL OPTIONS

- You are ready to complete your forms and place the order.
- Choose *Right* or *Left* leg
- Choose *Size*: If the greatest leg circumference is >90cm, order a Super.
- Choose *Flat* or *WAVEFOAM* (recommended)
- Choose the *Hip Attachment*, if desired (recommended)

FITTING TIPS

The ropey, thicker seam should be positioned laterally, and the tag is positioned at the top of the garment straps pull lateral to medial.

TAKE DIGITAL PHOTOS (REQUIRED) AFTER MEASUREMENT MARKINGS ARE MADE. EMAIL TO: us_orders@sigvaris.com

P.O. #	Account #	Contact	Date
--------	-----------	---------	------

Phone	Patient Name	Age	Height	Weight
-------	--------------	-----	--------	--------

Company Name

Dropship Name & Address

PRODUCT INFORMATION

Leg	Foam	Foot Options
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT	<input type="checkbox"/> FLAT <input type="checkbox"/> WAVEFOAM	<input type="checkbox"/> COMPREBOOT PLUS <input type="checkbox"/> CUSTOM MEDABOOT (additional charge)
Size*	Optional	
<input type="checkbox"/> REGULAR <input type="checkbox"/> SUPER *NOTE: If the greatest circumference measurement is >90cm, order a Super.	<input type="checkbox"/> HIP ATTACHMENT	

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

CIRCUMFERENCE MEASUREMENTS

TOP OF GARMENT

35 cm _____

30 cm _____

25 cm _____

20 cm _____

15 cm _____

10 cm _____

5 cm _____

Ø POINT TOP OF PATELLA

MID PATELLA

Ø POINT BOTTOM OF PATELLA

5 cm _____

10 cm _____

15 cm _____

20 cm _____

25 cm _____

30 cm _____

35 cm _____

40 cm _____

ANKLE BEND

● = Locations measured along Lateral aspect of leg

LENGTH MEASUREMENTS

Lateral

A1 _____

A2 _____

Medial

B1 _____

B2 _____

Posterior

C1 _____

C2 _____

Anterior

D1 _____

D2 _____

Knee

E _____

FOOT MEASUREMENTS

I _____

J _____

K _____

M _____

No straps provided over knee unless box is checked

LEGASSIST TH MEASURING INSTRUCTIONS

IMPORTANT

Measuring for this garment requires special measurements critical to proper fit. Before you begin, attend live training at a SIGVARIS Certified Compression Specialist course near you (info at sigvariseducation.com). Alternatively, call Customer Care at 800-322-7744, or e-mail us_orders@sigvaris.com, to receive a remote consultation/training.

SUPPLIES NEEDED

1. Cell phone with camera. Photos of the leg(s) with measurement markings MUST be e-mailed to: us_orders@sigvaris.com
2. Measuring Instructions and Forms
3. SIGVARIS Measuring tape and Body Pen (or eyeliner pencil)
4. Signed Custom Order Terms & Conditions Form (include with order)

INSTRUCTIONS

All measurements should be recorded in centimeters. Apply slight tension to hold the tape measure in place.

STEP 1: LOCATE LANDMARKS

1. Encircle leg with top edge of tape measure located at the Bottom of Patella. Mark leg on this edge of tape at the anterior, lateral, posterior and medial aspects, recording the circumference on the line labeled Bottom of Patella (Ø Point).
2. Encircle leg with the bottom edge of the tape measure aligned with the Ankle Bend. Mark leg on this edge of the tape at the anterior, lateral, posterior and medial aspects, recording the circumference on the line labeled Ankle Bend.
3. Encircle leg with bottom edge of the tape measure aligned with the Top of Patella. Mark leg on this edge of the tape at the anterior, lateral, posterior and medial aspects, recording the circumference on the line labeled Top of Patella.
4. You may choose to end the garment circumferentially at the gluteal crease, or extend the garment proximally on the lateral aspect of the thigh to the greater trochanter. Mark your preferred position for the top of the garment along the medial thigh, at least 4cm distal to the groin. Next, mark your preferred position for the top of the garment along the lateral thigh, anywhere from the level of the gluteal fold to the lateral greater trochanter. Encircle the leg with the proximal edge of the tape measure positioned just below your markings, trace a line around the leg along the top edge of the tape with your body pen, and record the circumference on Top of Garment line.

STEP 2: LOCATE CIRCUMFERENCE MEASUREMENT POSITIONS

1. Starting with "0" of the tape measure at your mark on the lateral aspect of the Bottom of Patella (Ø Point), mark the leg in 5cm increments down to your mark on the lateral aspect of the Ankle Bend.
2. Repeat step 1, only this time from the Top of the Patella (Ø Point) up to the lateral aspect of the Top of Garment.

STEP 3: MEASURE LENGTHS

LATERAL LENGTH (A1):

Measure upper leg from your mark on the lateral aspect of Top of Patella to your mark on lateral aspect of Top of Garment. Record length on Line A1.

MEDIAL LENGTH (B1):

Measure upper leg from your mark on medial aspect of Top of Patella to your mark on medial aspect of Top of Garment. Record length on Line B1.

POSTERIOR LENGTH (C1):

Measure upper leg from your mark on posterior aspect of Top of Patella to your mark on posterior aspect of Top of Garment. Record length on Line C1.

ANTERIOR LENGTH (D1):

Measure upper leg from your mark on anterior aspect of Top of Patella to your mark on anterior aspect of Top of Garment. Record length on Line D2.

KNEE SPACE:

Measure leg from your mark on the lateral aspect of Top of Patella to your mark on lateral aspect of Bottom of Patella. Record length on Line E.

LEGASSIST TH MEASURING INSTRUCTIONS

LATERAL LENGTH (A2):

Measure lower leg from your mark on the lateral aspect of Bottom of Patella to your mark on lateral aspect of Ankle Bend. Record length on Line A2.

MEDIAL LENGTH (B2):

Measure lower leg from your mark on medial aspect of Bottom of Patella to your mark on medial aspect of Ankle Bend. Record length on Line B2.

POSTERIOR LENGTH (C2):

Measure lower leg from your mark on posterior aspect of Bottom of Patella to your mark on posterior aspect of Ankle Bend. Record length on Line C2.

ANTERIOR LENGTH (D2):

Measure lower leg from your mark on anterior aspect of Bottom of Patella to your mark on anterior aspect of Ankle Bend. Record length on Line D2.

STEP 4: MEASURE CIRCUMFERENCES

UPPER: Encircle leg with tape measure at each mark located on lateral aspect between the Top of Patella and Top of Garment, and record each circumference on the corresponding line.

KNEE: Encircle leg with tape measure at Mid Patella and record each circumference on the corresponding line.

LOWER: Encircle leg with tape measure at each mark located on the lateral aspect between Bottom of Patella and Ankle Bend, and record each circumference on the corresponding line.

STEP 5: MEASURE FOOT

1. Measure foot from heel to 1st metatarsal head (or desired boot length) and record length on Line K.
2. Measure Top of Foot from the 3rd Metatarsal Head to Ankle Bend and record length on Line M.
3. Encircle the instep of foot and apex of Heel with a tape measure and record the circumference on Line I.
4. Encircle the foot across the metatarsal heads and record the circumference on Line J.

HIP ATTACHMENT (OPTIONAL)

1. Measure from lateral aspect of Gluteal Fold to lateral aspect of hip at the desired location for the bottom of the belt, and record length on Line H.
2. Measure circumference of torso at the desired location for the bottom of the belt, and record circumference on Line G.

TAKE DIGITAL PHOTOS (REQUIRED) AFTER MEASUREMENT MARKINGS ARE MADE. EMAIL TO: us_orders@sigvaris.com

NOTES: LEGASSIST THIGH HIGH is available in:

7099-TH Flat Foam, 7199-TH WAVEFOAM, 7299-TH Super Flat Foam, 7399-TH Super WAVEFOAM, TH with Hip Attachment is available in 7099-TA Flat Foam, 7199-TA Wave Foam, 7299-TA Super Flat Foam, 7399-TA Super WAVEFOAM (Includes MEDABOOT Custom 7099-MB)

SIGVARIS STANDARD MEASURING

Follow contour of the arm/leg by placing the measuring tape against the skin. Measure in centimeters.

ARM

With patient seated, place the arm extended and elbow bent slightly, with the palm down, on a flat surface.

CIRCUMFERENCES

Measure circumference at Palm and record on line E (*MEDAFIT ARM only*).

Measure at Wrist, and record on line A. Mark **dorsal** aspect at distal edge of tape.

Measure at Elbow, and record on line B.

Measure at mid Bicep, and record on line B1.

Measure at Axilla, and record on line C. Mark **dorsal** aspect at proximal edge of tape.

LENGTH

Measure length at **dorsal** aspect from mark at Point A to mark at Point C and record in box D.

BELOW KNEE

CIRCUMFERENCES

Measure circumferences where indicated by the black dots, beginning with the dot 5cm above the **Lateral Malleolus** (outer ankle bone), continue to the 15cm, 25cm and 30cm dots. Record in the corresponding lines A, B, C, C1(Left/Right).

LENGTH

Measure length from the **Lateral Malleolus** (outer ankle bone) to lateral aspect of the **Fibular Head**, and record in box G.

FOOT

Measure circumference of ankle bend and heel and record on line I. Measure circumference across metatarsal heads and record on line J. Measure length on medial side of foot from heel to 1st metatarsal head and record on line K.

THIGH COMPONENT

CIRCUMFERENCES

Measure circumferences where indicated by the black dots, beginning with the dot 5cm above the lateral aspect **Top of Patella** and continue to the 15cm, and 23cm dots. Record in the corresponding lines D, E, F (Left/Right).

LENGTH (*COMPREFLEX only*)

Measure length from the lateral aspect **Top of Patella** to lateral aspect **Gluteal Fold**, and record in box G1.

THIGH HIGH

CIRCUMFERENCES

Measure circumferences where indicated by the black dots, beginning with the dot 5cm above the **Lateral Malleolus** (outer ankle bone) and continue to the 15cm, 25cm and 30cm dots. Record in the corresponding lines A, B, C, C1(Left/Right). Measure circumferences at the 5cm dot from the lateral aspect **Top of Patella** and continue to the 15cm, and 23cm dots. Record in the corresponding lines D, E, F (Left/Right).

LENGTH

Measure length from the **Lateral Malleolus** (outer ankle bone) to lateral aspect of the **Gluteal Fold**, and record in box H.

FOOT

Measure circumference of ankle bend and heel and record on line I. Measure circumference across metatarsal heads and record on line J. Measure length on medial side of foot from heel to 1st metatarsal head and record on line K.

*For more information, please speak with your
SIGVARIS Sales Representative or Customer Care:
1-800-322-7744.*

Holland, MI:
FDA Registered
All products are 100% inspected.

Customer Care: Monday–Friday, 8am–7pm (EST)

Fax: 1-800-481-5488

Call: 1-800-322-7744, ext. 1

Online: **sigvaris-online.com**

ORDERING

ORDERING

Online: **sigvaris-online.com**

Call: 1-800-322-7744, ext. 1

Fax: 1-800-481-5488

Customer Care: Monday–Friday, 8am–7pm (EST)

QUALITY STANDARD

Peachtree City, GA:

Certified Quality System ISO 13485

FDA Registered

All products are 100% inspected.

sigvarisusa.com

SIGVARISUSA



SIGVARIS

The Worldwide Leader
in Compression Therapy



COMPRESSION SOLUTIONS 2018

INELASTIC WRAPS AND BANDAGING SYSTEMS
FOR LYMPHEDEMA, WOUND CARE &
COMPLEX EDEMAS