

WHAT YOU NEED

SafeStrip



FENTANYL TEST



Substance



SafeStrip Fentanyl Strip



Plastic bottlecap
1tbs (5ml)




10mg micro scoop
Helpful To Have


01 CRUSH AND SHAKE



02 MEASURE AND DILUTE



03 STIR AND TEST



Positive



FYL
FYL

Negative



FYL
FYL

FOR MORE INFORMATION AND INSTRUCTIONS VISIT WWW.SAFESTRIP.US

Tests are for forensic use.