

WHAT YOU NEED

SafeStrip




Substance



SafeStrip Xylazine Test Strip



Plastic bottlecap (5ml)




10mg micro scoop
Helpful To Have


01 CRUSH AND SHAKE



02 MEASURE AND DILUTE



03 STIR AND TEST



TESTS ARE FOR FORENSIC USE. THESE TESTS ARE DESIGNED TO BE USED WITH URINE. USING THIS TEST WITH ANYTHING ELSE IS TO BE CONSIDERED A PRECAUTION AND NOT A GUARANTEE