

blodem





D (mm)	L (mm)					
	6	8	10	12	14	16
3.8						
4.3						
4.8						

D (mm)	L (mm)					
	6	8	10	12	14	16
3.5						





2292

ISO:
13485:2012



Biodem GmbH

Speditionsstrasse 21 40221





Düsseldorf - **Germany**

www.bio-dem.com



www.bio-dem.com

SURGICAL KIT



Ratchet Hex Driver for Implant



Ratchet Hex Driver for Implant



Hand Hex Driver for Abutment long



Drill Extension



Ratchet Hex Driver for Implant



Ratchet Hex Driver for Implant



Hand Hex Driver for Abutment short



Deep Guide



Ratchet Hex Driver for Implant



Ratchet Hex Driver for Implant



Ratchet Hex Driver for Abutment short



2.2



3.2



3.5



3.9



Ratchet Hex Driver for Abutment long

CEMENTED REHABILITATION




CEMENTED REHABILITATION




SCREW-RETAINED REHABILITATION

MOVABLE PROSTHESIS



SCREW-RETAINED REHABILITATION

MOVABLE PROSTHESIS



TRANSFER

ANALOG

HEALING CAP




TRANSFER

ANALOG


HEALING CAP




PROSTHETIC PARTS




SFAT Standart Simante Abutment	NFAT Narrow Platform Simante Abutment
D(mm) = 4.5 L(mm) = 7,9,12,15	D(mm) = 4.0 L(mm) = 9,15




SEAT Straight Abutment	NEAT Narrow Straight Abutment
L(mm) = 1,2,3,5	L(mm) = 1,2,3,5




SBAT Wide Anatomic Abutment	NBAT Narrow Anatomic Abutment
D(mm) = 5.5 L(mm) = 1,2,3,5	D(mm) = 4.5 L(mm) = 1,2,3,5



AEAT Straight Abutment	NAET Narrow Straight Abutment
D(mm) = 4.5 L(mm) = 1,2,3,5	D(mm) = 4.5 L(mm) = 1,2,3,5




BAAT15 Standart Platform 15 Angular Abutment	NAAT15 Narrow Platform 25 Angular Abutment
L(mm) = 1,2,3,5	L(mm) = 1,2,3,5



SBA Ball Abutment	NBA Narrow Ball Abutment
D(mm) = 4.5 L(mm) = 1,2,3,5,7 D1(mm) 2.80	D(mm) = 4.5 L(mm) = 1,2,3,5 D1(mm) 2.20




BAAT25 Standart Platform 25 Angular Abutment	NAAT25 Narrow Platform 25 Angular Abutment
L(mm) = 1,2,3,5	L(mm) = 1,2,3,5



SDLC Locator	NDLC Narrow Locator
L(mm) = 1,2,3,5,7	L(mm) = 1,2,3,5,7

PROSTHETIC PARTS




BSDAT
Plastic Sleeve
for Casting

D(mm) = 4.5
L(mm) = 12
D1(mm) 2.47

BNDAT
Plastic Sleeve
for Casting

D(mm) = 4.0
L(mm) = 12
D1(mm) 2.03




PEEK
Temporary
Abutment

D(mm) = 4.5
L(mm) = 1, 2, 3, 5
D1(mm) 2.47

NPEEK
Temporary
Abutment

D(mm) = 4.5
L(mm) = 1, 2, 3, 5
D1(mm) 2.03




PEEK15
Temporary
Abutment

D(mm) = 4.5
L(mm) = 1, 2, 3, 5
D1(mm) 2.47

NPEEK15
Temporary
Abutment

D(mm) = 4.5
L(mm) = 1, 2, 3, 5
D1(mm) 2.03



GKV
Wide cover Screw

D(mm) = 6.5




BATS
Clips
Transfer

D(mm) = 4.5
L(mm) = 11
D1(mm) 2.47

BNTS
Clips
Transfer

D(mm) = 4.0
L(mm) = 11
D1(mm) 2.03




BSKT
Transfer

D(mm) = 4.5
L(mm) = 11, 15
D1(mm) 2.47

BNKT
Transfer

D(mm) = 4.0
L(mm) = 11, 15
D1(mm) 2.03



BSA
Standart
Analog

D(mm) = 4.0
L(mm) = 12

BSA
Narrow
Analog

D(mm) = 4.0
L(mm) = 12



Bsib
Healing Cap

D(mm) = 4.8
L(mm) = 1, 2, 3, 5, 7
D1(mm) 3.20

BNib
Healing Cap

D(mm) = 4.5
L(mm) = 1, 2, 3, 5
D1(mm) 2.60

biodem



Biodem GmbH

Speditionsstrasse 21 40221



Düsseldorf - **Germany**

www.bio-dem.com