

SPINE FIBO

Automatic Visual Inspection,
Counting and Bottle Filling of
Tablets & Capsules & Softgels



SENSUM



*Computer
Vision
Systems*

SPINE FIBO

**AUTOMATIC VISUAL INSPECTION, COUNTING AND BOTTLE FILLING
OF TABLETS & CAPSULES & SOFTGELS**



HIGHLIGHTS

- Multifunction machine for inspection, counting and bottle filling
- Exceptional 100 % visual inspection of tablets, capsules & softgels
- Product sorting before final counting and filling
- Easy adjustment for bottles of different size and type
- Filling bottles with count from 30 to 1000, and higher
- Capacity up to 165/100/10 bottles per minute at 30/100/1000 count
- Intuitive user-friendly graphical interface
- Ergonomic, compact design for simple cleaning and maintenance

REDUCE COSTS

SPINE FIBO combines the inspection performance of Sensum **SPINE** machine with high speed counting and bottle filling, that can be easily integrated into any new or existing bottle packaging line.




TWO PROCESSES COMBINED IN ONE

- 100% inspection and counting with bottle filling
- Smaller footprint than conventional counters
- Lower number of operators
- Less material handling
- Less material storage



RELIABLE BOTTLE FILLING

- System with no flaps
- Double counting with camera and optical system
- Less rejected bottles: only good products are filled into bottles because defective products are sorted out beforehand
- Bottle filling speed is not decreased due to defects in batch
- Reliable and controlled handling of each individual product
- Reliable filling of very small products (less than 6 mm)



LOWER RISKS



- Removed risk of product damage after inspection and before bottle filling due to material handling
- Reduced risk of human error after inspection
- Removed risk that two products count as one
- Removed risk that a defective product (e.g. laminated tablet) is filled in the bottle
- Reduced clogging at bottle neck: products are filled into a bottle one by one

IMPROVE YOUR QUALITY BY SORTING OUT DEFECTS

ROUND TABLETS



thickness



chips



shape



dark dots



bright dots



colour deviation



chipped coating



edge chips



edge cracks



edge dark dots



engraving defects



print defects

OVAL AND OBLONG TABLETS



lamination



chips



cracks



coating defects



engraving defects



print defects

CAPSULES



holes



dents



joints



length



empty transparent



print defects

SOFTGELS



size



damaged seam



dirt



hair



bubbles



broken


AND MANY MORE ...

INSPECTION OF TABLETS, CAPSULES OR SOFTGELS

INDICATIVE INSPECTION SPEEDS


ROUND TABLETS

Ø (DIAMETER) IN MM	TABLETS / HOUR
13 - 18	210.000
10 - 12	420.000
8 - 9	520.000
< 8	630.000




OVAL AND OBLONG TABLETS

L (LENGTH) IN MM	TABLETS / HOUR
> 18	150.000
14 - 18	230.000
11 - 13	310.000
8 - 10	390.000
< 8	470.000




CAPSULES

SIZE	CAPSULES / HOUR
0 & 1 & 2	200.000
3 & 4 & 5	300.000



SOFTGELS

L (LENGTH) IN MM	SOFTGELS / HOUR
> 15	160.000
14 - 15	280.000
11 - 13	380.000
8 - 10	470.000
< 8	570.000



ACTUAL INSPECTION SPEED depends on the shape, size and type of tablets, capsules and softgels, and on the number of inspected features. In special cases, even higher inspection speeds might be possible with dedicated solution.



BOTTLE FILLING AFTER INSPECTION AND SORTING



INDICATIVE BOTTLE FILLING SPEEDS

COUNT PER BOTTLE	BPM (BOTTLES PER MINUTE)
30	165
60	165
100	100
500	20
1000	10

Actual speeds depend on the inspection speed of the product.

YOUR FILLING SPEED FORMULA:

$$\text{BPM} = \frac{\text{INSPECTION SPEED (PR/H)}}{60 \text{ MIN} \times (\text{COUNT PER BOTTLE})}$$

(See inspection speed page)


Indicative bottle filling speeds are the upper limit for actual performance.

TECHNICAL DATA


APPLICATION	INSPECTION QUALITY	TYPICAL ELECTRICAL CONNECTION	COMPRESSED AIR	DIMENSIONS L x W x H	WEIGHT	CONSTRUCTION
Tablets, capsules and softgels	Same as Sensum SPINE machine	400/230V	7 Bar, 100 Nm ³ /h	L = 2560 mm W = 860 mm H = 2145 mm	1300 kg	CE, GMP, FDA acceptable materials


CONTINUOUS AND RELIABLE BOTTLE FILLING

Sensum **SPINE** FIBO automatically inspects the entire surface of products with 6 visual systems at high speed. The inspected products with defects are removed before bottle filling. The good products are **counted and filled** into bottles. Miscounted bottles are automatically **rejected** from the packaging line.



EFFICIENT BOTTLE FILLING WITHOUT FLAPS





SENSUM



Sensum

*Computer
Vision
Systems*

Tehnološki park 21
1000 Ljubljana
Slovenia, EU

www.sensum.eu
info@sensum.eu



FOLLOW US:

Follow us on LinkedIn for the latest updates on quality control in pharmaceutical and nutraceutical production.

ABOUT SENSUM:

We develop, produce and sell innovative automatic visual inspection solutions for the pharmaceutical industry. We have a strong R&D cooperation with the Imaging Technologies Lab, Department of Electrical Engineering, University of Ljubljana to keep in touch with the latest scientific developments.

