

# NEXT LEVEL PAIN MANAGEMENT



[WWW.SPRINGLIFE.COM](http://WWW.SPRINGLIFE.COM)

Easy-to-use and lightweight  
Pulsed Radio Frequency (PRF)-only  
generator for both transcutaneous  
and invasive treatments.



Clinical Guidelines



Spring2 Video



# NEXT LEVEL PAIN MANAGEMENT

## Quick overview

- PRF-only
- Compact in design and with a grip for easy handling
- Lightweight (only 2.6 kg)
- Intuitive user-friendly interface with touchscreen

## Transcutaneous PRF

- Automatically identified Springlife skin electrodes. For safety reasons the Spring2® recognizes the size of the attached electrodes and the maximum current output is adjusted to these sizes.
- Duty load: 15 msec (3 x 5ms/s)
- Up to 4 skin electrodes simultaneously
- Easy export of treatment details to USB memory stick

## Invasive multi-level PRF

- Treatments can be performed on 3 levels simultaneously
- Intuitive stimulation mode (sensory and motoric)
- Duty load: 15 msec (3 x 5ms/s)
- Easy export of treatment details to USB memory stick

