


DISPOSABLE TIPS FACT SHEET

Woorhi Mechatronics offers a wide range of disposable tips that you can rely on for your precise procedure every time.

Each Disposable Tips (Cartridges) are especially designed to maintain same quality throughout the entire procedure. Our disposable tips are designed to use only one time, however the points of needles are in the exact same position as at the beginning.

SHAPE



DESCRIPTION

Product Name	Description			
Massage Tip	Commonly used to massage face / Scalp / Various blood spot.			
1-Liner	Used for Fine hairlines / eyelines / lashes / small scars (acne) enhancement.			
3-Liner	Commonly used for eyeliner / eyebrows / lipliner / wider lines (fill in) /camouflage / shading.			
4-Flat	Primarily for shading certain areas but also used for hairlines / eyeliner / lip shading / brow foundation.			
5-Round	Extremely versatile needle used for soft color application / eyeliner / lip contour /shading / brow foundation / areola / scars.			
6-Flat	Covers larger areas than 4-flat hairlines / eyeliner / lip shading / brow foundation.			
9-Magnum	Used for Lip-shading / soft color application / eyeliner/ full lips areola designs /camouflage scars / lash enhancement.			
9-MTS 2.0	9-MTS shows effectiveness in Acne scar / Collagen induction / Cellulite treatment / Improving wrinkle / Skin pigmentation /			
9-MTS 0.2	Skin scar and Alopecia / Skin rejuvenation / Hyperpigmentation / Skin tone / Face lifting / Stretch mark / Skin pore.			

SPECIFICATION

Product Name	Weight	Length	Widht	Needle Gauge
Massage Tip	3.6g	70mm	Ø13mm	N/A
1-Liner	2.4g	62mm	Ø13mm	27ga
3-Liner	2.7g	62mm	Ø13mm	31ga
4-Flat	2.7g	62mm	Ø13mm	30ga
5-Round	2.7g	62mm	Ø13mm	31ga
6-Flat	2.7g	62mm	Ø13mm	30ga
9-Magnum	2.7g	62mm	Ø13mm	30ga
9-MTS 2.0	3.8g	76mm	Ø13mm	33ga
9-MTS 0.2	4.0g	76mm	Ø13mm	33ga



#301, Daekwang Plaza, 1027-15,Hogae-Dong,DongAn-Gu, Anyang-City, Kyunggi-Do, Korea #431-848 Tel: +8231-457-5032 / Fax: +8231-457-5031

E-mail: saleswrmc@gmail.com Http: www.woorhimecha.com