

NOBLE FACT SHEET

Woorhi Mecha's Auto MTS technology (NOBLE) is Most advance technology yet!

Noble is supplemental device that is ideal for non-surgical and non-ablative treatment of various skin conditions such as aging (wrinkles, stretching), scarring (acne, surgical), and hyperpigmentation.

Noble is also far more advance than old fashioned MTS (Roller Type). Clinical studies have shown AutoMTS (Automatic Micro-Needle Therapy System) to be more effective than ablative treatments like laser resurfacing, dermabrasion, and chemical peel.

Noble is also just as effective as non-ablative treatments like IPL, CO2 Laser, and Fraxel in stimulating collagen and elastin production to thicken the skin thereby erasing wrinkles and smoothing scars.

EFFECTS

- Improving wrinkle
- Skin rejuvenation
- Face lifting

- Skin pigmentation
- Skin tone
- Skin pore

BENEFITS

FEATURES	BENEFITS	
Automatic	Very easy to use, and also safe to use since the needle tips are	
Automatic	hidden inside the guide.	
	Ergonomically-designed Auto-MTS can help alleviate	
Ergonomic design	the discomfort caused by Repetitive Strain Injury (RSI),	
	Carpal Tunnel Syndrome (CTS)	
Safe using disposable	It is safe to share the hand-piece with other patient.	
needle tips	Be sure to change sterilized safety 9 needles every use.	
Economical	You do not have to buy a new device every time (like MTS).	
	Just replace the needle!	
Less Pain	The needle tips are inserted to your skin	
	perpendicularly. Therefore, pain is minimized.	
	It is more efficient & effective to treat narrow areas such as nose,	
Efficient & Effective	around eyes, or mouth, this hand-piece device allows you to treat	
	areas in whichever angles you want to.	

Improvements with NOBLE (AUTO-MTS)



SPEED & DEPTH


Level	Speed	Level	Speed
1	4200rpm (65~70Times/sec)	3	5400rpm (85~90Time/sec)
2	4800rpm (75-80Times/sec)	4	6000rpm (95~100Time/sec)

MTS ROLLER VS AUTO MTS NOBLE

	MTS Roller	Auto-MTS			
Areas	Difficult to operate in small or curved areas	Easy to maneuver around all facial contours Easy on nose near eyes			
Pain	Can be VERY painful, with some bleeding due to epidermal damage	Minimized epidermal damage means less pain and bleeding 90% less pain than roller			
Effects	Scarring, Enlarged pores, Wrinkles, Lack of skin firmness and elasticity, Acne and chicken pox scars, Under-eye bags, Sun damage and age spots, facial whitening, prevents hair loss.				
Application	Difficult to insert needles into scar tissue - needles can be bent	Because of vertical delivery, needles easily penetrate scar tissue without risk of bending			
Sanitary	Sterilize the roller every use and after use	Using sterilized disposable needles makes it easy and clean			



Auto-MTS



SPECIFICATION

PRODUCT	NOBLE	BELLA	RAFFINE	
INPUT	AC 100~250V			
OUTPUT	DC 5V			
WEIGHT	74.3g	145.3g	78.86g	
FREQUENCY	50~60HZ			
DIMENSION(WxHxD)	28x147x12mm	20X118X12mm	28X147X12mm	
OPERATION MODE	AUTO MTS	AUTO MTS /PERMANENT MAKE UP	AUTO MTS /PERMANENT MAKE UP	
NEEDLE DEPTH	0.2mm(FIXED)	0.1~2mm	0.1~2mm	



#301, Daekwang Plaza, 1027-15,Hogae-Dong,DongAn-Gu, Anyang-City, Kyunggi-Do, Korea #431-848 Tel: +8231-457-5032 / Fax: +8231-457-5031

E-mail: saleswrmc@gmail.com Http://www.woorhimecha.com