

Vantage *Plus* Scaler with Motor Pack



Experience the
Vantage *Plus*
Scaler with Motor Pack.

Focused LED Illumination
in both hand pieces. *pats. pending*

DELMARVA 2000™

**20+ Years of Distinction and
Superior Customer Service.**

MADE IN THE USA 

- **New** Scaler Handpiece with Focused LED Illumination *pat. pending*
- Automatic tuning recognizes and locks in on each insert.
- Powerful 25k magnetostrictive frequency.
- Programmed power control knob provides the optimum energy for all inserts.
- **Exclusive 30%** power increase setting for that “extra boost” to remove tough calculus.
- Specialty insert power setting for feline urethral catheters and other specialty inserts such as the periodontal insert.
- Unique low water pressure warning light.
- **New** Micro Motor handpiece with Focused LED Illumination *pat. pending*
- Micro Motor has variable speeds allowing for a wide RPM range for polishing, sectioning or drilling.
- Precision motor pack hand piece is a sealed unit...it generates no heat and requires **NO Lubrication**.
- Polish...Drill...Cutting...Rounding... Fine surface procedures.
- **New** Reversing Switch to untangle hair.
- Unit is self-contained-does not require a separate transformer pack for power. Small foot print.
- Indestructible foot pedal.
- Detachable hand pieces.
- Made in the USA with 3 year warranty.

“Lighting the Way”

Vantage Plus Scaler with Motor Pack

Scaling Inserts

Item#	Description
DE-5	TFI #3 Insert - 25k (Thru Flow Insert) Included with all dental units
DE-6	TFI #10 Insert - 25k (Thru Flow Insert)
DE-6A	SLI-10-S Insert - Periodontal - 25k (Focus Spray Slim Line Insert)
DE-7	Feline Urethral Catheter



Accessories

Item#	Description
DE-9	Prophy Angle
DE-15	Latch Type Contra Angle
DE-8	Prophy Cups (Snap on, bag 144)
DE-10	Prophy Paste (bag of 200)
DE-12	Dental Accessory Pack (6 burs, 1 disk, and contra angle)
DE-16B	Inline Water Filters (pack of 10)



Items pictured are not actual size.

DELMARVA 2000™

Milton, Delaware

302-645-2226

www.delmarva2000.com

email: customerservice@delmarva2000.com

“Lighting the Way”