

DELMARVA 2000

BiPolar Veterinary HFR

High Frequency Radiosurgery Device

The Delmarva 2000 Bipolar Veterinary HFR will easily meet the demanding needs of the discerning veterinarian. 125+ watts of power and 4mhz in radio frequency provides easy, fast, and safe surgery. The multi-function of bipolar and monopolar will afford any practice a wide range of procedures from removal a minor lesion to squamous cell carcinoma and anything in between.

- ✓ Totally isolated mono-polar configuration allowing for **Fine, Normal, and Coagulation** modes—with immediate **BiPolar** Switching.
- ✓ Isolated circuitry is compatible with all monitors
- ✓ 3 easy selectable modes of operation
 - ◆ **Cut** Fully rectified & filtered for smooth cutting and minimal hemostasis
 - ◆ **Blend** Fully rectified & not filtered for good cutting coupled with hemostasis
 - ◆ **Coag** Radio wave is modified for better coagulation controls bleeding and destroys unwanted tissue by desiccation
- ✓ Instant switching to **bipolar**, for pinpoint bleeder control or invert forceps for larger areas of coagulation
- ✓ Linear adjustable power control, assuring the correct power, whether patient is a pocket pet or equine.
- ✓ Autoclavable hand piece and cords
- ✓ Industrial grade hybrid electronics for stable, cool and safe operation
- ✓ Simple easy to use control panel
- ✓ HFR comes plug and play with all accessories and selection of electrodes.
- ✓ 3 year warranty

