



Radial Runway

Radial Wrist Positioning

The Radial Runway provides a radiolucent platform that keeps the patient's wrist in a deflected position for ideal radial access.

Disposable straps allow the wrist to be positioned in multiple neutral or extended angles, optimizing patient comfort and physician preference.

The foam pad and straps remain with the patient post-procedure, acting as an immobilizer for the patient's wrist during recovery.

Can be used for radial access on either wrist

Wrists may be positioned at multiple angles

Adjustable straps fit any size patient

Disposable straps, foam pad & brace

Keeps patient comfortable during the procedure & recovery

Latex Free

Use with Radial Board



Radial Runway Base
RRU-0008

Top Pads, Straps & Brace
RRU-0009



01925 387 121



Info@onetwomedical.com



Onetwomedical.com