



The Bedside Commode



This information is a general guideline. Please consult a healthcare provider for instructions specific to your situation.

A bedside commode can be used as a movable toilet or can be placed over the toilet to provide transfer support and a raised seat

Adjusting the leg height

- Press the snap button on the leg extension and slide the leg to the appropriate hole
- Make sure all buttons are fully engaged and adjusted to the same height
- After adjusting the height, test the stability of the commode by pushing down on the arm rests
- The user's feet should always be flat on the floor when seated on the commode. Do not adjust the legs too high for safe transfer