

Cardinal Health™ N250
Powder-free nitrile examination gloves



Standard protection nitrile exam glove

Cardinal Health™ N250 offers protection against chemicals and chemotherapy drugs, with more gloves per box.

Features and benefits:

- Standard level protection within our High Use portfolio
- Tested for use with chemotherapy drugs and chemicals^{1,2}
- USP 800 compliant*
- 25% more than FLEXAL, reducing packaging waste and frequency of box changes³



* This Cardinal Health product helps healthcare facilities comply with USP <800> standard for the safe handling of hazardous drugs.

High Use Nitrile:

non-sterile exam gloves designed to serve a broad range of everyday clinical needs.

Cardinal Health™ N250 Powder-free nitrile examination gloves



Cat. No.	Size	Pkg.
88LC01XS	X-Small	250/bx, 10 bx/cs
88LC02S	Small	250/bx, 10 bx/cs
88LC03M	Medium	250/bx, 10 bx/cs
88LC04L	Large	250/bx, 10 bx/cs
88LC05XL	X-Large	230/bx, 10 bx/cs
88LC06XXL	XX-Large	230/bx, 10 bx/cs

- Powder-free
- Not made with natural rubber latex
- Ambidextrous
- Color: Blue
- Nominal fingertip thickness: 3.7 mil (.1 mm)

For more information or to order, contact your Cardinal Health sales representative. Call 877.CARDINAL or visit cardinalhealth.com.

1 This glove has been tested for resistance to permeation of various chemotherapy drugs per ASTM D6978, "Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs." Caution: Gloves used for protection against chemotherapy drug and laboratory chemical exposure should be selected specifically for the type of chemicals being used.

Note: Review Safety Data Sheets (SDS) for the chemicals/drugs being used to determine required level of protection.

2 ASTM F739 - Standard Test Method for Permeation of Liquids and Gases through Protective Clothing Materials under Conditions of Continuous Contact

3 25% more N250 gloves per dispenser box when compared to FLEXAL™. 1000 gloves of N250 require 8 percent less storage space and use 9 percent less packaging material when compared to 1000 gloves of FLEXAL™.

