

SWISS PIEZOCLAST®

→ IN DEPTH WITH FOCUSED SHOCK WAVE THERAPY




EMS+

SWISS PIEZOCLAST®

> Treat deep musculoskeletal indications efficiently and comfortably with the EMS Focused ESWT solution.



**FOCUSED
ESWT**



ACCURATE →

> Deliver energy only where needed thanks to focused shock waves

DEEP →

> Treat pathologies down to 80 mm

PAIN FREE →

> Obtain comfortable treatments thanks to piezoceramic generated focused shock waves

PRODUCT REFERENCE

FT-204

TECHNOLOGY

Focused ESWT piezoelectric technology

POSITIVE ENERGY DENSITY (ED₊)

ED₊: 0.41 mJ/mm² - ED TOT: 0.82 mJ/mm²

CAVITATION LEVEL

Increases with intensity

FREQUENCY

Up to 8 Hz

STANDARD CONFIGURATION

PiezoClast® control unit, Therapy Source handpiece, gel pads, 3 pedal footswitch

EMS-MEDICAL.COM

For more information: welcome@ems-ch.com

**"I FEEL
GOOD"**