




WAHL CLIPPER CORPORATION | STERLING, IL 61081 | www.wahl.com



Insert comfort band end into position A and out position B as in the diagram above.

After assembly, place shield on head and measure for appropriate fit, just above the ear.

Once sized, insert shield strap tab in same method as used for the comfort band.

The shield should fit securely and comfortably with full face coverage. Adjust as desired.

Shield is intended for a single user, single use.

IMPORTANT INFORMATION

This face shield has not been FDA cleared or approved. This product has been authorized by FDA under an EUA for use by healthcare providers as personal protective equipment. It is only authorized for the duration of the declaration that circumstances justifying the authorization of emergency use under Section 564(b)(1) of the Act, 21 USC 360bbb-3(b)(1) unless the authorization is terminated or revoked sooner. The face shield is not intended for use in the presence of a high intensity heat source or flammable gas, use for antimicrobial or antiviral protection, use for infection prevention, or use for radiation protection.