




## Quickly and easily simulate “fuzz” and pill in carpet samples

The Pilfuz Carpet Test is used to simulate the “fuzz” and pill that can be associated with carpets under normal use.

The unit incorporates a velcro covered roller which interacts with the surface fibre of the sample for a pre-determined number of cycles. The samples are then visually assessed.



### Key Features

- ✓ *Quick and easy to use*
- ✓ *Instant results*
- ✓ *Can be used in the lab or in production facility*
- ✓ *No power required*



Consumables – PFZ:VEL – 5 meter roll of Velcro.

Order Code: PFZ:001