

### Cross-Reference Chart - HemaClear to Pneumatic Cuff Sizes

HemaClear Part Number	HemaClear Size - Adult	Pneumatic Cuff Size
PRH-028-PI-01A	HemaClear, Small Pink	8 inch Reusable or Disposable Cuff
PRH-028-PI-01A	HemaClear, Small Pink	12 inch Reusable or Disposable Cuff
PRH-040-YE-01A	HemaClear, Medium Yellow	18 inch Reusable or Disposable Cuff
PRH-060-OR-01A	HemaClear, Large Orange	24 inch Reusable or Disposable Cuff
PRH-090-BW-01A	HemaClear, Extra Large Black & White	30 inch Reusable or Disposable Cuff
PRH-090-BW-01A	HemaClear, Extra Large Black & White	34 inch Reusable or Disposable Cuff
PRH-090-BW-01A	HemaClear, Extra Large Black & White	42 inch Reusable or Disposable Cuff

HemaClear Part Number	HemaClear Size - Pediatric	Pneumatic Cuff Size
PRH-028-PI-01A	HemaClear, Small Pink	8 inch Reusable or Disposable Cuff
PRH-028-PI-01A	HemaClear, Small Pink	12 inch Reusable or Disposable Cuff
PRH-040-GR-01A	HemaClear, Medium Green	18 inch Reusable or Disposable Cuff
PRH-060-BL-01A	HemaClear, Large Blue	24 inch Reusable or Disposable Cuff
PRH-090-BW-01A	HemaClear, Extra Large Black & White	30 inch Reusable or Disposable Cuff
PRH-090-BW-01A	HemaClear, Extra Large Black & White	34 inch Reusable or Disposable Cuff
PRH-090-BW-01A	HemaClear, Extra Large Black & White	42 inch Reusable or Disposable Cuff

Note: HemaClear sizes may appear to be smaller than pneumatic cuff sizes; however the pneumatic cuff has a 4 to 8 inch Velcro overlap requirement to maintain pressure. The elastic HemaClear will stretch to accommodate the limb circumference size range shown above. Use a HemaClear sterile measuring tape to confirm proper size.