

HemaClear[®]
Sterile & Bloodless Surgical Field


Compare HemaClear Exsanguination Tourniquets to Pneumatic Cuffs

HemaClear[®] - a great alternative to pneumatic tourniquet systems. HemaClear significantly reduces tourniquet pain and post operative complications of tourniquet use.

All-in-One design, replaces Esmarch bandage exsanguination, limb stockinette, pneumatic cuff, hoses and pump. Always a dry sterile operating field.



Bilateral TKA with pneumatic tourniquet on right leg and HemaClear on left leg



Note the dryness of the HemaClear side

Cross Reference Charts		
Pneumatic Cuff Size	HemaClear Size - Adult	HemaClear Part Number
8 Inch Reusable or Disposable Cuff	HemaClear, Small, Pink	PRH-028-PI-01A
12 Inch Reusable or Disposable Cuff	HemaClear, Small, Pink	PRH-028-PI-01A
18 Inch Reusable or Disposable Cuff	HemaClear, Medium, Yellow	PRH-040-YE-01A
24 Inch Reusable or Disposable Cuff	HemaClear, Large, Orange	PRH-060-OR-01A
30 Inch Reusable or Disposable Cuff	HemaClear, Extra Large, Black & White	PRH-090-BW-01A
34 Inch Reusable or Disposable Cuff	HemaClear, Extra Large, Black & White	PRH-090-BW-01A
42 Inch Reusable or Disposable Cuff	HemaClear, Extra Large, Black & White	PRH-090-BW-01A
Pneumatic Cuff Size	HemaClear Size - Pediatric	HemaClear Part Number
8 Inch Reusable or Disposable Cuff	HemaClear, Small, Pink	PRH-028-PI-01A
12 Inch Reusable or Disposable Cuff	HemaClear, Small, Pink	PRH-028-PI-01A
18 Inch Reusable or Disposable Cuff	HemaClear, Medium, Green	PRH-040-GR-01A
24 Inch Reusable or Disposable Cuff	HemaClear, Large, Blue	PRH-060-BL-01A
30 Inch Reusable or Disposable Cuff	HemaClear, Extra Large, Black & White	PRH-090-BW-01A
34 Inch Reusable or Disposable Cuff	HemaClear, Extra Large, Black & White	PRH-090-BW-01A
42 Inch Reusable or Disposable Cuff	HemaClear, Extra Large, Black & White	PRH-090-BW-01A
NOTE 1: HemaClear sizes may appear to be smaller than pneumatic cuff sizes, however the pneumatic cuff has a 2 to 6 inch Velcro overlap requirement. The elastic HemaClear will stretch to accommodate the limb circumference size range shown above. Use the HemaClear sterile measuring tape to confirm proper size as		
NOTE 2: HemaClear All-in-One design includes exsanguination. No separate Esmarch Bandage purchase is required.		