











| Product | | Cat# | Units/ Case* | Blood Pressure | Skin Pressure | Circumference |
|-------------------------------------|---|----------------|-----------------|-------------------|---------------|---------------|
| SMALL Pink |  | PRH-028-PI-01A | 10 | < 130 mmHg | 197±24 mmHg | 14-28 cm |
| MEDIUM Yellow |  | PRH-040-YE-01A | 10 | < 190 mmHg | 254±43 mmHg | 24-40 cm |
| MEDIUM Red |  | PRH-040-RE-01A | 10 | < 160 mmHg | 246±42 mmHg | 24-40 cm |
| MEDIUM Green |  | PRH-040-GR-01A | 10 | < 130 mmHg | 234±42 mmHg | 24-40 cm |
| LARGE Brown |  | PRH-060-BR-01A | 10 | < 190 mmHg | 364±49 mmHg | 30-55 cm |
| LARGE Orange |  | PRH-060-OR-01A | 10 | < 160 mmHg | 329±47 mmHg | 30-55 cm |
| LARGE Blue |  | PRH-060-BL-01A | 10 | < 130 mmHg | 309±40 mmHg | 30-55 cm |
| X-LARGE- Black&White |  | PRH-090-BW-01A | 10 | < 160 mmHg | 320±19 mmHg | 50-85 cm |
| MODEL F |  | PRH-035-FA-01A | 25 | < 160 mmHg | 248±33 mmHg | 14-34 cm |
| MODEL A |  | PRH-032-MA-01A | 12 | < 160 mmHg | 176±9 mmHg | 22-32 cm |

*HemaClear® is sold per case