



Part Number	Application	Description	Limb Circumference Size	Blood Pressure	Unit of Measure	Pneumatic Tourniquet Cuff Size Equivalent
PRH-028-PI-01-A	Pediatric	HemaClear® SMALL - PINK	14-28 cm	<130 mmHg	1 Case (10 Units)	8 to 12 inch tourniquet cuff ¹
PRH-040-GR-01-A	Pediatric	HemaClear® MEDIUM - GREEN	24-40 cm	<130 mmHg	1 Case (10 Units)	10 to 18 inch tourniquet cuff
PRH-060-BL-01-A	Pediatric	HemaClear® LARGE - BLUE	30-60 cm	<130 mmHg	1 Case (10 Units)	12 to 24 inch tourniquet cuff
PRH-040-RE-01-A	Ped or Adult	HemaClear® MEDIUM - RED	24-40 cm	<160 mmHg	1 Case (10 Units)	10 to 18 inch tourniquet cuff
PRH-060-OR-01-A	Ped or Adult	HemaClear® LARGE - ORANGE	30-60 cm	<160 mmHg	1 Case (10 Units)	12 to 24 inch tourniquet cuff
PRH-090-BW-01-A	Ped or Adult	HemaClear® EX LARGE - BLACK&WHITE	50-85 cm	<160 mmHg	1 Case (10 Units)	20 to 34 inch tourniquet cuff
PRH-040-YE-01-A	Adult	HemaClear® MEDIUM - YELLOW	24-40 cm	<190 mmHg	1 Case (10 Units)	10 to 18 inch tourniquet cuff
PRH-060-BR-01A	Adult	HemaClear® LARGE - BROWN	30-60 cm	< 190 mmHg	1 Case (10 Units)	12 to 24 inch tourniquet cuff

1. HemaClear sizes may appear to be smaller than a pneumatic cuff size, however the pneumatic cuff has a 2 to 6 inch Velcro overlap requirement. The elastic HemaClear will stretch to accommodate larger limb circumferences. So there are fewer sizes of HemaClear. Use HemaClear sterile measuring tape to confirm the proper size.
2. When the sizes on the tape overlap, use the larger size.