



**Hydrogen Generator  
GGH301, 302**



## SALIENT FEATURES

- High Purity Hydrogen can be Produced by Turning On The Power
- Easy To Operate
- Automatic Display Pressure
- High Electrolysis Efficiency
- Purity & Flow Rate Will Not Decay
- Stable Output Pressure
- Digital Display Of Gas Flow Rate

## Hydrogen Generator GGH301, GGH302



HSN Code : 84051010

## SPECIFICATIONS

| MODEL                     | GGH301            | GGH302            |
|---------------------------|-------------------|-------------------|
| Type                      | Alkaline          | Pure Water        |
| Fluid Tank Capacity       | 1.2L              | 2L                |
| Make                      | LABMAN            |                   |
| Maximum Flow              | 300mL/min         |                   |
| Gas Purity                | 99.999%           |                   |
| Output Gas Pressure       | 0.3MPa            |                   |
| Operation Temperature     | 10°C-40°C         |                   |
| Output Flow Display       | LED Display       |                   |
| Power Wattage             | 150W              |                   |
| Phase                     | 1 Phase           |                   |
| Power supply              | AC 220V / 50Hz    |                   |
| Packing Dimension (LxWxH) | 450 x 340 x 500mm | 445 x 340 x 500mm |
| Gross Weight              | 17kg              | 18kg              |



★ Specifications are subject to change without notice.

**STANDARD ACCESSORIES** : Main Instrument, O-Ring, Power Cord, Instruction Manual, Test Certificate.  
Potassium Hydroxide(KOH)120g -(GGH301)

## LABMAN Product Range 2025-26

|                                                                       |                                                           |                                                      |                                                    |                                                        |
|-----------------------------------------------------------------------|-----------------------------------------------------------|------------------------------------------------------|----------------------------------------------------|--------------------------------------------------------|
| <p>LMFD61, 62, 63, 64</p> <p>Lyophilizer Freeze Dryer (-60°C)</p>     | <p>LMFD83, 84</p> <p>Lyophilizer Freeze Dryer (-80°C)</p> | <p>LMRE25</p> <p>Rotary Evaporator</p>               | <p>ADM3</p> <p>Automatic Density Meter</p>         | <p>MSOD6</p> <p>Microvolume Spectrophotometer</p>      |
| <p>NAE32</p> <p>Nucleic Acid Extractor</p>                            | <p>EMR4</p> <p>Elisa Microplate Reader</p>                | <p>GPTC96</p> <p>Gradient PCR Thermal Cycler</p>     | <p>HO30</p> <p>Hybridization Oven</p>              | <p>UVCL254</p> <p>UV - Crosslinker</p>                 |
| <p>TGH12, 24</p> <p>Tissue Grinder Homogenizer</p>                    | <p>PBH80T</p> <p>Paddle Blender Homogenizer</p>           | <p>TOC2NDir</p> <p>Total Organic Carbon Analyzer</p> | <p>TOC1CD</p> <p>Total Organic Carbon Analyzer</p> | <p>LMDS600, 1200</p> <p>Microwave Digestion System</p> |
| <p>HPLC1IS, HPLC2GS</p> <p>High Performance Liquid Chromatography</p> | <p>APT1800</p> <p>Automatic Potential Titrator</p>        | <p>KFT250VC</p> <p>Karl Fischer Titrator</p>         | <p>COD1500</p> <p>COD Analyzer with Digestor</p>   | <p>BOD4000</p> <p>BOD Analyzer</p>                     |
| <p>MTVM50</p> <p>Multi Tube Vortex Mixer</p>                          | <p>LMVM60D</p> <p>Vortex Mixer</p>                        | <p>OS20</p> <p>Orbital Shaker</p>                    | <p>TRM6</p> <p>Tube Roller Mixer</p>               | <p>DBNE12, 24</p> <p>Dry Bath Nitrogen Evaporator</p>  |

## LABMAN Product Range 2025-26




|                                                                                                                                               |                                                                                                                                               |                                                                                                                                           |                                                                                                                                                                                        |                                                                                                                                                                   |
|-----------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>DBI11T</b><br><br><b>Dry Bath Incubator</b>               | <b>DBI15SH, DBI16SC</b><br><br><b>Dry Bath Incubator</b>     | <b>IRS14, 35</b><br><br><b>Infrared Sterilizer</b>       | <b>GGA200</b><br><br><b>Air Generator</b>                                                           | <b>GGA301, 302</b><br><br><b>Hydrogen Generator</b>                            |
| <b>GGN303</b><br><br><b>Nitrogen Generator</b>               | <b>GHA304</b><br><br><b>Hydrogen Air Generator(2 in 1)</b>   | <b>LMVP10</b><br><br><b>Vacuum Pump</b>                  | <b>LMVP20</b><br><br><b>Vacuum Pump</b>                                                             | <b>LMCO1</b><br><br><b>Column Oven</b>                                         |
| <b>LMIO6</b><br><br><b>Mini Thermostatic Incubator</b>      | <b>GHP400</b><br><br><b>Graphite Hotplate</b>               | <b>LOBMS5</b><br><br><b>Oil Bath Magnetic Stirrer</b>   | <b>LMMS4LC, 20LC</b><br><br><b>Ceramic Hotplate Magnetic Stirrer</b>                               | <b>LMMS300</b><br><br><b>Hotplate Magnetic Stirrer</b>                        |
| <b>MMS15S</b><br><br><b>Multiposition Magnetic Stirrer</b> | <b>MMS12H</b><br><br><b>Multiposition Magnetic Stirrer</b> | <b>WBMS3</b><br><br><b>Water Bath Magnetic Stirrer</b> | <b>LMWB2H, 4H, 6H, 8H</b><br><br><b>Digital Water Bath</b>                                         | <b>HOS20L</b><br><br><b>Overhead Stirrer</b>                                 |
| <b>LMHH160</b><br><br><b>Hand Held Homogenizer</b>         | <b>PS1565</b><br><br><b>Peltier / Sipper System</b>        | <b>LUV6HG</b><br><br><b>PMT Spectrophotometer</b>      | <b>LUV2000T, 2100T</b><br><br><b>Single / Double Beam UV-Vis.<br/>Spectrophotometer with 21CFR</b> | <b>LMSPUV SERIES</b><br><br><b>Double Beam<br/>UV-Vis. Spectrophotometer</b> |



## LABMAN Product Range 2025-26

|                                                                                                                                                         |                                                                                                                                                    |                                                                                                                                                                |                                                                                                                                                        |                                                                                                                                                                |
|---------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>LMSPUV1200</b></p>  <p>UV-Vis. Single Beam Spectrophotometer</p> | <p><b>LMSPV320, 325</b></p>  <p>Visible Spectrophotometer</p>     | <p><b>ADP100CFR, 95CFR</b></p>  <p>Automatic Digital Polarimeter with CFR</p> | <p><b>ADP45, 50</b></p>  <p>Automatic Digital Polarimeter (±45°)</p> | <p><b>ADP75, 90</b></p>  <p>Automatic Digital Polarimeter (±90°)</p>        |
| <p><b>ADP40</b></p>  <p>Automatic Digital Polarimeter</p>              | <p><b>LMPM180</b></p>  <p>Manual Polarimeter</p>                  | <p><b>MAR100</b></p>  <p>Manual ABBE Refractometer</p>                        | <p><b>RFM100CFR</b></p>  <p>Digital Refractometer CFR</p>            | <p><b>RFM950, 970</b></p>  <p>Digital Refractometer</p>                     |
| <p><b>LMAR1317T</b></p>  <p>ABBE Refractometer</p>                    | <p><b>LMDV300H</b></p>  <p>Digital Viscometer with Thermosel</p> | <p><b>LMDV60, 100, 200</b></p>  <p>Digital Viscometer ( L Series )</p>       | <p><b>LMDV400R</b></p>  <p>Digital Viscometer ( R Series )</p>     | <p><b>LMKV2</b></p>  <p>Krebs Viscometer</p>                               |
| <p><b>LMDV111, 222</b></p>  <p>Digital Viscometer with SSA</p>       | <p><b>LVA30</b></p>  <p>Low Viscosity Adaptor</p>               | <p><b>LMCR80F</b></p>  <p>Visual Colorimeter</p>                            | <p><b>LMUC2D</b></p>  <p>Ultrasonic Cleaner</p>                    | <p><b>LMUC3D, 6D, 9D, 12D, 25D</b></p>  <p>Ultrasonic Cleaner - Degas</p> |
| <p><b>LMUC50, 100</b></p>  <p>Industrial Ultrasonic Cleaner</p>      | <p><b>LMPH15</b></p>  <p>pH Meter with Printer</p>              | <p><b>LMPH12, LMPH10H</b></p>  <p>pH Meter</p>                              | <p><b>LMPH9</b></p>  <p>pH Meter - Basic</p>                       | <p><b>LMCM20H</b></p>  <p>Conductivity / TDS Meter</p>                    |

## LABMAN Product Range 2025-26

|                                                                                                                                                         |                                                                                                                                                           |                                                                                                                                                   |                                                                                                                                                                      |                                                                                                                                                       |
|---------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>LMCM20H</b></p>  <p><b>Conductivity / TDS Meter</b></p>          | <p><b>LMMP30, 32</b></p>  <p><b>Multiparameter (7 in 1 / 9 in 1)</b></p> | <p><b>LMION40</b></p>  <p><b>Ion Meter</b></p>                   | <p><b>LMDO50</b></p>  <p><b>Dissolved Oxygen Meter</b></p>                         | <p><b>LMHD60</b></p>  <p><b>Water Hardness Meter</b></p>           |
| <p><b>OPH11</b></p>  <p><b>Online pH Controller</b></p>                | <p><b>OCM21</b></p>  <p><b>Online Conductivity Controller</b></p>        | <p><b>LMTB200</b></p>  <p><b>Turbidity Meter</b></p>             | <p><b>RCB1010T, 1020T, 1030T</b></p>  <p><b>Refrigerated Circulating Bath</b></p> | <p><b>RCB620T</b></p>  <p><b>Refrigerated Circulating Bath</b></p> |
| <p><b>RWC320</b></p>  <p><b>Recirculating Water Chiller</b></p>       | <p><b>LMDO4P</b></p>  <p><b>Portable DO Meter</b></p>                   | <p><b>LMPH1P</b></p>  <p><b>Portable pH / ORP Meter</b></p>     | <p><b>LMCM2P</b></p>  <p><b>Portable Conductivity / TDS Meter</b></p>             | <p><b>LMTB5P</b></p>  <p><b>Portable Turbidity Meter</b></p>      |
| <p><b>PRO656, 900, 1200, 1800</b></p>  <p><b>Probe Sonicator</b></p> | <p><b>PRO660R</b></p>  <p><b>Probe Sonicator with Remote</b></p>       | <p><b>PRO650</b></p>  <p><b>Probe Sonicator</b></p>            | <p><b>LMIF30T, 50T, 100T</b></p>  <p><b>Ice Flaker</b></p>                      | <p><b>LMIF200T</b></p>  <p><b>Ice Flaker</b></p>                 |
| <p><b>LMIM15</b></p>  <p><b>Ice Maker</b></p>                        | <p><b>LMIM20C</b></p>  <p><b>Ice Maker with Crusher</b></p>            | <p><b>ISE19, 62</b></p>  <p><b>Nitrate / Fluoride ISEs</b></p> | <p><b>ER SERIES</b></p>  <p><b>pH Electrode / DO Electrode</b></p>               | <p><b>CC1, CC2, CC3</b></p>  <p><b>Conductivity Cell</b></p>     |