




**Refrigerated Circulating Bath  
RCB1010T, 1020T, 1030T**





Touchscreen

## Refrigerated Circulating Bath

RCB1010T, 1020T, 1030T

### SALIENT FEATURES

- Touchscreen Display
- PT100 temperature sensor.
- Both sides folding handle for easy movement
- Overheat protection, PID temperature control system
- Stainless steel inner tank, Waterlevel alarm function
- Fast cooling, low noise, low vibration, stable and reliable performance
- Pump flow rate 8L/min. & Programmable temperature controller

• **Application** : To maintain the temperature, it can be attached in below instruments.

- Fermentation tank
- Spectrophotometer
- Rotary Evaporator
- Polarimeter
- Refractometer
- Viscometer
- Electrophoresis System
- Chromatographic Column



E - Catalogue

HSN Code : 84198990

### SPECIFICATIONS

| MODEL                          | RCB1010T                                    | RCB1020T          | RCB1030T          |
|--------------------------------|---------------------------------------------|-------------------|-------------------|
| Temperature Range              | -10°C to 100°C                              | -20°C to 100°C    | - 30°C to 100°C   |
| Make                           | LABMAN                                      |                   |                   |
| Temperature Display resolution | 0.01°C                                      |                   |                   |
| Temperature Accuracy           | ±0.05°C                                     |                   |                   |
| Temperature Control            | Microprocessor PID                          |                   |                   |
| Display                        | Touchscreen                                 |                   |                   |
| Refrigeration                  | Air Cooled refrigeration through compressor |                   |                   |
| Refrigerant                    | R-134A                                      |                   |                   |
| Pump circulation               | Internal & External                         |                   |                   |
| Temperature sensor             | PT100                                       |                   |                   |
| Pump flow rate                 | 8 L /min                                    |                   |                   |
| Bath Tank capacity             | 10L                                         |                   |                   |
| Tank size (LxWxD )             | 270 x 200 x 200mm                           |                   |                   |
| Bath Opening                   | 180 x 140 mm                                |                   |                   |
| Bath depth                     | 200 mm                                      |                   |                   |
| Power Requirements             | AC 220V / 50Hz                              |                   |                   |
| Packing Dimensions (LxWxH)     | 530 x 380 x 790mm                           | 535 x 380 x 800mm | 580 x 505 x 940mm |
| Gross Weight                   | 41kg                                        | 41kg              | 63kg              |

Liquid to be used : 5°C and below Alcohol, 5°C ~85°C, Pure Water, 85°C ~95°C - 15% Glycerin water solution


**STANDARD ACCESSORIES:** Main Instrument, SS-Lid, Circulating Tube, Thermal Sleeve for Circulating Tube, Fuse, Power Cord, Instructional Manual, Test Certificate

★ Specifications are subject to change without notice.


## LABMAN Product Range 2025-26

|                                                                       |                                                           |                                                      |                                                    |                                                        |
|-----------------------------------------------------------------------|-----------------------------------------------------------|------------------------------------------------------|----------------------------------------------------|--------------------------------------------------------|
| <p>LMFD61, 62, 63, 64</p> <p>Lyophilizer Freeze Dryer (-60°C)</p>     | <p>LMFD83, 84</p> <p>Lyophilizer Freeze Dryer (-80°C)</p> | <p>LMRE25</p> <p>Rotary Evaporator</p>               | <p>ADM3</p> <p>Automatic Density Meter</p>         | <p>MSOD6</p> <p>Microvolume Spectrophotometer</p>      |
| <p>NAE32</p> <p>Nucleic Acid Extractor</p>                            | <p>EMR4</p> <p>Elisa Microplate Reader</p>                | <p>GPTC96</p> <p>Gradient PCR Thermal Cycler</p>     | <p>HO30</p> <p>Hybridization Oven</p>              | <p>UVCL254</p> <p>UV - Crosslinker</p>                 |
| <p>TGH12, 24</p> <p>Tissue Grinder Homogenizer</p>                    | <p>PBH80T</p> <p>Paddle Blender Homogenizer</p>           | <p>TOC2NDir</p> <p>Total Organic Carbon Analyzer</p> | <p>TOC1CD</p> <p>Total Organic Carbon Analyzer</p> | <p>LMDS600, 1200</p> <p>Microwave Digestion System</p> |
| <p>HPLC1IS, HPLC2GS</p> <p>High Performance Liquid Chromatography</p> | <p>APT1800</p> <p>Automatic Potential Titrator</p>        | <p>KFT250VC</p> <p>Karl Fischer Titrator</p>         | <p>COD1500</p> <p>COD Analyzer with Digestor</p>   | <p>BOD4000</p> <p>BOD Analyzer</p>                     |
| <p>MTVM50</p> <p>Multi Tube Vortex Mixer</p>                          | <p>LMVM60D</p> <p>Vortex Mixer</p>                        | <p>OS20</p> <p>Orbital Shaker</p>                    | <p>TRM6</p> <p>Tube Roller Mixer</p>               | <p>DBNE12, 24</p> <p>Dry Bath Nitrogen Evaporator</p>  |





## LABMAN Product Range 2025-26

|                                                                                                                                               |                                                                                                                                               |                                                                                                                                           |                                                                                                                                                                                    |                                                                                                                                                               |
|-----------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>DBI11T</b><br><br><b>Dry Bath Incubator</b>               | <b>DBI15SH, DBI16SC</b><br><br><b>Dry Bath Incubator</b>     | <b>IRS14, 35</b><br><br><b>Infrared Sterilizer</b>       | <b>GGA200</b><br><br><b>Air Generator</b>                                                       | <b>GGA301, 302</b><br><br><b>Hydrogen Generator</b>                        |
| <b>GGN303</b><br><br><b>Nitrogen Generator</b>               | <b>GHA304</b><br><br><b>Hydrogen Air Generator(2 in 1)</b>   | <b>LMVP10</b><br><br><b>Vacuum Pump</b>                  | <b>LMVP20</b><br><br><b>Vacuum Pump</b>                                                         | <b>LMCO1</b><br><br><b>Column Oven</b>                                     |
| <b>LMIO6</b><br><br><b>Mini Thermostatic Incubator</b>      | <b>GHP400</b><br><br><b>Graphite Hotplate</b>               | <b>LOBMS5</b><br><br><b>Oil Bath Magnetic Stirrer</b>   | <b>LMMS4LC, 20LC</b><br><br><b>Ceramic Hotplate Magnetic Stirrer</b>                            | <b>LMMS300</b><br><br><b>Hotplate Magnetic Stirrer</b>                    |
| <b>MMS15S</b><br><br><b>Multiposition Magnetic Stirrer</b> | <b>MMS12H</b><br><br><b>Multiposition Magnetic Stirrer</b> | <b>WBMS3</b><br><br><b>Water Bath Magnetic Stirrer</b> | <b>LMWB2H, 4H, 6H, 8H</b><br><br><b>Digital Water Bath</b>                                     | <b>HOS20L</b><br><br><b>Overhead Stirrer</b>                             |
| <b>LMHH160</b><br><br><b>Hand Held Homogenizer</b>         | <b>PS1565</b><br><br><b>Peltier / Sipper System</b>        | <b>LUV6HG</b><br><br><b>PMT Spectrophotometer</b>      | <b>LUV2000T, 2100T</b><br><br><b>Single / Double Beam UV-Vis. Spectrophotometer with 21CFR</b> | <b>LMSPUV SERIES</b><br><br><b>Double Beam UV-Vis. Spectrophotometer</b> |

## LABMAN Product Range 2025-26

|                                                                                                                                                         |                                                                                                                                                    |                                                                                                                                                                |                                                                                                                                                        |                                                                                                                                                                |
|---------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>LMSPUV1200</b></p>  <p>UV-Vis. Single Beam Spectrophotometer</p> | <p><b>LMSPV320, 325</b></p>  <p>Visible Spectrophotometer</p>     | <p><b>ADP100CFR, 95CFR</b></p>  <p>Automatic Digital Polarimeter with CFR</p> | <p><b>ADP45, 50</b></p>  <p>Automatic Digital Polarimeter (±45°)</p> | <p><b>ADP75, 90</b></p>  <p>Automatic Digital Polarimeter (±90°)</p>        |
| <p><b>ADP40</b></p>  <p>Automatic Digital Polarimeter</p>              | <p><b>LMPM180</b></p>  <p>Manual Polarimeter</p>                  | <p><b>MAR100</b></p>  <p>Manual ABBE Refractometer</p>                        | <p><b>RFM100CFR</b></p>  <p>Digital Refractometer CFR</p>            | <p><b>RFM950, 970</b></p>  <p>Digital Refractometer</p>                     |
| <p><b>LMAR1317T</b></p>  <p>ABBE Refractometer</p>                    | <p><b>LMDV300H</b></p>  <p>Digital Viscometer with Thermosel</p> | <p><b>LMDV60, 100, 200</b></p>  <p>Digital Viscometer ( L Series )</p>       | <p><b>LMDV400R</b></p>  <p>Digital Viscometer ( R Series )</p>     | <p><b>LMKV2</b></p>  <p>Krebs Viscometer</p>                               |
| <p><b>LMDV111, 222</b></p>  <p>Digital Viscometer with SSA</p>       | <p><b>LVA30</b></p>  <p>Low Viscosity Adaptor</p>               | <p><b>LMCR80F</b></p>  <p>Visual Colorimeter</p>                            | <p><b>LMUC2D</b></p>  <p>Ultrasonic Cleaner</p>                    | <p><b>LMUC3D, 6D, 9D, 12D, 25D</b></p>  <p>Ultrasonic Cleaner - Degas</p> |
| <p><b>LMUC50, 100</b></p>  <p>Industrial Ultrasonic Cleaner</p>      | <p><b>LMPH15</b></p>  <p>pH Meter with Printer</p>              | <p><b>LMPH12, LMPH10H</b></p>  <p>pH Meter</p>                              | <p><b>LMPH9</b></p>  <p>pH Meter - Basic</p>                       | <p><b>LMCM20H</b></p>  <p>Conductivity / TDS Meter</p>                    |

## LABMAN Product Range 2025-26

|                                                                                                                                                         |                                                                                                                                                           |                                                                                                                                                   |                                                                                                                                                                      |                                                                                                                                                       |
|---------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>LMCM20H</b></p>  <p><b>Conductivity / TDS Meter</b></p>          | <p><b>LMMP30, 32</b></p>  <p><b>Multiparameter (7 in 1 / 9 in 1)</b></p> | <p><b>LMION40</b></p>  <p><b>Ion Meter</b></p>                   | <p><b>LMDO50</b></p>  <p><b>Dissolved Oxygen Meter</b></p>                         | <p><b>LMHD60</b></p>  <p><b>Water Hardness Meter</b></p>           |
| <p><b>OPH11</b></p>  <p><b>Online pH Controller</b></p>                | <p><b>OCM21</b></p>  <p><b>Online Conductivity Controller</b></p>        | <p><b>LMTB200</b></p>  <p><b>Turbidity Meter</b></p>             | <p><b>RCB1010T, 1020T, 1030T</b></p>  <p><b>Refrigerated Circulating Bath</b></p> | <p><b>RCB620T</b></p>  <p><b>Refrigerated Circulating Bath</b></p> |
| <p><b>RWC320</b></p>  <p><b>Recirculating Water Chiller</b></p>       | <p><b>LMDO4P</b></p>  <p><b>Portable DO Meter</b></p>                   | <p><b>LMPH1P</b></p>  <p><b>Portable pH / ORP Meter</b></p>     | <p><b>LMCM2P</b></p>  <p><b>Portable Conductivity / TDS Meter</b></p>             | <p><b>LMTB5P</b></p>  <p><b>Portable Turbidity Meter</b></p>      |
| <p><b>PRO656, 900, 1200, 1800</b></p>  <p><b>Probe Sonicator</b></p> | <p><b>PRO660R</b></p>  <p><b>Probe Sonicator with Remote</b></p>       | <p><b>PRO650</b></p>  <p><b>Probe Sonicator</b></p>            | <p><b>LMIF30T, 50T, 100T</b></p>  <p><b>Ice Flaker</b></p>                      | <p><b>LMIF200T</b></p>  <p><b>Ice Flaker</b></p>                 |
| <p><b>LMIM15</b></p>  <p><b>Ice Maker</b></p>                        | <p><b>LMIM20C</b></p>  <p><b>Ice Maker with Crusher</b></p>            | <p><b>ISE19, 62</b></p>  <p><b>Nitrate / Fluoride ISEs</b></p> | <p><b>ER SERIES</b></p>  <p><b>pH Electrode / DO Electrode</b></p>               | <p><b>CC1, CC2, CC3</b></p>  <p><b>Conductivity Cell</b></p>     |