



OXY CARE PLUS

5 LPM OXYGEN CONCENTRATOR

Feel the comfort of easy oxygen intake with Entrusted Oxygen Concentrators which delivers oxygen at a flow rate of 0 - 5 liters per minute. Our equipment delivers oxygen at a purity of 93 % \pm 3%. Our oxygen concentrator is designed with heavy duty compressors specially considering the voltage fluctuations in India. It comes with an added advantage that the premium quality zeolite from United States and France are being used to adsorb the nitrogen content in the atmospheric air. It comes with an innovative big LCD display and low noise which extends your comfort zone.

Anticipating your inevitable need of a nebulizer with an oxygen concentrator, we offer you built in nebulizer kit with our equipment.

Our LCD display parameters include switch times, operating pressure, present working time, accumulated working time, presetting time from 10 minutes to 40 hours.

It alerts you with an alarm during power failure, high and low pressure, low purity/low flow, high temperature, maintenance.



Specifications:

Flow rate: 0-5 LPM

Outlet pressure: 0.04-0.07 MPA

Weight: 16Kgs | Size: 348*280*510 mm

Sound level \leq 40 db Inbuilt nebulizer kit