



 **+91-73974 49991**

Work Mode	ST,S,T, CPAP
Therapy Data Display	AHI, SNI, Used Hours, etc for a period upto 1 year
Sound	<28 dB
Screen	2.4 inch Sharp LCD Display
Pressure Range	4-25 cm H2O
IPAP	4-25 cm H2O
EPAP	4-20 cm H2O
Data Storage	SD Card
Power Supply	Input 100-240V, 50-60Hz, 1-2A; Output +24V, 2.5 A
DC Voltage	24VDC, 2.5A Maximum
Cooling Down	Available
Leak Alert	Yes
Filter Change Alert	Yes
Backup Rate	5-40/min, adjustable
Mask Matching	Yes
Ramp Time	0-60 min adjustable 5min/step
Auto ON/OFF	Yes
Rise Time	100-900ms
Expiratory Sensitivity	Low, Medium, High
Inspiratory Sensitivity	Low, Medium, High
EPR	Yes, 0-3 Level
Auto - Altitude	Yes
Leak Compensation	Yes
Awake Detection	Yes
Breathe Event Detection	Yes
Water Capacity	250ml Max
Humidifier Level	Heated, 0-5 Level & detachable
Colour	Black
Dimensions(cm)	14Lx13Wx9H (main); 26Lx13Wx11H(with humidifier)
Weight	0.9Kg(Net weight); 1.5Kg (with humidifier)

- ✓ Compact, Light Weight & Ultra Portable (0.9Kg)
- ✓ User Friendly with Ergonomically Designed Rotary Dial for all inputs.
- ✓ Therapy Data - AHI, Snoring Index, Usage Hours, P90, Leak, etc., can be easily viewed from the device itself from previous night sleep upto an year.
- ✓ Includes Easily Detachable Heated Humidifier.
- ✓ Bio Compatible Easy to Clean Stainless Steel Humidifier Tank.
- ✓ Real Time Therapy Data Display.
- ✓ Includes S/T mode with Pressure Range of 4 - 25 cm h2o.
- ✓ Free SD card with Powerful Data Analysis Software.
- ✓ Free Comfo-Fit Full Face Mask with headgear.
- ✓ Easy to generate detailed printed reports for doctor's analysis.
- ✓ ISO, CE & IP 22 Certified Device.
- ✓ 2 Years Warranty