



 **+91-73974 49991**

- ✓ Compact, Light Weight & Ultra Portable (0.9Kg).
- ✓ User Friendly with Ergonomically Designed Rotary Dial for all inputs.
- ✓ Therapy Data - AHI, Snoring Index, Usage Hours, P90, Leak, etc., can be easily viewed from the device itself from previous night sleep upto an year.
- ✓ Includes Easily Detachable Heated Humidifier.
- ✓ Bio Compatible Easy to Clean Stainless Steel Humidifier Tank.
- ✓ Real Time Therapy Data Display.
- ✓ Includes Automatic & Manual CPAP modes.
- ✓ Expiratory Pressure Relief available.
- ✓ Free SD card with Powerful Data Analysis Software.
- ✓ Free Comfo-Fit Nasal Mask with headgear.
- ✓ Easy to generate detailed printed reports for doctor's analysis.
- ✓ ISO, CE & IP 22 Certified Device.
- ✓ 2 Years Warranty.

Work Mode	CPAP, Auto CPAP
Therapy Data Display	AHI, SNI, Used Hours, etc for a period upto 1 year
Sound	<28 dB
Screen	2.4 inch Sharp LCD Display
Pressure Range	4-20 cm H2O
Data Storage	SD Card
EPR	Yes, 0-3 Level
Screen Auto Sleep	Yes
Power Supply	Input 100-240V, 50-60Hz, 1-2A; Output +24V, 2.5 A
DC Voltage	24VDC, 2.5A Maximum
Cooling Down	Available
Leak Alert	Yes
Filter Change Alert	Yes
Backup Rate	5-40/min, adjustable
Mask Matching	Yes
Ramp Time	0-60 min adjustable 5min/step
Auto ON/OFF	Yes
Tubing	Heated Tubing
Mask	Comfo Fit Nasal Mask
Travel Bag	Yes
Humidifier tank	Stainless Steel
Auto - Altitude	Yes
Leak Compensation	Yes
Awake Detection	Yes
Breathe Event Detection	Yes
Water Capacity	250ml Max
Humidifier Level	Heated, 0-5 Level & detachable
Colour	White
Dimensions(cm)	14Lx13Wx9H (main); 26Lx13Wx11H(with humidifier)
Weight	0.9Kg(Net weight); 1.5Kg (with humidifier)