

Entrusted Oxy Care Digi 5 LPM Oxygen Concentrator with In Built Nebulizer & Digital Display indicating Flow Rate & Oxygen Purity Level



(ISO & CE Certified with 2 Years Warranty)

Technical Specifications of Entrusted Oxy Care Digi 5 LPM Model with Nebulizer



- 1) Flow rate: 1-5 LPM
- 2) 93% \pm 3% Oxygen Purity in all flow rate.
- 3) Outlet pressure: 0.02-0.07 MPA
- 4) Net Weight: 15.5 Kg
- 5) Gross Weight 18.5 Kg
- 6) Sound level \leq 54 db
- 7) Inbuilt nebulizer port
- 8) Includes one nasal cannula, one nebulizer, one power cord, user manual with one set spare filters
- 9) Has timer function in increment of 10 minutes upto 990 minutes (16.5 Hours)
- 10) Digital Display indicating accurate flow rate & oxygen purity level along with Hours & Minutes run.



- 11) Any error codes generated will display on the device screen for easy remedy & quick analysis.
- 12) Has wheels for easy maneuverability
- 13) Dimensions: 240x345x555 mm
- 14) Voltage requirement 220 V \pm 10 V.
- 15) Power Consumption = 400 Watts.
- 16) ISO & CE Certified
- 17) 2 Years Warranty

For more details call us at +91-7397 44 999 1 or login to www.entrusted.in