



Entrusted Oxy Care Digi 5 LPM Oxygen Concentrator with In Built Nebulizer & Digital Display indicating Flow Rate & Oxygen Purity Level (ISO & CE Certified with 2 Years Warranty)



low rate 5 LPM

Oxygen Purity in all flow rate

Outlet pressure MP

Net weight 5.5 kg

gross weight 8.5 kg

Sound level 55 db

inbuilt nebulizer port

includes one nasal cannula one nebulizer one power cord user manual with one set spare filters

timer function in increments of 5 minutes upto 5 hours

Digital Display indicating accurate flow rate oxygen purity level along with hours Minutes run

any error codes generated will display on the device screen for easy ready quick analysis

casters wheels for easy maneuverability

Dimensions 450 x 550 x 555

voltage requirement

Power Consumption 100 watts

ISO CE Certified

2 years warranty