




ENTRUSTED[®]
... Caring is our Calling

OXY CARE DIGI

5 LPM OXYGEN CONCENTRATOR

Oxy Care Digi comes equipped with an Inbuilt Oxygen Sensor & Flow Sensor to indicate live display of Output Oxygen Purity & Flow Rate set on its Crystal Clear Large LCD Display Screen. Easy & Noiseless Operation with Wheels for easy maneuverability. Has a separate port for Nebulization & electronic fuse to prevent it from electrical surge. Although promoted as a 5 LPM Oxygen Concentrator, same equipment can be used with reduced oxygen purity at 6 & 7 LPM during emergency.

- ✓ Flow Rate: 1-5 LPM
- ✓ 93% +/- 3% Oxygen Purity up to 5 LPM
- ✓ Net Weight : 15.5 Kgs
- ✓ Sound Level: <54 db
- ✓ Inbuilt Nebulizer Port
- ✓ Digital Display to Indicate Oxygen Flow Rate & Oxygen Purity on a sharp large LCD Display
- ✓ Dimensions: 240x345x555 mm
- ✓ Power Consumption: 400 Watts
- ✓ ISO & CE Certified



+91-73974 49991