

OXY CARE PRIME

10 LPM OXYGEN CONCENTRATOR



Heavy Duty Oxygen Concentrator with Flow Rate capacity of 10 LPM. Equipped with Durable Compressors & Imported Sieve Membrane to provide Oxygen at Highest Purity of 93% +/- 3%. Inbuilt Nebulizer. Easily portable with wheels. Innovative storage compartment to keep Nasal Cannula's & all disposables. Quiet Working & Preset Timer functions. Power Failure, Low Pressure / High Pressure, Low Oxygen Alarms for Easy Error Diagnosis. Display indicates Accumulation, Switch, Operation Time & Output Pressure.

- ✓ Flow Rate: 0-10 LPM
- ✓ 93% +/- 3% Oxygen Purity
- ✓ Net Weight : 28 Kgs
- ✓ Sound Level: <50 db
- ✓ Inbuilt Nebulizer Port & Storage Compartment
- ✓ Display Indicates Switch Time, Accumulation Time, Working Time, Outlet Pressure.
- ✓ Dimensions: 365x375x600 mm
- ✓ Power Consumption: 550 Watts
- ✓ ISO & CE Certified

Chennai +91-73974 49991 Bangalore: +91-97907 07896 Hyderabad +91-75502 79248 Pondycherry: +91-73388 84414

Email: entrustedmed@gmail.com Web: www.entrusted.in