



## Shoe Guide

### Acceptable Shoes

- Must be supportive and sturdy
- Must have a back or strap for support
- Must have sufficient material to attach hardware to
- Should have a low or sturdy heel

If you are unsure about the appropriate shoes for the MyQuest Soft AFO brace, contact us at [info@myquestorthosis.com](mailto:info@myquestorthosis.com).



## Shoe Guide

### **NOT Acceptable Shoes**

Must not be flimsy or have small straps

Must not be open at the back

Must not attach hardware to mesh, webbing or waterproof materials

Must not attach hardware to areas that have rhinestones or decorative features (buckles/velcro/etc.)

If you are unsure about the appropriate shoes for the MyQuest Soft AFO brace, contact us at [info@myquestorthosis.com](mailto:info@myquestorthosis.com).

