

CapN SHIELD™

**YOU SUPPLY THE CAP,
WE SUPPLY THE PROTECTION**

**100%
REUSABLE
AND
RECYCLABLE**

- INSTRUCTIONS -



SLIDE SHIELD ON TO BRIM OF YOUR CAP,
CENTER LEFT TO RIGHT ON BRIM.



2

ADJUST IF SHIELD TOUCHES YOUR NOSE:



TIP BRIM OF
CAP UP

OR



SLIGHTLY SLIDE
SHIELD ON BRIM
AWAY FROM
YOUR NOSE.

3

CLEAN INSIDE and OUT SIDE OF SHIELD AFTER
EACH USE WITH DISINFECTANT CLEANSER. ★



FOLLOW INSTRUCTIONS OF
CLEANSER FOR DISINFECTING
PLASTIC SURFACES.



WIPE SHIELD WITH
SOFT CLOTH. AVOID
ABRASIVE TOWELS
OR SCRUBBING.

WARNING: DO NOT STORE CapN SHIELD in a hot space for an extended period of time; you should be wearing it! If left in a hot car in the sun the shield may distort. It's like Kryptonite to the CapN Shield.

Reporting in the April 29, 2020 Journal of the American Medical Association, experts led by Dr. Eli Perencevich, of the university's department of internal medicine, and the Iowa City VA Health Care System... observed laboratory testing "suggests that cloth masks provide [only] some filtration of virus-sized aerosol particles." According to Perencevich's group, "face shields may provide a better option." The study noted, not only is "the use of a face shield a reminder to maintain social distancing, but it allows visibility of facial expressions and lip movements for speech perception."

★ **FDA Approved** - if cleaning procedure is followed
after each use, per EUA - Section 564(g)(2)(c)

MADE IN USA



Personality Protection

Protection for any crowded space where service should be with a smile!