

# HY101-01A

Simplicity is not simple. Compact yet able to perform an outstanding range of transfer, an economically priced HY101-01A devotes its ability in science to avoid physical strain from being placed on carers and offer greater independence for the patient. Curved design of lifting arm ensures the user free of contact with the patient lift. Detachable battery allows to have it charged anywhere. Manual emergency lowering mechanism. A simple use handset ensures the operator to control the movement simply by pushing a button. Free sling with purchase.



## TECHNICAL DATA

1200 mm  Length	620 mm  Width	860 mm  Max fork Distance
1375 mm  Lifting Distance	60 mm  Ground Clearance	24V 8000 n  Motor
60-80 Times  Battery Capacity	45 Kg  Net Weight	180 Kg  Max Loading Weight